




**WHAT IF YOUR BIGGEST
CHALLENGES WITH KERRISONS
COULD BE SOLVED WITH ONE PRODUCT?**



**INTRODUCING
THE RAPIDCLEAN® SYMMETRY SHARP KERRISON®**

WHAT IF...

**YOU COULD
HAVE A
PRODUCTIVE
BITE EVERY
TIME?**

**YOU NEVER
HAD TO USE A
DULL KERRISON
AGAIN?**

**YOU COULD
MINIMIZE THE
RISK OF DURAL
TEARS?**

**ALL OF THIS WAS
POSSIBLE WITH
ONE PRODUCT?**

**YOU COULD
REDUCE HAND
AND ARM
FATIGUE?**

**YOU COULD
EASILY OPEN THE
KERRISON FOR
CLEANING?**

**YOU NO LONGER
HAD TO OPEN
MULTIPLE TRAYS?**



SHARP TIP EVERY TIME


“ I choose to use the Symmetry Sharp Kerrisons for my ACD cases because they reduce the decompression time and improve the efficiency of the procedure.”

Neurosurgeon, Louisiana

“ If one has to tug and twist the Kerrison, the risk of a dural tear is much higher.”

Neurosurgeon, Virginia

Convertible handle for easier cleaning



Single-Use Tip

The single-use tips ensure a sharp tip every time, reducing the number of passes required and limiting surgeon hand fatigue.

Convenient hinge design opens easily for cleaning, inspection and lubrication while eliminating the risk of misplacing pieces during processing.

RAPIDCLEAN® SYMMETRY SHARP KERRISON®

RAPIDCLEAN® REUSABLE HANDLES

CERAMIC COATED

DESCRIPTION	CODE
5" (125 mm) Handle for 7" (178 mm) Kerrison	53-1674RC
6" (152 mm) Handle for 8" (203 mm) Kerrison	53-1675RC
7" (178 mm) Handle for 9" (229 mm) Kerrison	53-1676RC



SINGLE-USE TIPS

40° UP, 9 MM OPENING, 2" (50 MM) TIP LENGTH
PACKAGED STERILE, 3 PER BOX



DESCRIPTION	THIN FOOTPLATE	STANDARD FOOTPLATE
1 mm Tip	53-1661T	53-1661
2 mm Tip	53-1662T	53-1662
3 mm Tip	53-1663T	53-1663
4 mm Tip	53-1664T	53-1664

CALL YOUR SYMMETRY SURGICAL REPRESENTATIVE TODAY.
LEARN MORE AT SYMMETRYSURGICAL.COM OR CALL 1-800-251-3000.

