



15930 FluidGard® 120 Anti-Fog Surgical Mask with Extended Shield

Intended Use

Intended to be used in environments where low to moderate fluid exposure is a risk.



Complies with ASTM F2100

Fluid Resistance (mmHg): 120
BFE (%) @ 3.0 μm : $\geq 98\%$
PFE (%) @ 0.1 μm : $\geq 98\%$
 ΔP (mm H₂O/cm²): < 5
Flame Spread: Class 1

Key Features

- Complies with ASTM F2100 Level 2 requirements
- Coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Anti-fog foam inhibits fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Made in North America
- Not made with natural rubber latex

FluidGard® 120 Anti-Fog Surgical Mask with Extended Shield

Intended to be used in environments where low to moderate fluid exposure is a risk.



Technical Specifications

Part No.	15930
Description	FluidGard® 120 Anti-Fog Surgical Mask with Extended Shield
Color	Blue
Fluid Resistance	120 mm Hg
Filtration – BFE	≥99.0
Filtration – PFE	≥98.9
Delta P	≤3.3
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Qty	25/Box, 4/Boxes, 100/Case

Precept® and FluidGard® are a registered trademarks of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and model.

©2021 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1051547 rev 1 03-AUG-2021 ENG – US

Aspen Surgical Products, Inc.
phone +1 616.698.7100 toll-free +1 888.364.7004
aspensurgical.com

