


**Instructions for Use: REF# STP-SQUARE1**

1. Remove from packaging.
2. Use the STP-SQUARE1 as additional padding for patient bony prominences or in any situation where additional patient padding is required.
3. The pad may be cut or torn into smaller pieces if needed.
4. Discard after use.



Not made with natural rubber latex



Manufactured For:

  
**PRIME**  
MEDICAL  
1410 E Weisgarber  
Knoxville, TN 37909  
865-776-8994



QNET BV  
Kantstraat 19  
NL-5076NP Haaren  
The Netherlands