

Protek™ Probe Covers

For ultrasound scanning and biopsy

Intended Use

Intended to help minimize cross-contamination while protecting equipment, and reducing clean up time.

Key Features

- High-quality material is softer and more pliable for ultimate patient comfort
- Superior fold with easy grip cuff improves usability by creating a single point-of-entry for easier and faster application
- 3D end for seamless imaging
- Made in North America
- Cost savings


Protek™ Probe Covers
For ultrasound scanning and biopsy

Intended to help minimize cross-contamination while protecting equipment, and reducing clean up time.



Protek™ Probe Covers
Technical Specifications

Part No.	Description	Size	Sterility	Qty
1-519-2648	PU probe cover (telescopic-fold w/ 3D end & cuff) gel, 2 tapes, 2 bands	6" x 48"	Sterile	24/Box
1-519-2589	PU probe cover (telescopic-fold w/ 3D end & cuff) gel, 2 tapes, 2 bands	5" x 48"	Sterile	24/Box
1-519-2696	PU/PE probe cover (telescopic-fold w/ 3D end & cuff) gel, 2 tapes, 2 bands	6" x 96"	Sterile	24/Box
1-519-2596	PU/PE probe cover (telescopic-fold w/ 3D end & cuff) gel, 2 tapes, 2 bands	5" x 96"	Sterile	24/Box



*Arrow indicates probe placement (not hand).

Protek™ is a trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2023 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1088083 rev 2 20-Sept-2023 ENG - US

Aspen Surgical Products, Inc.
phone +1 616.698.7100 toll-free +1 888.364.7004
aspensurgical.com

