



# Product Catalog

Surgical Essentials and  
Infection Control Products

 **aspen** surgical®

We're Aspen Surgical, a supplier of surgical essentials. From the precision of the finest scalpel with our Bard-Parker® brand to our Precept® PPE and Protek™ ultrasound probe covers that help fight preventable infections and cross contamination. And our WriteSite® Plus skin markers that provide an accurate and reliable way to mark the incision site before surgery.

From one trusted supplier, you get a portfolio of essential, physician preferred medical disposables, sold through an integrated channel. With strong distributor relationships and a dedicated field service team, we ensure we stay grounded in our customer's needs every step of the way.

We help provide an easier set-up, clean up, and quicker turnover in the O.R. for surgical teams and facilities so you can focus on what matters the most – your patients.

From facilities across North America, our dedicated team designs, manufactures, and delivers high-quality products to healthcare facilities around the world.

Every day we strive to make a difference. Every day we hold ourselves to a higher standard, so whenever and wherever you need us, we'll be there.

Learn more about all that Aspen Surgical has to offer in the pages of this catalog, or online at [www.aspensurgical.com](http://www.aspensurgical.com)









Aspen Surgical Products, Inc. Corporate Headquarters, Caledonia, Michigan

The following brands are registered trademarks of Aspen Surgical Products, Inc:

Bard-Parker®, CinchMat®, ColorTrack®, Dr. Fog®, ErgoSupport®, EvacMat®, Hoffee®, Oxboro®, Precept®, Rib-Back®, Richard-Allan®, Saber®, SafeClean®, Sterion®, Strako®, Surgidyne®, SurgiSafe®, S-Vac®, WaterBug Quiet Floor Suction®, and WriteSite®.

The following brands are registered trademarks of Aspen Surgical Products, Inc by Colby Manufacturing:

AquaNot®, ViraSorb®, and WaterBoom®.

The following brands are trademarks of Aspen Surgical Products, Inc:

Colby™, Cord-Clamp™, DuoGuard™, Hermanator™, InnerVent™, Protek™, ProV-Access™, PullUp™, SafetyLock™, and SlimFit™.

Betadine® is a registered trademark of Avrio Health L.P.

Chloraprep® is a registered trademark of CareFusion 2200, Inc.

Duraprep® is a registered trademark of 3M Company.

Securline® is a registered trademark of Precision Dynamics, Corporation.

Sterrad® is a registered trademark of ASP Global Manufacturing GmbH.

Director™, Director XT™, Retrofit™, and Retrofit XT™ are trademarks of Protek Medical Products, Inc.

Trumpf Medical™ is a trademark of Trumpf Medizin Systeme GmbH & Co. KG and iLED® is a registered trademark of Trumpf Medizin Systeme GmbH & Co. KG.

TruLight™ is a trademark of Trumpf GmbH + Co. KG.

Please note, not all products are available in all regions.

©202 Aspen Surgical Products, Inc. All rights reserved.

# Table of Contents

## Blades & Scalpels... 7-13

Bard-Parker® Conventional Blades	7&8
Bard-Parker® Conventional Handles	9
Bard-Parker® Protected Blade System	10
Bard-Parker® Protected Disposable Scalpels	11
Bard-Parker® Conventional Disposable Scalpels	12

## Operating Room Essentials... 14-19

O.R. Camera Cover	14
Specialty Disposable Light Handle	14
Standard Disposable Light Handles	14
Bard-Parker® Blade Remover	15
Bard-Parker® Scalpel Holder	15
Sterion® Cautery Tip Cleaner	16
Surgical Headband	16
Precept® Cold Packs	16
Oxboro® Pin Covers	17
Richard-Allan® Snares	17
Richard-Allan® Safety Pins	17
Sterion® Clamp Covers	18
Sterion® Suture Aid Booties	18
Sterion® Vessel Loops	19

## Surgical Markers... 21-24

6" Aspen Surgical Ruler	21
Richard-Allan® Classic Regular Tip Skin Marker	21
Richard-Allan® Classic Fine Tip Skin Marker	21
Richard-Allan® Modern Regular Tip Skin Marker	21
Richard-Allan® Taper Tip Skin Marker	22
Securline® Taper Tip Skin Marker	22
WriteSite® Jr. Small Preoperative Skin Marker	22
WriteSite® Dual Skin Marker	22
WriteSite® Multi-Ink Marker	22
WriteSite® Plus Standard Surgical Skin Marker	23
WriteSite® Plus Jr. Surgical Skin Marker	23
WriteSite® Plus Multi-Ink Marker	23
Securline® Fine Tip Lab Marker	24
Richard-Allan® Regular Tip Utility Marker	24
Richard-Allan® Taper Tip Utility Marker	24

## Minimally Invasive... 25-43

Dr. Fog® Anti-Fog	26
Endoscopic Kittner	26
Protek™ Endocavity Needle Guide for use with Alpinion	27
Protek™ Endocavity Needle Guide for use with BK	27
Protek™ Endocavity Needle Guide for use with Carestream	27
Protek™ Endocavity Needle Guide for use with Chison	27
Protek™ Endocavity Needle Guide for use with ECM	27
Protek™ Endocavity Needle Guide for use with Esaote	27&28
Protek™ Endocavity Needle Guide for use with GE Healthcare	28
Protek™ Endocavity Needle Guide for use with Hitachi	28&29
Protek™ Endocavity Needle Guide for use with Mindray	29
Protek™ Endocavity Needle Guide for use with Philips	30
Protek™ Endocavity Needle Guide for use with Samsung	30&31
Protek™ Endocavity Needle Guide for use with Shimadzu	31
Protek™ Endocavity Needle Guide for use with Siemens	31&32
Protek™ Endocavity Needle Guide for use with SonoScape	32
Protek™ Endocavity Needle Guide for use with Sonosite	32
Protek™ Endocavity Needle Guide for use with Supersonic	32
Protek™ Endocavity Needle Guide for use with Terason	32&33
Protek™ Endocavity Needle Guide for use with Toshiba	33
Protek™ Endocavity Needle Guide for use with Ultrasonix	33
Protek™ Endocavity Needle Guide for use with Zonare	33
Protek™ Sterile Multi-Gauge Biopsy Needle Guide	34
Protek™ Sterile Multi-Gauge Biopsy Needle Guide Probe Cover Kit	34

Protek™ Elastic Bands . . . . .	34
Protek™ Sterile Single-Gauge Biopsy Needle Guide . . . . .	35
Protek™ Sterile Single-Gauge Biopsy Needle Guide and Probe Cover Kit . . . . .	35
Protek™ General Purpose Needle Guide Starter Kit for use with ECM . . . . .	36
Protek™ General Purpose Needle Guide Starter Kit for use with Esaote . . . . .	36
Protek™ General Purpose Needle Guide Starter Kit for use with GE Healthcare . . . . .	36&37
Protek™ General Purpose Needle Guide Starter Kit for use with Mindray . . . . .	37
Protek™ General Purpose Needle Guide Starter Kit for use with Philips . . . . .	37
Protek™ General Purpose Needle Guide Starter Kit for use with Samsung . . . . .	37, 38, 39
Protek™ General Purpose Needle Guide Starter Kit for use with Terason . . . . .	39
Protek™ General Purpose Needle Guide Starter Kit for use with Zonare . . . . .	39
Protek™ ProV-Access™ Single-Depth Needle Guides . . . . .	40
Protek™ ProV-Access™ 3-Depth Needle Guide Probe Cover Kits . . . . .	40
Protek™ Vascular Needle Guide Starter Kit for use with GE Healthcare . . . . .	41
Protek™ Vascular Needle Guide Starter Kit for use with Hitachi . . . . .	41
Protek™ Vascular Needle Guide Starter Kit for use with Mindray . . . . .	41
Protek™ Vascular Needle Guide Starter Kit for use with Sonosite . . . . .	42
Protek™ Vascular Needle Guide Starter Kit for use with Terason . . . . .	42
Protek™ Vascular Needle Guide Starter Kit for use with Toshiba (Canon) . . . . .	42
Protek™ Vascular Needle Guide Starter Kit for use with Zonare . . . . .	42
Hoffee® Blue Lasso Suture Retriever . . . . .	43
Green Genie Suture Passer . . . . .	43

## Orthopedic Patient Positioning... 44-58

Erin Beach Chair Face Mask . . . . .	45
Erin Beach Chair Face Mask Replacement Straps . . . . .	45
Universal Head Restraint . . . . .	45
Indy Beach Chair Face Mask . . . . .	46
Bethesda Beach Chair Face Mask . . . . .	46
Beach Chair Head and Chin Straps . . . . .	46
Marco Shoulder Stabilization Kit . . . . .	47
Berkette Beach Chair Limb Positioner Kit . . . . .	47
Jessie Arm Suspension Kit . . . . .	48
Norton Arm Restraint Kit . . . . .	48
Dr Don's® Arm Canoe . . . . .	49
Dr Don's® Arm Canoe-V . . . . .	49
Dr Don's® Arm Suspension Kit . . . . .	50
Dr Don's® Arm Suspension Kit – Montana . . . . .	50
Dr Don's® Arm Suspension Kit – Northern . . . . .	50
Carson 8" Perineal Post Pad . . . . .	51
Patient Care Boot Liner Pad . . . . .	51
Perineal Post Cover . . . . .	51
Lewiston Positioner Kit . . . . .	52
Lindy Traction Boot Liner . . . . .	52
Lateral Positioner Pad . . . . .	52
Racine Pelvic Positioning Pad . . . . .	53
Hermanator® Knee Positioning Pad . . . . .	53
XL Hermanator® Pink Knee Positioning Pad . . . . .	54
Strako® Knee Positioning Pad . . . . .	54
TLC Distraction Pad Set . . . . .	55
Leg Comfort Pad . . . . .	55
Smoot Ankle Distraction Strap . . . . .	56
Thigh Holder Pad . . . . .	56
Parent Post Knee Pad . . . . .	56
Prone Face Pillow . . . . .	57
Sedona Arm Positioner Set . . . . .	57
MD Traction Rope . . . . .	58
Brachytherapy Mattress . . . . .	58

## Fluid Management & Floor Safety... 59-67

SurgiSafe® Specialty Absorbent Floor Mats . . . . .	60
SurgiSafe® Standard Absorbent Floor Mats . . . . .	61
SurgiSafe® EvacMat® Absorbent Floor Mat . . . . .	62

SurgiSafe® CinchMat® Absorbent Floor Mat . . . . .	62
ErgoSupport® Anti-Fatigue Mat . . . . .	63
Entrance Mat (TacTec) . . . . .	63
ErgoSupport® Anti-Fatigue Step Stool Mat . . . . .	64
ErgoSupport® Anti-Fatigue Step Stool . . . . .	64
WaterBug Quiet Floor Suction® Device . . . . .	65
WaterBoom® Floor Suction Strip . . . . .	65
AquaNot® Floor Suction Device . . . . .	66
ViraSorb® Solidifier . . . . .	66
AspenBasin® Fluid Waste Container . . . . .	67

## Labor & Delivery... 68-72

AmniHook® Amniotic Membrane Perforator . . . . .	69
Double-Grip® Cord Clamp . . . . .	69
Cord-Clamp™ Clipper . . . . .	70
U-Bag® Urine Specimen Collector . . . . .	70
Ready-Rolled® FootPrinter. . . . .	71
Identification Form. . . . .	71
PlastiBell® Circumcision Device . . . . .	72

## Wound Care... 73-97

French Spring Eye ½ Circle Cutting . . . . .	75
French Spring Eye Kelly Intestinal ½ Circle Taper . . . . .	75
Kidney Blunt Straight . . . . .	75
Regular Surgeon's ¾ Circle Reverse Cutting . . . . .	76
Regular Surgeon's ½ Circle Reverse Cutting . . . . .	77
Adson Dura ½ Circle Taper Point . . . . .	78
Ferguson (Regular Surgeon's Eye) ½ Circle Taper Point . . . . .	78
Butler's Cervix ½ Circle Taper Point . . . . .	79
Eye Needle ½ Circle Taper Point . . . . .	79
Ferguson (Regular Eye) ½ Circle Taper Point . . . . .	79
Ferguson (Regular Eye) ½ Circle Taper Point . . . . .	80
Fistula Light ½ Circle Reverse Cutting . . . . .	80
Fistula Heavy ½ Circle Reverse Cutting . . . . .	81
French Spring Eye Intestinal ½ Circle Taper Point . . . . .	81
Kieth's Abdominal Straight Cutting . . . . .	82
Kelly's Cervix ½ Circle Reverse Cutting . . . . .	83
Kelly's Intestinal ½ Circle Taper Point . . . . .	83
King's ¾ Circle Reverse Cutting . . . . .	83
Mayo Intestinal Straight Taper Point . . . . .	84
Fine Intestinal Straight Taper Point . . . . .	84
Mayo Intestinal ½ Circle Taper Point . . . . .	85
Lane's Cleft Palate ¾ Circle Reverse Cutting. . . . .	85
Lane's Cleft Palate ½ Circle Reverse Cutting . . . . .	85
Lane's ½ Circle Reverse Cutting French Spring Eye. . . . .	86
Lane's ½ Circle Reverse Cutting. . . . .	86
Ligating ½ Circle Taper Point . . . . .	86
Martin's Uterine ½ Circle Reverse Cutting. . . . .	87
Mayo Catgut ½ Circle Taper Point. . . . .	87
Mayo Catgut ¾ Circle Taper Point (Heavy). . . . .	88
Mayo Catgut ¾ Circle Trocar Point. . . . .	88
Mayo Catgut ½ Circle Trocar Point . . . . .	88
Murphy's Trachelorrhaphy ½ Circle (Light Trocar) . . . . .	89
Murphy's Hemostatic ½ Circle Taper Point . . . . .	89
Murphy's Intestinal ½ Circle Taper Point . . . . .	89
Davis Tonsil . . . . .	89
Bunnels Straight . . . . .	90
½ Circle Improved Mayo Catgut Taper Point (Light). . . . .	90
Gallies Frascia ½ Circle Taper . . . . .	90
Mayo Catgut ½ Circle Trocar . . . . .	91
Fistula ½ Circle Reverse Cutting . . . . .	91
King Cutting ¾ Circle . . . . .	91
Ferguson's ½ Circle Taper Point . . . . .	92



Ferguson ½ Circle Taper . . . . .	92
Ferguson ½ Circle Taper . . . . .	92
S-Vac® Fluted Channel Drains . . . . .	93
Perforated Drains. . . . .	94&95
Surgidyne® S-Vac® Silicone Bulb Evacuator . . . . .	95
Sterion® Filament Reinforced Skin Closures. . . . .	96
Sterion® Flexible Skin Closures. . . . .	97

## Personal Protective Equipment... 98-140

Precept® Isolation Mask . . . . .	100
Precept® KittyKare™ Procedure Mask. . . . .	100
Precept® Children's Mask . . . . .	100
Precept® Lite Pouch Surgical Mask . . . . .	101
Precept® Lite & Cool Horizontal Tie Surgical Mask. . . . .	101
Precept® Fluid Resistant Procedure Mask . . . . .	102
Precept® Fluid Resistant Procedure Mask with Shield . . . . .	102
Precept® Fluid Resistant Procedure Face Masks with Anti-Glare Shield . . . . .	103
Precept® Lite & Cool Surgical Mask. . . . .	103
Precept® Standard Surgical Mask . . . . .	103
Precept® Comfort-Cool™ Surgical Mask. . . . .	104
Precept® Comfort-Plus™ Surgical Mask. . . . .	104
Precept® Lite & Cool Surgical Mask with Anti-Glare Shield. . . . .	104
Precept® Foam Anti-Fog Surgical Mask . . . . .	105
Precept® Film Anti-Fog Surgical Mask . . . . .	105
Precept® Tape Anti-Fog Surgical Mask . . . . .	105
Precept® Sensitive Skin Surgical Mask . . . . .	106
Precept® Sensitive Skin Foam Anti-Fog Surgical Mask. . . . .	106
Precept® Sensitive Skin Plus Surgical Mask . . . . .	106
Precept® Comfort-Plus™ Surgical Mask with Stretch Knit Ties. . . . .	107
Precept® Comfort-Plus™ Surgical Mask. . . . .	107
Precept® Comfort-Plus™ Surgical Mask with Stretch Knit Ties. . . . .	107
Precept® Foam Fog-Shield® Surgical Mask . . . . .	108
Precept® Tape Fog-Shield® Surgical Mask. . . . .	108
Precept® Tape Fog-Shield® Surgical Mask. . . . .	108
Precept® FluidGard® 120 Procedure Mask. . . . .	109
Precept® FluidGard® 120 Anti-Fog Procedure Mask with Extended Shield. . . . .	109
Precept® FluidGard® 120 Horizontal Tie Surgical Mask . . . . .	109
Precept® FluidGard® 120 FluidPouch™ Surgical Mask . . . . .	110
Precept® FluidGard® 120 Anti-Fog Surgical Mask. . . . .	110
Precept® FluidGard® 120 Horizontal Tie Surgical Mask with Anti-Glare Shield . . . . .	111
Precept® FluidGard® 120 FluidPouch™ Surgical Mask with Anti-Glare Shield. . . . .	111
Precept® FluidGard® 120 Anti-Fog Surgical Mask with Extended Shield . . . . .	112
Precept® FluidGard® 160 Procedure Mask . . . . .	112
Precept® FluidGard® 160 Anti-Fog Procedure Mask . . . . .	113
Precept® FluidGard® 160 Anti-Fog Procedure Mask with Extended Shield. . . . .	113
Precept® FluidGard® 160 Anti-Fog Procedure Mask with Anti-Glare Shield . . . . .	114
Precept® DualGard™ Foam Shield™ Anti-Fog Procedure Mask . . . . .	114
Precept® FluidGard® 160 Anti-Fog Surgical Mask. . . . .	115
Precept® DualGard™ High Filtration Surgical Mask. . . . .	115
Precept® FluidGard® 160 Anti-Fog Surgical Mask with Extended Shield . . . . .	116
Precept® FluidGard® 160 Anti-Fog Surgical Mask with Anti-Glare Shield. . . . .	116
Precept® DualGard™ Surgical Mask with Extended Shield . . . . .	117
Precept® DualGard™ Foam Shield™ Surgical Mask with Stretch Knit Ties . . . . .	117
Precept® High Filtration Surgical Mask with Stretch Knit Ties. . . . .	118
Precept® Face Shields. . . . .	119
Precept® Eye Shields . . . . .	119
Precept® Isolation Gown Full Coverage . . . . .	121
Precept® Isolation Gown Over-the-Head Thumb Loop Full Coverage . . . . .	122
Precept® Isolation Gown Full Coverage . . . . .	123
Precept® Isolation Gown Over-the-Head Thumb Loop Full Coverage. . . . .	124
Precept® Isolation Gown Full Coverage . . . . .	125
Precept® Protective Gown Over-the-Head Thumb Loop with Open Back . . . . .	126

Precept® Protective Gown Over-the-Head Thumb Loop with Open Back . . . . .	127
Precept® Blue Film Gown Over-the-Head Thumb Loop with Open Back. . . . .	128
Precept® Blue Film Gown Over-the-Head Thumb Loop with Open Back. . . . .	129
Precept® Chemotherapy Gown, Full Coverage . . . . .	130
Precept® Tyvek® Gown, Full Coverage . . . . .	131
Precept® Bouffant Cap . . . . .	132
Precept® Skull Cap. . . . .	132
Precept® Surgeon Cap . . . . .	132
Precept® Surgical Hood . . . . .	133
Precept® Legionnaire Surgical Hood. . . . .	133
Precept® Technique Cap . . . . .	133
Precept® Tyvek® Hood. . . . .	134
Precept® Beard Cover . . . . .	134
Precept® Wave™ Shoe Covers . . . . .	135
Precept® Tyvek® Shoe and Boot Covers . . . . .	135
Precept® Lab Coat and Jacket . . . . .	136
Precept® Coveralls, SMS and Tyvek® . . . . .	137
Precept® Patient Gown, Full Coverage . . . . .	138
Precept® Apron, Blue Film Over-the-Head Open Back. . . . .	138
Precept® Tyvek® Protective Sleeve . . . . .	139
Precept® Newborn Accessories . . . . .	139
Precept® Spill Pad and Entrance Mat . . . . .	140


## Equipment Covers... 141-154

Protek™ Sterile Soft PullUp™ General Purpose Probe Covers . . . . .	142
Protek™ Sterile Clear PullUp™ General Purpose Probe Covers. . . . .	142
Protek™ Sterile Soft Telescopically-Folded General Purpose Probe Covers . . . . .	143
Protek™ Sterile Clear Telescopically-Folded General Purpose Probe Covers . . . . .	143
Protek™ Sterile Clear Telescopically-Folded General Purpose Probe Covers . . . . .	144
Protek™ Sterile Soft Flat-Folded General Purpose Probe Covers. . . . .	144
Protek™ Sterile Clear Flat-Folded General Purpose Probe Covers. . . . .	145
Protek™ Sterile Wireless General Purpose Probe Cover Kits. . . . .	145
Protek™ Sterile Cuffed General Purpose Probe Cover Kits . . . . .	145
Protek™ SlimFit™ Endocavity Probe Cover Kits . . . . .	146
Protek™ Tapered Endocavity Probe Cover Kits. . . . .	146
Protek™ 3.5cm Rolled Latex Endocavity Probe Cover Kits . . . . .	147
Protek™ 2.0cm Rolled Latex Endocavity Probe Cover Kits . . . . .	147
Protek™ Sterile Accordion-Fold Surgical Probe Cover Kits. . . . .	148
Protek™ Sterile Telescopically-Folded Surgical Probe Cover Kits w/Arrow . . . . .	148
Protek™ Sterile Telescopically-Folded Tapered Surgical Probe Cover Kits. . . . .	149
Protek™ Sterile Telescopically-Folded Surgical Probe Cover Kits . . . . .	149
Protek™ Sterile Cuffed Surgical Probe Cover Kits . . . . .	149
Protek™ Sterile PullUp™ Tapered Surgical Probe Cover Kits . . . . .	150
Protek™ Sterile PullUp™ Surgical Probe Cover Kits . . . . .	150
Protek™ Sterile PullUp™ Surgical Probe Cover Kits w/ Form-Fit Tip . . . . .	151
Protek™ Transesophageal Probe Cover Kits . . . . .	151
Protek™ Sterile PullUp™ Instrument Cover for Vacuum Assisted Biopsy Systems. . . . .	152
Protek™ Sterile System Covers. . . . .	152
Protek™ Sterile ClingTouch™ Monitor Cover . . . . .	153
Protek™ Sterile Tablet Cover . . . . .	153
Protek™ Sterile Smartphone Cover. . . . .	154
Protek™ Probe Storage Covers. . . . .	154
Sanitary MRI Headset Covers. . . . .	154

## Central Sterile Essentials... 155-172

Osteotome Tip Protectors. . . . .	156
Standard Instrument Tip Protectors . . . . .	157&158
DuoGuard™ Instrument Tip Protectors . . . . .	159
Low Temperature Web Guards . . . . .	159
InnerVent™ Instrument Tip Protectors . . . . .	160
Gelpi Retractor Tip Protectors . . . . .	160
Weitlaner Instrument Tip Protectors. . . . .	161
Laparoscopic MIS Tip Protectors . . . . .	161

Endoscope Guards . . . . .	162
Round Tip Protectors . . . . .	162
ColorTrack® Instrument Identification Sheet Tape . . . . .	163&164
Roll Tape . . . . .	165
Sterion® Same Day Tape Remover . . . . .	165
Oxboro® Corner Protectors . . . . .	166
Sterion® Color-Coded Identification Locking Tags . . . . .	166
General Purpose Brushes . . . . .	167
Brush Cleaning Kits . . . . .	168
Suction Tube Cleaning Brushes . . . . .	169
Fan-Tip Brush . . . . .	169
Tracheostomy Brushes . . . . .	170
Toothbrush Style . . . . .	170
Double-Ended Nylon Cleaning Brushes . . . . .	170
Reusable Channel Cleaning Brushes . . . . .	171
Disposable Scope Channel Cleaning Brushes . . . . .	171
Hand/Nail Brushes . . . . .	171
Instrument Brush Storage . . . . .	172
Protek™ Needle Guide Cleaning Brushes . . . . .	172
Protek™ Disinfection Cup . . . . .	172






**Bard-Parker® Blades and Scalpels**

Physicians and hospitals have relied on the Bard-Parker® brand for more than 100 years. They choose Bard-Parker blades and scalpels for their exceptional precision and quality. Today, Aspen Surgical produces Bard-Parker surgical blades, handles, and disposable scalpels in a broad range of conventional and protected styles.

All Bard-Parker blades and scalpels are made from a proprietary carbon steel composition and a manufacturing process that result in an exceptionally strong and precise surgical instrument. To ensure high product quality and prompt delivery, we proudly produce Bard-Parker blades and scalpels in the United States, where our robust manufacturing capacity ensures a consistent, reliable product supply. You have good reason to depend on Bard-Parker.



**Bard-Parker® Conventional Blades**

**Product Details**

- Bard-Parker Conventional Blades are the blades that physicians rely on day in and day out.
- High-quality blades for exceptional control during incisions
  - Bard-Parker’s exclusive Carbon Steel Rib-Back® Blades are designed with added rigidity for better strength and control
  - Reusable handles available for every blade size
  - Foil pouch packaging with a unique dual inner liner for maximum protection of both the blade and the end user
  - Not made with natural rubber latex and intended for single-use

**Carbon Steel Rib-Back® Blades**

Individually packaged



Part No.	Size	Sterility	Qty
371110	10	Sterile	50/Box, 150/Case
371111	11	Sterile	50/Box, 150/Case
371112	12	Sterile	50/Box, 150/Case
371115	15	Sterile	50/Box, 150/Case
371120	20	Sterile	50/Box, 150/Case
371121	21	Sterile	50/Box, 150/Case
371122	22	Sterile	50/Box, 150/Case
371123	23	Sterile	50/Box, 150/Case



## Carbon Steel Rib-Back® Blades

Individually packaged (25 Strips of 6)

Part No.	Size	Sterility	Qty
371310	10	Non-sterile	150/Case
371311	11	Non-sterile	150/Case
371312	12	Non-sterile	50/Case
371315	15	Non-sterile	150/Case
371320	20	Non-sterile	150/Case
371321	21	Non-sterile	150/Case
371322	22	Non-sterile	150/Case
371323	23	Non-sterile	150/Case
371324	24	Non-sterile	150/Case
371325	25	Non-sterile	150/Case
371340*	60	Non-sterile	150/Case

\*Autopsy, #8 Handle only



## Stainless Steel Blades

Individually packaged

Part No.	Size	Sterility	Qty
371210	10	Sterile	50/Box, 150/Case
371211	11	Sterile	50/Box, 150/Case
371212	12	Sterile	50/Box, 150/Case
371215	15	Sterile	50/Box, 150/Case
371220	20	Sterile	50/Box, 150/Case
371221	21	Sterile	50/Box, 150/Case
371222	22	Sterile	50/Box, 150/Case
371223	23	Sterile	50/Box, 150/Case



## Stainless Steel, Specialty Surgeon's Blades

Individually packaged


Part No.	Size	Sterility	Qty
371700	10A	Sterile	50/Box, 150/Case
371712	12B	Sterile	50/Box, 150/Case
371716	15C	Sterile	50/Box, 150/Case



## SafetyLock™ Carbon Steel Rib-Back® Blades

Individually packaged

Part No.	Size	Sterility	Qty
371150	10	Sterile	50/Box, 150/Case
371151	11	Sterile	50/Box, 150/Case
371152	12	Sterile	50/Box, 150/Case
371153	15	Sterile	50/Box, 150/Case
371154	20	Sterile	50/Box, 150/Case
371155	21	Sterile	50/Box, 150/Case
371156	22	Sterile	50/Box, 150/Case
371157	23	Sterile	50/Box, 150/Case



**Bard-Parker® Conventional Handles**

**Product Details**

From our slender size 7 to our rugged size 4, there’s a Bard-Parker® handle suited to most applications and medical professionals.


- Reusable handles available for any blade size
- Not made with natural rubber latex



**Conventional Handles**

For use with Bard-Parker® Blades

Part No.	Size	Blade Size	Sterility	Qty
371030	3	10 – 15	Non-sterile	5/Box
371031	3L	10 – 15	Non-sterile	5/Box
371040	4	20 – 25	Non-sterile	5/Box
371050	5	10 – 15	Non-sterile	5/Box
371060	6	20 – 25	Non-sterile	5/Box
371070	7	10 – 15	Non-sterile	5/Box
371080	8	60	Non-sterile	5/Box
371090	9	10 – 15	Non-sterile	5/Box



## Bard-Parker® Protected Blade System

### Product Details

Sharps injuries can occur at any time: prior to, during, or after a procedure. The Bard-Parker Protected Blade System helps protect your staff before, during and after using the use of the scalpel. We designed the system with a similar weight and feel to our conventional products. It even uses the same stainless steel blades, but adds a locking retractable shield that features simple, one-handed activation.

- Similar weight and feel to our conventional scalpel
- Same stainless steel blade as our conventional scalpel
- Locking retractable shield with easy one-handed activation
- Long handle available with select blade sizes
- Can be used by a right or left handed surgeon
- Not made with natural rubber latex and intended for single-use

## Bard-Parker® Protected Stainless Steel Blades

Individually packaged

For use with Bard-Parker® Protected Handles

Part No.	Size	Sterility	Qty
373910	10	Sterile	50/Box, 150/Case
373911	11	Sterile	50/Box, 150/Case
373915	15	Sterile	50/Box, 150/Case
374305	15C	Sterile	50/Box, 150/Case
373920	20	Sterile	50/Box, 150/Case
373921	21	Sterile	50/Box, 150/Case
373922	22	Sterile	50/Box, 150/Case
373923	23	Sterile	50/Box, 150/Case

## Bard-Parker® Protected Handles

For use with Bard-Parker® Protected Surgical Blades

Part No.	Size	Blade Size	Sterility	Qty
374030	3	10, 11, 15, 15C	Non-sterile	5/Box
374031	3L	10, 11, 15, 15C	Non-sterile	5/Box
374040	4	20, 21, 22, 23	Non-sterile	5/Box
374041	4L	20, 21, 22, 23	Non-sterile	5/Box





## Bard-Parker® Protected Disposable Scalpels

### Product Details

Like our Protected Blade System, Bard-Parker® Protected Disposable Scalpels provide protection from sharps injuries, but are designed for single-use

- Similar weight and feel to our conventional disposable scalpel
- Same stainless steel blade as our conventional disposable scalpel
- Locking retractable shield with easy one-handed activation
- Long handle available with select blades sizes
- Can be used by a right or left handed surgeon
- Not made with natural rubber latex and intended for single-use

### Protected Disposable Scalpels, Stainless Steel Blade Individually packaged

Part No.	Size	Sterility	Qty
372610	10	Sterile	10/Box, 100/Case
372611	11	Sterile	10/Box, 100/Case
372612	12	Sterile	10/Box, 100/Case
372613	12B	Sterile	10/Box, 100/Case
372615	15	Sterile	10/Box, 100/Case
372620	20	Sterile	10/Box, 100/Case
372621	21	Sterile	10/Box, 100/Case
372622	22	Sterile	10/Box, 100/Case
372623	23	Sterile	10/Box, 100/Case
372710	Size 10, Long Handle	Sterile	10/Box, 100/Case
372711	Size 11, Long Handle	Sterile	10/Box, 100/Case
372715	Size 15, Long Handle	Sterile	10/Box, 100/Case

### Protected Disposable Scalpels, Stainless Steel Blade Individually packaged

Part No.	Size	Sterility	Qty
373810	10	Non-sterile	10 Strip, 50/Case
373811	11	Non-sterile	10 Strip, 50/Case
373815	15	Non-sterile	10 Strip, 50/Case
373820	20	Non-sterile	10 Strip, 50/Case
373822	22	Non-sterile	10 Strip, 50/Case



## Bard-Parker® Conventional Disposable Scalpels

Disposable Scalpels deliver the quality and precision of Bard-Parker® in a cost-effective scalpel designed for single-use.

- High quality stainless steel blade
- Available in most common blade sizes
- Not made with natural rubber latex and intended for single-use

## Disposable Scalpels, Stainless Steel Blade Individually packaged

Part No.	Size	Sterility	Qty
371610	10	Sterile	10/Box, 100/Case
371611	11	Sterile	10/Box, 100/Case
371612	12	Sterile	10/Box, 100/Case
371613	12B	Sterile	10/Box, 100/Case
371615	15	Sterile	10/Box, 100/Case
371620	20	Sterile	10/Box, 100/Case
371621	21	Sterile	10/Box, 100/Case
371623	23	Sterile	10/Box, 100/Case



## Disposable Mini Scalpel, Stainless Steel Blade Non-sterile

Part No.	Size	Sterility	Qty
371641	11	Non-sterile	100/Bag, 1000/Case





## O.R. Camera Cover

- Helps save time and reduce inconvenience associated with having to sterilize reusable camera covers and helps enable a quicker O.R. turnaround time
- The clear lens allows for complete visualization of the surgical field without distortion
- Adaptors are currently available for many common O.R. cameras
- Please contact your local Aspen Surgical representative for further adaptor information
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
LT-C02	Camera Cover	20



## Specialty Disposable Light Handle

- Designed to have a similar look and feel as the Trumpf Medical™ reusable handle
- Works with all Trumpf Medical™ ALC plus surgical lights, including the TruLight™ 5000 and iLED® 3/5
- Lens clarity is equivalent to the reusable handle, providing rapid response of the laser and sensor
- Ultra-durable lens and a snap-on feature that helps allow quick and easy installation and removal
- Light intensity can be adjusted with the touch sensors while the disposable handle is being used
- Intended for single-use and not made with natural rubber latex

Part No.	Description	Sterility	Qty/Box
LT-ALC01	Disposable Light Handle	Sterile	18



## Standard Disposable Light Handles

- Two systems offered for convenience: Rigid or Flexible Light Handle Cover
- Helps eliminate costly reprocessing and create uniformity throughout your hospital
- Adaptors are available for most of the popular lights
- Please contact your local Aspen Surgical representative for further adaptor information
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
LT-R01	Rigid Handle	24
LT-F01B	Flexible Cover (1 Each)	100
LT-F02	Flexible Cover (2 Each)	80
LT-F03	Flexible Cover (3 Each)	80





### Bard-Parker® Blade Remover

For use with Bard-Parker® Blades and Handles.

- Helps safely remove up to five blades from scalpels without the need to touch the blade
- Helps to protect users from sharps injuries with the ability to safely unload contaminated blades into a secure channel
- Clear lid allows easy counting of the blades that have been removed
- Locking tabs prevent blades from being removed from device and clearly indicate when device is full
- Adhesive backing secures device to many surfaces, including drapes
- Features an audible and tactile response to ensure that blade is safely secure in device before removing scalpel handle
- Helps to comply with OSHA Bloodborne Standard 1910.1030(d)(2)(vii)(B) and AORN Recommended Practices
- Intended for single-use and not made with natural rubber latex


Part No.	Description	Sterility	Color	Qty
ST-115	Blade Remover	Sterile	Red	24/Box
ST-115NS	Blade Remover	Non-sterile	Red	100/Box



### Bard-Parker® Scalpel Holder

- Provides a safe placement of sharps between handling
- Helps prevent accidental sharps injuries, while keeping scalpels ready
- Includes slots for three scalpels to be held in a safe position
- Includes an adhesive back to ensure a safe, secure, and convenient placement in the sterile field
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty
ST-100	Bard-Parker® Scalpel Holder	100/Box



### Sterion® Cautery Tip Cleaner

- Assists in the cleaning and removing of material from cautery tips during surgical procedures
- Consists of a foam pad with abrasive material on one side and strong adhesive on the other for attachment to both disposable and non-disposable drapes
- Cleans both mono-polar and bi-polar probes
- Radiopaque
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
032001PBX	Cautery Tip Cleaner	35
TC-100	Cautery Tip Cleaner	100
TC-100A	Cautery Tip Cleaner	25
TC-100B	Cautery Tip Cleaner	240



### Surgical Headband

- Used to wick sweat away from the clinician
- Non-sterile and not made with natural rubber latex

Part No.	Description	Qty/Box
400450	Surgical Headband	10



### Precept® Cold Packs

#### Product Details

- Premium Cold Packs available in large and small sizes
- Cold Packs available in universal size
- Wide mouth with opening aids that allow for extremely easy filling
- Clip sleeve provides a secure holding area for the clip
- Secure clip closure
- Sturdy and strong ties allow for secure attachment
- Three-layer design, soft outer layer for comfort, absorbent middle layer and an impervious inner plastic bladder which helps provide even cold therapy
- Not made with natural rubber latex

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
16100	Large	13"	8"	White	Premium Cold Pack	Rayon/Polyester	25/Box	100/Case
16101	Small	13"	6"	White	Premium Cold Pack	Rayon/Polyester	25/Box	100/Case
16105	Universal	13"	7½"	White	Cold Pack	Rayon/Polyester	25/Box	100/Case






## Oxboro® Pin Covers

- Help to protect patients, staff, and clothing from the sharp ends of Steinmann Pins and K-Wires
- Made of thermoplastic and are available in three color-coded sizes to fit all eleven sizes of Steinmann Pins and K-Wires
- Sterile, intended for single use, and not made with natural rubber latex


Part No.	Color	Fits This Size Range/Inches	Qty/Box
101001PBX	Yellow	0.028 – 0.062	1/ Pouch, 60 Pouches
101002PBX	Blue	0.078 – 0.125	2/Pouch, 30 Pouches
101003PBX	Red	0.141 – 0.188	2/Pouch, 30 Pouches



## Richard-Allan® Snares

- Formed in pre-cut, pre-bent sizes, ready for immediate use
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
2405	Nasal Snare Size 5	2/Pack, 50 Packs
2406	Nasal Snare Size 6	2/Pack, 50 Packs
2407	Tonsil Snare Size 7	2/Pack, 50 Packs
2408	Tonsil Snare Size 8	2/Pack, 50 Packs
2409	Tonsil Snare Size 9	2/Pack, 50 Packs



## Richard-Allan® Safety Pins

- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
2201	1"	2/Pack, 50 Packs
2202	1.5"	2/Pack, 50 Packs
2203	2"	2/Pack, 50 Packs



## Sterion® Clamp Covers

- Assist in the atraumatic clamping of delicate vessels during surgery
- Available in radiopaque silicone
- A variety of sizes to provide outstanding convenience and flexibility of use
- Sterile, intended for single-use and not made with natural rubber latex

Part No.	Color	Material	mm Length	Inches Length	Qty/Box
061011PBX	Red	Silicone	2.5 x 130	0.10 x 5	10 Pkgs of 2
061012PBX	White	Silicone	5.1 x 130	0.20 x 5	10 Pkgs of 2
061013PBX	Blue	Silicone	7.6 x 130	0.30 x 5	10 Pkgs of 2




## Sterion® Suture Aid Booties

- Protect, tag, and help locate sutures
- Vinyl suture booties are available in two, three, or five-pair configurations and are positioned in an adhesive-backed foam block for quick and easy usage
- Radiopaque
- Choose from a variety of sizes and colors as listed below
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Tray	Qty/Box
051003PBX	Standard Yellow	5 Pair	5 Trays
051011PBX	Standard Red	5 Pair	5 Trays
051012PBX	Standard Blue	5 Pair	5 Trays
051014PBX	Standard White	5 Pair	5 Trays
051015PBX	Standard Assorted	5 Pair	5 Trays
051004PBX	Large Pink	5 Pair	5 Trays
053003PBX	Standard Yellow	3 Pair	5 Trays
052004PBX	Large Pink	2 Pair	5 Trays
052003PBX	Standard Yellow	2 Pair	5 Trays
053030PBX	X-tra Mini Red	3 Pair	5 Pouches*

\* not in a tray



## Sterion® Vessel Loops

- Used in occluding, retracting, and identifying arteries, veins, tendons and nerves in surgical procedures
- Extruded from medical-grade silicone
- Radiopaque
- Intended for single-use and not made with natural rubber latex

## Sterion® Vessel Loops

2 Loops/Pkg; 10 Pkgs/Box/Double Pouched

Part No.	Color	Cross Section – Inches	Length – Inches	Cross Section – mm	Length – mm
011001PBX	Red Mini	0.06" x 0.03"	16"	1.5 x 0.88	406
011002PBX	Blue Mini	0.06" x 0.03"	16"	1.5 x 0.88	406
011003PBX	Yellow Mini	0.06" x 0.03"	16"	1.5 x 0.88	406
011004PBX	White Mini	0.06" x 0.03"	16"	1.5 x 0.88	406
011011PBX	Red Maxi	0.1" x .04"	16"	2.5 x 1.0	406
011012PBX	Blue Maxi	0.1" x .04"	16"	2.5 x 1.0	406
011013PBX	Yellow Maxi	0.1" x .04"	16"	2.5 x 1.0	406
011014PBX	White Maxi	0.1" x .04"	16"	2.5 x 1.0	406
011021PBX	White Supermaxi	0.2" x .05"	22"	5 x 1.3	559
011022PBX	Blue Supermaxi	0.2" x .05"	22"	5 x 1.3	559
011032PBX	Blue Micro	0.02" x 0.02"	16"	0.53 x 0.53	406
011301PBX	Red Mini	0.06" x 0.03"	30"	1.5 x 0.88	762
011302PBX	Blue Mini	0.06" x 0.03"	30"	1.5 x 0.88	762
011303PBX	Yellow Mini	0.06" x 0.03"	30"	1.5 x 0.88	762
011304PBX	White Mini	0.06" x 0.03"	30"	1.5 x 0.88	762
011311PBX	Red Maxi	0.1" x 0.04"	30"	2.5 x 1.0	762
011312PBX	Blue Maxi	0.1" x 0.04"	30"	2.5 x 1.0	762
011313PBX	Yellow Maxi	0.1" x 0.04"	30"	2.5 x 1.0	762
011314PBX	White Maxi	0.1" x 0.04"	30"	2.5 x 1.0	762

SEA  
SA  
LD

**WriteSite® Plus**  
Prep Resistant Skin Marker

## Richard-Allan®, Securline®, and WriteSite® Surgical Marking Pens

- Intended to provide an accurate and reliable way to mark the incision site prior to surgery
- Contains gentian violet ink
- A variety of barrel and tip options are available depending on preference
- Intended for single-use and not made with natural rubber latex

Select markers are available with Blank Labels, Rulers, or Medication Labels as noted.



9 Black and White Medication Labels



Blank Labels



Rulers



### 6" Aspen Surgical Ruler

Part No.	Description	Sterility	Qty
2554	6" Aspen Surgical Ruler	Sterile	25/Box



### Richard-Allan® Classic Regular Tip Skin Marker

Part No.	Description	Sterility	Qty
2650	Marker, Ruler, and Label Set	Sterile	50/Box



### Richard-Allan® Classic Fine Tip Skin Marker

Part No.	Description	Sterility	Qty
2630	Marker, Ruler, and Label Set	Sterile	50/Box



### Richard-Allan® Modern Regular Tip Skin Marker

Part No.	Description	Sterility	Qty
2750	Marker, Ruler, and Label Set	Sterile	50/Box



### Richard-Allan® Taper Tip Skin Marker

Part No.	Description	Sterility	Qty
2730	Marker, Ruler, and Label Set	Sterile	50/Box



### Securline® Taper Tip Skin Marker

Part No.	Description	Sterility	Qty
1000-00-PDG	Marker & Ruler Set	Sterile	10/Box, 10 Boxes/Case
100L-00-PDG	Marker, Ruler, & Label Set	Sterile	100/Box
1002-00-PDG	Marker & Ruler Set	Sterile	100/Box
1003-00-PDG	Marker & Ruler Set	Non-sterile	100/Box



### WriteSite® Jr. Small Preoperative Skin Marker

Part No.	Description	Sterility	Qty
2900BN	Marker Only	Non-sterile	100/Box



### WriteSite® Dual Skin Marker

Part No.	Description	Sterility	Qty
2777	Marker Only	Sterile	50/Box
2778BN	Marker Only	Non-sterile	100/Box



### WriteSite® Multi-Ink Marker

Part No.	Description	Sterility	Qty
2755	Marker Only	Sterile	12/Box
2755BN	Marker Only	Non-sterile	100/Box





## WriteSite®Plus Surgical Skin Markers

- Aspen Surgical carries a line of surgical marking pens that include an ink specifically designed for use with Chloraprep®, sold under the WriteSite® Plus brand
- Ink is also compatible with Duraprep® and Betadine® prep solutions
- Intended for single-use and not made with natural rubber latex



## WriteSite® Plus Standard Surgical Skin Marker

Part No.	Description	Sterility	Qty
2700	Marker Only	Sterile	50/Box
2701	Marker, Ruler, and Label Set	Sterile	50/Box
2700BN	Marker Only	Non-sterile	100/Box



## WriteSite® Plus Jr. Surgical Skin Marker

Part No.	Description	Sterility	Qty
2715	Marker Only	Sterile	50/Box
2716	Marker, Ruler, and Label Set	Sterile	50/Box
2715BN	Marker Only	Non-sterile	100/Jar



## WriteSite® Plus Multi-Ink Marker

**New!**

Part No.	Description	Sterility	Qty
2720	Marker, Ruler & Label Set	Sterile	50/Box
2720BN	Marker Only	Non-sterile	100/Box
2721BN	Marker Only	Non-sterile	500/Box



#### Securline® Fine Tip Lab Marker

- Used to make permanent markings on plastic, metal, and porcelain labware
- Can be used to write on most damp surfaces
- Available in four different colors
- Intended for single-use and not made with natural rubber latex

Part No.	Description	Color	Sterility	Qty/Box
1400-13-PDC	Marker Only	Blue	Non-sterile	10
1400-15-PDC	Marker Only	Green	Non-sterile	10
1400-16-PDC	Marker Only	Red	Non-sterile	10
1400-20-PDC	Marker Only	Black	Non-sterile	10



#### Richard-Allan® Utility Marker

- Intended to provide a permanent method of identifying instruments, specimens, medication labels, and other objects on and off the sterile field
- Contains a black, waterproof ink that can withstand autoclave sterilization and can be used on most surfaces
- Available in several barrel and tip options and with or without labels
- Intended for single-use and not made with natural rubber latex



#### Richard-Allan® Regular Tip Utility Marker

Part No.	Description	Sterility	Qty/Box
2850	Marker and 9 Label Set	Sterile	50
2850BN	Marker and 6 Label Set	Non-sterile	100
2851BN	Marker Only	Non-sterile	100



#### Richard-Allan® Taper Tip Utility Marker

Part No.	Description	Sterility	Qty/Box
2830BN	Marker and 6 Label Set	Non-sterile	100
2831BN	Marker Only	Non-sterile	100






#### Dr. Fog® Anti-Fog

- Offered with our easy squeeze bottle and sponge set or as a ready-to-use sponge, which is dry to touch and alcohol free
- Radiopaque
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
DF-3100	Original Solution with Sponge	20
DF-3120	Alcohol-Free Dry Treated Sponge	40



#### Endoscopic Kittner

- Our blunt dissecting instrument is 38cm long with a 5mm stainless steel shaft
- Available in both single and dual tip
- Radiopaque tip
- The brushed finish allows for safe use during laser surgery
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
KT-9101	Single Tip	20
KT-9105	Dual Tip	20

## Protek™ Endocavity Needle Guides

Endocavity needle guides ensure precise needle placement when performing a biopsy or aspiration

- Features an extra wide funnel for optimal needle entry and smooth profile to ensure ultimate patient comfort
- Please note, kits have both Sterile and Non-sterile components unless otherwise noted on product label

Below you will find single-use and reusable needle guides.



Sterile, single-use example



Non-sterile, reusable example

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-5250	Needle guide	16Ga	Alpinion E3-10	Sterile	24
1-535-5144	Needle guide	16Ga	Alpinion EC3-10, EC3-10H, EC3-10T, EC3-10X, EN3-10, EV3-10, EV3-10H, EV3-10T, EV3-10X	Sterile	24
1-535-5146	3.5cm Latex biopsy kit	16Ga	Alpinion EC3-10, EC3-10H, EC3-10T, EC3-10X, EN3-10, EV3-10, EV3-10H, EV3-10T, EV3-10X	Sterile	24
1-535-5149	12" Soft biopsy kit	16Ga	Alpinion EC3-10, EC3-10H, EC3-10T, EC3-10X, EN3-10, EV3-10, EV3-10H, EV3-10T, EV3-10X	Sterile	24
1-535-5183	48" SlimFit™ biopsy kit	16Ga	Alpinion EC3-10, EC3-10H, EC3-10T, EC3-10X, EN3-10, EV3-10, EV3-10H, EV3-10T, EV3-10X	Sterile	24
1-535-5142	Needle guide	16Ga	Alpinion VE3-10H	Sterile	24
1-535-6456	Stainless Steel Needle Guide	16Ga	Alpinion EC3-10, EN3-10, EV3-10, EC3-10T, EV3-10T, EC3-10H, EV3-10H, EC3-10X, EV3-10X	Non-sterile	1
1-535-5250	Needle guide	16Ga	BK Medical 8819, E10C4	Sterile	24
1-535-5250	Needle guide	16Ga	Carestream E10C4	Sterile	24
1-535-5585	Needle guide	16Ga	Chison D6C12L, V6-A	Sterile	24
1-535-6042	Stainless steel needle guide	16Ga	Chison D6C12L	Non-sterile	1
1-535-5595	Needle guide	16Ga	Chison D7C10L, V7-E	Sterile	24
1-535-5055	Needle guide	16Ga	ECM E610P	Sterile	24
1-535-5056	3.5cm Latex biopsy kit	16Ga	ECM E610P	Sterile	24
1-535-5058	12" Soft biopsy kit	16Ga	ECM E610P	Sterile	24
1-535-5059	24" SlimFit™ biopsy kit	16Ga	ECM E610P	Sterile	24
1-535-9000	Needle guide	16Ga	Esaote E8-5 R10	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Esaote E8-5 R10	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Esaote E8-5 R10	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Esaote E8-5 R10	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Esaote E8-5 R10	Sterile	24
1-535-5145	Needle guide	16Ga	Esaote EC123, EC1123, SE3123, SE3133, TRT12, IVT12	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	Esaote EC123, EC1123, SE3123, SE3133, TRT12, IVT12	Sterile	24

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-5151	24" Soft biopsy kit	16Ga	Esaote EC123, EC1123, SE3123, SE3133, TRT12, IVT12	Sterile	24
1-535-6390	Stainless steel needle guide	16Ga	Esaote EC123	Non-sterile	1
1-535-6360	Stainless steel needle guide	16Ga	Esaote EC1123, SE3133, E3-12	Non-sterile	1
1-535-5046	Needle guide	16Ga	GE Healthcare 618E	Sterile	25
1-535-6077	Stainless steel needle guide	16Ga	GE Healthcare 618E, TR7	Non-sterile	1
1-535-5014	24" Soft & 2.0cm Latex biopsy kit (no gel)	16Ga	GE Healthcare BE9C, BE9CS	Sterile	24
1-535-5017	24" Soft & 2.0cm Latex biopsy kit	16Ga	GE Healthcare BE9C, BE9CS	Sterile	24
1-535-6026	Stainless steel needle guide	16Ga	GE Healthcare BE9C	Non-sterile	1
1-535-6340	Stainless steel needle guide	14Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, MTZ, MZ	Non-sterile	1
1-535-6449	Stainless steel needle guide	16Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, MTZ, MZ	Non-sterile	1
1-535-5055	Needle guide	16Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, IC9B-RS, MTZ, MZ	Sterile	24
1-535-5056	3.5cm Latex biopsy kit	16Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, IC9B-RS, MTZ, MZ	Sterile	24
1-535-5058	12" Soft biopsy kit	16Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, IC9B-RS, MTZ, MZ	Sterile	24
1-535-5059	24" SlimFit™ biopsy kit	16Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, IC9B-RS, MTZ, MZ	Sterile	24
1-535-5145	Needle guide	16Ga	GE Healthcare IC5-9-RS, IC5-9W-RS	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	GE Healthcare IC5-9-RS, IC5-9W-RS	Sterile	24
1-535-5151	24" Soft biopsy kit	16Ga	GE Healthcare IC5-9-RS, IC5-9W-RS	Sterile	24
1-535-5570	Needle guide	16Ga	GE Healthcare IC9-RS	Sterile	24
1-535-6098	Stainless steel needle guide	16Ga	GE Healthcare IC9-RS	Non-sterile	1
1-535-5070	Needle guide	16Ga	GE Healthcare RIC5-9-D, RIC5-9H, RIC5-9W-D, RIC5-9W-RS, RIC6-12-D	Sterile	24
1-535-5071	3.5cm Latex biopsy kit	16Ga	GE Healthcare RIC5-9-D, RIC5-9H, RIC5-9W-D, RIC5-9W-RS, RIC6-12-D	Sterile	24
1-535-6097	Stainless steel needle guide	16Ga	GE Healthcare RIC5-9-D, RIC5-9A-RS, RIC5-9H, RIC5-9W-D, RIC5-9W-RS, RIC6-12-D	Non-sterile	1
1-535-5370	Needle guide	16Ga	Hitachi C41V1, EUP-V73W	Sterile	24
1-535-5371	24" SlimFit™ biopsy kit	16Ga	Hitachi C41V1, EUP-V73W	Sterile	24
1-535-5372	3.5cm Latex biopsy kit	16Ga	Hitachi C41V1, EUP-V73W	Sterile	24
1-535-5373	12" Soft biopsy kit	16Ga	Hitachi C41V1, EUP-V73W	Sterile	24
1-535-6370	Stainless steel needle guide	16Ga	Hitachi C41V1	Non-sterile	1
1-535-6381	Stainless steel needle guide (long)	14Ga	Hitachi EUP-V33, EUP-V33W, EUP-V53W, CC41R1 CC41R, CC41R1	Non-sterile	1
1-535-6382	Stainless steel needle guide (long)	16Ga	Hitachi EUP-V33, EUP-V33W, EUP-V53W, CC41R1 CC41R, CC41R1	Non-sterile	1
1-535-6384	Stainless steel needle guide (short)	16Ga	Hitachi EUP-V33, EUP-V33W, EUP-V53W, EUP-CC531, CC41R	Non-sterile	1



### Protek™ Endocavity Needle Guides

Endocavity needle guides ensure precise needle placement when performing a biopsy or aspiration

- Features an extra wide funnel for optimal needle entry and smooth profile to ensure ultimate patient comfort
- Please note, kits have both Sterile and Non-sterile components unless otherwise noted on product label

Below you will find single-use and reusable needle guides.



Sterile, single-use example



Non-sterile, reusable example

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-5380	Needle guide	16Ga	Hitachi EUP-V33, EUP-CC31, EUP-CC531, EUP-V12, EUP-V33W, EUP-V53W, EZU-PA5V, C41V, CC41R, CC41R1	Sterile	24
1-535-5381	12" Soft biopsy kit	16Ga	Hitachi EUP-V33, EUP-CC31, EUP-CC531, EUP-V12, EUP-V33W, EUP-V53W, EZU-PA5V, C41V, CC41R, CC41R1	Sterile	24
1-535-5382	3.5cm & 2.0cm Latex biopsy kit	16Ga	Hitachi EUP-V33, EUP-CC31, EUP-CC531, EUP-V12, EUP-V33W, EUP-V53W, EZU-PA5V, C41V, CC41R, CC41R1	Sterile	24
1-535-5250	Needle guide	16Ga	Hitachi UST-679	Sterile	24
1-535-6042	Stainless steel needle guide	16Ga	Mindray 65EC10HA, 65EC10EA		
1-535-9000	Needle guide	16Ga	Mindray EV6.5 MHZ/R10, 6CV1, 6CV1s, 65EC10EA, 65EC10EB, 65EC10HA, CB10-4, CB10-4e, CB10-4s, V10-4, V10-4B, V10-4BP, V10-4Bs, V10-4s, V11-3E, V11-3Hs, V11-3HU, V11-3WE, V11-3Ws	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Mindray EV6.5 MHZ/R10, 6CV1, 6CV1s, 65EC10EA, 65EC10EB, 65EC10HA, CB10-4, CB10-4e, CB10-4s, V10-4, V10-4B, V10-4BP, V10-4Bs, V10-4s, V11-3E, V11-3Hs, V11-3HU, V11-3WE, V11-3Ws	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Mindray EV6.5 MHZ/R10, 6CV1, 6CV1s, 65EC10EA, 65EC10EB, 65EC10HA, CB10-4, CB10-4e, CB10-4s, V10-4, V10-4B, V10-4BP, V10-4Bs, V10-4s, V11-3E, V11-3Hs, V11-3HU, V11-3WE, V11-3Ws	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Mindray EV6.5 MHZ/R10, 6CV1, 6CV1s, 65EC10EA, 65EC10EB, 65EC10HA, CB10-4, CB10-4e, CB10-4s, V10-4, V10-4B, V10-4BP, V10-4Bs, V10-4s, V11-3E, V11-3Hs, V11-3HU, V11-3WE, V11-3Ws	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Mindray EV6.5 MHZ/R10, 6CV1, 6CV1s, 65EC10EA, 65EC10EB, 65EC10HA, CB10-4, CB10-4e, CB10-4s, V10-4, V10-4B, V10-4BP, V10-4Bs, V10-4s, V11-3E, V11-3Hs, V11-3HU, V11-3WE, V11-3Ws	Sterile	24
1-535-5070	Needle guide	16Ga	Mindray DE10-3U, DE10-3WU	Sterile	24
1-535-5071	3.5cm Latex biopsy kit	16Ga	Mindray DE10-3U, DE10-3WU	Sterile	24
1-535-6097	Stainless steel needle guide	16Ga	Mindray DE10-3U, DE10-3WU	Non-sterile	1

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-5270	Needle guide	16Ga	Philips C8-4v, C9-3v, ClearVue™ C9-4v (on ClearVue 350™, 550™, 650™, and 850™ ultrasound systems), C10-3v	Sterile	24
1-535-5271	24" SlimFit™ biopsy kit	16Ga	Philips C8-4v, C9-3v, ClearVue™ C9-4v (on ClearVue 350™, 550™, 650™, and 850™ ultrasound systems), C10-3v	Sterile	24
1-535-5272	3.5cm Latex biopsy kit	16Ga	Philips C8-4v, C9-3v, ClearVue™ C9-4v (on ClearVue 350™, 550™, 650™, and 850™ ultrasound systems), C10-3v	Sterile	24
1-535-5144	Needle guide	16Ga	Philips C9-4ec	Sterile	24
1-535-5146	3.5cm Latex biopsy kit	16Ga	Philips C9-4ec	Sterile	24
1-535-5149	12" Soft biopsy kit	16Ga	Philips C9-4ec	Sterile	24
1-535-5183	48" SlimFit™ biopsy kit	16Ga	Philips C9-4ec	Sterile	24
1-535-5250	Needle guide	16Ga	Philips C9-4ec (on HD7 XE™ w/ 3.0.1 software and above, and HD9™ systems), BP10-5ec (on Affiniti 50™, Affiniti 70™, HD11™, HD115™, and HD7XE™ systems)	Sterile	24
1-535-5145	Needle guide	16Ga	Philips E6509	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	Philips E6509	Sterile	24
1-535-5151	24" SlimFit™ biopsy kit	16Ga	Philips E6509	Sterile	24
1-535-6390	Stainless steel needle guide	16Ga	Philips E6509	Non-sterile	1
1-535-5395	Needle guide	16Ga	Philips C9-4v (on Sparq® and Affiniti® Systems)	Sterile	24
1-535-5390	Needle guide	16Ga	Philips C9-5ec, C9-5ICT, C9-5	Sterile	24
1-535-6392	Stainless steel needle guide	16Ga	Philips C9-5ec, C9-5 ICT, 9-5 Endo 8	Non-sterile	1
1-535-5455	Needle guide	16Ga	Philips C10-4ec	Sterile	24
1-535-5154	Needle guide	16Ga	Philips V9-4v	Sterile	24
1-535-5154	Needle guide	16Ga	Samsung 3D4-9, 3DEC4-9, 3D4-9ES, 3DEC4-9ES	Sterile	24
1-535-5505	Needle guide	16Ga	Samsung 3D5-8EK, 3D5-9EK	Sterile	24
1-535-6190	Stainless steel needle guide	18Ga	Samsung 3D5-8EK, 3D5-9EK	Non-sterile	1
1-535-5145	Needle guide	16Ga	Samsung BCC4-9	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	Samsung BCC4-9	Sterile	24
1-535-5151	24" Soft biopsy kit	16Ga	Samsung BCC4-9	Sterile	24
1-535-5260	Needle guide	16Ga	Samsung E3-12A	Sterile	24
1-535-5261	24" SlimFit™ biopsy kit	16Ga	Samsung E3-12A	Sterile	24
1-535-5580	Needle guide (short)	16Ga	Samsung EA2-11AR, EA2-11AV	Sterile	24
1-535-5590	Needle guide (long)	16Ga	Samsung EA2-11AR, EA2-11AV	Sterile	24
1-535-6268	Stainless steel needle guide	18Ga	Samsung EA2-11AR, EA2-11AV	Non-sterile	1
1-535-5295	Needle guide	16Ga	Samsung EA2-11B	Sterile	24
1-535-5144	Needle guide	16Ga	Samsung EC4-9, EC4-9/10ED, EC4-9/10R, ER4-9, ER4-9/10ED, EV4-9, EV4-9/10ED, EVN4-9	Sterile	24
1-535-5146	3.5cm Latex biopsy kit	16Ga	Samsung EC4-9, EC4-9/10ED, EC4-9/10R, ER4-9, ER4-9/10ED, EV4-9, EV4-9/10ED, EVN4-9	Sterile	24
1-535-5149	12" Soft biopsy kit	16Ga	Samsung EC4-9, EC4-9/10ED, EC4-9/10R, ER4-9, ER4-9/10ED, EV4-9, EV4-9/10ED, EVN4-9	Sterile	24
1-535-5183	48" SlimFit™ biopsy kit	16Ga	Samsung EC4-9, EC4-9/10ED, EC4-9/10R, ER4-9, ER4-9/10ED, EV4-9, EV4-9/10ED, EVN4-9	Sterile	24

## Protek™ Endocavity Needle Guides

Endocavity needle guides ensure precise needle placement when performing a biopsy or aspiration

- Features an extra wide funnel for optimal needle entry and smooth profile to ensure ultimate patient comfort
- Please note, kits have both Sterile and Non-sterile components unless otherwise noted on product label

Below you will find single-use and reusable needle guides.



Sterile, single-use example



Non-sterile, reusable example

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-6042	Stainless steel needle guide	16Ga	Samsung EC4-9ES, EC4-9IS	Non-sterile	1
1-535-9000	Needle guide	16Ga	Samsung EC4-9ES, EC4-9IS, ER4-9ES, 10R EV/140D	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Samsung EC4-9ES, EC4-9IS, ER4-9ES, 10R EV/140D	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Samsung EC4-9ES, EC4-9IS, ER4-9ES, 10R EV/140D	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Samsung EC4-9ES, EC4-9IS, ER4-9ES, 10R EV/140D	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Samsung EC4-9ES, EC4-9IS, ER4-9ES, 10R EV/140D	Sterile	24
1-535-5275	Needle guide	16Ga	Samsung EV2-10A	Sterile	24
1-535-6267	Stainless steel needle guide	18Ga	Samsung EV2-10A	Non-sterile	1
1-535-6455	Stainless steel needle guide	16Ga	Samsung EV4-9/10ED, ER4-9/10ED, EC4-9, EVN4-9, EV4-9, ER4-9	Non-sterile	1
1-535-5250	Needle guide	16Ga	Samsung NER4-9ES, NEV4-9ES	Sterile	24
1-535-6097	Stainless steel needle guide	16Ga	Samsung V5-9	Non-sterile	1
1-535-5070	Needle guide	16Ga	Samsung V5-9, EV3-10B	Sterile	24
1-535-5071	3.5cm Latex biopsy kit	16Ga	Samsung V5-9, EV3-10B	Sterile	24
1-535-5055	Needle guide	16Ga	Samsung VR5-9	Sterile	24
1-535-5056	3.5cm Latex biopsy kit	16Ga	Samsung VR5-9	Sterile	24
1-535-5058	12" Soft biopsy kit	16Ga	Samsung VR5-9	Sterile	24
1-535-5059	24" SlimFit™ biopsy kit	16Ga	Samsung VR5-9	Sterile	24
1-535-6340	Stainless steel needle guide	14Ga	Samsung VR5-9	Non-sterile	1
1-535-6449	Stainless steel needle guide	16Ga	Samsung VR5-9	Non-sterile	1
1-535-5045	Needle guide	16Ga	Shimadzu TV11R-055U	Sterile	24
1-535-6005	Stainless steel needle guide	16Ga	Shimadzu TV11R-055U	Non-sterile	1
1-535-5145	Needle guide	16Ga	Siemens BE9-4	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	Siemens BE9-4	Sterile	24
1-535-5151	24" Soft biopsy kit	16Ga	Siemens BE9-4	Sterile	24
1-535-9000	Needle guide	16Ga	Siemens EC7, EC9-4 (on Sonoline G20™, G50S™, Acuson X500™, and X100™, systems) EV9-4 (on Acuson X500™, and X100™ systems), 6.5 EC10	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Siemens EC7, EC9-4 (on Sonoline G20™, G50S™, Acuson X500™, and X100™, systems) EV9-4 (on Acuson X500™, and X100™ systems), 6.5 EC10	Sterile	24

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-9002	12" Soft biopsy kit	16Ga	Siemens EC7, EC9-4 (on Sonoline G20™, G50S™, Acuson X500™, and X100™, systems) EV9-4 (on Acuson X500™, and X100™ systems), 6.5 EC10	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Siemens EC7, EC9-4 (on Sonoline G20™, G50S™, Acuson X500™, and X100™, systems) EV9-4 (on Acuson X500™, and X100™ systems), 6.5 EC10	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Siemens EC7, EC9-4 (on Sonoline G20™, G50S™, Acuson X500™, and X100™, systems) EV9-4 (on Acuson X500™, and X100™ systems), 6.5 EC10	Sterile	24
1-535-6042	Stainless steel needle guide	16Ga	Siemens EC7, EC9-4 (Acuson X500™ and Sonoline G20™, G50™, G60S™), EC-10C5	Non-sterile	1
1-535-6477	Stainless steel needle guide	14Ga	Siemens EC9-4 (on Acuson X150™, X300™, and Sonoline G40™)	Non-sterile	1
1-535-5480	Needle guide	16Ga	Siemens EC9-4 (on Acuson NX2™, NX3™, NX3 Elite™, P500™, X150™, X300™, and Sonoline G40™ systems), EV9-4 (on Acuson X150™, X300™, and Sonoline G40™ systems)	Sterile	24
1-535-5531	Needle guide	16Ga	Siemens EC9-4 (on Acuson S2000™ and Sonoline Antares™ systems)	Sterile	24
1-535-5535	24" SlimFit™ biopsy kit	16Ga	Siemens EC9-4 (on Acuson S2000™ and Sonoline Antares™ systems)	Sterile	24
1-535-5154	Needle guide	16Ga	Siemens EV9F3	Sterile	24
1-535-6482	Stainless steel needle guide	16Ga	Siemens P2	Non-sterile	1
1-535-5172	Needle guide	16Ga	SonoScape 6V1	Sterile	24
1-535-5178	3.5cm & 2.0cm Latex biopsy kit	16Ga	SonoScape 6V1	Sterile	24
1-535-5179	12" Soft biopsy kit	16Ga	SonoScape 6V1	Sterile	24
1-535-5380	Needle guide	16Ga	SonoScape 6V3 (2015 model and older)	Sterile	24
1-535-5381	12" Soft biopsy kit	16Ga	SonoScape 6V3 (2015 model and older)	Sterile	24
1-535-5382	3.5cm & 2.0cm Latex biopsy kit	16Ga	SonoScape 6V3 (2015 model and older)	Sterile	24
1-535-5665	Needle guide	16Ga	SonoScape 6V3 (2016 model and newer)	Sterile	24
1-535-5171	Needle guide	16Ga	SonoScape 6V4	Sterile	24
1-535-5650	Needle guide	16Ga	SonoScape EC9-5	Sterile	24
1-535-5651	24" SlimFit™ biopsy kit	16Ga	SonoScape EC9-5	Sterile	24
1-535-5652	3.5cm Latex biopsy kit	16Ga	SonoScape EC9-5	Sterile	24
1-535-6392	Stainless steel needle guide	16Ga	Sonosite C8 Series	Non-sterile	1
1-535-5575	Needle guide	16Ga	Sonosite ICT, ICTe, ICTx, ICTxp	Sterile	24
1-535-9000	Needle guide	16Ga	Supersonic Imagine SE12-3, SEV12-3	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Supersonic Imagine SE12-3, SEV12-3	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Supersonic Imagine SE12-3, SEV12-3	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Supersonic Imagine SE12-3, SEV12-3	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Supersonic Imagine SE12-3, SEV12-3	Sterile	24
1-535-5145	Needle guide	16Ga	Terason 8EC4	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	Terason 8EC4	Sterile	24
1-535-5151	24" Soft biopsy kit	16Ga	Terason 8EC4	Sterile	24
1-535-9000	Needle guide	16Ga	Terason 8EC4(V)	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Terason 8EC4(V)	Sterile	24

## Protek™ Endocavity Needle Guides

Endocavity needle guides ensure precise needle placement when performing a biopsy or aspiration

- Features an extra wide funnel for optimal needle entry and smooth profile to ensure ultimate patient comfort
- Please note, kits have both Sterile and Non-sterile components unless otherwise noted on product label

Below you will find single-use and reusable needle guides.



Sterile, single-use example



Non-sterile, reusable example

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-9002	12" Soft biopsy kit	16Ga	Terason 8EC4(V)	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Terason 8EC4(V)	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Terason 8EC4(V)	Sterile	24
1-535-5531	Needle guide	16Ga	Toshiba (Canon) PVF-651VT, PVF-641VT, PVF-621VT, PVG-601V, PVM-651VT, PVN-661VT, PVT-661VT, PVQ-641V	Sterile	24
1-535-5535	24" SlimFit™ biopsy kit	16Ga	Toshiba (Canon) PVF-651VT, PVF-641VT, PVF-621VT, PVG-601V, PVM-651VT, PVN-661VT, PVT-661VT, PVQ-641V	Sterile	24
1-535-5560	Needle guide	16Ga	Toshiba (Canon) PVT-781VT, PVU-781VT, 11C3	Sterile	24
1-535-5561	24" SlimFit™ biopsy kit	16Ga	Toshiba (Canon) PVT-781VT, PVU-781VT, 11C3	Sterile	24
1-535-5562	3.5cm Latex biopsy kit	16Ga	Toshiba (Canon) PVT-781VT, PVU-781VT, 11C3	Sterile	24
1-535-5144	Needle guide	16Ga	Ultrasonix EC9-5/10	Sterile	24
1-535-5146	3.5cm Latex biopsy kit	16Ga	Ultrasonix EC9-5/10	Sterile	24
1-535-5149	12" Soft biopsy kit	16Ga	Ultrasonix EC9-5/10	Sterile	24
1-535-5183	48" SlimFit™ biopsy kit	16Ga	Ultrasonix EC9-5/10	Sterile	24
1-535-6456	Stainless steel needle guide	16Ga	Ultrasonix EC9-5/10	Non-sterile	1
1-535-9000	Needle guide	16Ga	Zonare E9-3	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Zonare E9-3	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Zonare E9-3	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Zonare E9-3	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Zonare E9-3	Sterile	24
1-535-7411	12" Soft biopsy kit	16Ga	Zonare E9-4	Sterile	24



### Protek™ Sterile Multi-Gauge Biopsy Needle Guide

#### Product Details

- Easily snaps onto bracket for precise needle guidance
- Clear, wide funnel to visualize needle entry
- Quick-release for easy needle removal
- Sterile, single-use
- Retrofit XT™ fits GE Healthcare brackets, Director XT™ fits all other OEM brackets
- Set includes needle guide with 10 optional needle gauge inserts (8.5F, 14-18Ga, 20-23Ga)

Part No.	Description	Qty
1-535-4500	Director XT™ Needle Guide Set	24
1-535-5600	Retrofit XT™ Needle Guide Set	24



### Protek™ Sterile Multi-Gauge Biopsy Needle Guide Probe Cover Kit

#### Product Details

- Easily snaps onto bracket for precise needle guidance
- Clear, wide funnel to visualize needle entry
- Quick-release for easy needle removal
- 10 optional needle gauge inserts (8.5F, 14-18Ga, 20-23Ga)
- Retrofit XT™ fits GE Healthcare brackets, Director XT™ fits all other OEM brackets
- Sterile kit contains single-use needle guide with multi-gauge insert set, 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover, sterile gel, and 2 bands

Part No.	Description	Qty
1-535-4501	Director XT™ Needle Guide Probe Cover Kit	24
1-535-5601	Retrofit XT™ Needle Guide Probe Cover Kit	24



### Protek™ Elastic Bands

#### Product Details

- Use to secure ultrasound probe covers
- Not made with natural rubber latex
- Non-sterile

Part No.	Description	Sterility	Qty
1-519-0126	1.25" x 1/8" Blue Bands	Non-sterile	2150





#### Protek™ Sterile Single-Gauge Biopsy Needle Guide Product Details

- Easily snaps onto bracket for precise needle guidance
- Clear, wide funnel to visualize needle entry
- Quick-release for easy needle removal
- Color-coded gauges
- Sterile, single-use

#### Director™

Part No.	Gauge	Qty
1-535-4013	13Ga	24
1-535-4014	14Ga	24
1-535-4015	15Ga	24
1-535-4016	16Ga	24
1-535-4018	18Ga	24
1-535-4020	20Ga	24
1-535-4021	21Ga	24
1-535-4022	22Ga	24
1-535-4025	25Ga	24



#### Protek™ Sterile Single-Gauge Biopsy Needle Guide and Probe Cover Kit

##### Product Details

- Easily snaps onto bracket for precise needle guidance
- Clear, wide funnel to visualize needle entry
- Quick-release for easy needle removal
- Color-coded gauges
- Contains single-use needle guide, 5.5" x 48" (14cm x 122cm) Soft PullUp™ probe cover, sterile gel, and 2 bands

#### Director™

Part No.	Gauge	Qty
1-535-4216	16Ga	24
1-535-4218	18Ga	24
1-535-4222	22Ga	24

### Protek™ General Purpose Needle Guide Starter Kits

General purpose (in-plane) needle guide starter kits contain the necessary components for a safe and precise biopsy. The custom fit single or multi-angle bracket is designed to meet predetermined angles defined by ultrasound system software guidelines. Single-use biopsy needle guides easily snap onto the bracket providing a fixed path for the needle. PullUp™ probe cover kits help to maintain procedure sterility and reduce cross-contamination.

- Please note, kits have both Sterile and Non-sterile components unless otherwise noted on product label
- 1 Kit/Box

### Examples of Protek™ Needle Guide Starter Kits



Multi-angle example



Single-angle example

Part No.	For use with	Bracket Angles	Starter Kit Content
1-535-7138	ECM C360A, C360V	Single-angle	(1) Non-sterile reusable bracket; (6) Sterile Director™ needle guide/ PullUp™ probe cover kits (contains guide (2 ea. 16Ga, 18Ga, 22Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands)
1-535-7363	Esaote C1-8, AC2541	Multi-angle	(1) Non-sterile reusable bracket (5) Sterile Director XT™ multi-gauge needle guide/ PullUp™ probe cover kits (Contains: multi-gauge needle guide set (8.5F, 14-18,20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel and bands
1-535-7070	GE Healthcare 4C, 4C-A, 4C-D, 4C-RS, 4C-SC, 5CS	Multi-angle	(1) Non-sterile reusable bracket; (6) Sterile Retrofit™ single-use needle guides (2 ea. 18Ga, 20Ga, 22Ga); (6) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7068	GE Healthcare 9L, 9L-D, 9L-RS	Multi-angle	(1) Non-sterile reusable bracket; (6) Sterile Retrofit™ single-use needle guides (2 ea. 18Ga, 20Ga, 22Ga); (6) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7308	GE Healthcare C1-5-D, C1-5-RS	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Retrofit™ XT™ multi-gauge needle guide sets (8.5F, 14-18,20-23Ga); (5) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7309	GE Healthcare C1-6, C1-6-D	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Retrofit™ XT™ multi-gauge needle guide sets (8.5F, 14-18,20-23Ga); (5) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7306	GE Healthcare C2-9-D, C2-9VN-D	Multi-angle	(1) Non-sterile reusable bracket (5) Sterile Retrofit XT™ multi-gauge needle guide/ PullUp™ probe cover kits (Contains: multi-gauge needle guide set (8.5F, 14-18,20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel and bands

Part No.	For use with	Bracket Angles	Starter Kit Content
1-535-7107	GE Healthcare ML6-15, ML6-15-D, ML6-15-RS	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Retrofit XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands
1-535-7278	Mindray C5-1, C5-1E, C5-1s, C5-1U, SC6-1E, SC6-1s, SC6-1U, SC5-1E, SC5-1N, SC5-1Ns	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Retrofit XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands
1-535-7293	Mindray L14-5w, L14-5WE, L14-5WU,	Multi-Angle	(1) Non-sterile reusable bracket (5) Sterile Director XT™ multi-gauge needle guide/PullUp™ probe cover kits (Contains: multi-gauge needle guide set (8.5F, 14-18,20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel and bands
1-535-7294	Mindray L14-3WU, L14-3Ws	Multi-angle	(1) Non-sterile reusable bracket (5) Sterile Director XT™ multi-gauge needle guide/PullUp™ probe cover kits (Contains: multi-gauge needle guide set (8.5F, 14-18,20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel and bands
1-535-7385	Philips C5-2 (HD3 Ultrasound System)	Single-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ needle guide/PullUp™ probe cover kits (contains guide (8 ea. 16Ga, 18Ga, 22Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands)
1-535-7154	Samsung C2-5	Single-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 20Ga, 22Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7457	Samsung C2-8	Single-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ needle guide/PullUp™ probe cover kits (contains guide (8 ea. 16Ga, 18Ga, 22Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands)
1-535-7180	Samsung CA1-7A, CA1-7AD	Multi-angle	(1) Non-sterile reusable bracket; (6) Sterile Director™ single-use needle guides (2 ea. 16Ga, 18Ga, 22Ga); (6) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7167	Samsung CA1-7A, CA1-7AD	Multi-angle	(1) Non-sterile reusable bracket; (36) Sterile Director™ single-use needle guides (6 ea. 14Ga, 16Ga, 18Ga, 20Ga, 22Ga, 25Ga); (36) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7196	Samsung CA1-7S	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Director XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands
1-535-7179	Samsung CF2-8, CA2-8A, CA2-8AD, CA2-9A, CA2-9AD	Multi-angle	(1) Non-sterile reusable bracket; (6) Sterile Director™ single-use needle guides (2 ea. 16Ga, 18Ga, 22Ga); (6) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7166	Samsung CF2-8, CA2-8A, CA2-8AD, CA2-9A, CA2-9AD	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 22Ga, 25Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7265	Samsung CF4-9	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Director XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands

### Protek™ General Purpose Needle Guide Starter Kits

General purpose (in-plane) needle guide starter kits contain the necessary components for a safe and precise biopsy. The custom fit single or multi-angle bracket is designed to meet predetermined angles defined by ultrasound system software guidelines. Single-use biopsy needle guides easily snap onto the bracket providing a fixed path for the needle. PullUp™ probe cover kits help to maintain procedure sterility and reduce cross-contamination.

- Please note, kits have both Sterile and Non-sterile components unless otherwise noted on product label
- 1 Kit/Box

#### Examples of Protek™ Needle Guide Starter Kits



Multi-angle example



Single-angle example

Part No.	For use with	Bracket Angles	Starter Kit Content
1-535-7168	Samsung CV1-8A, CV1-8AD	Multi-angle	(1) Non-sterile reusable bracket; (36) Sterile Director™ single-use needle guides (6 ea. 14Ga, 16Ga, 18Ga, 20Ga, 22Ga, 25Ga); (36) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7162	Samsung L3-12A, LA3-14A, LA3-14AD, L5-13/50	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 20Ga, 22Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7181	Samsung LA2-9A	Multi-angle	(1) Non-sterile reusable bracket; (6) Sterile Retrofit™ single-use needle guides (2 ea. 18Ga, 20Ga, 22Ga); (6) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7285	Samsung LA4-18A, LA4-18AD	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Director XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands
1-535-7197	Samsung LA4-18B, LA4-18BD	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Director XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands
1-535-7165	Samsung LF5-13, LA3-16A	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 20Ga, 22Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands

Part No.	For use with	Bracket Angles	Starter Kit Content
1-535-7172	Samsung LF5-13, LA3-16A, LA3-16AD	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 20Ga, 22Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7461	Samsung LN5-12, LE5-12	Single-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ needle guide/ PullUp™ probe cover kits (contains guide (8 ea. 16Ga, 18Ga, 22Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands)
1-535-7164	Samsung SC1-6	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 22Ga, 25Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7250	Samsung VN4-8	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 20Ga, 21Ga, 22Ga, 25Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7138	Terason 4C2-V, 5C2, 5C2A	Single-angle	(1) Non-sterile reusable bracket; (6) Sterile Director™ needle guide/ PullUp™ probe cover kits (contains guide (2 ea. 16Ga, 18Ga, 22Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands)
1-535-7024	Zonare C5-2, C6-2	Single-angle	(1) Non-sterile reusable bracket; (6) Sterile Director™ single-use needle guides (2 ea. 16Ga, 18Ga, 22Ga); (6) Sterile 5.5" x 24" (14cm x 61cm) soft 3D ended probe cover kits with sterile gel and bands



### Protek™ ProV-Access™ Single-Depth Needle Guides

#### Product Details

- Attaches securely to bracket or Bard Site-Rite™ probe
- Clear pad to visualize needle entry
- Fits needle gauges 18 through 21
- Patented needle release for tilt out procedures
- Color-coded depths
- Sterile, single-use

Part No.	Depth	Sterility	Qty
1-535-5404	0.5cm (yellow)	Sterile	10
1-535-5403	1.0cm (blue)	Sterile	10
1-535-5402	1.5cm (white)	Sterile	10
1-535-5405	2.0cm (green)	Sterile	10
1-535-5406	2.5cm (pink)	Sterile	10
1-535-5401	3.5cm (red)	Sterile	10



### Protek™ ProV-Access™ 3-Depth Needle Guide Probe Cover Kits

#### Product Details

- Easily snaps onto bracket or probe
- Clear pad to visualize needle entry
- Fits needle gauges 18 through 21
- Patented needle release for tilt out procedures
- Color-coded depths
- Sterile, single-use
- Kit contains 3-depth needle guide set, 5" x 48" (13cm x 122cm) clear PullUp™ probe cover, sterile gel, and 2 bands

Part No.	Depths	Sterility	Qty
1-535-5414	2.0cm, 1.5cm, 1.0cm	Sterile	10
1-535-5410	3.5cm, 1.5cm, 1.0cm	Sterile	10
1-535-5416	3.5cm, 2.5cm, 1.5cm	Sterile	10



### Protek™ Vascular Needle Guide Starter Kit

Vascular (out-of-plane) needle guide starter kits contain the necessary components for a safe and precise biopsy. The custom fit bracket securely attaches to the probe with ease. Single-use, depth specific, ProV-Access™ needle guides snap onto the bracket providing a fixed path for the needle. PullUp™ probe cover kits help to maintain procedure sterility and reduce cross-contamination.

- Please note, kits have both Sterile and Non-sterile components unless otherwise noted on product label

#### Example of Protek™ Vascular Needle Guide Starter Kit



Part No.	For use with	Starter Kit Content
1-535-7311	GE Healthcare 4C, 4C-A, 4C-D, 4C-RS, 4C-SC, 5CS	(1) Non-sterile reusable vascular bracket; (3) Sterile 1.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7397	Hitachi L64	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7292	Mindray C4-1	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7276	Mindray C5-2E, C5-2s	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7277	Mindray L12-4s, 7L4s, L14-6Ns, L14-6NE, L12-3E	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7273	Mindray L14-5w, L14-5WE, L14-5WU	(1) Non-sterile reusable vascular bracket; (3) Sterile 1.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7274	Mindray L9-3U, L9-3s, and L9-3E	(1) Non-sterile reusable vascular bracket; (3) Sterile 1.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7279	Mindray SC8-2s, SC8-2U, SC8-2Us	(1) Non-sterile reusable vascular bracket; (3) Sterile 1.5cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands

Part No.	For use with	Starter Kit Content
1-535-7419	Sonosite C60, C60x	(1) Non-sterile reusable vascular bracket; (3) Sterile 1.5cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7492	Terason 15L4A, 15L4 Smart Mark™	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7536	Toshiba (Canon) PLI-1205BX (i18LX5)	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; Fits needle gauges 18 through 21; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7292	Zonare C4-1	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7296	Zonare L10-5	(1) Non-sterile reusable vascular bracket; (3) Sterile ProV-Access™ 18-21 Ga Disposable Needle Guide/PullUp™ probe cover kits (contains 3-depth needle guide set (2.0cm, 1.5cm, 1.0cm), 5" x 48" (13cm x 122cm) clear probe cover, sterile gel & bands)

## Suture Retrievers

The North American-made suture retrievers below provide high quality for the clinical staff, while at the same time offering cost savings to the medical facility. The knurled aluminum handles and distinctive coloring for ease of use and comfort in-hand. These suture retrievers are designed for single patient use. Intended to aid in the passing of sutures through bone and soft tissue during orthopedic surgery.



### Hoffee® Blue Lasso Suture Retriever

#### Product Details

- Comparable to Smith & Nephew® Hewson® Suture Retriever 7111-1579 and Arthrex® Transosseous Suture Passer AR-1000
- Hewson Suture Retriever replacement
- Aluminum knurled blue handle is weighted for ease of use
- Blue suture loop for increased visibility
- 10" length is suitable for most laparoscopic procedures

\*Please inquire about Non-sterile suture retrievers for custom packs

Part No.	Sterility	Qty
710000	Sterile, intended for single-use	5 Eaches/Box



### Green Genie Suture Passer

#### Product Details

- Comparable to SWANSON™ Suture Passer
- Aluminum knurled green handle is weighted for ease of use
- Wire suture loop
- 6.5" length
- Made in North America

Part No.	Sterility	Qty
710001	Sterile, intended for single-use	5 Eaches/Box



## Beach Chair Accessories

The North American-made head restraints are used with the following beach chairs: Allen®, Arthrex®, SchureMed™, Skytron®, Smith & Nephew® and Steris®. Each surgical positioning foam product is designed for single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.



### Erin Beach Chair Face Mask

#### Product Details

- For use with Smith & Nephew® T-MAX™ Beach Chair 7210996
- Comparable to Smith & Nephew® T-MAX™ Disposable Face Mask 7210559
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
711001	Non-sterile, intended for single-use	10/Box



### Erin Beach Chair Face Mask Replacement Straps

#### Product Details

- For use with Smith & Nephew® T-MAX™ Beach Chair 7210996
- Comparable to Smith & Nephew® T-MAX™ Replacement Strap 7210558
- Compatible with Erin Beach Chair Face Mask (711001)
- Reusable

Part No.	Sterility	Qty
711002	Non-sterile	1 Pair/Box



### Universal Head Restraint

#### Product Details

- For use with Allen® Beach Chair Positioner and Arthrex® Beach Chair Positioner
- Comparable to Allen® Universal™ Head Restraint A-90023 and Arthrex® Universal Head Restraint AR-1627-06
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
711004	Non-sterile, intended for single-use	10/Box



## Indy Beach Chair Face Mask

### Product Details

- For use with Allen® Flexiform™ Head Rest A-90017
- Comparable to Allen® Flexiform™ Head Rest Disposable A-90019
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
711003	Non-sterile, intended for single-use	10/Box



## Bethesda Beach Chair Face Mask

### Product Details

- For use with SchureMed™ Beach Chair 800-0004, Stryker® Shoulder Chair OZ 7519475 and OZ 7523075, and Skytron® Beach Chair
- Comparable to SchureMed™ Disposable Face Mask 800-0164, Stryker® Disposable Mask OZ 7547575, and Skytron® Disposable Face Mask 3-040-33
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
711000	Non-sterile, intended for single-use	10/Box



## Beach Chair Head and Chin Straps

### Product Details

- For use with Skytron® Beach Chair, SchureMed™ Beach Chair, Steris® Beach Chair, Stryker® Shoulder Chair OZ 7519475 and OZ 7523075
- Comparable to Skytron® Beach Chair Head & Chin Straps 3-040-22, SchureMed™ Disposable Head & Chin Strap 800-0039, Steris® Head Restraint Straps BF039, and Stryker® Chin Strap OZ 7518375
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
711005	Non-sterile, intended for single-use	20/Box



## Shoulder Positioner Kits

The North American-made sterile kits provide high quality for the clinical staff, while at the same time offering cost savings to the medical facility. For example, the lateral kits feature a snap clip to more securely attach the product to the frame and cost less than the brand name. Each kit is designed for single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.



### Marco Shoulder Stabilization Kit

#### Product Details

- For use with Smith & Nephew® SPIDER™ Limb Positioner 7210569
- Comparable to Smith & Nephew® SPIDER™ Shoulder Stabilization Kit 7210573
- Includes Drape, Arm & Hand Cuffs, Hand Grip, and 4" x 5yd 3M® Coban™
- Kit allows for quick and easy patient prep
- Soft, non-abrasive foam for patient comfort
- Not made of natural rubber latex

Part No.	Sterility	Qty
711100	Sterile, intended for single-use	12/Box



### Berkette Beach Chair Limb Positioner Kit

#### Product Details

- For use with Arthrex® Trimano® Shoulder Positioner AR-1640
- Comparable to Arthrex® Trimano® Beach Chair Shoulder Kit AR-1644
- Includes Foam Sleeve, Hand Grip, Drape, and 4" x 5yd 3M® Coban™
- Kit allows for quick and easy patient prep
- Soft, non-abrasive foam for patient comfort
- Not made of natural rubber latex

Part No.	Sterility	Qty
711200	Sterile, intended for single-use	6/Box



## Jessie Arm Suspension Kit

### Product Details

- For Use With McConnell Limb Positioner
- Comparable To McConnell Arm Support 12-401
- Includes Arm Wrap, Placement Strap, and Hand Hexapod
- Kit allows for quick and easy patient prep
- Soft, non-abrasive foam for patient comfort
- Not made of natural rubber latex

Part No.	Sterility	Qty
711300	Sterile, intended for single-use	10/Box



## Norton Arm Restraint Kit

### Product Details

- For Use With Allen® Beach Chair A-92000
- Comparable To Allen® Arm Positioner Restraint A-92001
- Includes Foam Arm Drape
- Kit allows for quick and easy patient prep
- Soft, non-abrasive foam for patient comfort
- Not made of natural rubber latex

Part No.	Sterility	Qty
711400	Sterile, intended for single-use	10/Box



## Dr Don's® Arm Canoe

### Product Details

- For Use With Arthrex® Lateral Decubitus Shoulder Traction Tower AR-1630 and Arthrex® 3-Point Shoulder Distraction System AR-1600M
- Comparable To Arthrex® STAR™ Sleeve AR-1606
- Includes Convolved Foam Sleeve, Arm Support Wrap, and 4" x 5yd 3M® Coban™
- Soft, non-abrasive foam for patient comfort
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure
- Not made of natural rubber latex

Part No.	Sterility	Qty
711500	Sterile, intended for single-use	5/Box



## Dr Don's® Arm Canoe-V

### Product Details

- For Use With Arthrex® Lateral Decubitus Shoulder Traction Tower AR-1630 and Arthrex® 3-Point Shoulder Distraction System AR-1600M
- Comparable To Arthrex® STAR™ Sleeve AR-1606V
- Includes Convolved Foam Sleeve and Arm Support Wrap
- Soft, non-abrasive foam for patient comfort
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure
- Not made of natural rubber latex

Part No.	Sterility	Qty
711501	Sterile, intended for single-use	5/Box



## Dr Don's® Arm Suspension Kit

### Product Details

- Comparable To Instrumed® Upper Extremity Drape IN-9300CLF and DeRoyal® Shoulder Suspension Kit 32-525
- Includes fluid resistant sleeve, rope and snap hook, and 4" x 5yd 3M® Coban™
- Fluid resistant stockinette allows for secure patient positioning
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure
- Not made of natural rubber latex

Part No.	Sterility	Qty
711502	Sterile, intended for single-use	10/Box



## Dr Don's® Arm Suspension Kit – Montana

### Product Details

- Comparable To Smith & Nephew® Shoulder Suspension Kit 72200195 and Zimmer Biomet® Shoulder Suspension Kit 905026
- Includes Reinforced sleeve, rope and snap hook and 4" x 5yd 3M® Coban™
- Reinforced, soft cloth lined stockinette allows for secure patient positioning and patient comfort
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure
- Not made of natural rubber latex

Part No.	Sterility	Qty
711503	Sterile, intended for single-use	10/Box



## Dr Don's® Arm Suspension Kit – Northern

### Product Details

- Comparable To Arthrex® Lateral Arm Traction Sleeve AR-1635
- Includes Jersey knit traction sleeve, grip lining, rope and snap hook, 4" x 5yd 3M® Coban™ and foam hand grip
- Jersey knit stockinette with interior grip lining allows for secure patient positioning
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure
- Not made of natural rubber latex

Part No.	Sterility	Qty
711504	Sterile, intended for single-use	10/Box

## Hip Positioners

The North American-made positioning foam products below are compatible with Innovative Medical Products® (IMP®), Mizuho OSI® and Smith & Nephew® hip positioning devices. These positioning foam products are designed for single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.



### Carson 8" Perineal Post Pad

#### Product Details

- For Use With Smith & Nephew® Supine Hip Positioning System 72200624 and Mizuho OSI® Hana®, ProFx®, and Hana SSXT® Tables
- Comparable To Smith & Nephew® Supine Perineal Post Pad 72200634
- Soft, non-abrasive foam for patient comfort
- Outside diameter: 9", Height: 11" Core Diameter: 1.5"

Part No.	Sterility	Qty
712200	Non-sterile, intended for single-use	1/Box



### Patient Care Boot Liner Pad

#### Product Details

- For Use With Smith & Nephew® Supine Hip Positioning System 72200624 and Smith & Nephew® SPIDER™ Boot 72201425
- Comparable To Smith & Nephew® Boot Insert Pads 72200635 and Smith & Nephew® SPIDER® Boot Liner 72201426
- Soft, non-abrasive foam for patient comfort
- IFU printed on product for quick reference

Part No.	Sterility	Qty
712300	Non-sterile, intended for single-use	1 Pair/Box



### Perineal Post Cover

#### Product Details

- For Use With Mizuho OSI® Hana®, ProFx®, and Hana SSXT® Tables
- Comparable To Mizuho OSI® Perineal Post Cover 6855-13
- Soft, non-abrasive foam for patient comfort
- Dimensions: 11.5" x 5.5"

Part No.	Sterility	Qty
712400	Non-sterile, intended for single-use	12/Box



## Lewiston Positioner Kit

### Product Details

- For Use With Arthrex® Hip Distraction System AR-6529S
- Comparable To Arthrex® Disposable Hip Kit AR-6529-11
- Includes Sueded perineal post pad cover and 1 pair of traction boot pads
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
712100	Non-sterile, intended for single-use	1/Box



## Lindy Traction Boot Liner

### Product Details

- For Use With Allen® Traction Boot A-93005 and A-93000 and Steris® OT 1000 Traction Boot
- Comparable To Allen® Traction Boot Disposable A-93001 and STERIS® Traction Boot Liner BF419
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
712900	Non-sterile, intended for single-use	6/Box



## Lateral Positioner Pad

### Product Details

- For Use With IMP® Universal Lateral Positioner® 1092
- Comparable To IMP® Universal Lateral Positioner® Economy Patient Protective Pad Set 1092-P
- Soft, non-abrasive foam for patient comfort
- Includes anterior support pads

Part No.	Sterility	Qty
712600	Non-sterile, intended for single-use	15/Box





## Racine Pelvic Positioning Pad

### Product Details

- For Use With IMP® McGuire™ Hip Positioner 1084
- Comparable To IMP® McGuire™ Hip Pads 1084-P
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
712800	Non-sterile, intended for single-use	15/Box

## Knee & Ankle Positioners

Each product below provides high quality for the clinical staff, while at the same time offering cost savings to the medical facility. These sterile kits and positioning foam products are designed for single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.



## Hermanator® Knee Positioning Pad

### Product Details

- For Use With IMP® De Mayo Knee Positioner® 803 AND 803-S, Zimmer® Alvarado Knee Holder 1320-000-00, SchureMed™ TKR Knee Positioner, and Innomed® Robb Knee Positioner
- Comparable To IMP® De Mayo Knee Positioner Pad 294-P-10 and 803-GP-10, Zimmer® Alvarado Disposable Liner 00-1320-015-10, Innomed® Sterile Protective Pad 2629-00, and SchureMed™ TKR Disposable Liner 800-0146
- Helps decrease device related mechanical pressure injury, especially to the high risk areas of the posterior, medial, and lateral malleolus
- Soft, high-density memory foam provides enhanced patient safety
- Not made of natural rubber latex

Part No.	Sterility	Qty
713100	Sterile, intended for single-use	10/Box



### XL Hermanator® Pink Knee Positioning Pad

#### Product Details

- For use with IMP® De Mayo Knee Positioner® 803 AND 803-S and Match Grade Ladderloc Knee Positioner
- Comparable to IMP® De Mayo Knee Positioner® 803-GP-10
- Helps decrease device related mechanical pressure injury, especially to the high risk areas of the posterior, medial, and later malleolus
- Soft, high-density memory foam provides enhanced patient safety
- Not made of natural rubber latex

Part No.	Sterility	Qty
713300	Sterile, intended for single-use	10/Box



### Sirako® Knee Positioning Pad

#### Product Details

- For use with the Stryker® Mako® Leg Positioner
- Comparable to Stryker® Mako® Leg Positioner Disposable Kit
- Secures patient's foot into the leg holder boot to help stabilize the knee during surgery
- Helps decrease device related mechanical pressure injury
- Sterile, intended for single-use and not made with natural rubber latex

Part No.	Sterility	Qty
713800	Sterile, intended for single-use	6/Box



### TLC Distraction Pad Set

#### Product Details

- For Use With IMP® De Mayo Universal Distractor® 907
- Comparable To IMP® De Mayo Pressure Protector Pad® Set 907-GP-10
- Includes Memory pad, 4" x 5yd 3M® Coban™ and Knee Distraction Guard
- Semi circle memory pad prevents device from rolling into the popliteal area while providing a safe view into the posterior capsule
- Helps protect the posterior thigh from mechanical pressure injury
- Soft, high-density memory foam provides enhanced patient safety
- Not made of natural rubber latex

Part No.	Sterility	Qty
713200	Sterile, intended for single-use	10/Box



### Leg Comfort Pad

#### Product Details

- For Use With Allen® Leg Holder R-734-33, Smith & Nephew® Acufex Knee Holder, SchureMed™ Schure Grip Leg Holder 800-0065, and Stryker® Arthroscopy Leg Holder OZ 7511975
- Comparable To Allen® Leg Holder Disposable R-S4133-V, Smith & Nephew® Disposable Leg Holder Pads 3757, SchureMed™ Schure Grip Leg Holder Pads 508-0103, and Stryker® Disposable Pad OZ 7512075
- Inner foam is soft and non-abrasive for patient comfort
- Outer foam is a dense material that helps provide appropriate clamping pressure to stabilize the leg

Part No.	Sterility	Qty
713600	Non-sterile, intended for single-use	12/Box



## Smoot Ankle Distraction Strap

### Product Details

- For Use With Smith & Nephew® GUHL™ Positioner and Arthrex® Ankle Distractor
- Comparable To Smith & Nephew® GUHL™ Ankle Distraction Strap 014407 and Arthrex® Ankle Distraction Strap AR-1712
- Soft, non-abrasive foam where ankle meets strap for patient comfort
- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure

Part No.	Sterility	Qty
713400	Sterile, intended for single-use	10/Box



## Thigh Holder Pad

### Product Details

- For Use With Smith & Nephew® FERKEL™ Thigh System 72201812
- Comparable To Smith & Nephew® FERKEL™ Thigh System 72201814
- Soft, non-abrasive foam for patient comfort

Part No.	Sterility	Qty
713500	Non-sterile, intended for single-use	6/Box



## Parent Post Knee Pad

### Product Details

- For Use With Arthrex® PivotPost® AR-1560
- Comparable To Arthrex® PivotPost® AR-1560-01
- Soft, non-abrasive foam for patient comfort
- Velcro® straps allow for efficient placement on PivotPost®

Part No.	Sterility	Qty
712500	Non-sterile, intended for single-use	5/Box

## Spine Positioning Foam

The North American-made positioning products are designed for single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.



### Prone Face Pillow

#### Product Details

- Comparable To Allen® Foam Head Block A-70862
- Consists of two specifically designed foam layers: bottom layer has supportive foam with pre-cut dual intubation slots with the top layer consisting of soft, non-abrasive foam designed for patient comfort
- Can be used with the C-Flex® Head Positioning System, Allen® Foam Block Stand, and the Allen® Advance Table Head Plate or any flat surface
- 5" Foam
- Memory Foam Topping
- Compressed Packaging

Part No.	Sterility	Qty
714100	Non-sterile, intended for single-use	12/Box



### Sedona Arm Positioner Set

#### Product Details

- Comparable To Mizuho OSI® Laminectomy Arm Cradle 5823-10
- Soft, non-abrasive foam pads the entire length of the patient's arm and has notches to allow for positioning of the patient's arm
- Dimensions: 24" x 5" x 3"
- Can be used in prone or supine position

Part No.	Sterility	Qty
714200	Non-sterile, intended for single-use	6 Pairs/Box

## Orthopedic Accessory Products

The North American-made sterile kits and positioning foam products are designed for specialty patient positioning needs and single patient use. Intended to support and maintain the position of an anesthetized patient during surgery.



\*Please inquire about Non-sterile MD Traction Rope for custom packs

### MD Traction Rope

#### Product Details

- Carabiner clip prevents accidental S-hook slips from occurring mid-procedure Snap Hook
- Rope Length Approx 89"

Part No.	Sterility	Qty
715000	Sterile, intended for single-use	50/Box



### Brachytherapy Mattress

#### Product Details

- Soft, non-abrasive foam for patient comfort
- 30" x 79" x 5"
- 18" x 6" Opening

Part No.	Sterility	Qty
715100	Non-sterile, intended for single-use	1/Box






## Colby™ Fluid Management Products

Our Colby™ line provides a convenient and effective means of managing infectious waste fluid generated during wet procedures. Our products help to reduce the likelihood of injuries associated with wet floors. They aid in protecting healthcare professionals and housekeeping personnel from contact with contaminated substances and help support your efforts to comply with OSHA's Bloodborne Pathogens Standard. Our products do not change the contents of medical waste, but rather help your team to more effectively contain it.

### SurgiSafe® Absorbent Floor Mats

- Help keep floors clean, dry, and slip-free, even during fluid intense procedures
- Contains fluid in one place to help make operating room turnover time faster and more efficient
- Absorb irrigates, blood, and bodily fluids through unique wicking, absorbency, and fluid retention properties, while conforming to uneven floor surfaces without bunching
- Help protect your staff from slips, trips, and falls in the operating room
- Help reduce staff exposure to infectious waste
- Helps eliminate using blankets and towels which can cause tripping hazards
- Non-sterile, intended for single-use, and not made with natural rubber latex



### SurgiSafe® Specialty Absorbent Floor Mats

#### Product Details

- Designed for strength and durability
- Features low profile design, which allows for foot pedals or stools to roll with ease
- Offered in large sizes to accommodate surgical procedures where more than one person needs to stand on the mat
- Blue mat features fluid barrier backing to help keep floor beneath dry and clean
- Yellow mat absorbs fluid from the top and bottom

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
82350	Blue	Medium	Yes	21" x 44"	50/Case
84625	Blue	Medium	Yes	32" x 46"	25/Case
84610	Blue	Medium	Yes	32" x 46"	10/Case
84672	Blue	Medium	Yes	32" x 72"	10/Case
81860-2	Yellow	Medium	No	16" x 45"	60/Case
83630-2	Yellow	Medium	No	32" x 45"	30/Case
83610-2	Yellow	Medium	No	32" x 45"	10/Case
83615-2L	Yellow	Medium	No	32" x 80"	15/Case
83625-2L	Yellow	Medium	No	32" x 80"	12/Case
83630-M	Yellow	Medium	No	32" x 45"	15/Case
	White	Medium	No	28" x 40"	15/Case



## SurgiSafe® Standard Absorbent Floor Mats

### Product Details

- Made with a minimum of 60% pre-consumer recycled material
- Available in three absorbency levels as well as with and without fluid barrier backing for convenience
- Pink – High level of fluid absorbency, will absorb 16x's its dry weight
- White – Medium level of fluid absorbency, will absorb 14x's its dry weight
- Green – Low level of fluid absorbency, will absorb 11x's its dry weight

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
83430	Pink	High	No	28" x 40"	30/Case
83472	Pink	High	No	28" x 72"	30/Case
83330	Pink	High	Yes	28" x 40"	30/Case
83372	Pink	High	Yes	28" x 72"	30/Case
81860	White	Medium	No	18" x 40"	60/Case
818100	White	Medium	No	18" x 21"	100/Case
81850	White	Medium	No	18" x 21"	50/Case
83610	White	Medium	No	28" x 40"	10/Case
83630	White	Medium	No	28" x 40"	30/Case
83625-L	White	Medium	No	28" x 72"	12/Case
83672	White	Medium	No	28" x 72"	30/Case
83530	White	Medium	Yes	28" x 40"	30/Case
83572	White	Medium	Yes	28" x 72"	30/Case
84030	Green	Low	No	28" x 40"	30/Case
83030	Green	Low	Yes	28" x 40"	30/Case
83072	Green	Low	Yes	28" x 72"	30/Case
830125	Green	Low	Yes	28" x 125"	24/Case



## SurgiSafe® EvacMat® Absorbent Floor Mat

### Product Details

- Combines the functionality of absorption and suction to aid in efficiently removing fluid waste from the O.R. floor
- Extremely quiet operation helps to decrease distraction of the O.R. staff– as quiet as an indoor conversation
- The top layer of the mat contains wicking to help prevent pooling
- Works well with wall suction
- Simple to use and clean up – red bag friendly
- Low profile allows for the O.R. staff to roll equipment over the mat without loss of suction
- Large enough for two people
- Provides cushion to help reduce fatigue during long procedures
- Aids in keeping the O.R. floor dry to reduce the risk of injuries elated to wet floors
- Non-sterile, intended for single-use, and not made with natural rubber latex

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
95005	Yellow	Medium	Yes	32" x 48"	5/Case



## SurgiSafe® CinchMat® Absorbent Floor Mat

### Product Details

- Combines the effective absorbency of the SurgiSafe® mat with an easier way to handle the mat during disposal
- Pull cord design helps make clean-up and O.R. turnaround time simple and fast
- The top layer has wicking abilities that quickly absorbs irrigants, blood, and bodily fluids from the O.R. floor to help prevent pooling and splashing
- Skid-resistant fluid barrier backing helps ensure that the floor beneath the mat remains dry and prevents staining
- Aids in keeping the O.R. floor dry to reduce the risk of injuries related to wet floors
- Non-sterile, intended for single-use, and not made with natural rubber latex

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
84625-C	Blue	Medium	Yes	32" x 46"	10/Case



## ErgoSupport® Anti-Fatigue Mat

- Intended to help reduce fatigue that healthcare personnel experience during long hours of standing in one location
- Uses a dual foam and gel technology to provide comfort under the foot
- The top surface provides a necessary friction level with damp or dry surgical booties
- The bottom surface is made of a high-traction material to provide slip resistance
- Available in four sizes
- Water-resistant
- Reusable for up to 6-12 months
- Non-sterile and not made with natural rubber latex

Part No.	Description	Color	Size	Qty/Case
82001	Anti-Fatigue Mat	Red	18" x 24"	1
82005	Anti-Fatigue Mat	Red	18" x 24"	5
83002	Anti-Fatigue Mat	Red	20" x 48"	2
84003	Anti-Fatigue Mat	Red	20" x 32"	3



## Entrance Mat (TacTec)

- Intended to help remove foot borne particulate before entering the cleanroom
- Full coverage adhesive backing for adhering to the floor for use
- Seamless one-piece construction
- Numbered tabs provide easy identification and sheet separation
- Easy mobility for carts and wheeled machinery
- Polyethelene film contains a special anti-microbial agent providing long lasting protection against bacterial growth
- Non-sterile and not made with natural rubber latex

Part No.	Description	Size	Qty/Case
22-8110	Entrance Mat (TacTec)	18" x 45"	4 (30 peel off layers)





## ErgoSupport® Anti-Fatigue Step Stool Mat

- Intended to help reduce fatigue that healthcare personnel experience during long hours of standing in one location
- Made with polyurethane foam to provide comfort under the foot
- Matte textured surface provides optimal friction to help prevent slips and falls
- Features a low-profile, durable edges, and slip-resistant bottom
- Sides are beveled to prevent tripping
- Small size fits on operating room step stool
- Fluid and dirt resistant with anti-microbial properties
- Non-sterile and not made with natural rubber latex

Part No.	Description	Sterility	Color	Qty/Case
85005	Anti-Fatigue Mat	Non-sterile	Black	5



## ErgoSupport® Anti-Fatigue Step Stool

- Stackable and interlocking cast aluminum step stool
- Protected with an anti-microbial coating
- Designed to eliminate sites for liquid pooling and debris accumulation to help reduce contamination
- High density rubber foot pegs
- Compatible with the ErgoSupport® Anti-Fatigue Step Stool Mat
- Anti-slip surface for better stability
- Non-sterile and not made with natural rubber latex

Part No.	Description	Sterility	Color	Qty/Case
ST-001	Step Stool	Non-sterile	Silver	4





## WaterBug Quiet Floor Suction® Device

### Product Details

- Efficiently suction fluid waste to help reduce risks and maintain a dry operating room floor
- Its raised design and 12' of flexible tubing allows for easy repositioning
- Suctions unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Tubing	Qty
90010	WaterBug Quiet Floor Suction® Device	Yes	10/Case
90020	WaterBug Quiet Floor Suction® Device	Yes	20/Case
90030	WaterBug Quiet Floor Suction® Device	Yes	30/Case



## WaterBoom® Floor Suction Strip

### Product Details

- Helps hold back, contain, and quietly suction an unlimited amount of waste fluid
- Provides strong suction, while conforming to your O.R. floor
- Low profile helps you step over it without suction loss
- 12' tubing attaches from the side to help reduce tripping hazards
- This unobtrusive device does not require repositioning and is a good choice for those who prefer nothing at their feet
- Non-sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Length	Tubing	Qty
94010	WaterBoom® Floor Suction Strip	32"	Yes	10/Case
94030	WaterBoom® Floor Suction Strip	32"	Yes	30/Case



### AquaNot® Floor Suction Device

#### Product Details

- Has similar characteristics as the WaterBug Quiet Floor Suction® Device, but has a lower profile
- Efficiently suctions fluid waste to help reduce risks and maintains a dry operating room floor
- Suctions unlimited amounts of irrigants, blood, and other bodily fluids from the floor
- For thorough pickup, fluid is captured from the top via holes as well as from the bottom through suction and wicking action
- Non-sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Tubing	Qty
97020	AquaNot® Floor Suction Device	Yes	20/Case
97060	AquaNot® Floor Suction Device	Yes	60/Case
97050-NT	AquaNot® Floor Suction Device	No	50/Case



### ViraSorb® Solidifier

#### Product Details

- Evenly gels canister contents in seconds without stirring or agitating
- Non-toxic, Glutaraldehyde-free, and can be safely incinerated
- It does not revert back to a fluid state, helping to eliminate spillage and leakage, qualifying it as Solid Regulated Medical Waste
- Pour nozzles are included with multi-use bottles
- Non-sterile and not made with natural rubber latex

Part No.	Description	Solidifies	Qty
1000	ViraSorb® Solidifier – Single-Use Bottle	1,200cc	96/Case
1001	ViraSorb® Solidifier – Single-Use Bottle	1,500cc	96/Case
1002	ViraSorb® Solidifier – Single-Use Bottle	2,000cc	36/Case
1003	ViraSorb® Solidifier – Single-Use Bottle	3,000cc	36/Case
1005	ViraSorb® Solidifier – Single-Use Pouch	500cc	200/Case
1010	ViraSorb® Solidifier – Multi-Use Bottle	10,000cc	12/Case
1020	ViraSorb® Solidifier – Multi-Use Bottle	20,000cc	12/Case



### AspenBasin® Fluid Waste Container

- Used for general disposal of waste fluids in the Catheterization Lab
- Helps reduce exposure to contaminated fluid and bloodborne pathogens
- Grooved disposal area allows for convenient placement of syringe during disposal
- Includes high fluid capacity exceeding 500ml
- Secure lid design allows for easy disposal
- Sterile, intended for single-use, and not made of natural rubber latex


Part No.	Description	Qty/Box
AB100	AspenBasin® Fluid Waste Container	10



## Stork® Labor and Delivery Products

The birth of a child is one of the riskiest and most costly procedures in a hospital. Sacrificing high quality for cost savings is not an option when liabilities are so high. Our traditional obstetric products have been trusted for well over 100 million births in consistently providing positive outcomes and the cost-effective risk mitigation.

All of our brands are produced to high manufacturing standards and exacting quality specifications. Our Labor and Delivery products must pass through our innovative, rigorous laser vision inspection systems to ensure that only products meeting 100% of our tough criteria will be the first and only product required for the procedure. We have been manufacturing these products and perfecting these processes in the same U.S. facility for decades since their inception. Our brands are built on this experience and our commitment to the highest quality, proven every day, every birth.



### AmniHook® Amniotic Membrane Perforator

#### Product Details

- 100% laser-vision inspected
- Rounded tip and protected point designed to rupture the amniotic membrane without injury to mother or fetus
- 10" long, with non-glare molded plastic and sides ensure a non-slip grip and precise control
- Disposable, pre-sterilized, individually packaged for convenience

Part No.	Description	Sterility	Qty
9601	Membrane Perforator	Sterile	100 ea/bx



### Double-Grip® Cord Clamp

#### Product Details

- Safe, disposable, easy-to-apply
- Dual gripping surface prevents slipping during application
- Maintains tight seal and applies constant pressure on any size cord as it shrinks
- Hinge guard confines cord between serrated jaws to eliminate movement
- No rough edges to snag clothing or get caught on bed linen

Part No.	Description	Sterility	Qty
9423	Cord Clamp	Sterile	125 ea/bx
9411	Cord Clamp	Non-sterile	250 ea/bx



### Cord-Clamp™ Clipper

#### Product Details

- Concealed cutting surface eliminates risk of sharp or pointed objects near baby
- One time single patient use only to prevent risk of cross-contamination
- Designed to work exclusively with our Double-Grip® Umbilical Cord Clamp, our Cord Clamp® Clipper safely cuts through the specially designed clamp hinge for easy removal from baby

Part No.	Description	Sterility	Qty
9441	Single Clipper	Non-sterile	1 ea/bx
9445	Bulk Clipper	Non-sterile	25 ea/bx



### U-Bag® Urine Specimen Collector

#### Product Details

- Transparent butterfly patch makes application and removal easy and accurate reducing valuable nursing time
- Compatible and comfortable for male and female anatomy
- Double-chambered construction with no-flow back valve keeps urine away from skin to help prevent irritation and chance of contaminating the urine specimen
- Adhesive patch can be folded in half to send sealed specimen to the lab
- Preemie size ideal for NICU use with low birth weight infants
- Designed with an infant's tender skin in mind

### U-Bag® with Porous Cloth Adhesive Single Specimen Style

Part No.	Description	Sterility	Qty/Ea/Box
7540	Preemie (10mL)	Sterile	100
7535	Newborn (100 mL)	Sterile	100
7531	Pediatric (200 mL)	Sterile	100

#### 24-hour Style

7537	Newborn (100 mL)	Sterile	50
7533	Pediatric (200 mL)	Sterile	50





### Ready-Rolled® FootPrinter

#### Product Details

- Disposable plastic plate fits easily in the palm of one hand
- 4½" x 2¾" less than ⅛" thick with stackable edges for easy, compact storage
- Non-toxic, black ink-easy to see and wash off of baby
- Great for our Heirloom Quality Birth Certificates, Newborn Identification Forms, or other hospital made mementos
- Designed to guard against over-inking

Part No.	Description	Qty
5402	Plastic FootPrinter	25 ea/pk – 4 pks/bx



### Identification Form

#### Product Details

- Disposable plastic plate fits easily in the palm of one hand
- Receives and retains every detail of a baby's tiny skin print
- Space for recording pertinent birth information and verifying signatures
- Enough space for two infant footprints and a mother's fingerprint
- Some forms allow for the attachment of infant's ID band
- Printed on high contrast, archival quality paper, perfect for both receiving and preserving infant foot prints.

Part No.	Description	Size	Sterility	Pouches/Box
5821	Newborn ID Form	8" x 5"	Non-sterile	50 sheets/pad, 5 pads/box
5823	Newborn ID Form	8" x 11"	Non-sterile	50 sheets/pad, 5 pads/box
5825	Newborn ID Form	8" x 11"	Non-sterile	50 sheets/pad, 5 pads/box
5826	Newborn ID Form	8" x 11"	Non-sterile	50 sheets/pad, 5 pads/box
8403	Patient Name Card	5¾"x¾"	Non-sterile	208 sheets/pack
8423	Patient Name Card	5¾"x1"	Non-sterile	144 sheets/pack



### PlastiBell® Circumcision Device

#### Product Details

- Disposable one piece device plus ligature, clear, plastic construction provides visual access at all times
- Produces a clean line of excision; helps to promote rapid healing
- No post operative care or special dressing required
- Sterile, individually packaged in six sizes and assortment packs
- Eliminates the need for repetitive and costly sterilization required of stainless steel clamps

Part No.	Description	Qty
9210	Hospital Assortment Pack	125 ea/pk
	1.1 cm (25 ea)	
	1.2 cm (35 ea)	
	1.3 cm (45 ea)	
	1.4 cm (10 ea)	
	1.5 cm (5 ea)	
	1.7 cm (5 ea)	
9230	Physicians Office Assortment Pack	25 ea/pk
	1.1 cm (5 ea)	
	1.2 cm (7 ea)	
	1.3 cm (8 ea)	
	1.4 cm (3 ea)	
	1.5 cm (1 ea)	
	1.7 cm (1 ea)	

#### Hospital (Bulk) Single Sizes

Part No.	Description	Qty
9211	1.1 cm	125 ea/bx
9212	1.2 cm	125 ea/bx
9213	1.3 cm	125 ea/bx
9214	1.4 cm	125 ea/bx
9215	1.5 cm	125 ea/bx
9217	1.7 cm	125 ea/bx

#### Physician Office Single Sizes

Part No.	Description	Qty
9231	1.1 cm	25 ea/bx
9232	1.2 cm	25 ea/bx
9233	1.3 cm	25 ea/bx
9234	1.4 cm	25 ea/bx
9235	1.5 cm	25 ea/bx
9237	1.7 cm	25 ea/bx



## Richard-Allan® Surgical Needles


Aspen Surgical carries many styles and sizes of Richard-Allan® Free-Eyed Disposable Needles used for suturing. This includes needles used in an operating room setting. Richard-Allan® Needles are known for their strength and are crafted to help guard against imperfections. This makes Richard-Allan® Surgical Needles a great choice for specialty procedures.


### Surgical Needle Reference Guide

2050 – Part Number

Description – Kidney Blunt ½ Circle Taper

Eye/Point Diagram


Wire Diameter  
Needle Length

205000 0.056"  
2.244"

Please note the  
needles drawn  
are not to scale.

Straight

Eye

Reverse Cutting

Cutting

Taper

Trocarn

Blunt

Regular Surgeon's

2134 VL-M

2146 L

None

2051 H

2147 VL

213413 M

2200 VL

¾ Circle

Eye

Reverse Cutting

Cutting

Taper

Trocarn

Blunt

Regular Surgeon's

2090 VL-H

2441 H

2168 H

2169 H

None

2124 VL

2195 M

2145 M

2153 VL-L

2158 VL

½ Circle

Eye

Reverse Cutting

Cutting

Taper

Trocarn

Blunt

Regular Surgeon's

2091 VL-M

None

2092 VL

2170 H

2050 H

2123 VL

2095 L-H

2179 M-H

2127 L-H

2116 M

2426 M-H

2128 L-H

2125 VL

2701 L

2135 M-H

2140 L-H

2154 VL-L

2150 VL

2157 VL

2151 L

2164 M-VH

2159 L-M

2432 L-M

2167 H

2180 M

2181 L-M

2222 M

2442 VL-L

2500 VL

2706 L

2707 L

Round

None

2126 VL-H

None

None

2223

VVH

Spring

2156

2000 VL

2129 VL

None

None


2002 M

## Scale

VL-Very Light Wire	L-Light Wire	M-Medium Wire	H-Heavy Wire	VH-Very Heavy Wire	VVH-Very, Very Heavy Wire
0.021"	0.027"	0.037"	0.050"	0.064"	0.080"
0.024"	0.030"	0.040"	0.056"		
	0.034"	0.043"			

## Part No. 2000

## French Spring Eye ½ Circle Cutting



Cutting

Spring Eye

200016 0.024"  
1.063"




200018 0.024"  
1.024"



## Part No. 2002

## French Spring Eye Kelly Intestinal ½ Circle Taper



Taper

Double  
Spring Eye

200202 0.037"  
1.496"



200204 0.037"  
1.496"




200203 0.037"  
1.260"



## Part No. 2051

## Kidney Blunt Straight



Blunt


Regular

205100 0.056"  
2.223"



Part No. 2090

Regular Surgeon's  $\frac{3}{8}$  Circle Reverse Cutting



Reverse Regular

209001 0.056"  
3.081"



209002 0.056"  
2.925"



209003 0.055"  
2.677"



209004 0.050"  
2.457"



209006 0.050"  
2.323"



209007 0.043"  
2.087"



209008 0.043"  
2.067"



209009 0.040"  
1.872"



209010 0.040"  
1.677"



209011 0.037"  
1.535"



209012 0.037"  
1.521"



209013 0.037"  
1.378"



209014 0.037"  
1.365"



209015 0.034"  
1.220"



209016 0.034"  
1.209"



209018 0.030"  
1.053"



209019 0.030"  
1.063"



209020 0.027"  
0.827"



209028 0.024"  
0.630"




209030 0.024"  
0.630"





Part No. 2091

Regular Surgeon's ½ Circle Reverse Cutting



Reverse    Regular

209101    0.056"  
3.081"



209102    0.055"  
2.953"



209103    0.056"  
2.652"



209104    0.050"  
2.480"



209106    0.050"  
2.301"



209107    0.043"  
2.087"



209108    0.043"  
2.067"



209109    0.040"  
1.872"



209110    0.040"  
1.677"



209112    0.037"  
1.521"



209113    0.037"  
1.378"



209114    0.037"  
1.365"



209115    0.034"  
1.220"



209116    0.034"  
1.209"



209117    0.030"  
1.063"



209118    0.030"  
1.053"



209119    0.030"  
1.063"



209120    0.027"  
0.819"



209124    0.024"  
0.748"



209128    0.024"  
0.630"




209130    0.024"  
0.630"



Part No. 2092

Adson Dura ½ Circle Taper Point



Taper

Regular

209201 0.024"  
0.741"



209203 0.024"  
0.507"




209202 0.024"  
0.624"



Part No. 2095

Ferguson (Regular Surgeon's Eye) ½ Circle Taper Point



Taper

Regular

209501 0.056"  
3.081"



209503 0.056"  
2.652"



209504 0.050"  
2.457"



209505 0.050"  
2.323"



209506 0.050"  
2.301"



209507 0.043"  
2.087"



209508 0.043"  
2.087"



209509 0.040"  
1.872"



209510 0.039"  
1.693"



209511 0.037"  
1.535"



209512 0.037"  
1.521"



209513 0.037"  
1.378"



209514 0.037"  
1.365"



209515 0.034"  
1.220"



209516 0.030"  
1.209"



209517 0.030"  
1.063"




209518 0.03"  
1.053"



Part No. 2116

Butler's Cervix ½ Circle Taper Point



Taper

Regular

211601 0.037"  
0.984"




211602 0.037"  
1.378"



Part No. 2125

Eye Needle ½ Circle Taper Point



Taper

Regular

212501 0.024"  
1.053"



212503 0.021"  
0.748"




212502 0.024"  
0.858"

212504 0.021"  
0.624"

Part No. 2126


Ferguson (Regular Eye) ½ Circle Taper Point



Taper

Regular

212601 0.056"  
3.801"



212602 0.056"  
2.925"




212603 0.050"  
2.613"



212604 0.050"  
2.457"



212605 0.043"  
2.340"



212609 0.037"  
1.482"



212610 0.034"  
1.287"



212611 0.034"  
1.299"



212612 0.030"  
1.014"



212613 0.027"  
0.858"




212614 0.024"  
0.819"



Part No. 2126

Ferguson (Regular Eye) ½ Circle Taper Point



Taper

Regular

212606    0.043"  
2.205"

212607    0.039"  
1.890"

212608    0.037"  
1.716"



212615    0.024"  
0.669"


212616    0.024"  
0.663"

212622    0.024"  
1.092"



Part No. 2127

Fistula Light ½ Circle Reverse Cutting



Reverse

Regular

212703    0.043"  
2.087"



212704    0.037"  
1.677"

212705    0.034"  
1.482"


212706    0.034"  
1.220"

212708    0.027"  
0.827"



Part No. 2128

Fistula Heavy ½ Circle Reverse Cutting



Reverse Regular

212801 0.056"  
2.340"



212802 0.056"  
2.106"



212803 0.050"  
1.833"



212804 0.043"  
1.638"



212805 0.037"  
1.378"




212806 0.034"  
1.092"



Part No. 2129

French Spring Eye Intestinal ½ Circle Taper Point



Reverse Regular

212901 0.024"  
1.260"



212902 0.024"  
1.102"



212903 0.024"  
0.984"



212904 0.024"  
0.866"



212905 0.024"  
0.748"




212906 0.024"  
0.630"



Part No. 2134

Kieth's Abdominal Straight Cutting




213401	0.040"	
	4.016"	
213402	0.040"	
	3.432"	
213404	0.034"	
	2.964"	
213405	0.034"	
	2.730"	
213406	0.034"	
	2.457"	
213407	0.030"	
	2.223"	
213408	0.030"	
	1.989"	
213409	0.030"	
	1.716"	
213410	0.027"	
	1.482"	
213412	0.024"	
	1.053"	



Part No. 2135

Kelly's Cervix ½ Circle Reverse Cutting



Reverse Regular

213501 0.056"  
2.126"



213503 0.043"  
1.654"



213502 0.050"  
1.693"




213504 0.037"  
1.378"



Part No. 2140

Kelly's Intestinal ½ Circle Taper Point



Taper Regular

214001 0.037"  
1.535"




214004 0.034"  
1.024"




Part No. 2145

King's ¾ Circle Reverse Cutting




Reverse Regular

214506 0.050"  
4.016"



214505 0.043"  
3.705"





214504 0.040"  
3.432"






Part No. 2146

Mayo Intestinal Straight Taper Point

		
Taper	Regular	
214601	0.034" 2.964"	
214602	0.034" 2.756"	
214603	0.030" 2.008"	
214604	0.030" 1.716"	
214605	0.027" 1.209"	


Part No. 2147

Fine Intestinal Straight Taper Point

		
Taper	Regular	
214701	0.024" 2.223"	
214702	0.024" 1.969"	
214703	0.024" 1.732"	
214704	0.024" 1.482"	
214706	0.024" 1.092"	

Part No. 2151

Mayo Intestinal ½ Circle Taper Point



Taper

Regular

215102    0.027"  
1.716"



215105    0.027"  
1.209"




215106    0.027"  
1.053"



Part No. 2153

Lane's Cleft Palate ¾ Circle Reverse Cutting



Reverse

Regular

215303    0.024"  
0.748"



215305    0.024"  
0.507"




215306    0.024"  
0.390"



Part No. 2154

Lane's Cleft Palate ½ Circle Reverse Cutting



Reverse

Regular

215402    0.027"  
0.827"



215405    0.024"  
0.512"



215403    0.024"  
0.748"



215406    0.024"  
0.390"




215404    0.024"  
0.630"



Part No. 2156

Lane's ½ Circle Reverse Cutting French Spring Eye



Reverse Double Spring Eye

215601 0.024"  
1.063"



215603 0.024"  
0.748"



215602 0.024"  
0.827"




215604 0.024"  
0.630"



Part No. 2157

Lane's ½ Circle Reverse Cutting



Reverse Regular

215701 0.024"  
1.053"



215703 0.024"  
0.748"



215702 0.024"  
0.858"




215704 0.024"  
0.630"



Part No. 2159

Ligating ½ Circle Taper Point



Taper Regular

215902 0.034"  
1.482"




215903 0.034"  
1.339"



## Part No. 2164

## Martin's Uterine ½ Circle Reverse Cutting



Reverse      Regular

216401    0.064"  
2.964"

216402    0.062"  
2.480"

216403    0.056"  
1.950"

216404    0.050"  
1.677"


216405    0.050"  
1.482"

216406    0.043"  
1.092"

216407    0.037"  
0.975"


## Part No. 2167

## Mayo Catgut ½ Circle Taper Point



Taper      Regular

216700    0.056"  
2.223"

216701    0.056"  
1.950"

216702    0.050"  
1.716"

216703    0.050"  
1.521"

216735    0.050"  
1.404"

216704    0.050"  
1.248"


216705    0.050"  
1.092"

216706    0.050"  
0.975"

216707    0.050"  
0.897"


Part No. 2168


Mayo Catgut  $\frac{3}{8}$  Circle Taper Point (Heavy)



Taper Regular

216801 0.056"  
1.950"

216802 0.050"  
1.716"




216803 0.050"  
1.521"

216804 0.050"  
1.248"



Part No. 2169

Mayo Catgut  $\frac{3}{8}$  Circle Trocar Point



Trocar Regular

216901 0.056"  
1.950"

216902 0.050"  
1.716"

216903 0.050"  
1.574"



216904 0.050"  
1.260"


216905 0.050"  
1.092"

216906 0.050"  
0.975"



Part No. 2170

Mayo Catgut  $\frac{1}{2}$  Circle Trocar Point



Trocar Regular

217001 0.056"  
1.950"

217002 0.050"  
1.716"

217003 0.050"  
1.521"



217004 0.050"  
1.248"

217005 0.050"  
1.092"


217006 0.050"  
0.975"





## Part No. 2179

## Murphy's Trachelorrhaphy ½ Circle (Light Trocar)



Trocar

Regular

217901 0.050"  
1.260"



217903 0.043"  
0.975"



217902 0.050"  
1.102"




217904 0.037"  
0.741"



## Part No. 2180

## Murphy's Hemostatic ½ Circle Taper Point



Trocar


Regular

218000 0.037"  
0.975"



## Part No. 2181

## Murphy's Intestinal ½ Circle Taper Point



Taper

Regular

218101 0.037"  
1.677"



218103 0.034"  
1.326"



218102 0.037"  
1.535"




218104 0.034"  
1.220"



## Part No. 2195

## Davis Tonsil



Taper

Regular

219500 0.043"  
0.975"



Part No. 2200

Bunnels Straight



220004 0.024"  
2.457"



220003 0.021"  
2.223"



Taper

Regular

220002 0.021"  
1.950"



Part No. 2222

½ Circle Improved Mayo Catgut Taper Point (Light)



222201 0.043"  
1.950"



222203 0.043"  
1.521"



Taper

Regular

222202 0.043"  
1.677"



222204 0.043"  
1.248"



222205 0.043"  
1.092"



Part No. 2223

Gallies Frascia ½ Circle Taper



222301 0.080"  
2.496"



222302 0.080"  
1.989"




Taper

Round

Part No. 2426

Mayo Catgut ½ Circle Trocar



Trocar Regular

242602 0.050"  
1.950"



242605 0.043"  
1.365"



242604 0.050"  
1.535"




242606 0.040"  
1.092"



Part No. 2432

Fistula ½ Circle Reverse Cutting



Reverse Regular

243204 0.043"  
1.654"



243207 0.037"  
0.984"



243205 0.043"  
1.365"



243208 0.034"  
0.866"



243206 0.040"  
1.248"




243209 0.030"  
0.741"



Part No. 2441







King Cutting ¾ Circle









Cutting Regular

244101 0.050"  
3.432"




Part No. 2442		Ferguson's ½ Circle Taper Point					
		<u>244216</u>			<u>244220</u>		
		0.027" 1.063"			0.024" 0.630"		
Taper	Regular	<u>244218</u>			<u>244222</u>		
		0.024" 0.748"			0.024" 0.512"		

Part No. 2706		Ferguson ½ Circle Taper	
		<div>270616</div> <div>0.030"</div> <div>1.063"</div>	
Taper	Regular		


Part No. 2707		Ferguson ½ Circle Taper	
		<div>270718</div> <div>0.027"</div> <div>0.780"</div>	
Taper	Regular		

### S-Vac® Fluted Channel Drains


- Designed to provide advanced drainage and increased patient comfort
- Fluted channel design helps minimize the chance of tissue in-growth creating the blockages that can cause premature drain removal when using standard drains
- 10mm Flat Fluted Drain provides 6-channels for optimal drainage, compared to the industry standard 4-channel design
- Soft, flat drain hubs offers enhanced patient comfort
- Exceptional flow rates provides improved drainage, helping reduce the chance of seroma formation
- Available with or without sharp stainless steel trocars
- Radiopaque
- Sterile, intended for single-use, and not made with natural rubber latex



6 Channel Fluted Flat Drain



4 Channel Fluted Flat Drain



4 Channel Fluted Round Drain

### Fluted Flat Drain

Part No.	Description	Qty/Case
322211	7mm Flat, Full Fluted, 4 Channel Drain	10
322214	10mm Flat, Full Fluted, 6 Channel Drain	10

### Fluted Flat Drain with Trocar

322212	7mm Flat, Full Fluted, 4 Channel Drain, Trocar	10
322215	10mm Flat, Full Fluted, 6 Channel Drain, Trocar	10

### Fluted Round Drain

332186	10 FR Round, Full Fluted, 4 Channel Drain	10
332188	15 FR Round, Full Fluted, 4 Channel Drain	10
332190	19 FR Round, Full Fluted, 4 Channel Drain	10

### Fluted Flat Drain with Trocar

322212	7mm Flat, Full Fluted, 4 Channel Drain, Trocar	10
322215	10mm Flat, Full Fluted, 6 Channel Drain, Trocar	10

### Fluted Round Drain with Trocar

332187	10 FR Round, Full Fluted, 4 Channel Drain, Trocar	10
332189	15 FR Round, Full Fluted, 4 Channel Drain, Trocar	10
332191	19 FR Round, Full Fluted, 4 Channel Drain, Trocar	10
332221	10 FR Round, 3/4 Fluted, 4 Channel Drain, Trocar	10
332223	15 FR Round, 3/4 Fluted, 4 Channel Drain, Trocar	10
332225	19 FR Round, 3/4 Fluted, 4 Channel Drain, Trocar	10

### Xtra Fluted Round Drain


322226	10 FR Round Xtra Fluted Drain	10
322227	10 FR Round Xtra Fluted Drain, With Trocar	10
322228	15 FR Round Xtra Fluted Drain	10
322229	15 FR Round Xtra Fluted Drain, With Trocar	10
322230	19 FR Round Xtra Fluted Drain	10
322231	19 FR Round Xtra Fluted Drain, With Trocar	10

### Drain System-Fluted with 100ml Bulb Evacuator


332263	10 FR Drain System Round Full Fluted, with Trocar and 100ml S-Vac® Bulb Evacuator	10
332264	15 FR Drain System Round Full Fluted, with Trocar and 100ml S-Vac® Bulb Evacuator	10

### Perforated Drains

- Discrete drain with our unique channel designed to help prevent collapsing or kinking
- Discrete drains contain internal ribs to help prevent collapse and occlusion of the drain
- 10mm discrete drain offers minimal transition from the drain hub to the drain tubing for ease of removal and added patient comfort
- Hubless styles also available
- Soft, silicone material
- Available with or without extra sharp stainless steel trocars
- Radiopaque
- Sterile, intended for single-use, and not made with natural rubber latex



Flat Perforated Drain



Round Perforated Drain

### Discrete Flat Perforated Drain

Part No.	Description	Qty/Case
370004	Discrete Drain, Flat 7mm x 20cm, $\frac{3}{4}$ Perf	10
370005	Discrete Drain, Flat 10mm x 20cm, $\frac{3}{4}$ Perf	10
370006	Discrete Drain, Flat 7mm x 20cm, Full Perf	10
370007	Discrete Drain, Flat 10mm x 20cm, Full Perf	10

### Discrete Flat Perforated Drain With Trocar

370000	Discrete Drain, Flat 7mm x 20cm, $\frac{3}{4}$ Perf, Trocar	10
370001	Discrete Drain, Flat 10mm x 20cm, $\frac{3}{4}$ Perf, Trocar	10
370002	Discrete Drain, Flat 7mm x 20cm, Full Perf, Trocar	10
370003	Discrete Drain, Flat 10mm x 20cm, Full Perf, Trocar	10

### Hubless Flat Perforated Drain With Trocar

360000	Hubless Drain, Flat 7mm x 20cm, $\frac{3}{4}$ Perf, Trocar	10
360002	Hubless Drain, Flat 7mm x 20cm, Full Perf, Trocar	10
360003	Hubless Drain, Flat 10mm x 20cm, Full Perf, Trocar	10

### Round Perforated Drain

370020	Discrete Drain, Round 7 FR, End-Perf	10
370021	Discrete Drain, Round 10 FR, End-Perf	10
370022	Discrete Drain, Round 15 FR, End-Perf	10
370023	Discrete Drain, Round 19 FR, End-Perf	10

### Round Perforated Drain With Trocar

370016	Discrete Drain, Round 7 FR, End-Perf, Trocar	10
370017	Discrete Drain, Round 10 FR, End-Perf, Trocar	10
370018	Discrete Drain, Round 15 FR, End-Perf, Trocar	10
370019	Discrete Drain, Round 19 FR, End-Perf, Trocar	10



### Drain Systems with S-Vac® 100ml Bulb Evacuator

Part No.	Description	Qty/Case
340001	10mm Drain System, Flat ¾ Perforated, with Trocar and 100ml S-Vac® Bulb Evacuator	10
340002	7mm Drain System, Flat Full Perforated, with Trocar and 100ml S-Vac® Bulb Evacuator	10
340003	10mm Drain System, Flat Full Perforated, with Trocar and 100ml S-Vac® Bulb Evacuator	10
340004	7mm Drain System, Flat ¾ Perforated, with 100ml S-Vac® Bulb Evacuator	10
340005	10mm Drain System, Flat ¾ Perforated, with 100ml S-Vac® Bulb Evacuator	10
340006	7mm Drain System, Flat Full Perforated, with 100ml S-Vac® Bulb Evacuator	10
340007	10mm Drain System, Flat Full Perforated, with 100ml S-Vac® Bulb Evacuator	10
340018	15 FR Drain System, Round End-Perforated, with Trocar and 100 ml S-Vac® Bulb Evacuator	10



### Surgidyne® S-Vac® Silicone Bulb Evacuator

- Can be used with all Surgidyne drains to help promote fluid drainage
- Preferred by some physician specialties
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Case
322160	S-VAC® Silicone 100ml Bulb Evacuator	10



#### Sterion® Filament Reinforced Skin Closures

- Reinforcing filaments offer added strength for large skin surface closures
- Filament Reinforced Skin Closure style is accepted as the standard by most healthcare professionals
- Breathable material designed to allow moisture transmission, to help reduce the chance of skin maceration
- Less painful to apply than suturing or stapling providing more comfort to the patient
- No adhesive enhancers required
- Strips are curl-resistant
- Sterile, intended for single-use, and not made with natural rubber latex




Part No.	Description	mm	inches	Strips/Pouch	Pouches/Box
751040PBX	Filament Reinforced Skin Closures	3 x 75	1/8" x 3"	5	50
751041PBX	Filament Reinforced Skin Closures	6 x 75	1/4" x 3"	3	50
751042PBX	Filament Reinforced Skin Closures	6 x 38	1/4" x 1 1/2"	6	50
751046PBX	Filament Reinforced Skin Closures	6 x 100	1/4" x 4"	10	50
751047PBX	Filament Reinforced Skin Closures	12 x 100	1/2" x 4"	6	50
751048PBX	Filament Reinforced Skin Closures	25.4 x 127	1" x 5"	4	25



#### Sterion® Flexible Skin Closures

- Non-woven flexible backing helps reduce the chance of skin shear that can occur when edema or flexing puts excessive pressure on the skin closure ends
- Less painful to apply than suturing or stapling, providing more comfort to the patient
- No adhesive enhancers required
- Sterile, intended for single-use, and not made with natural rubber latex








Part No.	Description	mm	inches	Strips/Pouch	Pouches/Box
753040PBX	Flexible Skin Closures	3 x 75	1/8" x 3"	5	50
753041PBX	Flexible Skin Closures	6 x 75	1/4" x 3"	3	50
753046PBX	Flexible Skin Closures	6 x 100	1/4" x 4"	10	50
753047PBX	Flexible Skin Closures	12 x 100	1/2" x 4"	6	50
753048PBX	Flexible Skin Closures	25 x 125	1" x 5"	4	25



## ASTM F2100 Face Standards

Three levels of protection are derived from standards that measure five performance characteristics:

-  **Fluid Penetration Resistance** *How well the mask protects against blood or other potentially infectious fluids. Measures resistance at one of three velocities corresponding to human blood pressures: 80, 120, and 160 mmHg.*
-  **Bacterial Filtration Efficiency (BFE)** *How well the mask protects against harmful bacteria. Effectiveness is measured using bacteria with a mean particle size of 3.0 µm (microns) ± 0.3 µm.*
-  **Sub-micron Particulate Filtration Efficiency (PFE)** *How well the mask protects against harmful aerosol particles. Effectiveness is measured using particulates down to the size of 0.1 µm (microns).*
-  **Differential Pressure (Delta P, or ΔP)** *How breathable is the mask. Measures resistance to air flow— the lower the number, the higher the air flow and, therefore, breathability.*
-  **Flame Resistance** *Measures the potential burn rate of the mask. Materials must meet Class 1 requirements indicating no unusual burning characteristics.*

## Differences between levels and how that affects protection & comfort



- Designed for lower risk fluid environments
- Provides good bacterial and particle filtration
- Greater comfort and breathability

FUNCTION	ASTM F2100 STANDARD	PERFORMANCE		
		GOOD	BETTER	BEST
Fluid Resistance	80 mmHg			
Bacterial Filtration	≥95%			
Particle Filtration	>95%			
Breathability	<4			
Flame Resistance	Class 1			



- Designed for moderate to high risk fluid environments
- Provides excellent bacterial and particle filtration
- Good balance of performance and comfort

FUNCTION	ASTM F2100 STANDARD	PERFORMANCE		
		GOOD	BETTER	BEST
Fluid Resistance	120 mmHg			
Bacterial Filtration	≥98%			
Particle Filtration	>98%			
Breathability	<5			
Flame Resistance	Class 1			



- Designed for high risk fluid environments
- Provides excellent bacterial and particle filtration
- Maintains breathability while providing the highest level of fluid protection

FUNCTION	ASTM F2100 STANDARD	PERFORMANCE		
		GOOD	BETTER	BEST
Fluid Resistance	160 mmHg			
Bacterial Filtration	≥98%			
Particle Filtration	>98%			
Breathability	<5			
Flame Resistance	Class 1			

## Precept® Comfort Masks

Comfort Isolation, Procedure and Surgical masks provide wear with soft comfortable protection as well as excellent bacterial (BFE) and particulate (PFE) filtration efficiencies. Integrated layers work together to provide masks that are soft, breathable, and cool. Comfort, Isolation, and Procedure masks have earloops that are non-irritating and exhibit excellent elasticity and durability. Comfort Surgical masks, available in pouch and horizontal tie styles, incorporate inner facing materials with excellent fiber pin-down, which helps prevent facial irritation. Child friendly, child sized earloop masks, and adult procedure mask have colorful fun prints to encourage compliance. Child size improves fit and protection for children.



### Precept® Isolation Mask Product Details

- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15100	Isolation Mask	Yellow	50/Box	500/Case
15101	Isolation Mask	Blue	50/Box	500/Case



### Precept® KittyKare™ Procedure Mask Product Details

- Child-friendly printed outer facing
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15130	Procedure Mask	KittyKare™	50/Box	500/Case



### Precept® Children's Mask<sup>①</sup> Product Details

- Child-sized to improve fit for children
- Child-friendly printed outer facing
- Soft comfortable earloops provide easier donning
- 75/Box, 10/Boxes, 750/Case
- Not made with natural rubber latex and non-sterile

<sup>①</sup>The Centers for Disease Control and Prevention (CDC) recommends the use of face masks in healthcare settings for all patients exhibiting symptoms of a respiratory infection. Precept's child sized face masks (15150) are designed specifically to be used by children ages 3-12, to decrease transmission of respiratory infections.

Part No.	Style	Color	Packaging	Qty
15150	Children's Mask	Smiley Face	75/Box	750/Case





## Precept® Lite Pouch Surgical Mask

### Product Details

- Pouch style mask provides off the face design giving more room to breathe
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15220	Surgical Mask	Blue	50/Box	300/Case



## Precept® Lite & Cool Horizontal Tie Surgical Mask

### Product Details

- Wider design provides more facial coverage than standard sized masks
- Soft longer horizontal ties promote a secure and comfortable fit
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- 50/Box, 12/Boxes, 600/Case
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15250	Surgical Mask	Blue	50/Box	600/Case

## Level 1 Fluid Resistant Face Masks

Precept ASTM F2100 Level 1 face masks are designed for lower risk fluid environments, providing good protection with filtration efficiencies (BFE and PFE) exceeding >97%, and all masks meet Class 1 flame resistance requirements. Precept Level 1 face masks are available in a variety of styles including Surgical and Procedure with and without extended shields. All masks are sonic welded and have and enclosed malleable nosepieces to ensure a snug, yet comfortable fit.

### Complies with ASTM F2100

Fluid Resistance (mmHg): 80  
BFE (%) @ 3.0  $\mu\text{m}$ :  $\geq 95\%$   
PFE (%) @ 0.1  $\mu\text{m}$ :  $\geq 95\%$   
 $\Delta P$  (mm H<sub>2</sub>O/cm<sup>2</sup>): < 4  
Flame Spread: Class 1



### Precept® Fluid Resistant Procedure Mask Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15110	Procedure Mask	Yellow	50/Box	500/Case
15111	Procedure Mask	Blue	50/Box	500/Case
15118	Procedure Mask	Blue	50/Box	500/Case Bulk
15120	Procedure Mask	Blue	50/Box	500/Case
15121	Procedure Mask	Yellow	50/Box	500/Case



### Precept® Fluid Resistant Procedure Mask with Shield Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15115	Procedure Mask w/Shield	Blue	25/Box	100/Case

## Complies with ASTM F2100

Fluid Resistance (mmHg): 80  
 BFE (%) @ 3.0  $\mu\text{m}$ :  $\geq 95\%$   
 PFE (%) @ 0.1  $\mu\text{m}$ :  $\geq 95\%$   
 $\Delta P$  (mm H<sub>2</sub>O/cm<sup>2</sup>): < 4  
 Flame Spread: Class 1



Level  
1

### Precept® Fluid Resistant Procedure Face Masks with Anti-Glare Shield

#### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Coated film shield with proprietary anti-glare, anti-fog treatment provides enhanced clarity
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15116	Procedure Mask w/Shield	Blue	25/Box	100/Case



Level  
1

### Precept® Lite & Cool Surgical Mask

#### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15200	Surgical Mask	Blue	50/Box	300/Case



Level  
1

### Precept® Standard Surgical Mask

#### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15201	Surgical Mask	Blue	50/Box	300/Case

## Complies with ASTM F2100

Fluid Resistance (mmHg): 80  
 BFE (%) @ 3.0  $\mu\text{m}$ :  $\geq 95\%$   
 PFE (%) @ 0.1  $\mu\text{m}$ :  $\geq 95\%$   
 $\Delta P$  (mm H<sub>2</sub>O/cm<sup>2</sup>):  $< 4$   
 Flame Spread: Class 1



Level  
1

### Precept® Comfort-Cool™ Surgical Mask Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft sewn ties and bindings promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3010	Surgical Mask	Blue	50/Box	300/Case



Level  
1

### Precept® Comfort-Plus™ Surgical Mask Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft sewn ties and bindings promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3110	Surgical Mask	Blue	50/Box	300/Case



Level  
1

### Precept® Lite & Cool Surgical Mask with Anti-Glare Shield Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Coated film shield with proprietary anti-glare, anti-fog treatment provides enhanced clarity
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15206	Surgical Mask	Blue	25/Box	100/Case

## Complies with ASTM F2100

Fluid Resistance (mmHg): 80  
 BFE (%) @ 3.0  $\mu\text{m}$ :  $\geq 95\%$   
 PFE (%) @ 0.1  $\mu\text{m}$ :  $\geq 95\%$   
 $\Delta P$  (mm H<sub>2</sub>O/cm<sup>2</sup>):  $< 4$   
 Flame Spread: Class 1



Level  
1

### Precept® Foam Anti-Fog Surgical Mask

#### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Fog-Shield® anti-fog foam inhibits eyewear fogging to help improve vision
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15210	Surgical Specialty Mask	Blue	50/Box	300/Case



Level  
1

### Precept® Film Anti-Fog Surgical Mask

#### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Anti-fog film inhibits eyewear fogging to help improve vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15211	Surgical Specialty Mask	Green	50/Box	300/Case



Level  
1

### Precept® Tape Anti-Fog Surgical Mask

#### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Anti-fog tape inhibits eyewear fogging to help improve vision
- Integrated layers work together to maximize breathability
- Soft sewn ties and bindings promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15212	Surgical Specialty Mask	Green	50/Box	300/Case

## Complies with ASTM F2100

Fluid Resistance (mmHg): 80

BFE (%) @ 3.0  $\mu\text{m}$ :  $\geq 95\%$

PFE (%) @ 0.1  $\mu\text{m}$ :  $\geq 95\%$

$\Delta P$  (mm H<sub>2</sub>O/cm<sup>2</sup>):  $< 4$

Flame Spread: Class 1



Level  
1

### Precept® Sensitive Skin Surgical Mask

#### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Color-free face mask with specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15215	Surgical Specialty Mask	White	50/Box	300/Case



Level  
1

### Precept® Sensitive Skin Foam Anti-Fog Surgical Mask

#### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Color-free face mask with specially designed inner facing improves comfort
- Anti-fog foam inhibits eyewear fogging to help improve vision
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15216	Surgical Specialty Mask	White	50/Box	300/Case



Level  
1

### Precept® Sensitive Skin Plus Surgical Mask

#### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Color-free face mask with specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft sewn ties and bindings promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3130	Surgical Specialty Mask	White	50/Box	300/Case



**Complies with ASTM F2100**  
 Fluid Resistance (mmHg): 80  
 BFE (%) @ 3.0 µm: ≥ 95%  
 PFE (%) @ 0.1 µm: ≥ 95%  
 ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 4  
 Flame Spread: Class 1



Level  
**1**

## Precept® Comfort-Plus™ Surgical Mask with Stretch Knit Ties

### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft sewn stretch knit ties provide a secure and comfortable fit
- 50/Box, 5/Boxes, 250/Case
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3112	Surgical Specialty Mask	Blue	50/Box	250/Case



Level  
**1**

## Precept® Comfort-Plus™ Surgical Mask

### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft sewn ties and bindings promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3120	Surgical Specialty Mask	Green	50/Box	300/Case



Level  
**1**

## Precept® Comfort-Plus™ Surgical Mask with Stretch Knit Ties

### Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft sewn stretch knit ties provide a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3122	Surgical Specialty Mask	Green	50/Box	250/Case

**Complies with ASTM F2100**  
 Fluid Resistance (mmHg): 80  
 BFE (%) @ 3.0  $\mu\text{m}$ :  $\geq 95\%$   
 PFE (%) @ 0.1  $\mu\text{m}$ :  $\geq 95\%$   
 $\Delta P$  (mm H<sub>2</sub>O/cm<sup>2</sup>):  $< 4$   
 Flame Spread: Class 1



Level  
1

## Precept® Foam Fog-Shield® Surgical Mask Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Fog-Shield® anti-fog foam inhibits eyewear fogging to help improve vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Strong polyester ties provide a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3320	Surgical Specialty Mask	Green	25/Box	150/Case



Level  
1

## Precept® Tape Fog-Shield® Surgical Mask Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Fog-Shield® anti-fog tape inhibits eyewear fogging to help improve vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft sewn ties and bindings promote secure comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3321	Surgical Specialty Mask	Green	50/Box	250/Case



Level  
1

## Precept® Tape Fog-Shield® Surgical Mask Product Details

- Complies with ASTM F2100 Level 1 Requirements
- Fog-Shield® anti-fog tape inhibits eyewear fogging to help improve vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Strong polyester ties provide a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3322	Surgical Specialty Mask	Green	50/Box	250/Case

## Level 2 Fluid Resistant Face Masks

Precept ASTM F2100 Level 2 face masks are designed for moderate to high risk fluid environments, providing good protection with filtration efficiencies (BFE and PFE) exceeding >98%, and all masks meet Class 1 flame resistance requirements. Precept Level 2 face masks are available in a variety of styles including Surgical and Procedure with and without extended shields. All masks are sonic welded and have enclosed malleable nosepieces to ensure a snug, yet comfortable fit.

### Complies with ASTM F2100

Fluid Resistance (mmHg): 80  
BFE (%) @ 3.0  $\mu\text{m}$ :  $\geq 98\%$   
PFE (%) @ 0.1  $\mu\text{m}$ :  $\geq 98\%$   
 $\Delta P$  (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5  
Flame Spread: Class 1



Level  
2

### Precept® FluidGuard® 120 Procedure Mask Product Details

- Complies with ASTM F2100 Level 2 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15900	Procedure Mask	Blue	50/Box	500/Case



Level  
2

### Precept® FluidGuard® 120 Anti-Fog Procedure Mask with Extended Shield Product Details

- Complies with ASTM F2100 Level 2 Requirements
- Coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Anti-fog foam inhibits fogging to help improve vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15910	Procedure Mask	Blue	25/Box	100/Case



Level  
2

### Precept® FluidGuard® 120 Horizontal Tie Surgical Mask Product Details

- Complies with ASTM F2100 Level 2 Requirements
- Wider design provides more facial coverage than standard sized masks
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer horizontal ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15251	Surgical Mask	Blue with Pattern	50/Box	600/Case

**Complies with ASTM F2100**  
 Fluid Resistance (mmHg): 80  
 BFE (%) @ 3.0 µm: ≥ 98%  
 PFE (%) @ 0.1 µm: ≥ 98%  
 ΔP (mm H2O/cm²): < 5  
 Flame Spread: Class 1



Level  
**2**

**Precept® FluidGard® 120 FluidPouch™ Surgical Mask**  
**Product Details**

- Complies with ASTM F2100 Level 2 Requirements
- Pouch style mask provides off the face design giving more room to breathe
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15820	Surgical Mask	Blue with Pattern	50/Box	300/Case



Level  
**2**

**Precept® FluidGard® 120 Anti-Fog Surgical Mask**  
**Product Details**

- Complies with ASTM F2100 Level 2 Requirements
- Anti-fog foam inhibits eyewear fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15920	Surgical Mask	Blue	50/Box	300/Case

**Complies with ASTM F2100**

Fluid Resistance (mmHg): 80  
 BFE (%) @ 3.0 µm: ≥ 98%  
 PFE (%) @ 0.1 µm: ≥ 98%  
 ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5  
 Flame Spread: Class 1



Level  
2

**Precept® FluidGuard® 120 Horizontal Tie Surgical Mask with Anti-Glare Shield**
**Product Details**

- Complies with ASTM F2100 Level 2 Requirements
- Coated film shield with proprietary anti-fog, anti-glare treatment provides enhanced clarity
- Wider design provides more facial coverage than standard sized masks
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer horizontal ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15255	Surgical Mask	Blue with Pattern	25/Box	100/Case



Level  
2

**Precept® FluidGuard® 120 FluidPouch™ Surgical Mask with Anti-Glare Shield**
**Product Details**

- Complies with ASTM F2100 Level 2 Requirements
- Coated film shield with proprietary anti-glare, anti-fog treatment provides enhanced clarity
- Pouch style mask provides off the face design giving more room to breathe
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15830	Surgical Mask	Blue with Pattern	25/Box	100/Case

**Complies with ASTM F2100**

Fluid Resistance (mmHg): 80  
 BFE (%) @ 3.0 µm: ≥ 98%  
 PFE (%) @ 0.1 µm: ≥ 98%  
 ΔP (mm H2O/cm²): < 5  
 Flame Spread: Class 1


Level  
2

**Precept® FluidGuard® 120 Anti-Fog Surgical Mask with Extended Shield**  
**Product Details**

- Complies with ASTM F2100 Level 2 Requirements
- Coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Anti-fog foam inhibits eyewear fogging to help improve vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15930	Surgical Mask	Blue	25/Box	100/Case

**Level 3 Fluid Resistant Face Masks**

ASTM F2100 Level 3 face masks are designed for high risk fluid environments, providing excellent fluid protection and filtration efficiencies (BFE and PFE) while maintaining excellent breathability. Aspen Level 3 face masks are available in a variety of styles including Surgical and Procedure as well as with and without extended shields. All masks are sonic welded and have an enclosed malleable nose piece to ensure a snug, yet comfortable fit.

**Complies with ASTM F2100**

Fluid Resistance (mmHg): 160  
 BFE (%) @ 3.0 µm: ≥ 98%  
 PFE (%) @ 0.1 µm: ≥ 98%  
 ΔP (mm H2O/cm²): < 5  
 Flame Spread: Class 1


Level  
3

**Precept® FluidGuard® 160 Procedure Mask**  
**Product Details**

- Complies with ASTM F2100 Level 3 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
14400	Procedure Mask	Blue	50/Box	500/Case
15300	Procedure Mask	Blue Diamond	50/Box	500/Case



**Complies with ASTM F2100**

Fluid Resistance (mmHg): 160

BFE (%) @ 3.0 µm: ≥ 98%

PFE (%) @ 0.1 µm: ≥ 98%

ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5

Flame Spread: Class 1


Level  
3

**Precept® FluidGard® 160 Anti-Fog Procedure Mask**  
**Product Details**

- Complies with ASTM F2100 Level 3 Requirements
- Anti-fog foam inhibits eyewear fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
14401	Procedure Mask	Blue	50/Box	500/Case
15301	Procedure Mask	Blue Diamond	50/Box	500/Case


Level  
3

**Precept® FluidGard® 160 Anti-Fog Procedure Mask with Extended Shield**  
**Product Details**

- Complies with ASTM F2100 Level 3 Requirements
- Coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Anti-fog foam inhibits fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
14410	Procedure Mask	Blue	25/Box	100/Case
15310	Procedure Mask	Blue Diamond	25/Box	100/Case



**Complies with ASTM F2100**  
 Fluid Resistance (mmHg): 160  
 BFE (%) @ 3.0  $\mu\text{m}$ :  $\geq 98\%$   
 PFE (%) @ 0.1  $\mu\text{m}$ :  $\geq 98\%$   
 $\Delta P$  (mm H<sub>2</sub>O/cm<sup>2</sup>):  $< 5$   
 Flame Spread: Class 1



Level  
**3**

## Precept® FluidGard® 160 Anti-Fog Procedure Mask with Anti-Glare Shield

### Product Details

- Complies with ASTM F2100 Level 3 Requirements
- Coated film shield with proprietary anti-glare, anti-fog treatment provides enhanced clarity
- Anti-fog foam inhibits eyewear fogging to help improve vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15311	Procedure Mask	Blue Diamond	25/Box	100/Case



Level  
**3**

## Precept® DualGard™ Foam Shield™ Anti-Fog Procedure Mask

### Product Details

- Complies with ASTM F2100 Level 3 Requirements
- Foam Shield™ provides added comfort with a foam strip that keeps shield off forehead
- Snap-on anti-fog Foam Shield™ can be removed from mask if needed
- Anti-fog foam inhibits fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3347	Procedure Mask	Blue	25/Box	100/Case

**Complies with ASTM F2100**

Fluid Resistance (mmHg): 160

BFE (%) @ 3.0 µm: ≥ 98%

PFE (%) @ 0.1 µm: ≥ 98%

ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5

Flame Spread: Class 1


Level  
3

**Precept® FluidGard® 160 Anti-Fog Surgical Mask**  
**Product Details**

- Complies with ASTM F2100 Level 3 Requirements
- Anti-fog foam inhibits fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- 25/Box, 4/Boxes, 100/Case
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
14420	Procedure Mask	Blue	50/Box	300/Case
15320	Procedure Mask	Blue Diamond	50/Box	300/Case


Level  
3

**Precept® DualGard™ High Filtration Surgical Mask**  
**Product Details**

- Complies with ASTM F2100 Level 3 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft sewn stretch knit ties provide a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3312	Surgical Mask	Blue Diamond	25/Box	100/Case

**Complies with ASTM F2100**

Fluid Resistance (mmHg): 160

BFE (%) @ 3.0 µm: ≥ 98%

PFE (%) @ 0.1 µm: ≥ 98%

ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5

Flame Spread: Class 1


Level  
3

**Precept® FluidGard® 160 Anti-Fog Surgical Mask with Extended Shield**
**Product Details**

- Complies with ASTM F2100 Level 3 Requirements
- Coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Anti-fog foam inhibits eyewear fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
14430	Surgical Mask	Blue	25/Box	100/Case
15330	Surgical Mask	Blue Diamond	25/Box	100/Case


Level  
3

**Precept® FluidGard® 160 Anti-Fog Surgical Mask with Anti-Glare Shield**
**Product Details**

- Complies with ASTM F2100 Level 3 Requirements
- Coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Anti-fog foam inhibits fogging to help improve vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
15331	Surgical Mask	Blue Diamond	25/Box	100/Case

**Complies with ASTM F2100**

Fluid Resistance (mmHg): 160

BFE (%) @ 3.0 µm: ≥ 98%

PFE (%) @ 0.1 µm: ≥ 98%

ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5

Flame Spread: Class 1


Level  
3

**Precept® DualGard™ Surgical Mask with Extended Shield**
**Product Details**

- Complies with ASTM F2100 Level 3 Requirements
- Snap-on coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Snap-on extended shield can be removed from mask if needed
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft sewn ties and bindings promote a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3344	Surgical Mask	Blue Diamond	25/Box	100/Case


Level  
3

**Precept® DualGard™ Foam Shield™ Surgical Mask with Stretch Knit Ties**
**Product Details**

- Complies with ASTM F2100 Level 3 Requirements
- Foam Shield™ coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft sewn stretch knit ties provide a secure and comfortable fit
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3348	Surgical Mask	Blue Diamond	25/Box	100/Case

**Complies with ASTM F2100**  
 Fluid Resistance (mmHg): 160  
 BFE (%) @ 3.0 µm: ≥ 98%  
 PFE (%) @ 0.1 µm: ≥ 98%  
 ΔP (mm H2O/cm²): < 5  
 Flame Spread: Class 1



**Precept® High Filtration Surgical Mask with Stretch  
 Knit Ties**  
**Product Details**

- Complies with ASTM F2100 Level 3 Requirements
- > 99% efficiency filtering laser plume particles of 0.1 microns or larger
- Soft sewn stretch knit ties provide a secure and comfortable fit
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Not made with natural rubber latex and non-sterile

Part No.	Style	Color	Packaging	Qty
65-3310	Surgical Mask	White	50/Box	250/Case

Full length and ¾ length face shields include a soft foam strip that cushions against the forehead, and the soft elastic headband helps provide a secure, comfortable fit. Each shield is treated with a specialty anti-fog coating to ensure optical clarity. Eyewear is lightweight and comfortable providing wrap around eye protection. Eye shield lenses are optically clear and anti-fog to provide enhanced clarity. The reusable frames are available in five attractive colors: Black, Blue, Green, Pink, and Purple. Users have the cost effective option of replacing frames and/or lenses. Face shields and eye shields are designed to be worn in conjunction with Aspen's various styles of face masks for added facial protection.



## Precept® Face Shields

### Product Details

- Two sizes of face shields: ¾ and full length
- Coated film shield with proprietary anti-fog treatment to provide enhanced clarity
- Soft, comfortable foam strip
- Stretchable head band provides a secure fit
- Convenient individual packaging
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Packaging	Qty
15990	7½"	Clear	Individual	50/Case
15995	9½"	Clear	Individual	50/Case



## Precept® Eye Shields

### Product Details

- Lightweight and comfortable providing wrap around eye protection
- Lenses are optically clear and anti-fog to provide enhanced clarity
- Reusable frames are available in five attractive colors
- Cost effective option of replaceable frames and/or lenses
- Not made with natural rubber latex and non-sterile

Part No.	Description	Color	Packaging	Qty
1599A	Eye shield assembled glasses	Variety pack (black, blue, green, pink, & purple)	10/Bag	50/Case
1599F	Eye shield reusable frames	Variety pack (black, blue, green, pink, & purple)	10/Bag	100/Case
1599L	Eye shield replacement lenses	Clear	25/Envelope	250/Case

## Isolation Gown Protection and Recommended Levels of Use

Protection is the most important feature to consider in making choices regarding the appropriate isolation gown for your needs. Aspen offers a variety of gowns and/or materials to meet the following AAMI<sup>1</sup> PB70<sup>2</sup> levels of fluid protection. Our intent is to allow each customer to identify and choose the product(s) that best meet their specific needs. The chart below provides basic information regarding the AAMI fluid resistance levels:



Designed for tasks or use in areas such as Med/Surg Units, Laundry, Housekeeping, Patient Transport, Basic Care

GOWN MATERIAL	TEST METHOD	AAMI LEVEL REQUIREMENT
Light-weight Multi-layer SMS	Water resistance: Impact Penetration AATCC 42	Spray Impact $\leq 4.5g$



Designed for tasks or use in areas such as ICU Med/Surg Units, Dialysis, Nursery, Pathology, Laboratories, Blood Draw, Suturing, IV Line Insertion

GOWN MATERIAL	TEST METHOD	AAMI LEVEL REQUIREMENT
Light-weight Multi-layer SMS	Water resistance: Impact Penetration AATCC 42	Spray Impact $\leq 1.0g$
	Water Resistance Hydrostatic Pressure AATCC 127	Hydrostatic Pressure $\geq 20cm$



Designed for tasks or use in areas such as ER, Trauma, Burn Unit, Critical Care Units, Arterial Blood Draw, IV Line Insertion

GOWN MATERIAL	TEST METHOD	AAMI LEVEL REQUIREMENT
Light-weight Multi-layer SMS	Water resistance: Impact Penetration AATCC 42	Spray Impact $\leq 1.0g$
	Water Resistance Hydrostatic Pressure AATCC 127	Hydrostatic Pressure $\geq 50cm$



Barrier Protection AAMI PB70:2012		
Test Method	Test Result	Anticipated Risk of Exposure
AATCC 42 Spray Impact	≤ 1.0g	Low
AATCC 127 Hydrostatic Pressure	≥ 20cm	Low



Level  
2

**Precept® Isolation Gown Full Coverage**  
**Product Details**

- Polycoated and multi-layer SMS material<sup>①</sup> delivers AAMI Level 2 protection and comfort
- Elastic wrists ensure proper glove donning
- Full coverage with overlapping back panel and secure neck ties
- Sonic-welded seams
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
51177	Large	45"	55"	Yellow	Secure Neck Tie with Elastic Wrist	SMS	10/Bag	100/Case
51178	X-Large	48"	60"	Yellow	Secure Neck Tie with Elastic Wrist	SMS	10/Bag	100/Case

①AATCC 42 – Water Resistance: Spray Impact Penetration Test / AATCC 127 – Water Resistance: Hydrostatic Pressure Test For more information on these test methods, please visit [aatcc.org](http://aatcc.org) and [aami.org](http://aami.org).

## Barrier Protection AAMI PB70:2012

Test Method	Test Result	Anticipated Risk of Exposure
AATCC 42 Spray Impact	≤ 1.0g	Low
AATCC 127 Hydrostatic Pressure	≥ 20cm	Low



Level  
2

### Precept® Isolation Gown Over-the-Head Thumb Loop Full Coverage

#### Product Details

- Multi-layer SMS or polycoated material<sup>①</sup> delivers AAMI Level 2 protection and comfort
- Thumb loops with elastic wrist ensure proper glove donning
- Full coverage with overlapping back panels
- Secure neck ties
- Sonic-welded seams
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
51193FC	X-Large	48"	61"	Yellow	OTH Thumb Loop w/ Elastic Wrist	SMS	10/Bag	100/Case
51193FCB	X-Large	48"	61"	Yellow	OTH Thumb Loop w/ Elastic Wrist	SMS	Individually Bagged	100/Case
51194FC	XX-Large	51"	62 ½"	Yellow	OTH Thumb Loop w/ Elastic Wrist	SMS	10/Bag	100/Case
51623FC	Universal	48"	61"	White	OTH Thumb Loop w/ Elastic Wrist	Polycoated	10/Bag	100/Case

①AATCC 42 – Water Resistance: Spray Impact Penetration Test / AATCC 127 – Water Resistance: Hydrostatic Pressure Test For more information on these test methods, please visit [aatcc.org](http://aatcc.org) and [aami.org](http://aami.org).

## Barrier Protection AAMI PB70:2012

Test Method	Test Result	Anticipated Risk of Exposure
AATCC 42 Spray Impact	≤ 1.0g	Low
AATCC 127 Hydrostatic Pressure	≥ 20cm	Low



Level  
2

### Precept® Isolation Gown Full Coverage Product Details

- Polycoated and multi-layer SMS material<sup>①</sup> delivers AAMI Level 2 protection and comfort
- Knit cuffs are comfortable and facilitate glove donning
- Full coverage with overlapping back panel
- Secure neck ties
- Sonic-welded seams
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
51112	XX-Large	52"	62"	White	Secure Neck Tie w/ Knit Cuffs	Polycoated	10/Bag	50/Case
51196	Universal	48"	60"	Blue	Secure Neck Tie w/ Knit Cuffs	SMS	10/Bag	100/Case
51198	Universal	48"	60"	Yellow	Secure Neck Tie w/ Knit Cuffs	SMS	10/Bag	100/Case

①AATCC 42 – Water Resistance: Spray Impact Penetration Test / AATCC 127 – Water Resistance: Hydrostatic Pressure Test For more information on these test methods, please visit [aatcc.org](http://aatcc.org) and [aami.org](http://aami.org).

Barrier Protection AAMI PB70:2012		
Test Method	Test Result	Anticipated Risk of Exposure
AATCC 42 Spray Impact	≤ 1.0g	Low
AATCC 127 Hydrostatic Pressure	≥ 20cm	Low



**Precept® Isolation Gown Over-the-Head Thumb Loop Full Coverage**  
**Product Details**

- Multi-layer SMS or polycoated material<sup>①</sup> delivers AAMI Level 2 protection and comfort
- Thumb loops with tapered wrist ensure proper glove donning
- Full coverage with overlapping back panels
- Secure neck ties
- Sonic-welded seams
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
51192FC	Large	48"	61"	Yellow	OTH Thumb Loop w/ Tapered Wrist	SMS	10/Bag	100/Case

①AATCC 42 – Water Resistance: Spray Impact Penetration Test / AATCC 127 – Water Resistance: Hydrostatic Pressure Test For more information on these test methods, please visit [aatcc.org](http://aatcc.org) and [aami.org](http://aami.org).

Barrier Protection AAMI PB70:2012		
Test Method	Test Result	Anticipated Risk of Exposure
AATCC 42 Spray Impact	≤ 1.0g	Low
AATCC 127 Hydrostatic Pressure	≥ 20cm	Low



Level  
2

**Precept® Isolation Gown Full Coverage**  
**Product Details**

- Multi-layer SMS material<sup>①</sup> delivers AAMI Level 2 protection and comfort
- Elastic wrist ensure proper glove donning
- Full coverage with overlapping back panels
- Tape tab neck closure
- Sonically sealed seams on sleeves, armholes, shoulders, and side ties
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
59979	Large	45"	55"	Yellow	Tape Tab Neck Closure w/ Elastic Wrist	SMS	10/Bag	100/Case
59989	X-Large	48"	60"	Yellow	Tape Tab Neck Closure w/ Elastic Wrist	SMS	10/Bag	100/Case

①AATCC 42 – Water Resistance: Spray Impact Penetration Test / AATCC 127 – Water Resistance: Hydrostatic Pressure Test For more information on these test methods, please visit [aatcc.org](http://aatcc.org) and [aami.org](http://aami.org).



## Precept® Protective Gown Over-the-Head Thumb Loop with Open Back

### Product Details

- Medium weight multi-layer SMS material<sup>①</sup> passes AAMI Level 2 protection and provides comfort
- Over-the-head donning helps improve compliance
- Thumb loops and elastic wrist ensure proper glove donning
- Sonic-welded seams
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
51192	Large	42"	41"	Yellow	OTH Thumb Loop w/ Elastic Wrist Open Back	SMS	10/Bag	100/Case
51193	X-Large	49"	43 ½"	Yellow	OTH Thumb Loop w/ Elastic Wrist Open Back	SMS	10/Bag	100/Case

①AATCC 42 – Water Resistance: Spray Impact Penetration Test / AATCC 127 – Water Resistance: Hydrostatic Pressure Test For more information on these test methods, please visit [aatcc.org](http://aatcc.org) and [aami.org](http://aami.org).



**Precept® Protective Gown Over-the-Head Thumb Loop with Open Back**

**Product Details**

- Medium weight multi-layer SMS material<sup>①</sup> passes AAMI Level 2 protection and provides comfort
- Over-the-head donning helps improve compliance
- Thumb loops ensure proper glove donning
- Sonic-welded seams
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
50095	Universal	42"	41"	Yellow	OTH Thumb Loop Open Back	SMS	10/Bag	100/Case

①AATCC 42 – Water Resistance: Spray Impact Penetration Test / AATCC 127 – Water Resistance: Hydrostatic Pressure Test For more information on these test methods, please visit [aatcc.org](http://aatcc.org) and [aami.org](http://aami.org).





## Precept® Blue Film Gown Over-the-Head Thumb Loop with Open Back

### Product Details

- Blue polyethylene film material<sup>®</sup> passes ASTM F1671 viral testing providing enhanced barrier protection
- Generous sizing improves comfort
- Over-the-head donning helps improve compliance
- Thumb loops with tapered wrist ensure proper glove donning
- Thermal impulse seams
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
51575	Large	45"	37"	Blue	OTH Thumb Loop w/ Tapered Wrist Open Back	PE Film	15/Box	75/Case
51576	X-Large	49"	44"	Blue	OTH Tapered Wrist w/ Thumb Loop Open Back	PE Film	15/Box	75/Case
8576B	X-Large	49"	44"	Blue	OTH Tapered Wrist w/ Thumb Loop Open Back	PE Film	Individually Bagged	100/Case

© ASTM F1671: Standard test method for resistance of material used in protective clothing to penetration by blood-borne pathogens using Phi-174 bacteriophage penetration as a test system. For more information on this test method, please visit [astm.org](http://astm.org).



## Precept® Blue Film Gown Over-the-Head Thumb Loop with Open Back

### Product Details

- Blue polyethylene film material<sup>®</sup> passes ASTM F1671 viral testing providing enhanced barrier protection
- Generous sizing improves comfort
- Over-the-head donning helps improve compliance
- Thumb loops with elastic wrist ensure proper glove donning
- Thermal impulse seams
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
51572	X-Large	49"	44"	Blue	OTH Thumb Loop w/ Elastic Wrist Open Back	PE Film	15/Box	75/Case
8572	Universal	46"	36 ½"	Blue	OTH Thumb Loop w/ Elastic Wrist Open Back	PE Film	15/Box	75/Case
8572B	Universal	46"	36 ½"	Blue	OTH Thumb Loop w/ Elastic Wrist Open Back	PE Film	Individually Bagged	100/Case

©ASTM F1671: Standard test method for resistance of material used in protective clothing to penetration by blood-borne pathogens using Phi-174 bacteriophage penetration as a test system. For more information on this test method, please visit [astm.org](http://astm.org).

## Barrier Protection AAMI PB70:2012

Test Method	Test Result	Anticipated Risk of Exposure
AATCC 42 Spray Impact	≤ 1.0g	Low
AATCC 127 Hydrostatic Pressure	≥ 20cm	Low



Level  
2

### Precept® Chemotherapy Gown, Full Coverage Product Details

- Tested for commonly used chemotherapy drugs per ASTM F739-20<sup>®</sup>
- Polycoated material delivers AAMI Level 2 protection and comfort
- Knit cuffs provide comfort
- Secure neck ties
- Full coverage with overlapping back panels
- Intended to be used in environments where there is a potential hazard in the handling of chemotherapy drugs
- Universal, one size fits most
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
50530	Universal	48"	60"	White	Secure Neck Tie w/ Knit Cuffs	Polycoated	10/Bag	100/Case

### Chemotherapy Commonly Tested Drugs

③The material and seams used in these gowns were tested against the following chemotherapy drugs per ASTM F739-20, and no breakthrough was observed during the test duration of 480 minutes based on average standard breakthrough time.

Carmustine (BCNU)	Cisplatin
Cyclophosphamide (Cytoxan)	Dacarbazine
Doxorubicin HCl	Etoposide
Fluorouracil	Methotrexate
Mitomycin C	Paclitaxel
ThioTepa	Vincristine Sulfate



**Precept® Tyvek® Gown, Full Coverage**

**Product Details**

- Lightweight Tyvek® material provides excellent protection against fluids and bacteria
- Tyvek® is noted for its strength and resistance to tears and punctures
- Full coverage, generous sizing
- Sewn seams on sleeves, armholes, and shoulders.
- Secure neck ties and knit cuff
- Low-linting, anti-static
- Intended to be used in environments where low linting is needed, and moderate protection from fluids and bacteria is needed
- Universal, one size fits most
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
35-5703	Universal	45"	56"	White	Secure Neck Tie with Knit Cuffs	Tyvek®	10/Bag	100/Case

Aspen Surgical's line of headwear is designed to meet the various needs of healthcare personnel. Headwear is available in a variety of styles, including bouffants, surgical caps, surgeon hoods, and beard covers.



## Precept® Bouffant Cap

### Product Details

- Designed for user comfort
- Elastic headband
- Made with polypropylene material
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Style	Material	Packaging	Qty
22133	24"	White	Bouffant Cap with Elastic Headband	Polypropylene	125/Box	500/Case



## Precept® Skull Cap

### Product Details

- Designed for user comfort
- Soft absorbent headband
- Elastic back keeps hair in place while providing a flexible fit
- Made with lightweight SMS material
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Style	Material	Packaging	Qty
1014P	Universal	Blue	Skull Cap with Soft Absorbent Headband	SMS	125/Box	500/Case



## Precept® Surgeon Cap

### Product Details

- Made with absorbent Kaycel fabric
- Designed for user comfort
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Style	Material	Packaging	Qty
5200	Universal	Blue	Surgeon Cap	Kaycel	125/Box	500/Case



## Precept® Surgical Hood

### Product Details

- Soft SMS material
- Tie closure
- Increased facial coverage
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Style	Material	Packaging	Qty
1011P	Universal	Blue	Surgical Hood with Tie Closure	SMS	100/Box	200/Case



## Precept® Legionnaire Surgical Hood

### Product Details

- Soft SMS material
- Tape closure
- Increased facial coverage
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Style	Material	Packaging	Qty
5246	Universal	Blue	Legionnaire Surgical Hood with Tape Closure	SMS	50/Bag	100/Case



## Precept® Technique Cap

### Product Details

- Soft SMS material
- Tape closure
- Universal sized
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Style	Material	Packaging	Qty
2210	Universal	Blue	Technique Cap with Tape Closure	SMS	100/Box	300/Case



## Precept® Tyvek® Hood

### Product Details

- Full face with snaps
- Large
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Style	Material	Packaging	Qty
35-60230L	Large	White	Tyvek® Hood with Snap Closure	Tyvek®	Bulk	100/Case



## Precept® Beard Cover

### Product Details

- Lightweight polypropylene material
- Soft comfortable elastic headband
- Excellent breathability
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Style	Material	Packaging	Qty
35-3015	Universal	White	Beard Cover with Elastic Headband	Polypropylene	500/Box	500/Case



Aspen offers shoe and boot covers with anti-skid embossed polyethylene sole, constructed of material such as polypropylene or Tyvek® by DuPont, and the convenience of dispenser box or bulk packaging.



## Precept® Wave™ Shoe Covers

### Product Details

- Anti-skid embossed polyethylene sole
- Triple-elastic runs that allow a smooth wrap around the foot to prevent tripping
- Available in Easy-Open bags, Easy Pull-Out or Econo dispensing boxes
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Style	Material	Packaging	Qty
SHS-0202-MF	Universal	Blue	Shoe cover anti-skid sole	Polypropylene	100/Bag, 3 Bags/Case	300/Case
SHS-0202-CF	Universal	Blue	Shoe cover anti-skid sole	Polypropylene	100/Bag, 6 Bags/Case	600/Case
SHS-0202-BF	Universal	Blue	Shoe cover anti-skid sole	Polypropylene	100/Bag, 10 Bags/Case	1,000/Case
SHS-0252-BF	Universal	Blue	Shoe cover anti-skid sole	Polypropylene	100/Bag, 10 Bags/Case	1,000/Case
SHS-0202-CB	Universal	Blue	Shoe cover anti-skid sole	Polypropylene	100/Box, 6 Dispenser/Case	600/Case
SHS-0202-MC	Universal	Blue	Shoe cover anti-skid sole	Polypropylene	150/Bag, 2 Bags/Case	300/Case
SHS-0102-MF	Universal	Blue	Shoe cover w/o anti-skid sole	Polypropylene	100/Bag, 3 Bags/Case	300/Case
SHS-0102-CF	Universal	Blue	Shoe cover w/o anti-skid sole	Polypropylene	100/Bag, 6 Bags/Case	600/Case
SHS-0102-BF	Universal	Blue	Shoe cover w/o anti-skid sole	Polypropylene	100/Bag, 10 Bags/Case	1,000/Case
SHS-0102-CB	Universal	Blue	Shoe cover w/o anti-skid sole	Polypropylene	100/Box, 6 Dispenser/Case	600/Case
SHS-6202-MF	X-Large	Blue	Shoe cover anti-skid sole	Polypropylene	100/Bag, 3 Bags/Case	300/Case
SHS-6202-CF	X-Large	Blue	Shoe cover anti-skid sole	Polypropylene	100/Bag, 4 Bags/Case	400/Case
SHS-6202-BF	X-Large	Blue	Shoe cover anti-skid sole	Polypropylene	100/Bag, 10 Bags/Case	1,000/Case
SHS-6202-CB	X-Large	Blue	Shoe cover anti-skid sole	Polypropylene	200/Box, 2 Dispenser/Case	400/Case
SHS-6102-MF	X-Large	Blue	Shoe cover w/o anti-skid sole	Polypropylene	100/Bag, 3 Bags/Case	300/Case
SHS-6102-CF	X-Large	Blue	Shoe cover w/o anti-skid sole	Polypropylene	100/Bag, 4 Bags/Case	400/Case
SHS-6102-BF	X-Large	Blue	Shoe cover w/o anti-skid sole	Polypropylene	100/Bag, 10 Bags/Case	1,000/Case
SHS-6102-CB	X-Large	Blue	Shoe cover w/o anti-skid sole	Polypropylene	200/Box, 2 Dispenser/Case	400/Case



## Precept® Tyvek® Shoe and Boot Covers

### Product Details

- Anti-skid embossed polyethylene sole
- Large, X-Large
- Not made with natural rubber latex and non-sterile

Part No.	Size	Color	Style	Material	Packaging	Qty
35-60420L	Large	White	Shoe cover anti-skid sole	Tyvek®	Bulk	100/Case
35-60420XL	X-Large	White	Shoe cover anti-skid sole	Tyvek®	Bulk	100/Case
35-60510L	Large	White	Boot cover 18" tall, anti-skid sole	Tyvek®	Bulk	100/Case
35-60510XL	X-Large	White	Boot cover 18" tall, anti-skid sole	Tyvek®	Bulk	100/Case




## Precept® Lab Coat and Jacket

### Product Details

- Soft knit cuffs and collar provide maximum comfort while enhancing fit
- Snap closures and 3 pockets provide efficiency and convenience
- Special finish provides durable tear resistant material that allows for ease in movement
- Sizing from petite to double extra large
- Intended to be used in environments where low fluid exposure is a risk
- Not made with natural rubber latex and non-sterile

### Sizing Chart — Approximation Only

Height in Inches Keep in mind these will cover your every day clothing. It is recommended to select one size larger than indicated in the chart.




Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
1239P	Petite	39 ½"	48"	Blue	Lab Coat, Knit Cuff and Collar	SMS	30/Bag	30/Case
1239M	Medium	40 ½"	50"	Blue	Lab Coat, Knit Cuff and Collar	SMS	30/Bag	30/Case
1239L	Large	41 ½"	52"	Blue	Lab Coat, Knit Cuff and Collar	SMS	30/Bag	30/Case
1239XL	X-Large	42 ½"	54"	Blue	Lab Coat, Knit Cuff and Collar	SMS	30/Bag	30/Case
1239XXL	XX-Large	43 ½"	56"	Blue	Lab Coat, Knit Cuff and Collar	SMS	30/Bag	30/Case
1725P	Petite	29"	48"	Blue	Lab Jacket, Knit Cuff and Collar	SMS	30/Bag	30/Case
1725M	Medium	29 ¾"	50"	Blue	Lab Jacket, Knit Cuff and Collar	SMS	30/Bag	30/Case
1725L	Large	31 ¾"	52"	Blue	Lab Jacket, Knit Cuff and Collar	SMS	30/Bag	30/Case
1725XL	X-Large	33"	54"	Blue	Lab Jacket, Knit Cuff and Collar	SMS	30/Bag	30/Case
1725XXL	XX-Large	31 ¾"	58"	Blue	Lab Jacket, Knit Cuff and Collar	SMS	30/Bag	30/Case



## Precept® Coveralls, SMS and Tyvek®

### Product Details

- Coveralls are available in a variety of styles and material such as multi-layer SMS or Tyvek®
- Special finish SMS material provides durable tear resistance that allows for ease in movement
- Low-linting, anti-static white Tyvek® material is noted for its strength and resistance to tears and punctures,
- Generous sizing
- Full collar with elastic or open wrist and ankles
- Not made with natural rubber latex and non-sterile



Part No.	Size	Color	Style	Material	Packaging	Qty
1547M	M	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1547L	L	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1547XL	XL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1547XXL	XXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1547XXXL	XXXL	Blue	Lightweight full collar with elastic wrists and ankles	SMS	25/Bag	25/Case
1548M	M	Blue	Lightweight full collar with open wrists and ankles	SMS	25/Bag	25/Case
1548L	L	Blue	Lightweight full collar with open wrists and ankles	SMS	25/Bag	25/Case
1548XL	XL	Blue	Lightweight full collar with open wrists and ankles	SMS	25/Bag	25/Case
1548XXL	XXL	Blue	Lightweight full collar with open wrists and ankles	SMS	25/Bag	25/Case
35-0429XXL	XX-Large	White	Full collar with elastic wrists and ankles	Tyvek®	25/Bag	25/Case
35-1412XL	X-Large	White	Full collar with elastic wrists and ankles	Tyvek®	25/Bag	25/Case



## Precept® Patient Gown, Full Coverage

### Product Details

- To maintain patient dignity the gown has  $\frac{3}{4}$  sleeves with ties at neck and waist
- Made with multi-layer SMS material
- Good opacity and tear resistance
- Disposable patient gowns help provide infection prevention and help eliminate the concern of missing buttons or ties
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
0456P	Universal	47"	57"	Dark Blue	Secure Neck and Waist Ties	SMS	100/Bag	100/Case



## Precept® Apron, Blue Film Over-the-Head Open Back

### Product Details

- Heavyweight blue film material
- Over-the-head donning
- Secure waist ties
- Die cut
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
8571	Universal	59"	43"	Blue	OTH Open Back	Blue Film	20/Bag	20/Case



## Precept® Tyvek® Protective Sleeve

### Product Details

- Lightweight Tyvek® material provides excellent protection against fluids and bacteria
- Tyvek® is noted for its strength and resistance to tears and punctures
- Low-linting, anti-static
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
35-2502	Universal	18"	NA	White	Sleeve with elastic	Lightweight Tyvek®	200/Bag	200/Case




## Precept® Newborn Accessories

### Product Details

- Disposable newborn accessories are designed with the comfort of infants in mind
- Materials are soft, yet durable, to help provide a nurturing environment
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
5180	Newborn	N/A	N/A	Pink & Blue Print	Infant Shirt with Tape Closure	Spunlace	100/Box	400/Case
5431	NA	19"	15"	Blue	Infant Scale Liner	SMS	250/Bag	1000/Case
60112	NA	30"	36"	White	Receiving Blanket	Rayon/ Polypropylene Blend	100/Bag	100/Case



**Precept® Spill Pad and Entrance Mat**

**Product Details**

- Intended to be used for the recovery of spilled liquids and clean-up
- Absorbent material
- Fluid impervious on one side
- Not made with natural rubber latex and non-sterile

Part No.	Size	Length	Width	Color	Style	Material	Packaging	Qty
99226	NA	13"	16"	Blue	Spill Pad	Absorbent blend spunlace laminated to polypropylene blend	200/Bag	200/Case









### Protek™ Sterile Soft PullUp™ General Purpose Probe Covers

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Telescopic-fold ensures easy gel and cover application
- Patented PullUp™ for even faster application
- Soft polyurethane (not made with natural rubber latex)
- Select covers offer a three-dimensional box end
- Kit includes probe cover, sterile gel and 2 bands

Part No.	Size	3D End	Qty
1-519-2231	4" x 24" (10cm x 61cm)	No	24
1-519-2044	5.5" x 24" (14cm x 61cm)	No	24
1-519-2099	5.5" x 24" (14cm x 61cm)	Yes	24
1-519-2048	5.5" x 48" (14cm x 122cm)	No	24
1-519-2097	5.5" x 48" (14cm x 122cm)	Yes	24



### Protek™ Sterile Clear PullUp™ General Purpose Probe Covers

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Telescopic-fold ensures easy gel and cover application
- Patented PullUp™ for even faster application
- Clear polyethylene (not made with natural rubber latex)
- Kit includes PullUp™ probe cover, sterile gel and 2 bands

Part No.	Size	3D End	Qty
1-519-1096	5" x 48" (13cm x 122cm)	No	24



### Protek™ Sterile Soft Telescopically-Folded General Purpose Probe Covers

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Telescopic-fold ensures easy gel and cover application
- Soft polyurethane (not made with natural rubber latex)
- Three-dimensional box end
- Kit includes probe cover, sterile gel and 2 bands

Part No.	Size	3D End	Qty
1-519-2320	5.5" x 36" (14cm x 91cm)	Yes	24
1-519-2450	5.5" x 48" (14cm x 122cm)	Yes	24
1-519-2210	6" x 58" (15cm x 147cm)	Yes	24



### Protek™ Sterile Clear Telescopically-Folded General Purpose Probe Covers

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Telescopic-fold ensures easy gel and cover application
- Clear polyethylene (not made with natural rubber latex)
- Select covers offer a three-dimensional box end
- Kit includes probe cover, sterile gel and 2 bands

Part No.	Size	3D End	Qty
1-519-1236	6" x 48" (15cm x 122cm)*	Yes	100
1-519-1460	6" x 48" (15cm x 122cm)	No	24

\*Kit does not contain sterile gel



### Protek™ Sterile Clear Telescopically-Folded General Purpose Probe Covers

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Telescopic-fold ensures easy gel and cover application
- Clear polyethylene (not made with natural rubber latex)
- Three-dimensional box end
- Kit includes probe cover with arrow, sterile gel, 2 tapes, and 2 bands

Part No.	Size	3D End	Qty
1-519-1090	6" x 48" (15cm x 122cm)	Yes	24



### Protek™ Sterile Soft Flat-Folded General Purpose Probe Covers

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Telescopic-fold ensures easy gel and cover application
- Soft polyurethane (not made with natural rubber latex)
- Select covers offer a three-dimensional box end
- Kit includes probe cover, sterile gel and 2 bands

Part No.	Size	3D End	Qty
1-519-2230	3.5" x 24" (9cm x 61cm)	No	24
1-519-2236	3.5" x 36" (9cm x 91cm)	No	24
1-519-2250	5.5" x 24" (14cm x 61cm)	Yes	24
1-519-2270	7" x 24" (18cm x 61cm)	No	24
1-519-2290	9" x 24" (23 x 61cm)	No	24



### Protek™ Sterile Clear Flat-Folded General Purpose Probe Covers

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Telescopic-fold ensures easy gel and cover application
- Clear polyethylene (not made with natural rubber latex)
- Select covers offer a three-dimensional box end
- Kit includes probe cover, sterile gel and 2 bands

Part No.	Size	3D End	Qty
1-519-1240	4" x 24" (10cm x 61cm)	Yes	24
1-519-1270	7" x 24" (18cm x 61cm)	No	24



### Protek™ Sterile Wireless General Purpose Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Telescopic-fold ensures easy gel and cover application
- Clear polyethylene (not made with natural rubber latex)
- Removable patented PullUp™ with adhesive seal ensures device is fully enclosed
- Kit includes probe cover, sterile gel and 2 bands

Part No.	Size	Qty
1-519-1512	5" x 12" (13cm x 30cm)	24



### Protek™ Sterile Cuffed General Purpose Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment and reduces post-procedure clean-up
- Superior fold with easy grip cuff improves usability by creating a single point-of-entry for easier and faster application
- Soft polyurethane (not made with natural rubber latex)
- 3D end for seamless imaging
- Sterile, intended for single-use

Part No.	Size	Qty
1-519-2648	6" x 48" (15cm x 122cm)	24



### Protek™ SlimFit™ Endocavity Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Intended for single-use
- Soft polyurethane (not made with natural rubber latex)
- Tapered end for a better fit
- Telescopic-fold ensures easy gel and cover application
- Patented PullUp™ for even faster application
- Non-sterile kits include SlimFit™ probe cover and bands.
- Sterile kit includes SlimFit™ probe cover, sterile gel and 2 bands

Part No.	Size	Sterility	Qty
1-519-2512	3.25" – 1.25" x 12" (8cm-3cm x 30cm)	Sterile	24
*1-519-2513	3.25-1.25" x 12" (8cm-3cm x 31cm)	Sterile	200
1-519-2517	3.25" – 1.25" x 12" (8cm-3cm x 30cm)	Non-sterile	24
1-519-2524	4" – 1.25" x 24" (10cm-3cm x 61cm)	Sterile	24
1-519-2548	4" – 1.25" x 48" (10cm-3cm x 122cm)	Sterile	24

\*Probe cover only



### Protek™ Tapered Endocavity Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Intended for single-use
- Soft polyurethane (not made with natural rubber latex)
- Tapered end for a better fit
- Telescopic-fold ensures easy gel and cover application
- Patented PullUp™ for even faster application
- Non-sterile kit includes probe cover and bands
- Sterile kit includes probe cover, sterile gel and 2 bands

Part No.	Size	Sterility	Qty
1-519-2220	5" – 2" x 24" (13cm-5cm x 61cm)	Sterile	24
1-519-2225	5" – 2" x 24" (13cm-5cm x 61cm)	Non-sterile	24




### Protek™ 3.5cm Rolled Latex Endocavity Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Made with natural rubber latex
- Intended for single-use
- Size: 3.5cm x 20cm (1.4" x 7.9")

Part No.	Kit Content	Sterility	Qty
1-519-3249	Probe cover, sterile gel, and 2 bands	Sterile	40
1-519-3230	Probe cover and 4 bands	Sterile	40
1-519-3235	Probe cover and 4 bands	Non-sterile	40
1-519-3241	Probe cover and 4 bands	Non-sterile	50
1-519-3238	Bag containing 100 probe covers	Non-sterile	1



### Protek™ 2.0cm Rolled Latex Endocavity Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Made with natural rubber latex
- Intended for single-use
- Size: 2.0cm x 20cm (0.8" x 7.9")

Part No.	Kit Content	Sterility	Qty
1-519-3233	Probe cover and 4 bands	Sterile	24
1-519-3232	Probe cover and 4 bands	Non-sterile	24



### Protek™ Sterile Accordion-Fold Surgical Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Clear polyethylene (not made with natural rubber latex)
- Size: 6" – 3" x 96" (15cm – 8cm x 244cm)

Part No.	Kit Content	Qty
1-519-1930	Probe cover, boot, sterile gel, and 4 bands	12
1-519-1931	Probe cover, boot, and 4 bands	12
1-519-1933	Probe cover and 4 bands	12



### Protek™ Sterile Telescopically-Folded Surgical Probe Cover Kits w/Arrow

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Clear polyethylene (not made with natural rubber latex)
- Three-dimensional box end
- Kit includes probe cover with arrow, sterile gel, 2 tapes, and 2 bands

Part No.	Size	Qty
1-519-1088	5" – 2¼" x 96" (13cm – 6cm x 244cm)	24
1-519-1092	6" – 3⅝" x 96" (15cm – 9cm x 244cm)	24





### Protek™ Sterile Telescopically-Folded Tapered Surgical Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Clear polyethylene (not made with natural rubber latex)
- Size: 6" – 3" x 96" (15cm-8cm x 244cm)

Part No.	Kit Content	Qty
1-519-1237	Probe cover and 2 bands	100
1-519-1934	Probe cover, boot, sterile gel, and 4 bands	24



### Protek™ Sterile Telescopically-Folded Surgical Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Clear polyethylene (not made with natural rubber latex)
- Kit includes probe cover, sterile gel, and 2 bands

Part No.	Size	Qty
1-519-1026	5.5" x 96" (14cm x 244cm)	24



### Protek™ Sterile Cuffed Surgical Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment and reduces post-procedure clean-up
- Superior fold with easy grip cuff improves usability by creating a single point-of-entry for easier and faster application
- Soft polyurethane (not made with natural rubber latex)
- 3D end for seamless imaging
- Sterile, intended for single-use

Part No.	Size	Qty
1-519-2696	6" – 3½" x 96" (15cm-9cm- x 244cm)	24



### Protek™ Sterile PullUp™ Tapered Surgical Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Patented PullUp™ for even faster application
- Clear polyethylene (not made with natural rubber latex)
- Kit includes probe cover, boot, sterile gel, and 2 tapes

Part No.	Size	Qty
1-519-1938	5" – 3" x 96" (13cm-8cm x 244cm)	24



### Protek™ Sterile PullUp™ Surgical Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Patented PullUp™ for even faster application
- Clear polyethylene and soft polyurethane (not made with natural rubber latex)
- Three-dimensional box end
- Kit includes probe cover, sterile gel, and 2 bands

Part No.	Size	Qty
1-519-2233	5" x 96" (13cm x 244cm)	24



### Protek™ Sterile PullUp™ Surgical Probe Cover Kits w/Form-Fit Tip

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Patented PullUp™ for even faster application
- Form-fit tip ensures a better fit
- Clear polyethylene with soft polyurethane tip (not made with natural rubber latex)

Part No.	Size	Kit Content	Qty
1-519-2016	5" - 0.39" x 96" (13cm-1cm x 244cm)	Probe cover, sterile gel, and 2 clips	10
1-519-2002	5" - 7/8" x 96" (13cm-2cm x 244cm)	Probe cover and 2 clips	10
1-519-2014	5" - 7/8" x 96" (13cm-2cm x 244cm)	Probe cover and 2 tapes	10
1-519-2087	5" - 1 3/4" x 96" (13cm-4cm x 244cm)	Probe cover, sterile gel, and 2 bands	10



### Protek™ Transesophageal Probe Cover Kits

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Sterile, intended for single-use
- Soft polyurethane (not made with natural rubber latex)
- Tapered end for a better fit
- Applicator tube provides access to base of cover for easy, mess-free application of gel
- Size: 1.75" - 1" x 48" (4cm-3cm x 122cm)

Part No.	Kit Content	Sterility	Qty
1-519-2143	Probe cover and 4 tapes	Non-sterile	12
1-519-2140	Probe cover, bite block, sterile gel, and 4 tapes	Sterile	12
1-519-2145	Probe cover, bite block, sterile gel, and 4 tapes	Sterile	12



### Protek™ Sterile PullUp™ Instrument Cover for Vacuum Assisted Biopsy Systems

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural clean up
- Patented needle isolator allows needle to pierce through cover while maintaining a sterile barrier
- Removable patented PullUp™ with self-adhesive ensures device is fully enclosed
- Telescopic-fold ensures easy gel and cover application
- Clear polyethylene (not made with natural rubber latex)
- Sterile, intended for single-use

Part No.	Size	Qty
1-519-1003	3.75" x 17" (10cm x 43cm)	24



### Protek™ Sterile System Covers

#### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Sterile, intended for single-use
- Clear polyethylene (not made with natural rubber latex)

Part No.	Size	Kit Content	Qty
1-519-1008	14.5" x 30" (37cm x 76cm)	CT control rail cover and 6 tapes	50
1-519-1010	6" x 10" (15cm x 25cm)	CT control cradle cover and 3 bands	50
1-519-1011	30" x 30" (76cm x 76cm)	System cover and 2 tapes	24
1-519-1212	12" x 24" (31cm x 61cm)	System cover and 4 tapes	24
1-519-1213	24" x 24" (61cm x 61cm)	System cover and 4 tapes	24
1-519-1624	24" W x 28" H (61cm x 71cm)	System cover w/10" deep pocket and 2 tapes	24
1-519-1628	48" W x 28" H (122cm x 71cm)	System cover w/10" deep pocket and 2 tapes	24
1-519-1217	28" x 17.5" (71cm x 45cm)	Portable system cover w/pocket & adhesive seal	24



## Protek™ Sterile ClingTouch™ Monitor Cover

### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Unrolls for easy application
- Top and bottom adhesive to secure cover
- Sterile, intended for single-use
- Clear polyethylene and vinyl (not made with natural rubber latex)

Part No.	Size	Qty
1-519-0180	17" x 17" (43cm x 43cm)	12
1-519-0181	17" x 12" (43cm x 30cm)	12



## Protek™ Sterile Tablet Cover

### Product Details

- Helps minimize cross-contamination
- Protects equipment
- Sterile, intended for single-use
- Clear polyethylene (not made with natural rubber latex)
- Self-adhesive strip ensures device is fully enclosed

Part No.	Size	Qty
1-519-1028	8" x 10" (20cm x 25cm)	24



### Protek™ Sterile Smartphone Cover

#### Product Details

- Helps minimize cross-contamination
  - Protects equipment
  - Sterile, intended for single-use
  - Clear polyethylene (not made with natural rubber latex)
  - Self-adhesive strip ensures device is fully enclosed
- \*Fits iPhone® 8+ and smaller

Part No.	Size	Qty
1-519-1039	4" x 6" (10cm x 15cm)	24



### Protek™ Probe Storage Covers

#### Product Details

- Protects entire probe including handle after cleaning
- Designated space to apply disinfection label
- Clear polyethylene (not made with natural rubber latex)
- Non-sterile, intended for single-use
- Cover size 4" x 16" (10cm x 41cm)

Part No.	Description	Qty
1-519-1160	Dispenser box containing 100 probe covers	4




### Sanitary MRI Headset Covers

#### Product Details

- Minimizes cross-contamination during MRI scan
- Protects and extends equipment life
- Elastic opening ensures easy application and secure placement
- Comfortable, soft material (not made with natural rubber latex)
- Non-sterile, intended for single-use

Part No.	Description	Size	Color	Qty
1-719-0601	Sanitary MRI Headset Covers	Large	White	1,000/Box







### Oxboro® Instrument Tip Protectors


- Designed to protect valuable instruments from breaking and chipping during sterilization and handling
- Help eliminate reprocessing costs caused by torn pouches or wrap and protect staff from sharp sticks during handling



#### Osteotome Tip Protectors


- Help protect Osteotome instrument edges from chipping with these tip protectors
- Available tinted or clear and vented or non-vented
- All types allow visual instrument inspection and are gas and steam permeable
- Not made with natural rubber latex, Non-sterile, and intended for single-use

#### Osteotome Tip Protectors – Clear




Part No. Vented	Part No. Non-Vented	Size mm	Qty/Bag
094031BBG	094021BBG	6.4 x 19	50
094032BBG	094022BBG	9.5 x 25	50
094033BBG	094023BBG	12.7 x 25	50
094034BBG	094024BBG	19 x 25.4	50
094035BBG	094025BBG	25 x 25	50
094036BBG	094026BBG	32 x 32	50
094037BBG	094027BBG	38 x 38	50
094038BBG	094028BBG	44.5 x 32	50
094040BBG		Assorted	50

#### Osteotome Tip Protectors – Blue



Part No. Vented	Part No. Non-Vented	Size mm	Qty/Bag
094011BBG	094001BBG	6.4 x 19	50
094012BBG	094002BBG	9.5 x 25	50
094013BBG	094003BBG	12.7 x 25	50
094014BBG	094004BBG	19 x 25.4	50
094015BBG	094005BBG	25 x 25	50
094016BBG	094006BBG	32 x 32	50
094017BBG	094007BBG	38 x 38	50
094018BBG	094008BBG	44.5 x 32	50
094020BBG		Assorted	50



### Standard Instrument Tip Protectors

- Available solid, tinted or clear and vented or non-vented
- Gas or steam permeable and can be used in flash or conventional autoclave processes
- Solid colors are radiopaque
- Not made with natural rubber latex, Non-sterile, and intended for single-use

### Standard Tip Protector Sizing Guide

Use this guide to determine the standard tip protector that will best fit your instrument.



1.6mm x 19mm	Part No.	Qty/Bag
Solid Vented	091011BBG	100
Solid Non-Vented	091001BBG	100
Tinted Vented	093011BBG	100
Tinted Non-Vented	093001BBG	100
Clear Vented	092011BBG	100
Clear Non-Vented	092001BBG	100



2mm x 19mm	Part No.	Qty/Bag
Solid Vented	091012BBG	100
Solid Non-Vented	091002BBG	100
Tinted Vented	093012BBG	100
Tinted Non-Vented	093002BBG	100
Clear Vented	092012BBG	100
Clear Non-Vented	092002BBG	100



2.8mm x 19mm	Part No.	Qty/Bag
Solid Vented	091013BBG	100
Solid Non-Vented	091003BBG	100
Tinted Vented	093013BBG	100
Tinted Non-Vented	093003BBG	100
Clear Vented	092013BBG	100
Clear Non-Vented	092003BBG	100



3.2mm x 25.4mm	Part No.	Qty/Bag
Solid Vented	091014BBG	100
Solid Non-Vented	091004BBG	100
Tinted Vented	093014BBG	100
Tinted Non-Vented	093004BBG	100
Clear Vented	092014BBG	100
Clear Non-Vented	092004BBG	100



5mm x 25mm	Part No.	Qty/Bag
Solid Vented	091015BBG	100
Solid Non-Vented	091005BBG	100
Tinted Vented	093015BBG	100
Tinted Non-Vented	093005BBG	100
Clear Vented	092015BBG	100
Clear Non-Vented	092005BBG	100



10mm x 19mm	Part No.	Qty/Bag
Solid Vented	091016BBG	100
Solid Non-Vented	091006BBG	100
Tinted Vented	093016BBG	100
Tinted Non-Vented	093006BBG	100
Clear Vented	092016BBG	100
Clear Non-Vented	092006BBG	100



### Standard Instrument Tip Protectors

- Available solid, tinted or clear and vented or non-vented
- Gas or steam permeable and can be used in flash or conventional autoclave processes
- Solid colors are radiopaque
- Not made with natural rubber latex, Non-sterile, and intended for single-use

### Standard Tip Protector Sizing Guide

Use this guide to determine the standard tip protector that will best fit your instrument.



2mm x 9mm x 25mm Part No.	Qty/Bag
Solid Vented 091017BBG	100
Solid Non-Vented 091007BBG	100
Tinted Vented 093017BBG	100
Tinted Non-Vented 093007BBG	100
Clear Vented 092017BBG	100
Clear Non-Vented 092007BBG	100

Assorted - Solid	Part No.	Qty/Bag
Solid Vented	091020BBG	100
Solid Non-Vented	091010BBG	100



2mm x 16mm x 25mm Part No.	Qty/Bag
Solid Vented 091018BBG	100
Solid Non-Vented 091008BBG	100
Tinted Vented 093018BBG	100
Tinted Non-Vented 093008BBG	100
Clear Vented 092018BBG	100
Clear Non-Vented 092008BBG	100

Assorted - Tinted	Part No.	Qty/Bag
Tinted Vented	093020BBG	100
Tinted Non-Vented	093010BBG	100



3mm x 25mm x 25mm Part No.	Qty/Bag
Solid Vented 091019BBG	100
Solid Non-Vented 091009BBG	100
Tinted Vented 093019BBG	100
Tinted Non-Vented 093009BBG	100
Clear Vented 092019BBG	100
Clear Non-Vented 092009BBG	100

Assorted - Clear	Part No.	Qty/Bag
Clear Vented	092020BBG	100
Clear Non-Vented	092010BBG	100



### DuoGuard™ Instrument Tip Protectors

#### Product Details

- Ensure that your instruments stay in the proper open position during processing with DuoGuard™ Instrument Tip Protectors
- Available vented and non-vented, they are tinted for easy viewing of instrument tips
- Not made with natural rubber latex, Non-sterile, and intended for single-use

Tinted Part No. Vented	Part No. Non-Vented	Color	Size mm	Qty
093031BBG	093021BBG	White	1.6 x 20.3	50
093032BBG	093022BBG	Blue	2 x 20.3	50
093033BBG	093023BBG	Green	2.8 x 20.3	50
093034BBG	093024BBG	Red	3.2 x 25.4	50
093035BBG	093025BBG	Yellow	4.7 x 25.4	50
093037BBG	093027BBG	Orange	9.5 x 26.4	50
093040BBG		Assorted	Assorted	50



### Low Temperature Web Guards

#### Product Details

- An Aspen exclusive
- These unique guards with a flow through, open weave are designed to be used only with low heat, plasma sterilization processes, such as the Sterrad® Sterilization System
- Great for helping to protect endoscopes and their accessories
- Not made with natural rubber latex, Non-sterile, and intended for single-use

Part No.	Color	mm	Qty/Bag
099001BBG	Orange Mesh	10 x 76	24
099002BBG	Green Mesh	15 x 76	24
099003BBG	Blue Mesh	20 x 76	24
099004BBG	Red Mesh	40 x 76	24
099005BBG	Black Mesh	60 x 76	16



Solid/Radiopaque



Clear



Tinted

### InnerVent™ Instrument Tip Protectors

#### Product Details

- Offer the performance of vented tip protectors with the safety of non-vented tip protectors
- Raised inner nubs allow improved flow through of steam and gas while preventing the tip protectors from "vacuuming" onto the instrument during drying
- Not made with natural rubber latex, Non-sterile, and intended for single-use

Part No.	Color	Size/mm	Qty/Bag
091117BBG	Orange, Solid	2 x 9 x 25	100
091118BBG	Brown, Solid	2 x 16 x 25	100
091119BBG	Black, Solid	3 x 25 x 25	100
092117BBG	Clear	2 x 9 x 25	100
092118BBG	Clear	2 x 16 x 25	100
092119BBG	Clear	3 x 25 x 25	100
093117BBG	Orange, Tinted	2 x 9 x 25	100
093118BBG	Brown, Tinted	2 x 16 x 25	100
093119BBG	Black, Tinted	3 x 25 x 25	100



### Gelpi Retractor Tip Protectors

#### Product Details

- Protect expensive Gelpi retractors with Aspen Surgical's exclusive Gelpi Retractor Tip Protectors
- Available in three sizes to fit the most commonly used Gelpi retractors
- Gas and steam sterilization compatible
- Not made with natural rubber latex, Non-sterile, and intended for single-use

Part No.	Color	Size/mm	Qty/Bag
096109BBG	Clear	38 Length	55
096129BBG	Clear	51 Length	45
096139BBG	Clear	76 Length	40



### Weitlaner Instrument Tip Protectors

#### Product Details

- Help protect your expensive Weitlaner instruments and reduce reprocessing costs associated with compromised wrap
- Specifically designed for use with up to a 4" Weitlaner retractor and offer easy one-handed application
- Not made with natural rubber latex, Non-sterile, and intended for single-use

Vented Part No.	Size	Color	Qty/Bag
096409BBG	Small	Magenta, solid	50
096410BBG	Small	Clear	50
096419BBG	Medium	Green, solid	50
096420BBG	Medium	Clear	50



### Laparoscopic MIS Tip Protectors

#### Product Details

- Venting of these specialty tip protectors provides extra flow through to Trocar probes and Veress needles during high level disinfection or sterilization
- Not made with natural rubber latex, Non-sterile, and intended for single-use

Vented and Clear Part No.	Size/mm	Applications	Qty/Bag
097011BBG	19	Use for Veress Needle	50
097012BBG	20 x 5	Use for Trocar Probe	25
097013BBG	43 x 5	Use for Trocar Probe	20
097014BBG	20 x 10	Use for Trocar Probe	25
097015BBG	57 x 10	Use for Trocar Probe	20



### Endoscope Guards

#### Product Details

- Help safeguard the eye pieces of delicate endoscopic equipment during storage with these soft, conformable guards
- Not made with natural rubber latex, Non-sterile, and intended for single-use

#### Blue, Non-Vented

Part No.	Description	Qty/Bag
095016BBG	Small Eye Piece	6



### Round Tip Protectors

#### Product Details

- Designed to help protect valuable instruments, ranging in sizes from the smallest Veress Needle to the largest endoscope, during handling and storage
- Tinted for easy visual inspection and vented for improved steam and gas flow
- Not made with natural rubber latex, Non-sterile, and intended for single-use


#### Vented and Tinted

Part No. 38.1mm	Part No. 50.8mm	Color	Fits Tip Sizes Size mm	Qty/Bag
096001BBG	096021BBG	White	2-3	70
096002BBG	096022BBG	Blue	3-4	70
096003BBG	096023BBG	Green	4-5	70
096004BBG	096024BBG	Red	5-6	70
096005BBG	096025BBG	Yellow	6-7	70
096006BBG	096026BBG	Maroon	7-8	70
096007BBG	096027BBG	Orange	8-9	70



## ColorTrack® Instrument Identification Tape

- Helps make organizing and tracking valuable instruments easy
- Flexible, adheres securely to instruments, and will not chip, flake, peel, or crack
- With 75+ colors, symbols, and character options, the color-coding options are nearly endless
- If a change is needed, simply remove the ColorTrack® I.D. Tape and replace with your new color choice
- Not made with natural rubber latex and Non-sterile



### Sheet Tape

#### Product Details

- Each 8½" x 11" sheet has 374 total inches of tape and contains:
- 44 each of 2" strips
- 44 each of 1½" strips
- 220 each of 1" strips

### Solid Color

¼" Height / Various widths

Part No.	Description
152001EEA	White
152002EEA	Royal Blue
152003EEA	Kelly Green
152004EEA	Red
152005EEA	Yellow
152006EEA	Purple
152007EEA	Orange
152008EEA	Brown
152009EEA	Black
152010EEA	Gray
152011EEA	Lavender
152012EEA	Beige
152013EEA	Light Blue
152014EEA	Lime Green
152015EEA	Pink
152017EEA	Plum
152018EEA	Mint Green

### 4-Color

¼" Height / Various widths

Note: Assorted color package will contain separate sheets each with its own color.

Part No.	Description
152031EEA	Red/Yellow/ Purple/Orange
152032EEA	Blue/White/ Black/Green

### White Stripe

¼" Height / Various widths

Part No.	Description
154002EEA	White/Royal Blue
154003EEA	White/Kelly Green
154004EEA	White/Red
154005EEA	White/Yellow
154007EEA	White/Orange
154008EEA	White/Brown
154009EEA	White/Black
154011EEA	White/Lavender
154013EEA	White/Light Blue
154014EEA	White/Lime Green
154015EEA	White/Pink
154017EEA	White/Plum

### Specialty Stripe

¼" Height / Various widths

Part No.	Description
156003EEA	Blue/Green
156004EEA	Blue/Red
157002EEA	Yellow/Blue
157003EEA	Yellow/Green
157004EEA	Yellow/Red

### Black Stripe

¼" Height / Various widths





Part No.	Description
153001EEA	Black/White
153003EEA	Black/Kelly Green
153004EEA	Black/Red
153005EEA	Black/Yellow
153007EEA	Black/Orange
153010EEA	Black/Gray
153011EEA	Black/Lavender
153012EEA	Black/Beige
153013EEA	Black/Light Blue
153014EEA	Black/Lime Green
153015EEA	Black/Pink
153017EEA	Black/Plum

### Tri-Color Stripe


¼" Height / Various widths

Part No.	Description
155051EEA	Red/White/Blue
155052EEA	Red/White/Black
155053EEA	Red/Yellow/Black
155054EEA	Red/Black/Green



## Dots ¼" Height / Various widths

Part No.	Description
 158040EEA	Fuchsia Dots
 158041EEA	Medical Blue Dots
 158042EEA	Emerald Green Dots
 158043EEA	Violet Dots

## Stars ¼" Height / Various widths





Part No.	Description
 158030EEA	Fuchsia Stars
 158031EEA	Medical Blue Stars
 158032EEA	Emerald Green Stars
 158033EEA	Violet Stars

## Triangles ¼" Height / Various widths



Part No.	Description
 158050EEA	Fuchsia Triangles
 158051EEA	Medical Blue Triangles
 158052EEA	Emerald Green Triangles
 158053EEA	Violet Triangles

## Numbers ¼" Height / Various widths

Each sheet contains numbers 0-8

Part No.	Description
 158020EEA	Fuchsia Numbers
 158021EEA	Medical Blue Numbers
 158022EEA	Emerald Green Numbers
 158023EEA	Violet Numbers








## Assorted Symbols ¼" Height / Various widths


Part No.	Description
 158060EEA	Red Hearts on White
 158080EEA	Alphabet, Fuchsia/White, A-E

## Diagonal Stripes ¼" Height / Various widths

Part No.	Description
 159001EEA	Purple/Black Diagonal
 159002EEA	Orange/Black Diagonal
 159003EEA	Blue/Black Diagonal
 159004EEA	Rubine/Black Diagonal
 159005EEA	Yellow/Black Diagonal
 159006EEA	Green/Black Diagonal
 159007EEA	White/Black Diagonal

## Code to Specialty ¼" Height / Various widths

Part No.	Description
 158090EEA	Mint/White ENT
 158091EEA	Purple/White ER
 158092EEA	Royal Blue/White GEN
 158093EEA	Lt Blue/White OB/GYN
 158094EEA	Kelly Green/White ORTHO
 158095EEA	Orange/White NEURO
 158096EEA	Red/ White CARDIAC



### Roll Tape

#### Application Instructions:

For secure placement through multiple processing cycles, ColorTrack® I.D. Tape should be wrapped 1½ times around the instrument and excess length should be trimmed. Instrument sterility under ColorTrack® I.D. Tape is assured after processing with properly applied tape.\*

\* Test results available on request size and color of tapes shown may vary from actual products.

### Roll Tape

#### Roll Tape ¼" x 300"

Part No.	Description
151001EEA	White
151002EEA	Royal Blue
151003EEA	Kelly Green
151004EEA	Red
151005EEA	Yellow
151006EEA	Purple
151007EEA	Orange
151008EEA	Brown
151009EEA	Black

### Roll Tape

#### Roll Tape ⅛" x 250"

Part No.	Description
151011EEA	White
151012EEA	Royal Blue
151013EEA	Kelly Green
151014EEA	Red
151015EEA	Yellow
151016EEA	Purple
151017EEA	Orange
151018EEA	Brown
151019EEA	Black

Helps  
to Easily  
Remove  
Our Tape!



### Sterion® Same Day Tape Remover

#### Product Details

- An easy, instrument-friendly solution to remove ColorTrack® Tape
- Helps eliminate the risk of cuts and sticks common with scalpel or utility knife tape removal
- Neutral scent helps make Same Day Tape Remover comfortable to use

Part No.	Description
151066EEA	1 Gallon Container



### Oxboro® Corner Protectors

#### Product Details

- Protect your instrument trays from rips, holes, and tears in sterilization wrap due to sharp tray corners, improper handling, or heavy surgical trays
- Aspen Surgical's Corner Protectors are designed to create a flexible barrier between the tray's sharp corners and the sterilization wrap
- They can be used with steam, gas, ozone, and low temperature plasma sterilization and are reusable
- Not made with natural rubber latex



Part No.	Description	Qty/Bag
CP-80NS	Corner Protectors	80



### Sterion® Color-Coded Identification Locking Tags


#### Product Details

- These tamper evident locking tags can be customized to a variety of uses such as identifying instruments, carts, trays, containers or cabinets
- Available pre-printed with "Repair," "Sharpen," or "Good" on one side and writable on the other or blank on one side and writable on the other
- Not made with natural rubber latex

Part No.	Description	Qty/Bag
251001BBG	Red, Imprint: Repair, Write-On	100
251002BBG	Yellow, Imprint: Sharpen, Write-On	100
251003BBG	Green, Imprint: Good, Write-On	100
251004BBG	Blue, Write-On	100
251005BBG	White, Write-On	100
251006BBG	Green, Write-On	100
251007BBG	Black, Write-On	100
251008BBG	Purple, Write-On	100
251009BBG	Orange, Write-On	100
251010BBG	Assorted Colors, Write-On	150

### SafeClean® Instrument and Channel Cleaning Brushes

- Designed to clean surgical/medical instruments and endoscopic devices
- The complete line of SafeClean® Cleaning Brushes are strong enough to effectively remove surgical residue, yet gentle enough to use on the most delicate instruments and devices
- Not made with natural rubber latex



#### General Purpose Brushes

##### Product Details

- Can be used throughout all of your instrument and device reprocessing areas to clean uniquely shaped and larger bore surgical, medical, and endoscopic instruments and devices
- Not made with natural rubber latex

Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
241007BBG	1mm (0.039")	18"	0.5"	2
241010BBG	1.6mm (0.063")	18"	0.5"	2
241015BBG	1.8mm (0.071")	18"	0.5"	2
241019BBG	2mm (0.079")	12"	1"	2
241020BBG	2.2mm (0.086")	18"	0.5"	2
241024BBG	2.5mm (0.098")	24"	2.5"	2
241025BBG	2.7mm (0.106")	18"	1"	2
241029BBG	3mm (0.118")	12"	1"	2
241030BBG	3.2mm (0.126")	18"	2"	2
241032BBG	4mm (0.158")	12"	1"	2
241033BBG	5mm (0.197")	12"	2.5"	2
241034BBG	5mm (0.197")	24"	2.5"	2
241035BBG	5mm (0.197")	16"	2.5"	2
241040BBG	6mm (0.236")	16"	2.5"	2
241041BBG	6mm (0.236")	24"	2.5"	2
241042BBG	7mm (0.275")	16"	2.5"	2
241043BBG	7mm (0.275")	24"	2.5"	2
241044BBG	10mm (0.393")	16"	2.5"	2
241045BBG	11mm (0.433")	16"	2.5"	2
241046BBG	11mm (0.433")	24"	2.5"	2
241047BBG	10mm (0.393")	24"	2.5"	2
241050BBG	12mm (0.472")	16"	2.5"	2
241055BBG	15mm (0.590")	16"	2.5"	2
241060BBG	20mm (0.787")	16"	2.5"	2

### Brush Cleaning Kits

Help make ordering easier with all of the brushes you need in one kit.

#### 241100BBG – Small Cleaning Brush Kit

- Includes various small general purpose cleaning brushes
- Not made with natural rubber latex

This kit contains:

Part No.	Qty/Each
241010BBG	1
241015BBG	1
241020BBG	1
241025BBG	1
241030BBG	1

#### 241200BBG – Large Cleaning Brush Kit

- Includes various large general purpose cleaning brushes
- Not made with natural rubber latex

This kit contains:

Part No.	Qty/Each
241035BBG	1
241040BBG	1
241045BBG	1
241050BBG	1
241055BBG	1
241060BBG	1

#### 241300BBG – "Best Sellers" Cleaning Brush Kit

- Includes a selection of our most popular general purpose cleaning brushes
- Not made with natural rubber latex

This kit contains:

Part No.	Qty/Each
241010BBG	1
241015BBG	1
241020BBG	1
241025BBG	1
241030BBG	1
241035BBG	1
241040BBG	1
241045BBG	1
241050BBG	1
241055BBG	1
241060BBG	1

#### 241700BBG – Micro/Ophthalmic Brush Kit

- Contains a selection of our smallest general purpose cleaning brushes
- Not made with natural rubber latex

This kit contains:

Part No.	Qty/Each
241007BBG	3
241010BBG	3
241015BBG	3
241020BBG	3
241025BBG	3
241030BBG	3

#### 241600BBG – General Cleaning Brush Kit

- Contains an extensive selection of general purpose brushes, as well as toothbrush and fan tip brushes
- Not made with natural rubber latex

This kit contains:


Part No.	Qty/Each
241025BBG	3
241030BBG	3
241035BBG	3
241040BBG	3
241041BBG	3
241045BBG	3
241046BBG	3
241050BBG	3
241055BBG	3
241060BBG	3
242001BBG	3
242002BBG	3
242003BBG	3
243001BBG	3
243004BBG	3

#### 243100BBG – Tooth Brush Kit

- Not made with natural rubber latex

This kit contains:

Part No.	Qty/Each
243001BBG	1
243002BBG	1
243003BBG	1



### Suction Tube Cleaning Brushes

#### Product Details

- Brushes are designed to fit into and effectively clean very small bore devices and instruments used in surgery, MIS, and therapy procedures today
- Not made with natural rubber latex

Part No.	Brush Head Diameter	Brush Total Total	Bristle Length	Qty/ Bag
241003BBG	3 French (1mm, 0.039")	6"	0.5"	2
241005BBG	5 French (1.67mm, 0.066")	6"	0.5"	2
241006BBG	6 French (2mm, 0.079")	8"	0.5"	2
241008BBG	7 French (2.3mm, 0.092")	8"	0.5"	2
241009BBG	8 French (2.7mm, 0.105")	8"	0.5"	2
241011BBG	9 French (3mm, 0.118")	8"	0.5"	2
241012BBG	10 French (3.3mm, 0.131")	8"	0.5"	2
241013BBG	11 French (3.7mm, 0.144")	8"	0.5"	2
241014BBG	12 French (4mm, 0.158")	8"	0.5"	2



### Fan-Tip Brush

#### Product Details

- Features bristles that wrap around the end of the brush to clean bottles, jars, and other large instruments
- Not made with natural rubber latex

Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
242001BBG	4.8mm-12.7mm Tapered (0.189"-0.5")	7"	3"	2
242002BBG	13mm (0.512")	6"	3"	2
242003BBG	16mm (0.623")	5"	2.5"	2
242046BBG	9.6mm (0.384")	24"	2.5"	2





### Tracheostomy Brushes

#### Product Details

- Uniquely designed to clean tubes
- Not made with natural rubber latex

Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
247001BBG	Fits Trach Sizes 00–5FR 9mm (0.354")	4.75"	1.85"	2
247002BBG	Fits Trach Sizes 6–12FR 13mm (0.512")	8"	3.66"	2



### Toothbrush Style

#### Product Details

- Variety of bristle types available
- Designed to clean all types of instruments and stains including nylon, brass, and steel
- Not made with natural rubber latex

Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/ Bag
243001BBG	13mm (0.512"), Nylon	7.25"	1.5"	2
243002BBG	13mm (0.512"), Brass	7.25"	1.5"	2
243003BBG	13mm (0.512"), Steel	7.25"	1.5"	2

Available  
Sterile and  
Non-sterile



### Double-Ended Nylon Cleaning Brushes

#### Product Details

- Features one end with an extra large brush head for general cleaning and one end with a single row of bristles for more intricate cleaning
- Not made with natural rubber latex

Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Qty
243004BBG	Non-sterile Short End, 2mm (0.079") Long End, 8mm (0.315")	7"	0.867" 1.38"	2/Bag
243014PBX	Sterile Short End, 2mm (0.079") Long End, 8mm (0.315")	7"	0.867" 1.38"	1/Pack, 12/Packs/Box



### Reusable Channel Cleaning Brushes

#### Product Details

- Channel Brushes with tapered brush heads, elevator sheath, and comfortable round, plastic handles
- Available in a variety of lengths and materials to fit the range of endoscopic devices
- Not made with natural rubber latex

Part No.	Description Diameter	Total Length	Sheath	Qty/ Bag
248107BBG	3-5mm	70cm	1.6mm	1
248115BBG	3-5mm	115cm	1.6mm	1
248160BBG	3-5mm	160cm	1.6mm	1
248180BBG	3-5mm	230cm	1.8mm	1
248190BBG	3-5mm	120cm	1.8mm	1




### Disposable Scope Channel Cleaning Brushes

#### Product Details

- Brush heads that range from to 2.5mm - 5mm and will fit scopes as small as 1.1mm
- Features more than double the bristle length than competing low-cost brushes
- Strong, sturdy handle making it easy to pull through the scope, while still cleaning the scope effectively
- Intended for single-use and not made with natural rubber latex

Part No.	Description	Color	Total Length	Sheath Diameter	Qty/ Box
248102BBG	2.5mm	Yellow	230cm	1mm	50
248103BBG	3mm	Purple	230cm	1.55mm	50
248104BBG	5mm – Double Ended	Pink	230cm	1.85mm	50
248105BBG	5mm	Pink	230cm	1.85mm	50



### Hand/Nail Brushes

#### Product Details

- Used to clean dirt on hands and nails
- Can also be used for general instrument cleaning
- Not made with natural rubber latex

Part No.	Brush Head Area	Qty/Bag
244001BBG	35mm x 93mm (1.375" x 1.75")	2



## Instrument Brush Storage

### Product Details

- Convenient storage with our acrylic racks to organize your cleaning brushes
- Helps ensure you will have the right brush for the job
- Not made with natural rubber latex

Part No.	Description	Qty/Box
246001BBG	Brush Rack, Medium (without drip tray)	1



## Protek™ Needle Guide Cleaning Brushes

### Product Details

- Use with disinfection cup to clean reusable needle guides and brackets
- 7" long (16 Ga)
- Non-sterile

Part No.	Description	Qty/Box
1-235-5105	Cleaning Brushes	5



## Protek™ Disinfection Cup

### Product Details

- A convenient reservoir for disinfecting ultrasound probes, reusable needle guides, and brackets
- Extended length for maximum probe coverage
- Easy wipe down cleaning
- Slide proof table grips
- Hooks for cord and connection management
- Wall mounting hardware

Part No.	Description	Qty
1-577-6050	Disinfection Cup Kit	1
1-577-6051	Replacement Cups	3

Notes:

### Notes:

[illegible]

### Notes:

[illegible]

Allen, C-Flex, Flexiform, and Universal are either registered trademarks or trademarks of Allen Medical Systems, Inc. ("Allen Medical"). Arthrex, PivotPost, and S.T.a.R. are either registered trademarks or trademarks of Arthrex, Inc. ("Arthrex"). Trimano is a registered trademark of MAQUET GmbH ("MAQUET"). DeRoyal is a registered trademark of DeRoyal Industries, Inc. ("DeRoyal"). IMP, Innovative Medical Products, Inc., McGuire, Universal Lateral Positioner, Knee Positioning Triangles, De Mayo Knee Positioner, De Mayo Universal Distractor, and Pressure Protector Pad, are either registered trademarks or trademarks of Innovative Medical Products, Inc. ("IMP"). InstraMed is a trademark of InstraMed, Inc. ("InstraMed"). Affiniti, Affiniti 50, Affiniti 70, ClearVue, ClearVue C9-4v, ClearVue 350, ClearVue 550, ClearVue 650, ClearVue 850, HD7 XE, HD9, HD11, HD115, and Sparq are either registered trademarks or trademarks of Koninklijke Philips N.V. ("Philips"). Mizuho OSI, Hana, ProFX, and SSXT are either registered trademarks or trademarks of Mizuho OSI, Inc. ("Mizuho OSI"). SchureMed is a trademark of Schuerch Corporation ("Schuerch"). ACUSON, ACUSON X500, ACUSON X100, ACUSON X150, ACUSON X300, ACUSON NX, ACUSON NX2, ACUSON NX3, ACUSON NX3 Elite, ACUSON P500, ACUSON S2000, Sonoline G20, Sonoline G40, Sonoline G50S, and Sonoline Antares are either registered trademarks or trademarks of Siemens Medical Solutions USA, Inc. ("Siemens"). Skytron is registered trademark of Skytron, LLC. ("Skytron"). Smith & Nephew, GUHL, T-MAX, SPIDER, SPIDER2, FERKEL, and SWANSON SUTURE PASSER are either registered trademarks or trademarks of Smith & Nephew PLC ("Smith & Nephew"). HEWSON is a registered trademark of Smith & Nephew, Inc. ("Smith & Nephew"). Steris is a registered trademark of Steris Inc. ("Steris"). Stryker is a trademark of Stryker Corporation ("Stryker"). Velcro is a registered trademark of Velcro IP Holdings LLC. Zimmer is a registered trademark of Zimmer, Inc. ("Zimmer").

Aspen Surgical Products, Inc. is not affiliated with, endorsed by, or sponsored by Allen Medical, Arthrex, DeRoyal, IMP, InstraMed, Mizuho OSI, Philips, Schuerch, Siemens, Skytron, Smith & Nephew, Steris, Stryker, or Zimmer.

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2024 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1072409 rev 3 02-JULY-2024 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

USA 616-698-7100  
Canada 888-364-7004  
[www.aspensurgical.com](http://www.aspensurgical.com)

