


# FastClamp<sup>™</sup>

## Laparoscopic Clamping System

The FastClamp<sup>™</sup> clamping system is a high specification laparoscopic clamping solution. It ensures quick, stable, unobtrusive instrument and scope fixation.

### FastClamp<sup>™</sup> clamping arm:

- Has multiple degrees of freedom, yet can be securely locked to position instruments accurately and quickly
- Securely grips either 5mm or 10mm instruments with 'soft jaws', which also minimize the risk of instrument damage
- Uniquely designed to allow one or two arms to be used with each table clamp
- Two arm lengths available: standard and bariatric, ideally complementing the PretzelFlex<sup>™</sup> retractor and other retraction instruments



### FastClamp<sup>™</sup> table clamp:

- Easily positioned with a wide range of options of adjustment for secure, yet rapid clamping to any standard table rail
- The clamp handle's toggle action ensures a rigid connection to the operating table
- Entire system can be released in seconds

Visit [symmetrysurgical.com](http://symmetrysurgical.com) or contact your local Sales Representative for more information.

### FastClamp<sup>™</sup> Laparoscopic Clamping System

#### FastClamp<sup>™</sup> Clamping Arm Options

CODE	DESCRIPTION
A-1000	FastClamp <sup>™</sup> System Includes: 1 Table Clamp - TC-1000-01 1 Instrument Clamp - IC-1000-01 2 Vertical Bars Standard Length - BC-1030-01 & BC-1050-01
A-1001	Extended FastClamp <sup>™</sup> System Includes: 1 Table Clamp - TC-1000-01 1 Extended Instrument Clamp - IC-1001-01 2 Vertical Bars - BC-1030-01 & BC-1090-01

#### FastClamp<sup>™</sup> Table Clamp Options

CODE	DESCRIPTION
TC-1000-01	Table Clamp
IC-1000-01	Instrument Clamp
IC-1001-01	Extended Instrument Clamp

#### Bar Options

CODE	DESCRIPTION
BC-1030-01	30mm Bar
BC-1050-01	50mm Bar
BC-1070-01	70 mm Bar
BC-1090-01	90mm Bar



3034 Owen Drive | Antioch, TN 37013

800.251.3000 | 615.964.5532

[customerservice@symmetrysurgical.com](mailto:customerservice@symmetrysurgical.com)

[symmetrysurgical.com](https://www.symmetrysurgical.com)

©2020 Symmetry Surgical, Inc. All rights reserved.