



## Specialty Light Handle and Camera Covers designed to work with STERIS® Light Handles & Cameras

### Intended Use

Intended to provide a sterile interface between the non-sterile operating room camera & surgical light from the sterile field of the nurse or surgeon.

### Key Features for Light Handle & Camera Cover

- Creates uniformity throughout hospital by utilizing one light handle and one camera cover product for multiple STERIS® Models
- Helps to save time and money associated with reprocessing reusable handles & covers
- Helps enable quicker O.R. turnaround times
- A sturdy opening with a large lip helps for an easier install
- Sterile, intended for single use, and not made with natural rubber latex

### Additional Camera Cover Features

- Clear lens helps allow for complete visualization of surgical field
- No distortion of image through lens
- Ultra-durable clear lens to help avoid breakage during surgery

# Specialty Light Handle and Camera Covers designed to work with STERIS® Light Handles and Cameras

Intended to provide a sterile interface between the non-sterile operating room camera & surgical light from the sterile field of the nurse or surgeon.



## Specialty Light Handle and Camera Covers designed to work with STERIS® Light Handles & Cameras Technical Specifications

Part No.	Description	Color	Sterility	Qty
LT-53	Disposable O.R. Light Handle Cover, comparable to STERIS® LB53 (Compatible with HarmonyAIR® E-Series, HarmonyAIR® M-Series, HarmonyAIR® G-Series, Harmony LED, Harmony vLED, Harmony LA)	Blue	Sterile	120/Box
LT-54	Disposable O.R. Camera Cover, comparable to STERIS® LB54 (Compatible with HarmonyAIR® E-Series, HarmonyAIR® M-Series, HarmonyAIR® G-Series, FREE5® Camera System, Other legacy STERIS® Harmony Surgical Lighting Systems)	Blue	Sterile	30/Box

HarmonyAIR® and FREE5® are registered trademarks of STERIS® Corporation. STERIS® is a registered trademark of STERIS® Inc.

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2023 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1090033 rev 2 17-May-2023 ENG – US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
aspensurgical.com

