

 symmetry surgical® | **Bovie**®
FAMILY PRACTICE PRODUCTS



HIGH FREQUENCY DESICCATORS AND ELECTROSURGICAL GENERATORS



Bovie DERM 102

Feature: Easy to use and comes with a starter accessory package

Benefit: Ready to use to maximize your ROI

Feature: Adjust power from 0.1 watt-10 watts

Benefit: No grounding pad is needed

Feature: 3-Button Handpiece (included)

Benefit: Allows the surgeon to adjust the power output up and down as well as activate the unit, even when in foot control mode

Feature: Large display screen

Benefit: Allows the user to see power setting from across the room



Bovie DERM 942

Feature: Designed to require fewer accessories for full functionality

Benefit: Saves time during procedures and requires the purchase of fewer components

Feature: Adjust power from 0.1 watt – 40 watts

Benefit: Higher resolution in bipolar mode

Feature: 3-Button Handpiece (included)

Benefit: Allows the surgeon to adjust the power output up and down as well as activate the unit, even when in foot control mode



Bovie Bantam | PRO A952

Feature: Designed and engineered to be the most reliable and durable generator on the market

Benefit: The unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts

Feature: 50 watts of power

Benefit: No grounding pad is needed when in desiccation mode

Feature: High frequency desiccator with cutting ability

Benefit: Provides a two-in-one generator system offering 50 watts of cutting power as well as coag capabilities

Feature: 3-Button Handpiece (included)

Benefit: Allows the surgeon to adjust the power output up and down as well as activate the unit, even when in foot control mode

Feature: Split Pad Bovie NEM Technology

Benefit: Complete patient safety to prevent burns during procedures



UNSURPASSED 4-YEAR WARRANTY ON ALL BOVIE ELECTROSURGERY GENERATORS!



GENERATOR COMPARISON

FEATURES	DERM 102	A942	A952
Cut	N/A	N/A	50 W
Blend	N/A	N/A	50 W
Coag	10 W	40 W	50 W
Fulguration/Spray Coag	10 W	40 W	40 W
Standard Bipolar	10 W	40 W	40 W
Pad Alarm	N/A	N/A	Included
Pad Sensing (tissue impedance)	N/A	N/A	Bovie NEM

FOOT CONTROL/HANDPIECE ACCESSORIES	DERM 102	A942	A952
Foot Control	A803	A803	A803
Pencils/Handpiece	A902	A902	A902

ILLUMINATING THE PATH TO BETTER HEALTHCARE

MEDICAL LIGHTING

Bovie® LED medical lights are energy efficient to provide 50,000 hours of cool, clear visibility, making them ideal for patient exam and procedure rooms. Bovie procedure lights (available in floor-, wall- and ceiling-mount options) include a sterilizable handle for dimming and illumination field size control with one hand, keeping the user in the sterile field.

Bovie continues to prioritize safety with the UV24 Air Purification systems. These units function as ceiling lights while also providing 24/7 protection to your patients and staff against airborne bacteria, viruses, and mold in all areas of your facility.



RECOMMENDED COMPLEMENTARY PRODUCTS

Bovie® Smoke Shark® III Smoke Evacuator

The Bovie® Smoke Shark® III incorporates powerful and effective smoke evacuation to reduce the biological, physical, and chemical hazards present in the surgical smoke generated from the use of any energy devices including electrosurgical, laser, and ultrasonic products. It is designed as a smaller device to uniquely address the limited space and shorter procedures performed in the physician's office.



STYLE # SE03-120

▼ Products Included

- SE03-120** (1) Smoke Shark® III
- SF06** (1) 6-hour filter
- SS95-1** (1) SharkSkin adaptor
- SERS2** (1) Smart Shark Remote Activator Switch
- SE03-M** (1) Wall mount kit
- 110 VAC Hospital-grade power cord*
- User's guide (**P000027200**)*

DIMENSIONS

- Height:** 18cm (7.1")
- Width:** 23.9cm (9.4")
- Depth:** 15.5cm (6.1")
- Weight:** 2.8kg (6.25 lbs.)

LINE VOLTAGE

120 VAC

LINE FREQUENCY

50 – 60 Hz

WARRANTY

Two years

PRODUCTS INCLUDED

			
<p>SF06 Smoke Shark® III Multi-procedure 6-hour filter</p>	<p>SS95 SharkSkin SharkSkin Electrosurgery Pencil Adaptor with 3/8" x 10' tubing</p>	<p>SERS2 Bovie® Smart Shark Remote Activator Switch</p>	<p>SE03-M Wall mount kit</p>



RECOMMENDED COMPLEMENTARY PRODUCTS (CONTINUED)

Diagnose and Treat

	
<p>DermaScout™ II The DermaScout II is an efficient and dependable instrument for visualization of dermal lesions.</p>	<p>High-Temperature Cauteries Bovie® is the world leader in battery-operated cauteries.</p>

Pencils, Sheaths and Penlights

			
<p>A902 Three-button power control & autoclavable handpiece</p>	<p>A910ST Disposable handpiece and cord sheath, sterile</p>	<p>A910 Disposable handpiece sheath, non-sterile</p>	<p>6666 Physician's Disposable Penlights</p>

Bipolar Accessories

		
<p>A824 Reusable bipolar forceps, 4.5", curved, micro (0.5mm) tips</p>	<p>A827V Reusable bipolar forceps power cord</p>	<p>A803 Footswitch for A952, A942, and DERM 102</p>

RECOMMENDED COMPLEMENTARY PRODUCTS (CONTINUED)

Assortment of Electrode Styles

Mobile Stand

					
ES01 Standard Blade Electrode	ES02 Standard Needle Electrode	ES07 Extended Ball Electrode	ES22-8 1/4" Disposable Loop Electrode	ES24-8 3/8" Disposable Loop Electrode	A812-C Mobile stand for A952, A942, and DERM 102 (47" high) <i>*Does not accommodate the SE03-120 Smoke Shark® III</i>

Dermal Tip Disposable Electrodes

				
A804 (non-sterile) & A805 (sterile) Sharp dermal tip	A806 (non-sterile) & A807 (sterile) Blunt dermal tip	A806DE (non-sterile) & A807DE (sterile) Derm-Elite blunt dermal tip	H10012 (non-sterile) & H10008 (sterile) H-type sharp tip electrode	H10112 (non-sterile) & H10108 (sterile) H-type blunt tip electrode





3034 Owen Drive | Antioch, TN 37013 | 800.251.3000 (US) | +41 (0) 52 557 59 00
customerservice@symmetrysurgical.com | [symmetrysurgical.com](https://www.symmetrysurgical.com)
©2021 Symmetry Surgical Inc. All rights reserved.
BOVIE FAMILY PRACTICE BROCHURE 05.25.2022 v1