

BOOKWALTER® III RETRACTOR SYSTEM



"Good exposure is the key to good surgery."

John Bookwalter, M.D.




Dr. John Bookwalter designed and developed the **Bookwalter® III Retractor System** to enhance function and visualization during general surgery. It is designed with a natural toe-in motion that mimics hand-held retractors.

The system is constructed with lightweight components to reduce set weight and includes an assortment of blades to support numerous surgical procedures.

BOOKWALTER® III RETRACTOR SYSTEM

The **Bookwalter® III Retractor System** is included in the Bookwalter patented portfolio representing more than 40 years of innovation, performance and reliability.



The **Bookwalter® III Retractor System** includes the following:

CODE	DESCRIPTION	QTY
50-5710	Bookwalter III Retractor System	
50-4688	Bookwalter 12" Lightweight Table Post	1
50-4554	Bookwalter Lightweight Table Post Coupling	1
50-4582	Bookwalter Lightweight Horizontal Flexbar	1
50-4583	Bookwalter Medium Segmented Ring	1
50-4636	Bookwalter 4" Straight Ring Segment	1
50-4580	Bookwalter Tilting Ratchet	4
50-4696	Bookwalter Low Profile Self-Retaining RoTilt Ratchet	4
50-4568	Bookwalter 2x3 Kelly Retractor Blade	2
50-4569	Bookwalter 2x4 Kelly Retractor Blade	2
50-4570	Bookwalter 2x5 Kelly Retractor Blade	1
50-4571	Bookwalter 2x6 Kelly Retractor Blade	1
50-4566	Bookwalter Harrington Retractor Blade	1
50-4576	Bookwalter Balfour Retractor	1
50-4672	Bookwalter 1 1/2 x6 Malleable Retractor Blade	1
50-4673	Bookwalter 2x6 Malleable Retractor Blade	2
50-4674	Bookwalter 3x6 Malleable Retractor Blade	1
BW0338015	Bookwalter Twister Ratchet	1