



Symmetry Surgical Karlin Magnum, PG Magnum and XL Magnum instruments

- Lightweight
- Ribbed, satin-black handle to assure non-slip grip
- Optimally weighted for surgeon comfort and control.



Karlin Magnum Instruments

Karlin Magnum Instruments

- Strengthened and reinforced handles for mallet impact durability
- Ribbed, satin-black handles facilitate non-slip grip and reduce glare
- Hardened, razor sharp steel cutting edges facilitate precise osteotomies
- Optimal weight and balance for surgeon comfort and control



Total Length: 10.7" (27cm) · Handle Length: 5.9" (15cm) · Handle Diameter: 1.25" (3cm) · Working Distance: 4.7" (12cm)

Straight/Curved Magnum Chisels

	30mm	25mm	20mm	15mm	10mm	6mm
Straight	46-3106	46-3105	46-3104	46-3103	46-3102	46-3101
Curved	46-3116	46-3115	46-3114	46-3113	46-3112	46-3111



Karlin PG Magnum and XL Magnum Instruments

- Two-handed grip delivers controlled power
- Dual wrist action produces fine motor precision
- Extra length enhances visualization and addresses the challenge of large patients
- Razor sharp, lightweight and durable blades reduce surgeon effort and minimize operative time
- Patented* ribbed, satin-black handles assure non-slip grip and reduce glare

Karlin PG: Total Length: 12.6" (32cm) · Handle Length: 7.8" (20cm) · Handle Diameter: 1.12" (3cm) · Working Distance: 4.7" (12cm)

Karlin XL: Total Length: 13.25" (33.7cm) · Handle Length: 8.5" (22cm) · Handle Diameter: 1.25" (3.2cm)

Karlin PG Magnum and XL Magnum Ring Curettes

Unique, broad-neck fulcrum provides atraumatic leverage on opposing vertebral end plates. Dual-coned, precision sharpened rings facilitate rapid, efficient bone sculpting and removal. Each instrument is engraved with depth indicator to minimize potential for deep penetration.



PG Magnum	46-3406	46-3405	46-3404	46-3403	46-3402	46-3401
XL Magnum	46-3221	46-3220	46-3219	46-3218	46-3217	46-3216

Karlin PG Magnum and XL Magnum Cup Currettes

Three tip configurations: Forward Straight (FS), Forward Angled (FA), and Down Pushing (DP.) Large capacity cup sizes permit rapid tissue removal.

Forward
Straight

19mm



17mm



15mm



13mm



11mm



9mm



7mm



PG Magnum
XL Magnum

46-3417
46-3237

46-3416
46-3236

46-3415
46-3230

46-3414
46-3229

46-3413
46-3228

46-3412
46-3227

46-3411
46-3226

Forward
Angled

11mm



9mm



7mm



PG Magnum
XL Magnum

46-3420
46-3233

46-3419
46-3232

46-3418
46-3231

Down Pushing (Backward Down)

9mm



7mm



46-3422
46-3235

46-3421
46-3234

Karlin PG Magnum and XL Magnum Cobb Elevators

Broad, razor sharp blades accelerate periosteal stripping while reducing potential for spinal canal penetration.

35mm



30mm



25mm



20mm



15mm



10mm



PG Magnum
XL Magnum

46-3431
46-3246

46-3430
46-3245

46-3429
46-3244

46-3428
46-3243

46-3427
46-3242

46-3426
46-3241

Visit symmetrysurgical.com or contact your local Sales Representative for more information.

CODE	DESCRIPTION	QTY
46-3100	Magnum Chisels, Straight, Set of 6 & Case Includes	
46-3101	Chisel, Straight, 6mm	1
46-3102	Chisel, Straight, 10mm	1
46-3103	Chisel, Straight, 15mm	1
46-3104	Chisel, Straight, 20mm	1
46-3105	Chisel, Straight, 25mm	1
46-3106	Chisel, Straight, 30mm	1
46-3107	Sterilizing Case for any 6 Magnum Instruments	1

46-3425	PG Magnum Cobb Elevators, Set of 6 & Case Includes	
46-3426	Cobb Elevator, Straight, 10mm	1
46-3427	Cobb Elevator, Straight, 15mm	1
46-3428	Cobb Elevator, Straight, 20mm	1
46-3429	Cobb Elevator, Straight, 25mm	1
46-3430	Cobb Elevator, Straight, 30mm	1
46-3431	Cobb Elevator, Straight, 35mm	1
46-3407	Sterilizing Case for any 12PG Magnum Instruments	1

46-3400	PG Magnum Ring Curettes, Set of 6 & Case Includes	
46-3401	Ring Curette, Straight, 6mm	1
46-3402	Ring Curette, Straight, 8mm	1
46-3403	Ring Curette, Straight, 10mm	1
46-3404	Ring Curette, Angled, 6mm	1
46-3405	Ring Curette, Angled, 8mm	1
46-3406	Ring Curette, Angled, 10mm	1
46-3407	Sterilizing Case for any 12PG Magnum Instruments	1

46-3410	PG Magnum Cup Curettes, Set of 12 & Case Includes	
46-3411	Cup Curette, Forward, Straight, 7mm	1
46-3412	Cup Curette, Forward, Straight, 9mm	1
46-3413	Cup Curette, Forward, Straight, 11mm	1
46-3414	Cup Curette, Forward, Straight, 13mm	1
46-3415	Cup Curette, Forward, Straight, 15mm	1
46-3416	Cup Curette, Forward, Straight, 17mm	1
46-3417	Cup Curette, Forward, Straight, 19mm	1
46-3418	Cup Curette, Forward, Angled, 7mm	1
46-3419	Cup Curette, Forward, Angled, 9mm	1
46-3420	Cup Curette, Forward, Angled, 11mm	1
46-3421	Cup Curette, Down Pushing, 7mm	1
46-3422	Cup Curette, Down Pushing, 9mm	1
46-3407	Sterilizing Case for any 12PG Magnum Instruments	1

CODE	DESCRIPTION	QTY
46-3110	Magnum Chisels, Curved, Set of 6 & Case Includes	
46-3111	Chisel, Curved, 6mm	1
46-3112	Chisel, Curved, 10mm	1
46-3113	Chisel, Curved, 15mm	1
46-3114	Chisel, Curved, 20mm	1
46-3115	Chisel, Curved, 25mm	1
46-3116	Chisel, Curved, 30mm	1
46-3107	Sterilizing Case for any 6 Magnum Instruments	1

46-3240	XL Magnum Cobb Elevators, Set of 6 & Case Includes	
46-3241	Cobb Elevator, Straight, 10mm	1
46-3242	Cobb Elevator, Straight, 15mm	1
46-3243	Cobb Elevator, Straight, 20mm	1
46-3244	Cobb Elevator, Straight, 25mm	1
46-3245	Cobb Elevator, Straight, 30mm	1
46-3246	Cobb Elevator, Straight, 35mm	1
46-3222	Sterilizing Case for any 12XL Magnum Instruments	1

46-3215	XL Magnum Ring Curettes, Set of 6 & Case Includes	
46-3216	Ring Curette, Straight, 6mm	1
46-3217	Ring Curette, Straight, 8mm	1
46-3218	Ring Curette, Straight, 10mm	1
46-3219	Ring Curette, Angled, 6mm	1
46-3220	Ring Curette, Angled, 8mm	1
46-3221	Ring Curette, Angled, 10mm	1
46-3222	Sterilizing Case for any 12XL Magnum Instruments	1

46-3225	XL Magnum Cup Curettes, Set of 12 & Case Includes	
46-3226	Cup Curette, Forward, Straight, 7mm	1
46-3227	Cup Curette, Forward, Straight, 9mm	1
46-3228	Cup Curette, Forward, Straight, 11mm	1
46-3229	Cup Curette, Forward, Straight, 13mm	1
46-3230	Cup Curette, Forward, Straight, 15mm	1
46-3236	Cup Curette, Forward, Straight, 17mm	1
46-3237	Cup Curette, Forward, Straight, 19mm	1
46-3231	Cup Curette, Forward, Angled, 7mm	1
46-3232	Cup Curette, Forward, Angled, 9mm	1
46-3233	Cup Curette, Forward, Angled, 11mm	1
46-3234	Cup Curette, Down Pushing, 7mm	1
46-3235	Cup Curette, Down Pushing, 9mm	1
46-3222	Sterilizing Case for any 12XL Magnum Instruments	1

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

*Patent No. D318,118

©2014 Symmetry Surgical, Inc. All rights reserved. Cat No. SM 0039



3034 Owen Drive | Antioch, TN 37013 | 800.251.3000
 customerservice@symmetrysurgical.com | symmetrysurgical.com