



# Aspen™ EvacMat®

 SURGICAL AND MEDICAL ESSENTIALS

## INTENDED USE

Combines the functionality of absorption and suction to aid in efficiently removing fluid waste from the O.R. floor.

# Aspen™

## EvacMat®

 SURGICAL AND MEDICAL ESSENTIALS

Combines the functionality of absorption and suction to aid in efficiently removing fluid waste from the O.R. floor.

### KEY FEATURES

- Extremely quiet operation helps to decrease distraction of the O.R. staff – as quiet as an indoor conversation
- The top layer of the mat contains wicking to help prevent pooling
- Works well with wall suction
- Simple to use and clean up – red bag friendly
- Low profile allows for the O.R. staff to roll equipment over the mat without loss of suction
- Large enough for two people
- Provides cushion to help reduce fatigue during long procedures
- Aids in keeping the O.R. floor dry to reduce the risk of injuries related to wet floors
- Non-Sterile, intended for single-use, and not made with natural rubber latex

### Aspen™ EvacMat®

Technical Specifications

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
95005	Yellow	Medium	Yes	32" x 48"	5/Case

EvacMat® is a registered trademark of Aspen Surgical Products, Inc. Aspen™ is a trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

PHONE / +1 616.698.7100 TOLL-FREE / +1 888.364.7004 WEBSITE / [aspensurgical.com](http://aspensurgical.com)