



## SurgiSafe® EvacMat® Absorbent Floor Mat

### Intended Use

Combines the functionality of absorption and suction to aid in efficiently removing fluid waste from the O.R. floor.

### Key Features

- Extremely quiet operation helps to decrease distraction of the O.R. staff – as quiet as an indoor conversation
- The top layer of the mat contains wicking to help prevent pooling
- Works well with wall suction
- Simple to use and clean up – red bag friendly
- Low profile allows for the O.R. staff to roll equipment over the mat without loss of suction
- Large enough for two people
- Provides cushion to help reduce fatigue during long procedures
- Aids in keeping the O.R. floor dry to reduce the risk of injuries related to wet floors
- Non-Sterile, intended for single-use, and not made with natural rubber latex

# SurgiSafe® EvacMat® Absorbent Floor Mat

Combines the functionality of absorption and suction to aid in efficiently removing fluid waste from the O.R. floor.



## SurgiSafe® EvacMat® Absorbent Floor Mat Technical Specifications

Part No.	Color	Absorbency Level	Fluid Barrier Backing	Size	Qty
95005	Yellow	Medium	Yes	36" x 48"	5/Case

EvacMat® and SurgiSafe® are registered trademarks of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1014391 rev 3 9-Dec-2020 ENG - US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
[aspensurgical.com](http://aspensurgical.com)

