

Symmetry Surgical[®] Secto[®] Dissectors & Tonsil Sponges




Symmetry Surgical's **Secto[®] Dissectors** are designed to provide a gentle, atraumatic way to separate tissue during surgical procedures. These sterile, cotton-constructed dissectors are x-ray detectable and available in four patterns: Cherry, Triangle, Peanut and Kittner.

- Available in a dispenser box including 96 packages of 5 each; 4 dispenser boxes in each case

Symmetry Surgical's **Secto[®] Tonsil Sponges** are sterile, cotton, non-linting, x-ray detectable gauzes with identification tape and string for use during tonsillectomies.

- Available 5 each per sterile pack, 25 packs per dispenser carton, 2 cartons per case, for a total of 250 per case.



Secto[®] Cherry
Dissectors
81-1002




Secto[®] Triangle
Dissectors
81-1006



Secto[®] Peanut
Dissectors
81-1000



Secto[®] Kittner
Dissectors
81-1001



Secto[®] Tonsil
Sponges
81-1012

Visit symmetrysurgical.com or contact your local Sales Representative for more information.

Secto® Cherry Dissectors

Dispenser box includes 96 packages of 5 each; 4 dispenser boxes to a case.

CODE	DESCRIPTION	DIAMETER
81-1002	Sterile, X-Ray Detectable	3/8 in 10 mm

Secto® Kittner Dissectors

Dispenser box includes 96 packages of 5 each; 4 dispenser boxes to a case.

CODE	DESCRIPTION	DIAMETER
81-1001	Sterile, X-Ray Detectable	5/16 in 8 mm

Secto® Peanut Dissectors

Dispenser box includes 96 packages of 5 each; 4 dispenser boxes to a case.

CODE	DESCRIPTION	DIAMETER
81-1000	Sterile, X-Ray Detectable	3/8 in 10 mm

Secto® Triangle Dissectors

Dispenser box includes 96 packages of 5 each; 4 dispenser boxes to a case.

CODE	DESCRIPTION	DIAMETER
81-1006	Sterile, X-Ray Detectable	1/2 in 13 mm

Secto® Tonsil Sponges

X-ray detectable, sterile, 5 per card/pouch, 25 pouches per dispenser, 2 dispenser boxes per case.

CODE	CORD LENGTH	DIAMETER
81-1010	10in 254mm	5/8 in 16 mm
81-1011	10in 254mm	1 in 25 mm
81-1012	10in 254mm	1-1/4 in 32 mm

