



 symmetry surgical® | **Bovie**®  
**AMBULATORY SURGERY CENTER PRODUCTS**

# THE NEXT GENERATION OF ELECTROSURGERY GENERATORS



**A3350 OR | PRO 300**



**A2350 Surgi-Center | PRO**

**Feature:** The dual monopolar ports allow two surgeons to work simultaneously in the spray coagulation setting.

**Benefit:** Two surgeons can work in separate areas, allowing the procedure to be completed faster.

**Feature:** The Auto Bipolar setting allows the surgeon to activate the bipolar forceps without a foot pedal.

**Benefit:** The surgeon does not have to look down to locate the foot pedal, allowing for quicker activation of the bipolar forcep, saving time and motion during the procedure.

**Feature:** The Bovie® Bipolar (BBp mode) has more power than the Macro, Micro, and Standard bipolar settings.

**Benefit:** This mode provides enough power for the most demanding surgical procedures.

**Feature:** Gentle Coag Mode

**Benefit:** This mode provides gentler coagulation when a Gold Probe™\* is used during a colonoscopy to prevent perforation of the intestine.

**Feature:** Split Pad

**Benefit:** Complete patient safety to prevent burns during procedures.

**Feature:** The dual monopolar ports allow two surgeons to work simultaneously in the spray coagulation setting.

**Benefit:** Two surgeons can work in separate areas, allowing the procedure to be completed faster.

**Feature:** The surgeon can move smoothly from monopolar to bipolar.

**Benefit:** Saves the surgeon time and motion during the procedures.

**Feature:** The A2350 Surgi-Center | PRO generator has four blend settings.

**Benefit:** The four blend settings allow the surgeon to customize how much coag they want in the cut setting.

**Feature:** The CUT II setting allows the surgeon to do laparoscopic procedures.

**Benefit:** Performing laparoscopic procedures in an outpatient environment is usually more convenient for the patient and provides more revenue to the surgeon. This mode will help prevent direct coupling during a laparoscopic procedure.

**Feature:** Split Pad

**Benefit:** Complete patient safety to prevent burns during procedures.

\*Gold Probe is a product of Boston Scientific



**UNSURPASSED 4-YEAR WARRANTY  
ON ALL BOVIE ELECTROSURGERY  
GENERATORS!**



## GENERATOR COMPARISON

FEATURES	A3350	A2350
Cut (# of modes)	300 W (2)	200 W (2)
Blend (# of steps)	200 W (4)	200 W (4)
Pinpoint Coag	120 W	120 W
Fulguration/Spray Coag	120 W	120 W
Gentle Coag	120 W	120 W
Macro Bipolar	150 W	80 W
Standard Bipolar	50 W	50 W
Micro Bipolar	80 W	80 W
Bovie Bipolar	225 W	N/A
Pad Alarm	YES	YES
Pad Sensing (tissue impedance)	YES	YES
Presets (# of presets)	10	10

FOOT CONTROL / HANDPIECE ACCESSORIES	A3350	A2350
Foot Control	YES	YES
Standard 3-Pin	YES	YES
ESPR2 Reusable Pencil	YES	YES
Two Pencil Configuration	YES	YES
US, Canada	YES	YES
CE	YES	YES

# ILLUMINATING THE PATH TO BETTER HEALTHCARE

## MEDICAL LIGHTING


The System Two LED series presents an unmatched value for surgery centers seeking lighting options. The energy-efficient design consumes half the electricity of its halogen counterpart to offer significant long-term savings for the healthcare facility. The long-life 50,000-hour rated LEDs eliminate procedure room downtime and are output-rated at 130,000 Lux with a Kelvin color temperature of 4300°. The System Two LED series is designed with the unique spring arm system which can be adjusted for the life of the product providing optimal surgical site visibility and positioning. This lighting system is available in solo, duo, and tri mounting styles; with or without camera and monitors.

Bovie continues to prioritize safety with the UV24 Air Purification systems. These units function as ceiling lights while also providing 24/7 protection to your patients and staff against airborne bacteria, viruses and mold in all areas of your facility.




## RECOMMENDED COMPLEMENTARY PRODUCTS

### Electrosurgical Pencils



### Return Electrodes and Cord



## RECOMMENDED COMPLEMENTARY PRODUCTS (continued)

### Footswitches



**BV-1254B**

Bipolar footswitch for A3350 and A2350



**BV-1253B**

Monopolar footswitch for A3350, A2350 and A1250S

### Mobile cart



**BV-IDS-CS2**

Mobile stand for A3350, A2350 and A1250S; 16" wide x 20" deep x 35" high

### Smoke Evacuation System



**Smoke Shark II®**

Bovie® Smoke Shark® II is designed to be easy to use and is effective for controlling surgical smoke generated from the use of electro-surgical, laser, and ultrasonic energy products.



**ISEP1000**

Bovie® Orca™, Smoke Evacuation Pencil, 10ft Tubing, Sterile



**SEPAT**

Smoke Evacuation Pencil Adaptor and Tubing, 3/8" x 4' Hose with 7/8" to 6' Tube



**SERS2**









Bovie® Smart Shark Remote Activator Switch




















## RECOMMENDED COMPLEMENTARY PRODUCTS (continued)




### Disposable Electrodes\*

 <p><b>ES581</b> Resistick™ II 6.35cm coated blade</p>	 <p><b>ES57T</b> Resistick™ II 10.2cm coated needle</p>	 <p><b>ES21T</b> Resistick™ II 4mm coated ball</p>	 <p><b>ES54T</b> Resistick™ II 10.2cm coated modified blade</p>
 <p><b>ES60</b> SuperCut™ 3cm tungsten needle</p>	 <p><b>ES63</b> SuperCut™ 4.5cm tungsten needle</p>	 <p><b>AR00</b> Arthroscopic hook, 90°</p>	 <p><b>AR01</b> Arthroscopic hook, 45°</p>

### Laparoscopic Electrodes\*

 <p><b>LJ07</b> J-hook wire</p>	 <p><b>LH08</b> L-hook wire</p>	 <p><b>LN03</b> Needle</p>	 <p><b>LC04</b> Curved blade</p>		
 <p><b>ES09</b> 10mm x 10mm loop</p>	 <p><b>ES11</b> 15mm x 10mm loop</p>	 <p><b>ES12</b> 20mm x 8mm loop</p>	 <p><b>ES49</b> 20mm x 10mm loop</p>	 <p><b>ES07</b> Extended 5mm ball</p>	
 <p><b>ES01</b> Standard 7cm blade</p>	 <p><b>ES02</b> Standard 7cm needle</p>	 <p><b>ES39</b> Extended 15.2cm modified blade</p>	 <p><b>ES40</b> Extended 15.2cm modified needle</p>	 <p><b>ES37</b> Modified 7cm blade</p>	 <p><b>ES38</b> Modified 7cm needle</p>



























### Reusable Bipolar Forcep & Cord\*

 <p><b>A820</b> Bayonet bipolar forceps, 18.4cm long, 1.0mm fine tip</p>	 <p><b>A843</b> Tenzel bipolar forceps, 10.2cm long, 0.5mm micro tip</p>	 <p><b>A827V</b> Bipolar forceps cord, fits A3350, A2350, A1250S, A952, A942 and DERM 102</p>
---	---	--


\*More electrode styles available

# RECOMMENDED COMPLEMENTARY PRODUCTS (continued)

## Microdissection Needles\*

 <p><b>97110BX</b> Needle Electrode, Insulated, Sharp Tip, Tungsten Wire, 3.8cm (1-1/2in)</p>	 <p><b>97111BX</b> Needle Electrode, Insulated, Sharp Tip, Tungsten Wire, 5.1cm (2in)</p>	 <p><b>97112BX</b> Needle Electrode, Insulated, Sharp Tip, Tungsten Wire, 6cm (2-3/8in)</p>
 <p><b>97114BX</b> Needle Electrode, Insulated, Sharp Tip, Tungsten Wire, 15.2cm (6in)</p>	 <p><b>97115BX</b> Needle Electrode, Insulated, Sharp Tip, Tungsten Wire, 45° Angled Shaft, 6cm (2-3/8in)</p>	 <p><b>97116BX</b> Needle Electrode, Insulated, Sharp Tip, Tungsten Wire, 8.9cm (3-1/2in)</p>
 <p><b>97117BX</b> Needle Electrode, Insulated, Sharp Tip, Tungsten Wire, 11.4cm (4-1/2in)</p>	 <p><b>97125BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Bayonet Shaft, 15.2cm (6in)</p>	 <p><b>97130BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 3.2cm (1-1/4in)</p>
 <p><b>97131BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 4.1cm (1-5/8in)</p>	 <p><b>97132BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 5.1cm (2in)</p>	 <p><b>97133BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 6.4cm (2-1/2in)</p>
 <p><b>97134BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 11.7cm (4-5/8in)</p>	 <p><b>97135BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 14cm (5-1/2in)</p>	 <p><b>97136BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 16.8cm (6-5/8in)</p>
 <p><b>97141BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 3mm Angled Tip, 4.1cm (1-5/8in)</p>	 <p><b>97143BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 3mm Angled Tip, 5.1cm (2in)</p>	 <p><b>97144BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 3mm Angled Tip, 10.2cm (4in)</p>
 <p><b>97145BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 3mm Angled Tip, 11.7cm (4-5/8in)</p>	 <p><b>97147BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 3mm Angled Tip, 16.8cm (6-5/8in)</p>	 <p><b>97148BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 3mm Angled Tip, 17.8cm (7in)</p>
 <p><b>97150BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 10mm Angled Shaft, 5.1cm (2in)</p>	 <p><b>97152BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 10mm Angled Shaft, 12.4cm (4-7/8in)</p>	 <p><b>97154BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, 10mm Angled Shaft, 17.5cm (6-7/8in)</p>
 <p><b>97170BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, Curved Shaft, 10.2cm (4in)</p>	 <p><b>97171BX</b> Needle Electrode, Insulated, Ultra Sharp Tip, Tungsten Wire, Straight, Curved Shaft, 15cm (5-7/8in)</p>	

# AMBULATORY SURGERY CENTER INSTRUMENTS



Item Code	Description	Item Code	Description
01-2703	Symmetry® Knife Handle, #3, Standard, 5in, 127mm	10-1213	Symmetry® Needle Holder, Halsey, Serrated, 5-3/8in, 127mm
01-2804	Symmetry® Knife Handle, #4, Standard, 5-1/4 in, 133mm	10-1813	Symmetry® Needle Holder, Mayo-Hegar, 5-1/4in, 133mm
01-2814	Symmetry® Knife Handle, #4L, Long, 8-3/8in, 212mm	10-1815	Symmetry® Needle Holder, Mayo-Hegar, 6-2/3in, 160mm
02-2011	Symmetry® Scissors, Iris Standard Pattern, Straight, 4-1/2in	10-3121	Symmetry® Needle Holder, Heaney, 8-1/4in, 209mm
02-2111	Symmetry® Scissors, Iris Standard Pattern, Curved, 4-1/2in	14-0609	Symmetry® Scissors, Spencer Stitch, 3-1/2in
02-5002	Symmetry® Scissors, Tenotomy, Straight, Sharp, 4-1/4in	15-3209	Symmetry® Forceps, Hartman Baby Mosquito, Straight, 3-7/8in
02-5003	Symmetry® Scissors, Tenotomy, Curved, Sharp, 4-1/8in	15-3309	Symmetry® Forceps, Hartman Baby Mosquito, Curved, 4in
03-0415	Symmetry® Scissors, Operating, Straight, Sharp/Sharp, 6in	15-4412	Symmetry® Forceps, Halsted Mosquito, Straight, 5in
03-0613	Symmetry® Scissors, Operating, Straight, Sharp/Blunt, 5-1/8in	15-4512*	Symmetry® Forceps, Halsted Mosquito, Curved, 5in
03-5017	Symmetry® Scissors, Mayo, Straight, Beveled Blades, 6-3/4in	15-7016	Symmetry® Forceps, Kelly Artery, Straight, 6-1/4in
03-5117	Symmetry® Scissors, Mayo, Curved, Beveled Blades, 6-3/4in	15-7114	Symmetry® Forceps, Kelly Artery, Curved, 5-1/2in
04-3214	Symmetry® Scissors, Metzenbaum Dissecting, Straight, Standard Pattern, 5-3/4in	15-7116	Symmetry® Forceps, Kelly Artery, Curved, 6-1/4in
04-3314	Symmetry® Scissors, Metzenbaum, Curved, Standard Pattern, 5-3/4in	15-7316	Symmetry® Forceps, Crile Artery, Curved, 6-1/4in
04-3318	Symmetry® Scissors, Metzenbaum, Curved, Standard Pattern, 7in	16-1016	Symmetry® Forceps, Rochester-Pean Hemostatic, Straight, 6-5/8in
06-0014	Symmetry® Forceps, Dressing, Serrated, 5-1/2in	16-1116	Symmetry® Forceps, Rochester-Pean Hemostatic, Curved, 6-3/5in
06-0020	Symmetry® Forceps, Dressing, Serrated, 8in	17-5508	Symmetry® Clamp, Backhaus Towel, 3-1/2in
06-0514	Symmetry® Forceps, Tissue, 1x2 Teeth, 5-3/4in	17-5513	Symmetry® Clamp, Backhaus Towel, 5-1/4in, 133mm
06-0520	Symmetry® Forceps, Tissue, 1x2 Teeth, 8in	19-5024	Symmetry® Forceps, Foerster, Straight, Serrated, 9-3/4in
06-2012	Symmetry® Forceps, Adson Dressing, Serrated Jaws, 4-3/4in	19-5579	Classic® Clamp, Gemini, 5-1/2in
06-2112	Symmetry® Forceps, Adson Tissue, 1x2 Teeth, 4-3/4in	20-0508*	Symmetry® Retractor, Senn, Double-Ended, Sharp, 6-1/4in
06-2113	Symmetry® Forceps, Adson Tissue, Serrated, 1x2 Teeth, 4-3/4in	20-5001	Symmetry® Hook, Joseph Skin, Single-Prong, 6.0mm, 6-1/4in, 158mm
06-2512	Symmetry® Forceps, Adson-Brown Tissue, 4-3/4in	20-5005	Symmetry® Hook, Joseph Skin, Double-Prong, 5mm Spread, 4.5mm, 6-1/4in, 158mm
06-2590	Symmetry® Forceps, Callison Tying Tissue, 1.5mm Teeth, Tying Platform, 1x2 Teeth, 4-7/8in	22-1418	Symmetry® Retractor, Parker, Large, 7in, (Set of 2)
06-8798	Forceps, Heaney Tissue, 2x3 Teeth, 5-1/2in	23-5534	Symmetry® Curette, Sims Uterine, Size 1, Sharp, Malleable Shaft, 11in
09-8309	Symmetry® Suction Tube, Frazier, Medium Angle, 9 Fr, 4-1/3in, 110mm	23-7503	Symmetry® Retractor, Deaver, 1in Wide, Standard Handle, 9-1/4in
10-1012	Symmetry® Needle Holder, Webster, Straight, Blunt, Smooth Jaws, 5-1/4in, 133mm	24-2211	Symmetry® Retractor, Weitlaner, Blunt, 2x3 Prong, 4-1/2in
10-1212	Symmetry® Needle Holder, Derf, 5in, 127mm	30-5561	Symmetry® Forceps, Bozeman Dressing, Curved, 260mm

\*Image shown above



## AMBULATORY SURGERY CENTER INSTRUMENTS (continued)

Item Code	Description	Item Code	Description
30-6556	Symmetry® Forceps, Bishop Harmon Iris, Delicate	52-2793	Speculum, Auvarud Vaginal, 1 3/4 in x 4 in Blades, 2.5 lbs
32-110	Symmetry® Needle Holder, Baumgartner, Tungsten Carbide, 5 in	52-4001	Symmetry® Sound, Sims Uterine, Graduated In Centimeters, 13 1/4 in, 335 mm
32-405	Symmetry® Needle Holder, Webster, Tungsten Carbide, Large Ring Handle, 4 3/4 in	52-4426	Symmetry® Forceps, Braun Tenaculum, Straight, 9 1/2 in
32-430	Symmetry® Needle Holder, Crile-Wood, 6 in, 152 mm	54-2555	Symmetry® Scissors, Shortbent Stitch
32-500*	Symmetry® Forceps, Adson Dressing, Tungsten Carbide, 4 3/4 in	54-4030	Symmetry® Scissors, Operating, Straight, 5 1/2 in
50-0215	Symmetry® Forceps, Allis Tissue, 4x5 Teeth, 6 in	60-1813	Symmetry® Scissors, Stitch, Curved, Delicate, Sharp, Large Ribbon Ring Handles w/Flat Shanks, 4 in
50-0219	Symmetry® Forceps, Allis Tissue, 5x6 Teeth, 7 1/2 in	60-1856	Symmetry® Scissors, Strabismus, Straight, Ribbon Handles, 4 1/4 in
50-5540	Symmetry® Retractor, Jackson Vaginal, Small	60-1924*	Symmetry® Needle Holder, Castroviejo, Straight, w/Lock, 12 mm, 5 1/2 in
50-5541	Symmetry® Retractor, Jackson Vaginal, Medium	60-1925	Symmetry® Needle Holder, Castroviejo, Curved, w/Lock, 12 mm, 5 1/2 in
50-5542	Symmetry® Retractor, Jackson Vaginal, Large	67-5507	Scissors, Ragnell, Curved, Beveled Blades, Dissecting Tips, 5 in
52-0610	Symmetry® Speculum, Graves Vaginal, 1 1/4 in x 4 in Blades, Medium	86-5013	Symmetry® Scissors, Stevens Tenotomy, Curved, Serrated, Sc, 4 in

\*Image shown in page 8



# MINIMALLY INVASIVE SURGERY PROCEDURAL DEVICES



Item Code	Description
<b>DUMI-350A</b>	DUMI ManipulatOR™, 350mm adjustable diagnostic uterine manipulator injector
<b>H-250</b>	HolstOR™, Small, 250mm instrument holder used in open surgery, yellow
<b>H-330</b>	HolstOR™, Large, 330mm instrument holder used in laparoscopic surgery, pink
<b>H-410</b>	HolstOR™, X-large, 410mm instrument holder used in bariatric surgery, blue
<b>H-Pro-250</b>	HolstOR PRO™, Small, 250mm instrument holder used in open surgery, yellow
<b>H-Pro-330</b>	HolstOR PRO™, Large, 330mm instrument holder used in laparoscopic surgery, pink
<b>H-Pro-410</b>	HolstOR PRO™, X-large, 410mm instrument holder used in bariatric surgery, blue
<b>IN-12-4</b>	InsufflatOR Needle™, 120mm needle (Gauge: 14GA)
<b>IN-15-4</b>	InsufflatOR Needle™, 150mm needle (Gauge: 14GA)
<b>LSH-01</b>	Liquid Scope HelpOR™, Liquid Scope WarmOR™, two trocar swabs (5-8mm and 10-12mm), one microfiber cloth
<b>PC-120 *</b>	Port ClosOR PRO™, 150mm Port ClosOR (Gauge:14GA) and two trocar swabs (5-8mm and 10-12mm)
<b>PFITF-10</b>	InsufflatOR Tubing Set™, with dual luer lock fitting, 0.1µm filter
<b>PS0012</b>	InsufflatOR Tubing Set™, w/dual connector 0.1µm filter (Storz)
<b>PS003</b>	AntifogOR™, 6ml anti-fog solution with radiopaque sponge
<b>PS-500</b>	Liquid Scope WarmOR™, Long version (three buttons)
<b>PS-500-M *</b>	Liquid Scope WarmOR™, Short version (two buttons)
<b>RSG-14F</b>	Suture PassOR PRO™, 150mm Suture PassOR, one 10-12mm guide and two trocar swabs (5-8mm and 10-12mm)
<b>RSG-14/18F *</b>	Suture PassOR PRO™, 150mm Suture PassOR, two guides (10-12mm and 15mm) and two trocar swabs (5-8mm and 10-12mm)
<b>RSG-18F-XL</b>	Suture PassOR PRO™, 180mm Suture PassOR, two bariatric (XL) length guides (10-12mm and 15mm) and two trocar swabs (5-8mm and 10-12mm)
<b>SC-01</b>	Scope CleanOR™, Two trocar swabs (5-8mm and 10-12mm), one microfiber cloth
<b>SHK-01</b>	Scope HelpOR™, Lap Scope WarmOR™ warming sheath, two trocar swabs (5-8mm and 10-12mm) one microfiber cloth
<b>SP102</b>	InsufflatOR Tubing Set™, w/quick disconnect (Smith&Nephew)
<b>SW-PA-6V</b>	Lap Scope WarmOR™, Medical grade, reusable 6V AC transformer
<b>TM6000</b>	InsufflatOR Tubing Set™, w/funnel connector (Stryker)

\*Image shown above





3034 Owen Drive | Antioch, TN 37013 | 800.251.3000 (US) | +41 (0) 52 557 59 00  
customerservice@symmetrysurgical.com | [symmetrysurgical.com](https://www.symmetrysurgical.com)  
©2021 Symmetry Surgical Inc. All rights reserved.  
BOVIE BROCHURE AMBULATORY 06.15 .2022 v1