

 symmetry surgical® | **Bovie**®
DERMATOLOGY PRACTICE PRODUCTS



THE NEXT GENERATION OF HIGH FREQUENCY DESSICATORS AND ELECTROSURGICAL GENERATORS



A942 Bovie DERM 942

Feature: Designed to require fewer accessories for full functionality

Benefit: Saves time during procedures and requires the purchase of fewer components

Feature: Adjust power from 0 – 10 watts in 0.1 watt increments

Benefit: Higher resolution in bipolar mode

Feature: 3-Button Handpiece (included)

Benefit: Allows the surgeon to adjust the power output up and down as well as activate the unit, even when in foot control mode

Feature: Digital Error Detection

Benefit: Continuously monitors every aspect of the output and instantly disables the unit at the sign of any problem



UNSURPASSED 4-YEAR WARRANTY ON ALL BOVIE ELECTROSURGERY GENERATORS!



A952 Bovie Bantam | PRO

Feature: Designed and engineered to be the most reliable and durable generator on the market

Benefit: Unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts

Feature: 50 watts of power

Benefit: No grounding pad needed when in desiccation mode

Feature: High frequency desiccator with cutting ability

Benefit: Provides a 2-in-one generator system offering 50 watts of cutting power as well as coag capabilities

Feature: 3-Button Handpiece (included)

Benefit: Allows the surgeon to adjust the power output up and down as well as activate the unit, even when in foot control mode

Feature: Split Pad Bovie NEM Technology

Benefit: Complete patient safety to prevent burns during procedures

Perfect option for MOHS procedures



A1250S Specialist | PRO

Feature: The affordable solution for most in-office electrosurgery procedures

Benefit: The surgeon receives a much faster return on their investment

Feature: Cut and coag activation is on the handpiece

Benefit: The surgeon can move seamlessly from cut to coag without leaving the sterile field

Feature: The push-button controls allow easier use of the generator

Benefit: The push-button front panel saves time and motion in adjusting the generator. The smoothly embossed push-button front panel allows the staff to clean the generator more quickly since there are no knobs

Feature: Split Pad Bovie NEM Technology

Benefit: Complete patient safety to prevent burns during procedures

Perfect option for MOHS procedures



GENERATOR COMPARISON

FEATURES	A942	A952	A1250S
Cut (# of modes)	N/A	50W	120W (1)
Blend (# of steps)	N/A	50W	90W (1)
Coag	40W	50W	80W (1)
Fulguration/Spray Coag	40W	40W	40W
Standard Bipolar	40W	40W	30W
Pad Alarm	N/A	YES	YES
Pad Sensing (tissue impedance)	N/A	BOVIE NEM	BOVIE NEM

FOOT CONTROL / HANDPIECE ACCESSORIES	A942	A952	A1250S
Footswitch	A803	A803	BV-1253B
Pencils/handpiece	A902	A902	ESPR2

ILLUMINATING THE PATH TO BETTER HEALTHCARE

MEDICAL LIGHTING

Bovie® LED medical lights are energy efficient to provide 50,000 hours of cool, clear visibility, making them ideal for patient exam and procedure rooms. Bovie procedure lights (available in floor-, wall- and ceiling-mount options) include a sterilizable handle for dimming and illumination field size control with one hand, keeping the user in the sterile field.


Bovie continues to prioritize safety with the UV24 Air Purification systems. These units function as ceiling lights while also providing 24/7 protection to your patients and staff against airborne bacteria, viruses and mold in all areas of your facility.




*Incorporates office lighting options only.

RECOMMENDED COMPLEMENTARY PRODUCTS

Smoke Evacuation System



Diagnose and Treat



RECOMMENDED COMPLEMENTARY PRODUCTS (continued)

Pens and Sheaths

 <p>A902 Three-button power control & autoclavable handpiece</p>	 <p>ESP1 Sterile push-button pencil with blade electrode</p>	 <p>A910ST Disposable handpiece and cord sheath, sterile</p>	 <p>A910 Disposable handpiece sheath, non-sterile</p>
--	--	---	---

Pens and Sheaths

 <p>A824 Reusable bipolar forceps, 4.5", curved, micro (0.5mm) tips</p>	 <p>A827V Bipolar forceps cord, fits A3350, A2350, A1250S, A952, A942 and DERM 102</p>	 <p>A803 Footswitch for A952, A942, and DERM 102</p>	 <p>BV-1253B Monopolar footswitch for A3350, A2350 and A1250S</p>
--	--	---	---





RECOMMENDED COMPLEMENTARY PRODUCTS (continued)

Specialty Electrodes

 <p>Tungsten wire loop electrodes</p>	 <p>Microdissection needle electrodes</p>	 <p>ES61 SuperCut™ tungsten needle, modified super fine, 3cm, sterile</p>	 <p>ES20 Standard 3/16" (4.8mm) ball electrode</p>	 <p>ES02 Needle electrode, 2.75"</p>
--	---	---	--	--

Dermal Tip Electrodes

 <p>A804 (non-sterile) & A805 (sterile) Sharp dermal tip</p>	 <p>A806 (non-sterile) & A807 (sterile) Blunt dermal tip</p>	 <p>A806DE (non-sterile) & A807DE (sterile) Derm-Elite blunt dermal tip</p>	 <p>H10012 (non-sterile) & H10008 (sterile) H-type sharp tip electrode</p>	 <p>H10112 (non-sterile) & H10108 (sterile) H-type blunt tip electrode</p>
--	--	---	--	--

DERMATOLOGY INSTRUMENTS



Item Code	Procedure Instrumentation Description	Qty
01-2703	Symmetry® Knife Handle, #3, Standard, 5 in, 127 mm	1/each
01-2804	Symmetry® Knife Handle, #4, Standard, 5-1/4 in, 133 mm	1/each
01-2814	Symmetry® Knife Handle, #4L, Long, 8-3/8 in, 212 mm	1/each
02-2011	Symmetry® Scissors, Iris Standard Pattern, Straight, 4-1/2 in	1/each
02-2111	Symmetry® Scissors, Iris Standard Pattern, Curved, 4-1/2 in	1/each
02-5002	Symmetry® Scissors, Tenotomy, Straight, Sharp, 4-1/4 in	1/each
02-5003	Symmetry® Scissors, Tenotomy, Curved, Sharp, 4-1/8 in	1/each
03-0415	Symmetry® Scissors, Operating, Straight, Sharp/Sharp, 6 in	1/each
03-0613	Symmetry® Scissors, Operating, Straight, Sharp/Blunt, 5-1/8 in	1/each
03-5017	Symmetry® Scissors, Mayo, Straight, Beveled Blades, 6-3/4 in (173 mm)	1/each
03-5117	Symmetry® Scissors, Mayo, Curved, Beveled Blades, 6-3/4 in	1/each
04-3214	Symmetry® Scissors, Metzenbaum Dissecting, Straight, Standard Pattern, 5-3/4 in	1/each
04-3314	Symmetry® Scissors, Metzenbaum, Curved, Standard Pattern, 5-3/4 in	1/each
04-3318	Symmetry® Scissors, Metzenbaum, Curved, Standard Pattern, 7 in (178 mm)	1/each
06-0014	Symmetry® Forceps, Dressing, Serrated, 5-1/2 in	1/each
06-0020	Symmetry® Forceps, Dressing, Serrated, 8 in	1/each
06-0514	Symmetry® Forceps, Tissue, 1x2 Teeth, 5-3/4 in (140 mm)	1/each
06-2012	Symmetry® Forceps, Adson Dressing, Serrated Jaws, 4-3/4 in	1/each
06-2112	Symmetry® Forceps, Adson Tissue, 1x2 Teeth, 4-3/4 in	1/each
06-2113	Symmetry® Forceps, Adson Tissue, Serrated, 1x2 Teeth, 4-3/4 in	1/each
06-2512	Symmetry® Forceps, Adson-Brown Tissue, 4-3/4 in (121 mm), 7 x 7 Teeth	1/each
06-2590	Symmetry® Forceps, Callison Tying Tissue, 1.5 mm Teeth, Tying Platform, 1x2 Teeth, 4-7/8 in	1/each
06-8798	Forceps, Heaney Tissue, 2x3 Teeth, 5-1/2 in	1/each

*Image shown above

DERMATOLOGY INSTRUMENTS (continued)

Item Code	Procedure Instrumentation Description	Qty
09-8309	Symmetry [®] Suction Tube, Frazier, Medium Angle, 9 Fr, 4-1/3 in, 110 mm	1/each
10-1012	Symmetry [®] Needle Holder, Webster, Straight, Blunt, Smooth Jaws, 5-1/4 in, 133 mm	1/each
10-1212	Symmetry [®] Needle Holder, Derf, 5 in, 127 mm	1/each
10-1213*	Symmetry [®] Needle Holder, Halsey, Serrated, 5-3/8 in, 127 mm	1/each
10-1813	Symmetry [®] Needle Holder, Mayo-Hegar, 5-1/4 in, 133 mm	1/each
10-1815	Symmetry [®] Needle Holder, Mayo-Hegar, 6-2/3 in, 160 mm	1/each
10-3121	Symmetry [®] Needle Holder, Heaney, 8-1/4 in, 209 mm	1/each
14-0609	Symmetry [®] Scissors, Spencer Stitch, 3-1/2 in	1/each
15-3209	Symmetry [®] Forceps, Hartman Baby Mosquito, Straight, 3-7/8 in	1/each
15-3309	Symmetry [®] Forceps, Hartman Baby Mosquito, Curved, 4 in	1/each
15-4412	Symmetry [®] Forceps, Halsted Mosquito, Straight, 5 in	1/each
15-4512	Symmetry [®] Forceps, Halsted Mosquito, Curved, 5 in (122 mm)	1/each
15-7016	Symmetry [®] Forceps, Kelly Artery, Straight, 6-1/4 in	1/each
15-7114	Symmetry [®] Forceps, Kelly Artery, Curved, 5-1/2 in (140 mm)	1/each
15-7116	Symmetry [®] Forceps, Kelly Artery, Curved, 6-1/4 in	1/each
15-7316	Symmetry [®] Forceps, Crile Artery, Curved, 6-1/4 in	1/each
16-1016	Symmetry [®] Forceps, Rochester-Pean Hemostatic, Straight, 6-5/8 in	1/each
16-1116	Symmetry [®] Forceps, Rochester-Pean Hemostatic, Curved, 6-3/5 in	1/each
17-5508	Symmetry [®] Clamp, Backhaus Towel, 3-1/2 in (89 mm)	1/each
17-5513	Symmetry [®] Clamp, Backhaus Towel, 5-1/4 in (133 mm)	1/each
19-5024	Symmetry [®] Forceps, Foerster, Straight, Serrated, 9-3/4 in	1/each
19-5579	Classic [®] Clamp, Gemini, 5-1/2 in	1/each
20-0508	Symmetry [®] Retractor, Senn, Double-Ended, Sharp, 6-1/4 in	1/each
20-5001	Symmetry [®] Hook, Joseph Skin, Single-Prong, 6.0 mm, 6-1/4 in, 158 mm	1/each
20-5005*	Symmetry [®] Hook, Joseph Skin, Double-Prong, 5 mm Spread, 4.5 mm, 6-1/4 in, 158 mm	1/each
24-2211	Symmetry [®] Retractor, Weitlaner, Blunt, 2x3 Prong, 4-1/2 in	1/each
30-5561	Symmetry [®] Forceps, Bozeman Dressing, Curved, 260 mm	1/each
30-6556*	Symmetry [®] Forceps, Bishop Harmon Iris, Delicate	1/each
32-110	Symmetry [®] Needle Holder, Baumgartner, Tungsten Carbide, 5 in	1/each
32-405	Symmetry [®] Needle Holder, Webster, Tungsten Carbide, Large Ring Handle, 4-3/4 in	1/each
32-430	Symmetry [®] Needle Holder, Crile-Wood, 6 in, 152 mm	1/each
32-500	Symmetry [®] Forceps, Adson Dressing, Tungsten Carbide, 4-3/4 in	1/each
50-0215	Symmetry [®] Forceps, Allis Tissue, 4x5 Teeth, 6 in (152 mm)	1/each
50-0219	Symmetry [®] Forceps, Allis Tissue, 5x6 Teeth, 7-1/2 in	1/each
52-2793	Speculum, Auvard Vaginal, 1-3/4 in x 4 in Blades, 2.5 lbs	1/each
52-4426	Symmetry [®] Forceps, Braun Tenaculum, Straight, 9-1/2 in	1/each
54-2555	Symmetry [®] Scissors, Shortbent Stitch	1/each
54-4030	Symmetry [®] Scissors, Operating, Straight, 5-1/2 in	1/each
60-1813	Symmetry [®] Scissors, Stitch, Curved, Delicate, Sharp, Large Ribbon Ring Handles w/Flat Shanks, 4 in	1/each
60-1856	Symmetry [®] Scissors, Strabismus, Straight, Ribbon Handles, 4-1/4 in	1/each
60-1924	Symmetry [®] Needle Holder, Castroviejo, Straight, w/Lock, 12 mm, 5-1/2 in	1/each
60-1925	Symmetry [®] Needle Holder, Castroviejo, Curved, w/Lock, 12 mm, 5-1/2 in	1/each
67-5507*	Scissors, Ragnell, Curved, Beveled Blades, Dissecting Tips, 5 in	1/each
86-5013*	Symmetry [®] Scissors, Stevens Tenotomy, Curved, Serrated, Sc, 4 in	1/each

*Image shown in page 7





3034 Owen Drive | Antioch, TN 37013 | 800.251.3000 (US) | +41 (0) 52 557 59 00
customerservice@symmetrysurgical.com | [symmetrysurgical.com](https://www.symmetrysurgical.com)
©2021 Symmetry Surgical Inc. All rights reserved.
BOVIE BROCHURE DERM 05.05.2022 v1