



symmetry surgical®

MODULAR LAPAROSCOPIC INSTRUMENTS



CONTENTS

LAPAROSCOPIC INSTRUMENTS PORTFOLIO OVERVIEW 4 – 5

Features and Benefits 6

Modular Jaw Inserts – 5mm 7–16

- Traumatic Graspers 7
- Marylands and Right Angles 9
- Glassman and DeBakeys 10
- Fenestrated Graspers 11
- Dolphin Dissectors 12
- Reusable Scissors 12
- Blunt Graspers 13
- Alligator Forceps 13
- Babcocks 14
- Atraumatic Graspers 15
- Needle Holder 16
- Spoon 16
- Biopsy 16

Modular Jaw Inserts – 10mm 17-26

- Babcocks 17
- Right Angle Dissectors with 5mm Jaws 17
- Right Angles 18
- Glassman 18
- Anvil 18
- Traumatic Forceps 19
- Spoon 19

Free Rotation 20

Modular Handle Guide – 5mm Modular Handles 21-24

- Axial Handles 25

Modular Handle Guide – 10mm Modular Handles 26-27

CONTENTS

Modular Handle Guide – Dedicated Scissor Handles 28


- Single-use Scissor Inserts 29

Accessories 30-33

- Guide – Knot Pushers 30
- Guide – Coagulators 31-32
- Laparoscopic Probe 33

Cleaning & Instrument Care 33-34

Customer Support 35




LAPAROSCOPIC PORTFOLIO

OVERVIEW

Our portfolio of laparoscopic hand instruments provides ergonomics, versatility and durability. Our laparoscopic hand instruments help you safely and efficiently perform laparoscopic surgery. These instruments help meet your minimally invasive surgery needs with comfort and precision.

Our modular laparoscopic instrumentation offers more than 8,000 combinations to ensure you have the right instrument for every procedure. These two-piece instruments are designed to come completely apart for cleaning and can be sterilized fully assembled which removes complexity from both sterile processing and the OR.


Symmetry Surgical is proud to offer the quality you know and trust in this opportunity to simplify your laparoscopic instrument trays. We strive to enable clinicians to achieve the best surgical outcome for their patients while minimizing cost.




Symmetry Surgical
U.S. Phone 1-800-251-3000
symmetrysurgical.com
customerservice@symmetrysurgical.com

LAPAROSCOPIC PORTFOLIO

Laparoscopic Instrument Assembly



Helping to remove complexity from SPD and the Operating Room and provide some of the best Class II laparoscopic instruments.



The above is showing the instrument separate.

SYMMETRY LAPAROSCOPIC


Features & Benefits



- Made in the USA.
- Can be completely disassembled for easy cleaning.
- More than 200 inserts for custom solutions.
- Durable and ergonomic composite and Stainless-Steel handle options.

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 5mm Traumatic Graspers



5mm Allis Forcep	
Part #	Description
99-80002	Standard Length
99-80002-45	Bariatric Length



5mm Cobra Jaw Grasping Forcep	
Part #	Description
99-80034	Standard Length
99-80034-45	Bariatric Length



5mm Kocher Grasping Forcep	
Part #	Description
99-80042	Standard Length
99-80042-45	Bariatric Length



5mm Tenaculum	
Part #	Description
99-80058	Standard Length
99-80058-45	Bariatric Length




5mm D/A Traumatic Grasper, 2/3 Teeth	
Part #	Description
99-80004	Standard Length
99-80004-45	Bariatric Length



5mm S/A Traumatic Grasper, 2x3 Teeth	
Part #	Description
99-80047	Standard Length
99-80047-45	Bariatric Length



5mm Ripple Claw Forcep	
Part #	Description
99-80050	Standard Length
99-80050-45	Bariatric Length



5mm Tooth Grasper Insert with Cup	
Part #	Description
99-80061	Standard Length
99-80061-45	Bariatric Length



5mm Claw Forcep with 2x3 Teeth	
Part #	Description
99-80023	Standard Length
99-80023-45	Bariatric Length



5mm Long Allis Forcep	
Part #	Description
99-80003	Standard Length
99-80003-45	Bariatric Length



5mm Murray Grasper	
Part #	Description
99-80073	Standard Length
99-80073-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts




SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 5mm Marylands & Right Angles

MODULAR JAW INSERTS



5mm Short Curved Dissecting Forcep	
Part #	Description
99-80014	Standard Length
99-80014-45	Bariatric Length



5mm Petrovich Curved Dissector	
Part #	Description
99-80069	Standard Length
99-80069-45	Bariatric Length



5mm Curved Maryland Dissector	
Part #	Description
99-80043	Standard Length
99-80043-45	Bariatric Length



5mm Medium Curved Dissector	
Part #	Description
99-80059	Standard Length
99-80059-45	Bariatric Length



5mm Curved Dissecting Forcep	
Part #	Description
99-80013	Standard Length
99-80013-45	Bariatric Length



5mm Maryland with Cross-Serrations	
Part #	Description
99-80071	Standard Length
99-80071-45	Bariatric Length




5mm Credi Curved Dissector	
Part #	Description
99-80065	Standard Length
99-80065-45	Bariatric Length




5mm Mini Right Angle Dissector	
Part #	Description
99-80019	Standard Length
99-80019-45	Bariatric Length



5mm Mixer Dissector	
Part #	Description
99-80015	Standard Length
99-80015-45	Bariatric Length



5mm Heavy Curved Dissector	
Part #	Description
99-80056	Standard Length
99-80056-45	Bariatric Length



5mm Graves Grasper	
Part #	Description
99-80079	Standard Length
99-80079-45	Bariatric Length



5mm Angeled Dissector	
Part #	Description
99-80076	Standard Length
99-80076-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 5mm Glassman & DeBakeys



5mm Glassman 20mm Forcep	
Part #	Description
99-80055	Standard Length
99-80055-45	Bariatric Length



5mm S/A DeBakey Forcep	
Part #	Description
99-80081	Standard Length
99-80081-45	Bariatric Length



5mm Glassman 30mm Forcep	
Part #	Description
99-80062	Standard Length
99-80062-45	Bariatric Length



5mm Glassman Forcep	
Part #	Description
99-80006	Standard Length
99-80006-45	Bariatric Length




5mm Stevy Grasper	
Part #	Description
99-80001	Standard Length
99-80001-45	Bariatric Length



5mm Curved DeBakey Tissue Forcep	
Part #	Description
99-80068	Standard Length
99-80068-45	Bariatric Length




5mm DeBakey Tissue Forcep	
Part #	Description
99-80027	Standard Length
99-80027-45	Bariatric Length



5mm DeBakey 50mm Forcep	
Part #	Description
99-80080	Standard Length
99-80080-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 5mm Fenestrated Graspers



5mm Mini-Micro Dissecting Forcep	
Part #	Description
99-80051	Standard Length
99-80051-45	Bariatric Length



5mm St. Joseph Dissecting Forcep	
Part #	Description
99-80066	Standard Length
99-80066-45	Bariatric Length




5mm Strong Fenestrated Grasping Forcep	
Part #	Description
99-80008	Standard Length
99-80008-45	Bariatric Length



5mm Micro Dissecting Forcep	
Part #	Description
99-80031	Standard Length
99-80031-45	Bariatric Length




5mm Meah Grasper	
Part #	Description
99-80077	Standard Length
99-80077-45	Bariatric Length



5mm Maxi Fenestrated Grasping Forcep	
Part #	Description
99-80040	Standard Length
99-80040-45	Bariatric Length



5mm Fenestrated Grasping Forcep	
Part #	Description
99-80030	Standard Length
99-80030-45	Bariatric Length



5mm Johan Grasper	
Part #	Description
99-80053	Standard Length
99-80053-45	Bariatric Length



5mm S/A Johan Grasper	
Part #	Description
99-80054	Standard Length
99-80054-45	Bariatric Length



5mm Dorsey Fenestrated Grasper	
Part #	Description
99-80025	Standard Length
99-80025-45	Bariatric Length




5mm Dorsey 40mm Grasper	
Part #	Description
99-80057	Standard Length
99-80057-45	Bariatric Length



5mm Fenestrated Bowel Grasper	
Part #	Description
99-80063	Standard Length
99-80063-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 5mm Dolphin Dissectors




5mm Straight Mini Dissector	
Part #	Description
99-80018	Standard Length
99-80018-45	Bariatric Length




5mm Duckbill Dissector	
Part #	Description
99-80016	Standard Length
99-80016-45	Bariatric Length




5mm Dolphin Nose Grasping Forcep	
Part #	Description
99-80044	Standard Length
99-80044-45	Bariatric Length




5mm Delicate Cone Dissector	
Part #	Description
99-80021	Standard Length
99-80021-45	Bariatric Length




5mm Micro Dolphin Dissecting Forcep	
Part #	Description
99-80022	Standard Length
99-80022-45	Bariatric Length



5mm Micro Dolphin Dissector with Cup	
Part #	Description
99-80045	Standard Length
99-80045-45	Bariatric Length



5mm Dolphin Nose Dissector with Cup	
Part #	Description
99-80033	Standard Length
99-80033-45	Bariatric Length




5mm Pointed Dissector with Cup	
Part #	Description
99-80011	Standard Length
99-80011-45	Bariatric Length



5mm Straight Tapered Dissector	
Part #	Description
99-80017	Standard Length
99-80017-45	Bariatric Length


Modular Jaw Inserts – 5mm Reusable Scissors



3/4" Curved Metz	
Part #	Description
99-70001	Standard Length
99-70001-45	Bariatric Length



1/2" Delicate Curved Metz	
Part #	Description
99-70002	Standard Length
99-70002-45	Bariatric Length



Hook Scissors	
Part #	Description
99-70004	Standard Length
99-70004-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 5mm Blunt Graspers



5mm Standard Grasping Forcep	
Part #	Description
99-80029	Standard Length
99-80029-45	Bariatric Length



5mm Flat Dissector	
Part #	Description
99-80052	Standard Length
99-80052-45	Bariatric Length



5mm Wave Jaw Grasping Forcep	
Part #	Description
99-80036	Standard Length
99-80036-45	Bariatric Length



5mm Grasping Forcep with Cup	
Part #	Description
99-80049	Standard Length
99-80049-45	Bariatric Length



5mm Atraumatic Cross-Serrated Dissector	
Part #	Description
99-80067	Standard Length
99-80067-45	Bariatric Length




5mm Long Grasping Forcep	
Part #	Description
99-80039	Standard Length
99-80039-45	Bariatric Length

5mm Alligator Forceps



5mm Extracting Forcep with Reverse Teeth	
Part #	Description
99-80028	Standard Length
99-80028-45	Bariatric Length



5mm Alligator Forcep	
Part #	Description
99-80010	Standard Length
99-80010-45	Bariatric Length



5mm Retrograde Grasping Forcep	
Part #	Description
99-80041	Standard Length
99-80041-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 5mm Babcocks



5mm Pennington Forcep	
Part #	Description
99-80007	Standard Length
99-80007-45	Bariatric Length



5mm Andrew Babcock Forcep	
Part #	Description
99-80070	Standard Length
99-80070-45	Bariatric Length



5mm Babcock Grasping Forcep	
Part #	Description
99-80032	Standard Length
99-80032-45	Bariatric Length



5mm Paddle Babcock	
Part #	Description
99-80037	Standard Length
99-80037-45	Bariatric Length



5mm Endo Allis Forcep	
Part #	Description
99-80005	Standard Length
99-80005-45	Bariatric Length



5mm Round Trip Allis Grasping Forcep	
Part #	Description
99-80012	Standard Length
99-80012-45	Bariatric Length



5mm Endo Babcock Forcep	
Part #	Description
99-80026	Standard Length
99-80026-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 5mm Atraumatic Graspers



5mm Atraumatic Grasper	
Part #	Description
99-80035	Standard Length
99-80035-45	Bariatric Length



5mm S/A Atraumatic Grasper	
Part #	Description
99-80048	Standard Length
99-80048-45	Bariatric Length



5mm Short Atraumatic Grasper	
Part #	Description
99-80060	Standard Length
99-80060-45	Bariatric Length



5mm Tri-Grasper	
Part #	Description
99-80038	Standard Length
99-80038-45	Bariatric Length



5mm McKernan Tri-Grasping Fenestrated Forcep	
Part #	Description
99-80009	Standard Length
99-80009-45	Bariatric Length



5mm Atraumatic Fundus Grasper	
Part #	Description
99-80082	Standard Length
99-80082-45	Bariatric Length



5mm Oviduct Atraumatic Forcep	
Part #	Description
99-80046	Standard Length
99-80046-45	Bariatric Length



5mm Atraumatic Tube Forcep	
Part #	Description
99-80020	Standard Length
99-80020-45	Bariatric Length



5mm Holmes Atraumatic Grasper	
Part #	Description
99-80078	Standard Length
99-80078-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 5mm Needle Holder



5mm S/A Straight Needle Holder	
Part #	Description
99-80083	Standard Length
99-80083-45	Bariatric Length



5mm S/A Left Curve Needle Holder	
Part #	Description
99-80084	Standard Length
99-80084-45	Bariatric Length




5mm S/A Right Curve Needle Holder	
Part #	Description
99-80085	Standard Length
99-80085-45	Bariatric Length



5mm S/A Self-Righting Needle Holder	
Part #	Description
99-80086	Standard Length
99-80086-45	Bariatric Length



5mm S/A Super Self-Righting Needle Holder	
Part #	Description
99-80087	Standard Length
99-80087-45	Bariatric Length



5mm S/A Needle Holder	
Part #	Description
*99-80064	Standard Length
*99-80064-45	Bariatric Length

*Needs a 10mm trocar

5mm Spoon



5mm D/A Spoon Forcep	
Part #	Description
99-80074	Standard Length
99-80074-45	Bariatric Length

5mm Biopsy



5mm Biopsy Forcep with 2 Teeth	
Part #	Description
99-80024	Standard Length
99-80024-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 10mm Babcocks




10mm Babcock with Pyramid Serations	
Part #	Description
99-80104	Standard Length
99-80117-45	Bariatric Length



10mm Babcock with DeBaKey Teeth	
Part #	Description
99-80099	Standard Length
99-80113-45	Bariatric Length



10mm Babcock with Straight Serations	
Part #	Description
99-80100	Standard Length
99-80114-45	Bariatric Length



10mm Paddle Babcock	
Part #	Description
99-80101	Standard Length
99-80115-45	Bariatric Length



10mm Pennington Forcep	
Part #	Description
99-80096	Standard Length
99-80110-45	Bariatric Length


10mm Right Angle Dissectors with 5mm Jaws



5mm Right Angle Dissector	
Part #	Description
99-80094	Standard Length
99-80108-45	Bariatric Length

SYMMETRY LAPAROSCOPIC


Modular Jaw Inserts – 10mm Right Angles



10mm Mixer Dissector	
Part #	Description
99-80095	Standard Length
99-80109-45	Bariatric Length




10mm Right Angle Dissector	
Part #	Description
99-80093	Standard Length
99-80107-45	Bariatric Length




10mm Maryland	
Part #	Description
99-80103	Standard Length

10mm Glassman



10mm Glassman Forcep	
Part #	Description
99-80092	Standard Length
99-80106-45	Bariatric Length

10mm Anvil




10mm Angled Circle Grasper	
Part #	Description
99-80091	Standard Length
99-80105-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Jaw Inserts – 10mm Traumatic Forceps



10mm S/A Claw Forcep	
Part #	Description
99-80097	Standard Length
99-80111-45	Bariatric Length



10mm Tenaculum	
Part #	Description
99-80102	Standard Length
99-80116-45	Bariatric Length

10mm Spoon




10mm S/A Spoon Forcep	
Part #	Description
99-80098	Standard Length
99-80112-45	Bariatric Length



SYMMETRY LAPAROSCOPIC Modular Handle Guide – 5mm Modular Handles


5mm Modular Slide-Lock Handle, with Rotation

Part #	Description
99-90000	Standard Length
99-90000-45	Bariatric Length




5mm Modular Slide-Lock Handle, with Free Rotation

Part #	Description
99-90014	Standard Length
99-90014-45	Bariatric Length



MODULAR HANDLE GUIDE


What is Free Rotation?

Rotations  will NOT TURN when grasping tissue unless handles are marked with  indicating non-locking Free Rotation.

 means the  WILL TURN when grasping tissue.


SYMMETRY LAPAROSCOPIC

Modular Handle Guide – 5mm Modular Handles




5mm Modular Steel Handle, with Rotation

Part #	Description
99-90010	Standard Length
99-90010-45	Bariatric Length




5mm Modular Steel Handle, with Free Rotation

Part #	Description
99-90001	Standard Length
99-90001-45	Bariatric Length



5mm Modular Ratcheting Steel Handle, with Rotation

Part #	Description
99-90011	Standard Length
99-90011-45	Bariatric Length




5mm Modular Ratcheting Steel Handle, with Free Rotation

Part #	Description
99-90012	Standard Length
99-90012-45	Bariatric Length


SYMMETRY LAPAROSCOPIC

Modular Handle Guide – 5mm Modular Handles



5mm Modular Monopolar Handle, with Free Rotation

Part #	Description
99-90002	Standard Length
99-90002-45	Bariatric Length




5mm Modular Monopolar Handle, with Rotation

Part #	Description
99-90004	Standard Length
99-90004-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Handle Guide – 5mm Modular Handles



5mm Modular Monopolar Ratcheting Handle, with Free Rotation

Part #	Description
99-90003	Standard Length
99-90003-45	Bariatric Length




5mm Modular Monopolar Ratcheting Handle, with Rotation

Part #	Description
99-90007	Standard Length
99-90007-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Handle Guide – 5mm Modular Axial Handles



5mm Axial Handle with Rotation


Part #	Description
99-90017	Standard Length
99-90017-45	Bariatric Length

5mm Axial Ratcheting Handle with Rotation


Part #	Description
99-90016	Standard Length
99-90016-45	Bariatric Length

SYMMETRY LAPAROSCOPIC


Modular Handle Guide – 10mm Modular Handles



10mm Modular Slide-Lock Handle, with Rotation	
Part #	Description
99-90022	Standard Length
99-90022-45	Bariatric Length




10mm Modular Steel Handle, with Rotation	
Part #	Description
99-90020	Standard Length
99-90020-45	Bariatric Length




10mm Modular Ratcheting Steel Handle, with Rotation	
Part #	Description
99-90021	Standard Length
99-90021-45	Bariatric Length

SYMMETRY LAPAROSCOPIC

Modular Handle Guide – 10mm Modular Handles




10mm Modular Monopolar Handle, with Rotation	
Part #	Description
99-90018	Standard Length
99-90018-45	Bariatric Length





10mm Modular Monopolar Ratcheting Handle, with Rotation	
Part #	Description
99-90019	Standard Length
99-90019-45	Bariatric Length

SYMMETRY LAPAROSCOPIC


Modular Handle Guide – Dedicated Scissor Handles




5mm Modular Steel Scissors Handle (Free Rotation)	
Part #	Description
99-70010	Standard Length
99-70010-45	Bariatric Length

5mm Modular Monopolar Scissors Handle (Free Rotation)	
Part #	Description
99-70009	Standard Length
99-70009-45	Bariatric Length



Knob Color Matters




Blue handles will **ONLY** accept Symmetry Scissor Inserts




Black handles will **NOT** accept Symmetry Scissor Inserts

SYMMETRY LAPAROSCOPIC


Single-use Scissor Inserts



3/4" Curved Metz	
Part #	Description
99-70005K	Standard Length (single-use 5mm 3/4" curved metz scissors insert, sold in packages of 10)
99-70005K-45	Bariatric Length (single-use 5mm 3/4" curved metz scissors insert, sold in packages of 10)




1/2" Delicate Curved Metz	
Part #	Description
99-70006K	Standard Length (single-use 5mm 1/2" curved metz scissors insert, sold in packages of 10)
99-70006K-45	Bariatric Length (single-use 5mm 1/2" curved metz scissors insert, sold in packages of 10)



Hook Scissors	
Part #	Description
99-70008K	Standard Length (single-use 5mm hook scissors insert, sold in packages of 10)
99-70008K-45	Bariatric Length (single-use 5mm hook scissors insert, sold in packages of 10)


SYMMETRY LAPAROSCOPIC

Accessories Guide – Knot Pushers




4mm Open Ended Knot Pusher with Knurled Handle

Part #	Description
99-60002	Standard Length
99-60002L	Bariatric Length



4mm Closed End Knot Pusher

Part #	Description
99-60003	Standard Length
99-60003L	Bariatric Length




Appel Type Knot Pusher with Knurled Handle

Part #	Description
99-60004	Standard Length
99-60004L	Bariatric Length

SYMMETRY LAPAROSCOPIC

Accessories Guide – Coagulators




Blunt Needle Coagulator

Suction Part #	Non-Suction Part #	Description
99-60006	99-60012	Standard Length
99-60006L	99-60012L	Bariatric Length



Right Angle Coagulator

Suction Part #	Non-Suction Part #	Description
99-60007	99-60013	Standard Length
99-60007L	99-60013L	Bariatric Length



V Hook Coagulator


Suction Part #	Non-Suction Part #	Description
99-60008	99-60014	Standard Length
99-60008L	99-60014L	Bariatric Length

SYMMETRY LAPAROSCOPIC


Accessories Guide – Coagulators



J-Hook Coagulator		
Suction Part #	Non-Suction Part #	Description
99-60009	99-60015	Standard Length
99-60009L	99-60015L	Bariatric Length



Spatula Coagulator		
Suction Part #	Non-Suction Part #	Description
99-60010	99-60016	Standard Length
99-60010L	99-60016L	Bariatric Length



Blunt Nose Coagulator		
Suction Part #	Non-Suction Part #	Description
99-60011	99-60017	Standard Length
99-60011L	99-60017L	Bariatric Length

SYMMETRY LAPAROSCOPIC


Accessories Guide – Laparoscopic Probe




Laparoscopic Probe	
Suction Part #	Description
99-60001	Standard Length
99-60001L	Bariatric Length

Cleaning and Instrument Care

Reusable Channel Cleaning Brushes				
Part #	Description Diameter	Total Length	Sheath	Qty/Bag
248107BBG	3-5mm	70cm	1.6mm	1
248115BBG	3-5mm	115cm	1.6mm	1
248160BBG	3-5mm	160cm	1.6mm	1
248180BBG	3-5mm	230cm	1.8mm	1
248190BBG	3-5mm	120cm	1.8mm	1




Disposable Scope Channel Cleaning Brushes					
Part #	Description	Color	Total Length	Sheath Diameter	Qty/Box
248102BBG	2.5mm	Yellow	230cm	1mm	50
248103BBG	3mm	Purple	230cm	1.55mm	50
248104BBG	5mm – Double Ended	Pink	230cm	1.85mm	50
248105BBG	5mm	Pink	230cm	1.85mm	50



SYMMETRY LAPAROSCOPIC

Cleaning and Instrument Care

General Purpose Brushes				
Part #	Brush Head Diameter	Brush Total Length	Bristle Length	Qty/Bag
241007BBG	1mm (0.039")	18"	0.5"	2
241010BBG	1.6mm (0.063")	18"	0.5"	2
241015BBG	1.8mm (0.071")	18"	0.5"	2
241019BBG	2mm (0.079")	12"	1"	2
241020BBG	2.2mm (0.086")	18"	0.5"	2
241024BBG	2.5mm (0.098")	24"	2.5"	2
241025BBG	2.7mm (0.106")	18"	1"	2
241029BBG	3mm (0.118")	12"	1"	2
241030BBG	3.2mm (0.126")	18"	2"	2
241032BBG	4mm (0.158")	12"	1"	2
241033BBG	5mm (0.197")	12"	2.5"	2
241034BBG	5mm (0.197")	24"	2.5"	2
241035BBG	5mm (0.197")	16"	2.5"	2
241040BBG	6mm (0.236")	16"	2.5"	2
241041BBG	6mm (0.236")	24"	2.5"	2
241042BBG	7mm (0.275")	16"	2.5"	2
241043BBG	7mm (0.275")	24"	2.5"	2
241044BBG	10mm (0.393")	16"	2.5"	2
241045BBG	11mm (0.433")	16"	2.5"	2
241046BBG	11mm (0.433")	24"	2.5"	2
241047BBG	10mm (0.393")	24"	2.5"	2
241050BBG	12mm (0.472")	16"	2.5"	2
241055BBG	15mm (0.590")	16"	2.5"	2
241060BBG	20mm (0.787")	16"	2.5"	2



Our highly-trained customer support team is available to take your calls from 7:00 am to 5:30 pm (CST), Monday - Friday. This group is dedicated to providing a full range of support services, including:

- Order Placement
- Order Inquiries
- Delivery Information
- Product Information
- Pricing Inquiries

If we can be of assistance, please contact us using the methods below:

PHONE

US & Canada: +1 800-251-3000
International: +1 615-964-5532

EMAIL

customerservice@symmetrysurgical.com
orders@symmetrysurgical.com

WEB

www.symmetrysurgical.com

Live Chat on symmetrysurgical.com:

Click the "Can We Help You" window in the bottom corner of any page on symmetrysurgical.com to be connected directly to a member of our sales or customer service team.

We can chat in real time and answer your questions about our company and surgical instrument portfolio.

(Available Monday through Friday, 7:30 am - 5:30 pm CST)

MAILING ADDRESS

Symmetry Surgical Inc.
Attn: Customer Service
3034 Owen Drive
Antioch, TN 37013



Symmetry Surgical is now part of Aspen Surgical®



Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

M-SM-CAT-1119835 rev 2 21-SEPT-2023 ENG – US

For further information about this product or a service, please contact your local Symmetry Surgical or Aspen Surgical representative or visit one of our webpages:

Symmetry Surgical: +1 615-964-5532 (INTL) | +1 800-251-3000 (US)
customerservice@symmetrysurgical.com | [symmetrysurgical.com](https://www.symmetrysurgical.com)

Aspen Surgical: USA +1 616.698.7100 | Canada +1 888.364.7004
customerservice@aspensurgical.com | [aspensurgical.com](https://www.aspensurgical.com)

©2023 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.

©2023 Symmetry Surgical Inc. All rights reserved.