



## 15220 Lite Pouch Surgical Mask

### Intended Use

Intended to be used in environments where exposure to fluids is not a risk.

### Key Features

- Pouch style mask provides off the face design giving more room to breathe
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote secure comfortable fit
- Made in North America
- Not made with natural rubber latex

# Lite Pouch Surgical Mask

Intended to be used in environments where exposure to fluids is not a risk.



## Technical Specifications

Part No.	15220
Description	Lite Pouch Surgical Mask
Color	Blue
Fluid Resistance	N/A
Filtration - BFE	$\geq 98.8$
Filtration - PFE	$\geq 99.1$
Delta P	$\leq 2.9$
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Qty	50/Box, 6 Boxes, 300/Case

Precept® is a registered trademark of Aspen Surgical Products, Inc.

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2021 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1051529 rev 1 13-APR-2021 ENG – US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
[aspensurgical.com](http://aspensurgical.com)

