

Probe Covers, Equipment Covers, Needle Guides, and Cleaning Accessories



We're Aspen Surgical, a supplier of surgical essentials. From the precision of the finest scalpel with our Bard-Parker® brand to our Precept® PPE and Protek™ ultrasound probe covers that help fight preventable infections and cross contamination. And our WriteSite® Plus skin markers that provide an accurate and reliable way to mark the incision site before surgery.

From one trusted supplier, you get a portfolio of essential, physician preferred medical disposables, sold through an integrated channel. With strong distributor relationships and a dedicated field service team, we ensure we stay grounded in our customer's needs every step of the way.

We help provide an easier set-up, clean up, and quicker turnover in the O.R. for surgical teams and facilities so you can focus on what matters the most – your patients.

From facilities across North America, our dedicated team designs, manufactures, and delivers high-quality products to healthcare facilities around the world.

Every day we strive to make a difference. Every day we hold ourselves to a higher standard, so whenever and wherever you need us, we'll be there.

Learn more about all that Aspen Surgical has to offer in the pages of this catalog, or online at www.aspensurgical.com





Aspen Surgical Products, Inc. Corporate Headquarters located in Caledonia, Michigan.

The following brands are registered trademarks of Aspen Surgical Products, Inc:
Precept® and WriteSite®.

The following brands are trademarks of Aspen Surgical Products, Inc:
Protek™, ProV-Access™, PullUp™, and SlimFit™.

Director™, Director XT™, Retrofit™, and Retrofit XT™ are trademarks of Protek Medical Products, Inc.

Please note, not all products are available in all regions.

©2023 Aspen Surgical Products, Inc. All rights reserved.

Table of Contents

Needle Guides... 9-24

Endocavity Needle Guides	
Alpinion	9
BK	9
Carestream	9
Chison	9
ECM	9
Esaote	9&10
GE Healthcare	10&11
Hitachi	11
Mindray	11&12
Philips	12
Samsung	12&13
Shimadzu	13
Siemens	13&14
SonoScape	14&15
Sonosite	15
Supersonic	15
Terason	15
Toshiba	15
Ultrasonix	15
Zonare	15
General Purpose Needle Guides Starter Kits	
ECM	16
Esaote	16
GE Healthcare	16&17
Mindray	17
Philips	17
Samsung	17, 18, 19
Terason	19
Zonare	19
General Purpose Needle Guides	
Multi-Gauge Director Needle Guides	20
Single-Gauge Needle Guides	21
Vascular Needle Guides Starter Kits	
GE Healthcare	22
Hitachi	22
Mindray	22
Sonosite	23
Terason	23
Toshiba (Canon)	23
Zonare	23
Vascular Needle Guides	
Single-Depth Needle Guides	24
3-Depth Needle Guide Kits	24

Probe Covers... 26-33


General Purpose Probe Covers	
PullUp™ Probe Covers	26
Telescopically-Folded Probe Covers	26
Flat-Folded Probe Covers	27
Wireless Probe Cover	27
Endocavity Probe Covers	
SlimFit™ Probe Covers	28
Tapered Probe Covers	28
Rolled Latex Probe Covers	29
Surgical Probe Covers	
PullUp™ Probe Covers	30
PullUp™ Probe Cover Kits w/ Form-Fit Tip	30
Telescopically-Folded Probe Covers	31
Accordion-Folded Probe Covers	31
Transesophageal Probe Covers	
Probe Storage Covers	
Elastic Bands	33

Equipment Covers... 34-37

ClingTouch™ Monitor Covers	34
Ultrasound System Covers	34
Tablet Cover	35
Smartphone Cover	35
MRI Headset Covers	36
Biopsy Instrument Cover	36

Cleaning Accessories... 37

Disinfection Cup	37
Cleaning Brushes	37



Protek™ Endocavity Needle Guides


Endocavity needle guides ensure precise needle placement when performing a biopsy or aspiration

- Features an extra wide funnel for optimal needle entry and smooth profile to ensure ultimate patient comfort

Below you will find single-use and reusable needle guides.



Sterile, single-use example



Non-sterile, reusable example

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-5250	Needle guide	16Ga	Alpinion E3-10	Sterile	24
1-535-5144	Needle guide	16Ga	Alpinion EC3-10, EC3-10H, EC3-10T, EC3-10X, EN3-10, EV3-10, EV3-10H, EV3-10T, EV3-10X	Sterile	24
1-535-5146	3.5cm Latex biopsy kit	16Ga	Alpinion EC3-10, EC3-10H, EC3-10T, EC3-10X, EN3-10, EV3-10, EV3-10H, EV3-10T, EV3-10X	Sterile	24
1-535-5149	12" Soft biopsy kit	16Ga	Alpinion EC3-10, EC3-10H, EC3-10T, EC3-10X, EN3-10, EV3-10, EV3-10H, EV3-10T, EV3-10X	Sterile	24
1-535-5183	48" SlimFit™ biopsy kit	16Ga	Alpinion EC3-10, EC3-10H, EC3-10T, EC3-10X, EN3-10, EV3-10, EV3-10H, EV3-10T, EV3-10X	Sterile	24
1-535-5142	Needle guide	16Ga	Alpinion VE3-10H	Sterile	24
1-535-6456	Stainless Steel Needle Guide	16Ga	Alpinion EC3-10, EN3-10, EV3-10, EC3-10T, EV3-10T, EC3-10H, EV3-10H, EC3-10X, EV3-10X	Non-sterile	1
1-535-5250	Needle guide	16Ga	BK Medical 8819, E10C4	Sterile	24
1-535-5250	Needle guide	16Ga	Carestream E10C4	Sterile	24
1-535-5585	Needle guide	16Ga	Chison D6C12L, V6-A	Sterile	24
1-535-6042	Stainless steel needle guide	16Ga	Chison D6C12L	Non-sterile	1
1-535-5595	Needle guide	16Ga	Chison D7C10L, V7-E	Sterile	24
1-535-5055	Needle guide	16Ga	ECM E610P	Sterile	24
1-535-5056	3.5cm Latex biopsy kit	16Ga	ECM E610P	Sterile	24
1-535-5058	12" Soft biopsy kit	16Ga	ECM E610P	Sterile	24
1-535-5059	24" SlimFit™ biopsy kit	16Ga	ECM E610P	Sterile	24
1-535-9000	Needle guide	16Ga	Esaote E8-5 R10	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Esaote E8-5 R10	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Esaote E8-5 R10	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Esaote E8-5 R10	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Esaote E8-5 R10	Sterile	24
1-535-5145	Needle guide	16Ga	Esaote EC123, EC1123, SE3123, SE3133, TRT12, IVT12	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	Esaote EC123, EC1123, SE3123, SE3133, TRT12, IVT12	Sterile	24

Protek™ Endocavity Needle Guides

Endocavity needle guides ensure precise needle placement when performing a biopsy or aspiration

- Features an extra wide funnel for optimal needle entry and smooth profile to ensure ultimate patient comfort

Below you will find single-use and reusable needle guides.



Sterile, single-use example



Non-sterile, reusable example

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-5151	24" Soft biopsy kit	16Ga	Esate EC123, EC1123, SE3123, SE3133, TRT12, IVT12	Sterile	24
1-535-6390	Stainless steel needle guide	16Ga	Esate EC123	Non-sterile	1
1-535-6360	Stainless steel needle guide	16Ga	Esate EC1123, SE3133, E3-12	Non-sterile	1
1-535-5046	Needle guide	16Ga	GE Healthcare 618E	Sterile	25
1-535-6077	Stainless steel needle guide	16Ga	GE Healthcare 618E, TR7	Non-sterile	1
1-535-5014	24" Soft & 2.0cm Latex biopsy kit (no gel)	16Ga	GE Healthcare BE9C, BE9CS	Sterile	24
1-535-5017	24" Soft & 2.0cm Latex biopsy kit	16Ga	GE Healthcare BE9C, BE9CS	Sterile	24
1-535-6026	Stainless steel needle guide	16Ga	GE Healthcare BE9C	Non-sterile	1
1-535-6340	Stainless steel needle guide	14Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, MTZ, MZ	Non-sterile	1
1-535-6449	Stainless steel needle guide	16Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, MTZ, MZ	Non-sterile	1
1-535-5055	Needle guide	16Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, IC9B-RS, MTZ, MZ	Sterile	24
1-535-5056	3.5cm Latex biopsy kit	16Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, IC9B-RS, MTZ, MZ	Sterile	24
1-535-5058	12" Soft biopsy kit	16Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, IC9B-RS, MTZ, MZ	Sterile	24
1-535-5059	24" SlimFit™ biopsy kit	16Ga	GE Healthcare E8C, E721, E8C-RS, E72, IC5-9-D, IC5-9H, IC9B-RS, MTZ, MZ	Sterile	24
1-535-5145	Needle guide	16Ga	GE Healthcare IC5-9-RS, IC5-9W-RS	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	GE Healthcare IC5-9-RS, IC5-9W-RS	Sterile	24
1-535-5151	24" Soft biopsy kit	16Ga	GE Healthcare IC5-9-RS, IC5-9W-RS	Sterile	24
1-535-5570	Needle guide	16Ga	GE Healthcare IC9-RS	Sterile	24
1-535-6098	Stainless steel needle guide	16Ga	GE Healthcare IC9-RS	Non-sterile	1
1-535-5070	Needle guide	16Ga	GE Healthcare RIC5-9-D, RIC5-9H, RIC5-9W-D, RIC5-9W-RS, RIC6-12-D	Sterile	24

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-5071	3.5cm Latex biopsy kit	16Ga	GE Healthcare RIC5-9-D, RIC5-9H, RIC5-9W-D, RIC5-9W-RS, RIC6-12-D	Sterile	24
1-535-6097	Stainless steel needle guide	16Ga	GE Healthcare RIC5-9-D, RIC5-9A-RS, RIC5-9H, RIC5-9W-D, RIC5-9W-RS, RIC6-12-D	Non-sterile	1
1-535-5370	Needle guide	16Ga	Hitachi C41V1, EUP-V73W	Sterile	24
1-535-5371	24" SlimFit™ biopsy kit	16Ga	Hitachi C41V1, EUP-V73W	Sterile	24
1-535-5372	3.5cm Latex biopsy kit	16Ga	Hitachi C41V1, EUP-V73W	Sterile	24
1-535-5373	12" Soft biopsy kit	16Ga	Hitachi C41V1, EUP-V73W	Sterile	24
1-535-6370	Stainless steel needle guide	16Ga	Hitachi C41V1	Non-sterile	1
1-535-6381	Stainless steel needle guide (long)	14Ga	Hitachi EUP-V33, EUP-V33W, EUP-V53W, CC41R1 CC41R, CC41R1	Non-sterile	1
1-535-6382	Stainless steel needle guide (long)	16Ga	Hitachi EUP-V33, EUP-V33W, EUP-V53W, CC41R1 CC41R, CC41R1	Non-sterile	1
1-535-6384	Stainless steel needle guide (short)	16Ga	Hitachi EUP-V33, EUP-V33W, EUP-V53W, EUP-CC531, CC41R	Non-sterile	1
1-535-5380	Needle guide	16Ga	Hitachi EUP-V33, EUP-CC31, EUP-CC531, EUP-V12, EUP-V33W, EUP-V53W, EZU-PA5V, C41V, CC41R, CC41R1	Sterile	24
1-535-5381	12" Soft biopsy kit	16Ga	Hitachi EUP-V33, EUP-CC31, EUP-CC531, EUP-V12, EUP-V33W, EUP-V53W, EZU-PA5V, C41V, CC41R, CC41R1	Sterile	24
1-535-5382	3.5cm & 2.0cm Latex biopsy kit	16Ga	Hitachi EUP-V33, EUP-CC31, EUP-CC531, EUP-V12, EUP-V33W, EUP-V53W, EZU-PA5V, C41V, CC41R, CC41R1	Sterile	24
1-535-5250	Needle guide	16Ga	Hitachi UST-679	Sterile	24
1-535-6042	Stainless steel needle guide	16Ga	Mindray 65EC10HA, 65EC10EA		
1-535-9000	Needle guide	16Ga	Mindray EV6.5 MHZ/R10, 6CV1, 6CV1s, 65EC10EA, 65EC10EB, 65EC10HA, CB10-4, CB10-4e, CB10-4s, V10-4, V10-4B, V10-4BP, V10-4Bs, V10-4s, V11-3E, V11-3Hs, V11-3HU, V11-3WE, V11-3Ws	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Mindray EV6.5 MHZ/R10, 6CV1, 6CV1s, 65EC10EA, 65EC10EB, 65EC10HA, CB10-4, CB10-4e, CB10-4s, V10-4, V10-4B, V10-4BP, V10-4Bs, V10-4s, V11-3E, V11-3Hs, V11-3HU, V11-3WE, V11-3Ws	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Mindray EV6.5 MHZ/R10, 6CV1, 6CV1s, 65EC10EA, 65EC10EB, 65EC10HA, CB10-4, CB10-4e, CB10-4s, V10-4, V10-4B, V10-4BP, V10-4Bs, V10-4s, V11-3E, V11-3Hs, V11-3HU, V11-3WE, V11-3Ws	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Mindray EV6.5 MHZ/R10, 6CV1, 6CV1s, 65EC10EA, 65EC10EB, 65EC10HA, CB10-4, CB10-4e, CB10-4s, V10-4, V10-4B, V10-4BP, V10-4Bs, V10-4s, V11-3E, V11-3Hs, V11-3HU, V11-3WE, V11-3Ws	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Mindray EV6.5 MHZ/R10, 6CV1, 6CV1s, 65EC10EA, 65EC10EB, 65EC10HA, CB10-4, CB10-4e, CB10-4s, V10-4, V10-4B, V10-4BP, V10-4Bs, V10-4s, V11-3E, V11-3Hs, V11-3HU, V11-3WE, V11-3Ws	Sterile	24
1-535-5070	Needle guide	16Ga	Mindray DE10-3U, DE10-3WU	Sterile	24
1-535-5071	3.5cm Latex biopsy kit	16Ga	Mindray DE10-3U, DE10-3WU	Sterile	24

Protek™ Endocavity Needle Guides

Endocavity needle guides ensure precise needle placement when performing a biopsy or aspiration

- Features an extra wide funnel for optimal needle entry and smooth profile to ensure ultimate patient comfort

Below you will find single-use and reusable needle guides.



Sterile, single-use example



Non-sterile, reusable example

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-6097	Stainless steel needle guide	16Ga	Mindray DE10-3U, DE10-3WU	Non-sterile	1
1-535-5270	Needle guide	16Ga	Philips C8-4v, C9-3v, ClearVue™ C9-4v (on ClearVue 350™, 550™, 650™, and 850™ ultrasound systems), C10-3v	Sterile	24
1-535-5271	24" SlimFit™ biopsy kit	16Ga	Philips C8-4v, C9-3v, ClearVue™ C9-4v (on ClearVue 350™, 550™, 650™, and 850™ ultrasound systems), C10-3v	Sterile	24
1-535-5272	3.5cm Latex biopsy kit	16Ga	Philips C8-4v, C9-3v, ClearVue™ C9-4v (on ClearVue 350™, 550™, 650™, and 850™ ultrasound systems), C10-3v	Sterile	24
1-535-5144	Needle guide	16Ga	Philips C9-4ec	Sterile	24
1-535-5146	3.5cm Latex biopsy kit	16Ga	Philips C9-4ec	Sterile	24
1-535-5149	12" Soft biopsy kit	16Ga	Philips C9-4ec	Sterile	24
1-535-5183	48" SlimFit™ biopsy kit	16Ga	Philips C9-4ec	Sterile	24
1-535-5250	Needle guide	16Ga	Philips C9-4ec (on HD7 XE™ w/ 3.0.1 software and above, and HD9™ systems), BP10-5ec (on Affiniti 50™, Affiniti 70™, HD11™, HD115™, and HD7XE™ systems)	Sterile	24
1-535-5145	Needle guide	16Ga	Philips E6509	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	Philips E6509	Sterile	24
1-535-5151	24" SlimFit™ biopsy kit	16Ga	Philips E6509	Sterile	24
1-535-6390	Stainless steel needle guide	16Ga	Philips E6509	Non-sterile	1
1-535-5395	Needle guide	16Ga	Philips C9-4v (on Sparq® and Affiniti® Systems)	Sterile	24
1-535-5390	Needle guide	16Ga	Philips C9-5ec, C9-5ICT, C9-5	Sterile	24
1-535-6392	Stainless steel needle guide	16Ga	Philips C9-5ec, C9-5 ICT, 9-5 Endo 8	Non-sterile	1
1-535-5455	Needle guide	16Ga	Philips C10-4ec	Sterile	24
1-535-5154	Needle guide	16Ga	Philips V9-4v	Sterile	24
1-535-5154	Needle guide	16Ga	Samsung 3D4-9, 3DEC4-9, 3D4-9ES, 3DEC4-9ES	Sterile	24
1-535-5505	Needle guide	16Ga	Samsung 3D5-8EK, 3D5-9EK	Sterile	24
1-535-6190	Stainless steel needle guide	18Ga	Samsung 3D5-8EK, 3D5-9EK	Non-sterile	1

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-5145	Needle guide	16Ga	Samsung BCC4-9	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	Samsung BCC4-9	Sterile	24
1-535-5151	24" Soft biopsy kit	16Ga	Samsung BCC4-9	Sterile	24
1-535-5260	Needle guide	16Ga	Samsung E3-12A	Sterile	24
1-535-5261	24" SlimFit™ biopsy kit	16Ga	Samsung E3-12A	Sterile	24
1-535-5580	Needle guide (short)	16Ga	Samsung EA2-11AR, EA2-11AV	Sterile	24
1-535-5590	Needle guide (long)	16Ga	Samsung EA2-11AR, EA2-11AV	Sterile	24
1-535-6268	Stainless steel needle guide	18Ga	Samsung EA2-11AR, EA2-11AV	Non-sterile	1
1-535-5295	Needle guide	16Ga	Samsung EA2-11B	Sterile	24
1-535-5460	Needle guide	16Ga	Samsung miniER7	Sterile	24
1-535-5144	Needle guide	16Ga	Samsung EC4-9, EC4-9/10ED, EC4-9/10R, ER4-9, ER4-9/10ED, EV4-9, EV4-9/10ED, EVN4-9	Sterile	24
1-535-5146	3.5cm Latex biopsy kit	16Ga	Samsung EC4-9, EC4-9/10ED, EC4-9/10R, ER4-9, ER4-9/10ED, EV4-9, EV4-9/10ED, EVN4-9	Sterile	24
1-535-5149	12" Soft biopsy kit	16Ga	Samsung EC4-9, EC4-9/10ED, EC4-9/10R, ER4-9, ER4-9/10ED, EV4-9, EV4-9/10ED, EVN4-9	Sterile	24
1-535-5183	48" SlimFit™ biopsy kit	16Ga	Samsung EC4-9, EC4-9/10ED, EC4-9/10R, ER4-9, ER4-9/10ED, EV4-9, EV4-9/10ED, EVN4-9	Sterile	24
1-535-6042	Stainless steel needle guide	16Ga	Samsung EC4-9ES, EC4-9IS	Non-sterile	1
1-535-9000	Needle guide	16Ga	Samsung EC4-9ES, EC4-9IS, ER4-9ES, 10R EV/140D	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Samsung EC4-9ES, EC4-9IS, ER4-9ES, 10R EV/140D	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Samsung EC4-9ES, EC4-9IS, ER4-9ES, 10R EV/140D	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Samsung EC4-9ES, EC4-9IS, ER4-9ES, 10R EV/140D	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Samsung EC4-9ES, EC4-9IS, ER4-9ES, 10R EV/140D	Sterile	24
1-535-5275	Needle guide	16Ga	Samsung EV2-10A	Sterile	24
1-535-6267	Stainless steel needle guide	18Ga	Samsung EV2-10A	Non-sterile	1
1-535-6455	Stainless steel needle guide	16Ga	Samsung EV4-9/10ED, ER4-9/10ED, EC4-9, EVN4-9, EV4-9, ER4-9	Non-sterile	1
1-535-5250	Needle guide	16Ga	Samsung NER4-9ES, NEV4-9ES	Sterile	24
1-535-6097	Stainless steel needle guide	16Ga	Samsung V5-9	Non-sterile	1
1-535-5070	Needle guide	16Ga	Samsung V5-9, EV3-10B	Sterile	24
1-535-5071	3.5cm Latex biopsy kit	16Ga	Samsung V5-9, EV3-10B	Sterile	24
1-535-5055	Needle guide	16Ga	Samsung VR5-9	Sterile	24
1-535-5056	3.5cm Latex biopsy kit	16Ga	Samsung VR5-9	Sterile	24
1-535-5058	12" Soft biopsy kit	16Ga	Samsung VR5-9	Sterile	24
1-535-5059	24" SlimFit™ biopsy kit	16Ga	Samsung VR5-9	Sterile	24
1-535-6340	Stainless steel needle guide	14Ga	Samsung VR5-9	Non-sterile	1
1-535-6449	Stainless steel needle guide	16Ga	Samsung VR5-9	Non-sterile	1
1-535-5045	Needle guide	16Ga	Shimadzu TV11R-055U	Sterile	24
1-535-6005	Stainless steel needle guide	16Ga	Shimadzu TV11R-055U	Non-sterile	1
1-535-5145	Needle guide	16Ga	Siemens BE9-4	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	Siemens BE9-4	Sterile	24
1-535-5151	24" Soft biopsy kit	16Ga	Siemens BE9-4	Sterile	24

Protek™ Endocavity Needle Guides

Endocavity needle guides ensure precise needle placement when performing a biopsy or aspiration

- Features an extra wide funnel for optimal needle entry and smooth profile to ensure ultimate patient comfort

Below you will find single-use and reusable needle guides.



Sterile, single-use example



Non-sterile, reusable example

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-9000	Needle guide	16Ga	Siemens EC7, EC9-4 (on Sonoline G20™, G50S™, Acuson X500™, and X100™, systems) EV9-4 (on Acuson X500™, and X100™ systems), 6.5 EC10	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Siemens EC7, EC9-4 (on Sonoline G20™, G50S™, Acuson X500™, and X100™, systems) EV9-4 (on Acuson X500™, and X100™ systems), 6.5 EC10	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Siemens EC7, EC9-4 (on Sonoline G20™, G50S™, Acuson X500™, and X100™, systems) EV9-4 (on Acuson X500™, and X100™ systems), 6.5 EC10	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Siemens EC7, EC9-4 (on Sonoline G20™, G50S™, Acuson X500™, and X100™, systems) EV9-4 (on Acuson X500™, and X100™ systems), 6.5 EC10	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Siemens EC7, EC9-4 (on Sonoline G20™, G50S™, Acuson X500™, and X100™, systems) EV9-4 (on Acuson X500™, and X100™ systems), 6.5 EC10	Sterile	24
1-535-6042	Stainless steel needle guide	16Ga	Siemens EC7, EC9-4 (Acuson X500™ and Sonoline G20™, G50™, G60S™), EC-10C5	Non-sterile	1
1-535-6477	Stainless steel needle guide	14Ga	Siemens EC9-4 (on Acuson X150™, X300™, and Sonoline G40™)	Non-sterile	1
1-535-5480	Needle guide	16Ga	Siemens EC9-4 (on Acuson NX2™, NX3™, NX3 Elite™, P500™, X150™, X300™, and Sonoline G40™ systems), EV9-4 (on Acuson X150™, X300™, and Sonoline G40™ systems)	Sterile	24
1-535-5531	Needle guide	16Ga	Siemens EC9-4 (on Acuson S2000™ and Sonoline Antares™ systems)	Sterile	24
1-535-5535	24" SlimFit™ biopsy kit	16Ga	Siemens EC9-4 (on Acuson S2000™ and Sonoline Antares™ systems)	Sterile	24
1-535-5154	Needle guide	16Ga	Siemens EV9F3	Sterile	24
1-535-6482	Stainless steel needle guide	16Ga	Siemens P2	Non-sterile	1
1-535-5172	Needle guide	16Ga	SonoScape 6V1	Sterile	24
1-535-5178	3.5cm & 2.0cm Latex biopsy kit	16Ga	SonoScape 6V1	Sterile	24
1-535-5179	12" Soft biopsy kit	16Ga	SonoScape 6V1	Sterile	24

Part No.	Description	Gauge	For use with	Sterility	Qty
1-535-5380	Needle guide	16Ga	SonoScape 6V3 (2015 model and older)	Sterile	24
1-535-5381	12" Soft biopsy kit	16Ga	SonoScape 6V3 (2015 model and older)	Sterile	24
1-535-5382	3.5cm & 2.0cm Latex biopsy kit	16Ga	SonoScape 6V3 (2015 model and older)	Sterile	24
1-535-5665	Needle guide	16Ga	SonoScape 6V3 (2016 model and newer)	Sterile	24
1-535-5171	Needle guide	16Ga	SonoScape 6V4	Sterile	24
1-535-5650	Needle guide	16Ga	SonoScape EC9-5	Sterile	24
1-535-5651	24" SlimFit™ biopsy kit	16Ga	SonoScape EC9-5	Sterile	24
1-535-5652	3.5cm Latex biopsy kit	16Ga	SonoScape EC9-5	Sterile	24
1-535-6392	Stainless steel needle guide	16Ga	Sonosite C8 Series	Non-sterile	1
1-535-5575	Needle guide	16Ga	Sonosite ICT, ICTe, ICTx, ICTxp	Sterile	24
1-535-9000	Needle guide	16Ga	Supersonic Imagine SE12-3, SEV12-3	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Supersonic Imagine SE12-3, SEV12-3	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Supersonic Imagine SE12-3, SEV12-3	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Supersonic Imagine SE12-3, SEV12-3	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Supersonic Imagine SE12-3, SEV12-3	Sterile	24
1-535-5145	Needle guide	16Ga	Terason 8EC4	Sterile	24
1-535-5147	2.0cm Latex biopsy kit	16Ga	Terason 8EC4	Sterile	24
1-535-5151	24" Soft biopsy kit	16Ga	Terason 8EC4	Sterile	24
1-535-9000	Needle guide	16Ga	Terason 8EC4(V)	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Terason 8EC4(V)	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Terason 8EC4(V)	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Terason 8EC4(V)	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Terason 8EC4(V)	Sterile	24
1-535-5531	Needle guide	16Ga	Toshiba (Canon) PVF-651VT, PVF-641VT, PVF-621VT, PVG-601V, PVM-651VT, PVN-661VT, PVT-661VT, PVQ-641V	Sterile	24
1-535-5535	24" SlimFit™ biopsy kit	16Ga	Toshiba (Canon) PVF-651VT, PVF-641VT, PVF-621VT, PVG-601V, PVM-651VT, PVN-661VT, PVT-661VT, PVQ-641V	Sterile	24
1-535-5560	Needle guide	16Ga	Toshiba (Canon) PVT-781VT, PVU-781VT, 11C3	Sterile	24
1-535-5561	24" SlimFit™ biopsy kit	16Ga	Toshiba (Canon) PVT-781VT, PVU-781VT, 11C3	Sterile	24
1-535-5562	3.5cm Latex biopsy kit	16Ga	Toshiba (Canon) PVT-781VT, PVU-781VT, 11C3	Sterile	24
1-535-5144	Needle guide	16Ga	Ultrasonix EC9-5/10	Sterile	24
1-535-5146	3.5cm Latex biopsy kit	16Ga	Ultrasonix EC9-5/10	Sterile	24
1-535-5149	12" Soft biopsy kit	16Ga	Ultrasonix EC9-5/10	Sterile	24
1-535-5183	48" SlimFit™ biopsy kit	16Ga	Ultrasonix EC9-5/10	Sterile	24
1-535-6456	Stainless steel needle guide	16Ga	Ultrasonix EC9-5/10	Non-sterile	1
1-535-9000	Needle guide	16Ga	Zonare E9-3	Sterile	24
1-535-9001	3.5cm Latex biopsy kit	16Ga	Zonare E9-3	Sterile	24
1-535-9002	12" Soft biopsy kit	16Ga	Zonare E9-3	Sterile	24
1-535-9003	2.6cm Latex biopsy kit	16Ga	Zonare E9-3	Sterile	24
1-535-9005	24" SlimFit™ biopsy kit	16Ga	Zonare E9-3	Sterile	24
1-535-7411	12" Soft biopsy kit	16Ga	Zonare E9-4	Sterile	24

Protek™ General Purpose Needle Guide Starter Kits

General purpose (in-plane) needle guide starter kits contain the necessary components for a safe and precise biopsy. The custom fit single or multi-angle bracket is designed to meet predetermined angles defined by ultrasound system software guidelines. Single-use biopsy needle guides easily snap onto the bracket providing a fixed path for the needle. PullUp™ probe cover kits help to maintain procedure sterility and reduce cross-contamination.

- Please note, kits have both Sterile and Non-sterile components unless otherwise noted on product label

Examples of Protek™ Needle Guide Starter Kits



Multi-angle example



Single-angle example

Part No.	For use with	Bracket Angles	Starter Kit Content
1-535-7138	ECM C360A, C360V	Single-angle	(1) Non-sterile reusable bracket; (6) Sterile Director™ needle guide/ Pullup™ probe cover kits (contains guide (2 ea. 16Ga, 18Ga, 22Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands)
1-535-7363	Esaote C1-8, AC2541	Multi-angle	(1) Non-sterile reusable bracket (5) Sterile Director XT™ multi-gauge needle guide/PullUp™ probe cover kits (Contains: multi-gauge needle guide set (8.5F, 14-18,20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel and bands
1-535-7070	GE Healthcare 4C, 4C-A, 4C-D, 4C-RS, 4C-SC, 5CS	Multi-angle	(1) Non-sterile reusable bracket; (6) Sterile Retrofit™ single-use needle guides (2 ea. 18Ga, 20Ga, 22Ga); (6) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7068	GE Healthcare 9L, 9L-D, 9L-RS	Multi-angle	(1) Non-sterile reusable bracket; (6) Sterile Retrofit™ single-use needle guides (2 ea. 18Ga, 20Ga, 22Ga); (6) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7308	GE Healthcare C1-5-D, C1-5-RS	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Retrofit™ XT™ multi-gauge needle guide sets (8.5F, 14-18,20-23Ga); (5) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7309	GE Healthcare C1-6, C1-6-D	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Retrofit™ XT™ multi-gauge needle guide sets (8.5F, 14-18,20-23Ga); (5) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7306	GE Healthcare C2-9-D, C2-9VN-D	Multi-angle	(1) Non-sterile reusable bracket (5) Sterile Retrofit XT™ multi-gauge needle guide/PullUp™ probe cover kits (Contains: multi-gauge needle guide set (8.5F, 14-18,20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel and bands

Part No.	For use with	Bracket Angles	Starter Kit Content
1-535-7107	GE Healthcare ML6-15, ML6-15-D, ML6-15-RS	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Retrofit XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands
1-535-7278	Mindray C5-1, C5-1E, C5-1s, C5-1U, SC6-1E, SC6-1s, SC6-1U, SC5-1E, SC5-1N, SC5-1Ns	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Retrofit XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands
1-535-7293	Mindray L14-5w, L14-5WE, L14-5WU,	Multi-Angle	(1) Non-sterile reusable bracket (5) Sterile Director XT™ multi-gauge needle guide/PullUp™ probe cover kits (Contains: multi-gauge needle guide set (8.5F, 14-18,20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel and bands
1-535-7294	Mindray L14-3WU, L14-3Ws	Multi-angle	(1) Non-sterile reusable bracket (5) Sterile Director XT™ multi-gauge needle guide/PullUp™ probe cover kits (Contains: multi-gauge needle guide set (8.5F, 14-18,20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel and bands
1-535-7385	Philips C5-2 (HD3 Ultrasound System)	Single-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ needle guide/ Pullup™ probe cover kits (contains guide (8 ea. 16Ga, 18Ga, 22Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands)
1-535-7154	Samsung C2-5	Single-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 20Ga, 22Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7457	Samsung C2-8	Single-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ needle guide/ Pullup™ probe cover kits (contains guide (8 ea. 16Ga, 18Ga, 22Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands)
1-535-7180	Samsung CA1-7A, CA1-7AD	Multi-angle	(1) Non-sterile reusable bracket; (6) Sterile Director™ single-use needle guides (2 ea. 16Ga, 18Ga, 22Ga); (6) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7167	Samsung CA1-7A, CA1-7AD	Multi-angle	(1) Non-sterile reusable bracket; (36) Sterile Director™ single-use needle guides (6 ea. 14Ga, 16Ga, 18Ga, 20Ga, 22Ga, 25Ga); (36) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7196	Samsung CA1-7S	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Director XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands
1-535-7179	Samsung CF2-8, CA2-8A, CA2-8AD, CA2-9A, CA2-9AD	Multi-angle	(1) Non-sterile reusable bracket; (6) Sterile Director™ single-use needle guides (2 ea. 16Ga, 18Ga, 22Ga); (6) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7166	Samsung CF2-8, CA2-8A, CA2-8AD, CA2-9A, CA2-9AD	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 22Ga, 25Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7265	Samsung CF4-9	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Director XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands

Protek™ General Purpose Needle Guide Starter Kits

General purpose (in-plane) needle guide starter kits contain the necessary components for a safe and precise biopsy. The custom fit single or multi-angle bracket is designed to meet predetermined angles defined by ultrasound system software guidelines. Single-use biopsy needle guides easily snap onto the bracket providing a fixed path for the needle. PullUp™ probe cover kits help to maintain procedure sterility and reduce cross-contamination.

- Please note, kits have both Sterile and Non-sterile components unless otherwise noted on product label

Examples of Protek™ Needle Guide Starter Kits




Multi-angle example



Single-angle example

Part No.	For use with	Bracket Angles	Starter Kit Content
1-535-7168	Samsung CV1-8A, CV1-8AD	Multi-angle	(1) Non-sterile reusable bracket; (36) Sterile Director™ single-use needle guides (6 ea. 14Ga, 16Ga, 18Ga, 20Ga, 22Ga, 25Ga); (36) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7162	Samsung L3-12A, LA3-14A, LA3-14AD, L5-13/50	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 20Ga, 22Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7181	Samsung LA2-9A	Multi-angle	(1) Non-sterile reusable bracket; (6) Sterile Retrofit™ single-use needle guides (2 ea. 18Ga, 20Ga, 22Ga); (6) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7285	Samsung LA4-18A, LA4-18AD	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Director XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands
1-535-7197	Samsung LA4-18B, LA4-18BD	Multi-angle	(1) Non-sterile reusable bracket; (5) Sterile Director XT™ needle guide/PullUp™ probe cover kits (contains multi-gauge needle guide set (8.5F, 14-18Ga, 20-23Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands
1-535-7165	Samsung LF5-13, LA3-16A	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 20Ga, 22Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands

Part No.	For use with	Bracket Angles	Starter Kit Content
1-535-7172	Samsung LF5-13, LA3-16A, LA3-16AD	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 20Ga, 22Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7461	Samsung LN5-12, LE5-12	Single-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ needle guide/ Pullup™ probe cover kits (contains guide (8 ea. 16Ga, 18Ga, 22Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands)
1-535-7164	Samsung SC1-6	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 16Ga, 18Ga, 22Ga, 25Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7250	Samsung VN4-8	Multi-angle	(1) Non-sterile reusable bracket; (24) Sterile Director™ single-use needle guides (6 ea. 20Ga, 21Ga, 22Ga, 25Ga); (24) Sterile 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover kits with sterile gel and bands
1-535-7138	Terason 4C2-V, 5C2, 5C2A	Single-angle	(1) Non-sterile reusable bracket; (6) Sterile Director™ needle guide/ PullUp™ probe cover kits (contains guide (2 ea. 16Ga, 18Ga, 22Ga), 5.5" x 48" (14cm x 122cm) soft cover, sterile gel, and bands)
1-535-7024	Zonare C5-2, C6-2	Single-angle	(1) Non-sterile reusable bracket; (6) Sterile Director™ single-use needle guides (2 ea. 16Ga, 18Ga, 22Ga); (6) Sterile 5.5" x 24" (14cm x 61cm) soft 3D ended probe cover kits with sterile gel and bands



Protek™ Multi-Gauge Needle Guides

Product Details

- Easily snaps onto bracket for precise needle guidance
- Clear, wide funnel to visualize needle entry
- Quick-release for easy needle removal
- Sterile, single-use
- Retrofit XT™ fits GE Healthcare brackets, Director XT™ fits all other OEM brackets
- Set includes needle guide with 10 optional needle gauge inserts (8.5F, 14-18Ga, 20-23Ga)

Part No.	Description	Qty
1-535-4500	Director XT™ Needle Guide Set	24
1-535-5600	Retrofit XT™ Needle Guide Set	24



Protek™ Multi-Gauge Needle Guide Probe Cover Kits

Product Details

- Easily snaps onto bracket for precise needle guidance
- Clear, wide funnel to visualize needle entry
- Quick-release for easy needle removal
- 10 optional needle gauge inserts (8.5F, 14-18Ga, 20-23Ga)
- Retrofit XT™ fits GE Healthcare brackets, Director XT™ fits all other OEM brackets
- Sterile kit contains single-use needle guide with multi-gauge insert set, 5.5" x 48" (14cm x 122cm) soft PullUp™ probe cover, sterile gel, and 2 bands

Part No.	Description	Qty
1-535-4501	Director XT™ Needle Guide Probe Cover Kit	24
1-535-5601	Retrofit XT™ Needle Guide Probe Cover Kit	24



Protek™ Single-Gauge Needle Guides

Product Details

- Easily snaps onto bracket for precise needle guidance
- Clear, wide funnel to visualize needle entry
- Quick-release for easy needle removal
- Color-coded gauges
- Sterile, single-use

Director™		
Part No.	Gauge	Qty
1-535-4013	13Ga	24
1-535-4014	14Ga	24
1-535-4015	15Ga	24
1-535-4016	16Ga	24
1-535-4018	18Ga	24
1-535-4020	20Ga	24
1-535-4021	21Ga	24
1-535-4022	22Ga	24
1-535-4025	25Ga	24



Protek™ Single-Gauge Needle Guide and Probe Cover Kits

Product Details

- Easily snaps onto bracket for precise needle guidance
- Clear, wide funnel to visualize needle entry
- Quick-release for easy needle removal
- Color-coded gauges
- Contains single-use needle guide, 5.5" x 48" (14cm x 122cm) Soft PullUp™ probe cover, sterile gel, and 2 bands

Director™		
Part No.	Gauge	Qty
1-535-4216	16Ga	24
1-535-4218	18Ga	24
1-535-4222	22Ga	24

Protek™ Vascular Needle Guide Starter Kit

Vascular (out-of-plane) needle guide starter kits contain the necessary components for a safe and precise biopsy. The custom fit bracket securely attaches to the probe with ease. Single-use, depth specific, ProV-Access™ needle guides snap onto the bracket providing a fixed path for the needle. PullUp™ probe cover kits help to maintain procedure sterility and reduce cross-contamination.


- Please note, kits have both Sterile and Non-sterile components unless otherwise noted on product label

Example of Protek™ Vascular Needle Guide Starter Kit



Part No.	For use with	Starter Kit Content
1-535-7311	GE Healthcare 4C, 4C-A, 4C-D, 4C-RS, 4C-SC, 5CS	(1) Non-sterile reusable vascular bracket; (3) Sterile 1.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7397	Hitachi L64	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7292	Mindray C4-1	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7276	Mindray C5-2E, C5-2s	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7277	Mindray L12-4s, 7L4s, L14-6Ns, L14-6NE, L12-3E	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7273	Mindray L14-5w, L14-5WE, L14-5WU	(1) Non-sterile reusable vascular bracket; (3) Sterile 1.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7274	Mindray L9-3U, L9-3s, and L9-3E	(1) Non-sterile reusable vascular bracket; (3) Sterile 1.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7279	Mindray SC8-2s, SC8-2U, SC8-2Us	(1) Non-sterile reusable vascular bracket; (3) Sterile 1.5cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands

Part No.	For use with	Starter Kit Content
1-535-7419	Sonosite C60, C60x	(1) Non-sterile reusable vascular bracket; (3) Sterile 1.5cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7492	Terason 15L4A, 15L4 Smart Mark™	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7536	Toshiba (Canon) PLI-1205BX (i18LX5)	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; Fits needle gauges 18 through 21; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7292	Zonare C4-1	(1) Non-sterile reusable vascular bracket; (3) Sterile 2.0cm ProV-Access™ single-use needle guides; (3) Sterile 5" x 48" (13cm x 122cm) clear PullUp™ probe cover kits with sterile gel and bands
1-535-7296	Zonare L10-5	(1) Non-sterile reusable vascular bracket; (3) Sterile ProV-Access™ 18-21 Ga Disposable Needle Guide/PullUp™ probe cover kits (contains 3-depth needle guide set (2.0cm, 1.5cm, 1.0cm), 5" x 48" (13cm x 122cm) clear probe cover, sterile gel & bands)



- Attaches securely to bracket or Bard Site-Rite™ probe
- Clear pad to visualize needle entry
- Fits needle gauges 18 through 21
- Patented needle release for tilt out procedures
- Color-coded depths
- Sterile, single-use

Part No.	Depth	Sterility	Qty
1-535-5404	0.5cm (yellow)	Sterile	10
1-535-5403	1.0cm (blue)	Sterile	10
1-535-5402	1.5cm (white)	Sterile	10
1-535-5405	2.0cm (green)	Sterile	10
1-535-5406	2.5cm (pink)	Sterile	10
1-535-5401	3.5cm (red)	Sterile	10



- ### Product Details
- Easily snaps onto bracket or probe
 - Clear pad to visualize needle entry
 - Fits needle gauges 18 through 21
 - Patented needle release for tilt out procedures
 - Color-coded depths
 - Sterile, single-use
 - Kit contains 3-depth needle guide set, 5" x 48" (13cm x 122cm) clear PullUp™ probe cover, sterile gel, and 2 bands

Part No.	Depths	Sterility	Qty
1-535-5414	2.0cm, 1.5cm, 1.0cm	Sterile	10
1-535-5410	3.5cm, 1.5cm, 1.0cm	Sterile	10
1-535-5416	3.5cm, 2.5cm, 1.5cm	Sterile	10





Protek™ PullUp™ Probe Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Patented PullUp™ for even faster application
- Telescopic-fold ensures easy gel and cover application
- Available in clear polyethylene or soft polyurethane (not made with natural rubber latex)
- Select covers offer a 3D box end for seamless imaging
- Kit includes PullUp™ probe cover, sterile gel and 2 bands
- Sterile, intended for single-use

Part No.	Size	3D End	Material	Qty
1-519-2231	4" x 24" (10cm x 61cm)	No	Soft	24
1-519-2044	5.5" x 24" (14cm x 61cm)	No	Soft	24
1-519-2099	5.5" x 24" (14cm x 61cm)	Yes	Soft	24
1-519-2048	5.5" x 48" (14cm x 122cm)	No	Soft	24
1-519-2097	5.5" x 48" (14cm x 122cm)	Yes	Soft	24
1-519-1096	5" x 48" (13cm x 122cm)	No	Clear	24



Protek™ Telescopically-Folded Probe Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Telescopic-fold ensures easy gel and cover application
- Available in clear polyethylene or soft polyurethane (not made with natural rubber latex)
- Select covers offer a 3D box end for seamless imaging
- Kit includes probe cover, sterile gel and 2 bands
- Sterile, intended for single-use

Part No.	Size	3D End	Material	Qty
1-519-2320	5.5" x 36" (14cm x 91cm)	Yes	Soft	24
1-519-2450	5.5" x 48" (14cm x 122cm)	Yes	Soft	24
1-519-2210	6" x 58" (15cm x 147cm)	Yes	Soft	24
1-519-1236	6" x 48" (15cm x 122cm)*	Yes	Clear	100
1-519-1460	6" x 48" (15cm x 122cm)	No	Clear	24
1-519-1090	6" x 48" (15cm x 122cm)**	Yes	Clear	24
1-519-2648	6" x 48" (15cm x 122cm)**	Yes	Soft	24
1-519-2458	4" x 58" (10cm x 147cm)	Yes	Soft	24
1-519-2589	5" x 48" (13cm x 122cm)**	Yes	Soft	24
1-519-2558	5.5" x 58" (14cm x 147cm)**	Yes	Soft	24

*Kit does not contain sterile gel

**Kit includes 2 tapes and cover with arrow



Protek™ Flat-Folded Probe Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Available in clear polyethylene or soft polyurethane (not made with natural rubber latex)
- Select covers offer a 3D box end for seamless imaging
- Kit includes probe cover, sterile gel and 2 bands
- Sterile, intended for single-use

Part No.	Size	3D End	Material	Qty
1-519-2230	3.5" x 24" (9cm x 61cm)	No	Soft	24
1-519-2236	3.5" x 36" (9cm x 91cm)	No	Soft	24
1-519-2250	5.5" x 24" (14cm x 61cm)	Yes	Soft	24
1-519-2270	7" x 24" (18cm x 61cm)	No	Soft	24
1-519-2290	9" x 24" (23 x 61cm)	No	Soft	24
1-519-1240	4" x 24" (10cm x 61cm)	Yes	Clear	24
1-519-1270	7" x 24" (18cm x 61cm)	No	Clear	24



Protek™ Wireless Probe Cover

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Removable patented PullUp™ with adhesive seal ensures device is fully enclosed
- Telescopic-fold ensures easy gel and cover application
- Clear polyethylene (not made with natural rubber latex)
- Kit includes probe cover, sterile gel and 2 bands
- Sterile, intended for single-use

Part No.	Size	Qty
1-519-1512	5" x 12" (13cm x 30cm)	24



Protek™ SlimFit™ Probe Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Patented PullUp™ for even faster application
- Telescopic-fold ensures easy gel and cover application
- Soft polyurethane (not made with natural rubber latex)
- Tapered end for a better fit
- Sterile kit includes SlimFit™ probe cover, sterile gel and 2 bands
- Non-sterile kits include SlimFit™ probe cover and bands.
- Intended for single-use

Part No.	Size	Sterility	Qty
1-519-2512	3.25" – 1.25" x 12" (8cm-3cm x 30cm)	Sterile	24
1-519-2513	3.25-1.25" x 12" (8cm-3cm x 31cm)*	Sterile	200
1-519-2517	3.25" – 1.25" x 12" (8cm-3cm x 30cm)	Non-sterile	24
1-519-2524	4" – 1.25" x 24" (10cm-3cm x 61cm)	Sterile	24
1-519-2548	4" – 1.25" x 48" (10cm-3cm x 122cm)	Sterile	24

*Probe cover only




Protek™ Tapered Probe Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Soft polyurethane (not made with natural rubber latex)
- Tapered end for a better fit
- Sterile kit includes probe cover, sterile gel and 2 bands
- Non-sterile kit includes probe cover and bands
- Intended for single-use

Part No.	Size	Sterility	Qty
1-519-2220	5" – 2" x 24" (13cm-5cm x 61cm)	Sterile	24
1-519-2225	5" – 2" x 24" (13cm-5cm x 61cm)	Non-sterile	24



Protek™ Rolled Latex Probe Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Made with natural rubber latex
- Length 20cm (7.9")
- Intended for single-use

Part No.	Width	Kit Content	Sterility	Qty
1-519-3249	3.5cm (1.4")	Probe cover, sterile gel, and 2 bands	Sterile	40
1-519-3230	3.5cm (1.4")	Probe cover and 4 bands	Sterile	40
1-519-3235	3.5cm (1.4")	Probe cover and 4 bands	Non-sterile	40
1-519-3241	3.5cm (1.4")	Probe cover and 4 bands	Non-sterile	50
1-519-3238	3.5cm (1.4")	Bag containing 100 probe covers	Non-sterile	1
1-519-3233	2cm (0.8")	Probe cover and 4 bands	Sterile	24
1-519-3232	2cm (0.8")	Probe cover and 4 bands	Non-sterile	24



Protek™ PullUp™ Probe Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Patented PullUp™ for even faster application
- Telescopic-fold ensures easy gel and cover application
- Available in clear polyethylene or soft polyurethane (not made with natural rubber latex)
- Select covers offer a 3D box end for seamless imaging
- Kit includes probe cover, sterile gel, bands, and/or tape
- Sterile, intended for single-use

Part No.	Size	3D End	Material	Qty
1-519-1938	5" – 3" x 96" (13cm-8cm x 244cm)*	No	Clear	24
1-519-2233	5" x 96" (13cm x 244cm)	Yes	Soft	24

*Kit includes 5.5" x 12" polyurethane boot



Protek™ PullUp™ Probe Covers w/Form-Fit Tip

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Patented PullUp™ for even faster application
- Telescopic-fold ensures easy gel and cover application
- Clear polyethylene with soft polyurethane tip (not made with natural rubber latex)
- Form-fit tip ensures a better fit
- Sterile, intended for single-use

Part No.	Size	Kit Content	Qty
1-519-2016	5" – 0.39" x 96" (13cm-1cm x 244cm)	Probe cover, sterile gel, and 2 clips	10
1-519-2002	5" – 7⁄8" x 96" (13cm-2cm x 244cm)	Probe cover and 2 clips	10
1-519-2014	5" – 7⁄8" x 96" (13cm-2cm x 244cm)	Probe cover and 2 tapes	10
1-519-2087	5" – 1¾" x 96" (13cm-4cm x 244cm)	Probe cover, sterile gel, and 2 bands	10



Protek™ Telescopically-Folded Probe Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Telescopic-fold ensures easy gel and cover application
- Available in clear polyethylene or soft polyurethane (not made with natural rubber latex)
- Select covers offer a 3D box end for seamless imaging
- Kit includes probe cover, sterile gel, bands, and/or tape
- Sterile, intended for single-use

Part No.	Size	3D End	Material	Qty
1-519-1026	5.5" x 96" (14cm x 244cm)	No	Clear	24
1-519-1088	5" – 2¼" x 96" (13cm – 6cm x 244cm)	Yes	Clear	24
1-519-1092	6" – 3⅝" x 96" (15cm – 9cm x 244cm)	Yes	Clear	24
1-519-1237	6" – 3" x 96" (15cm-8cm x 244cm)**	No	Clear	100
1-519-1934	6" – 3" x 96" (15cm-8cm x 244cm)*	No	Clear	24
1-519-2696	6" x 96" (15cm x 244cm)***	Yes	Soft	24
1-519-2596	5" x 96" (13cm x 244cm)***	Yes	Soft	24

*Kit includes 5.5" x 12" polyurethane boot

**Gel not included

***Kit includes 2 tapes and cover with arrow



Protek™ Accordion-Folded Probe Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Clear polyethylene (not made with natural rubber latex)
- Size: 6" – 3" x 96" (15cm – 8cm x 244cm)
- Sterile, intended for single-use

Part No.	Kit Content	Qty
1-519-1930	Probe cover, boot, sterile gel, and 4 bands	12
1-519-1931	Probe cover, boot, and 4 bands	12
1-519-1933	Probe cover and 4 bands	12



Protek™ Transesophageal Probe Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Soft polyurethane (not made with natural rubber latex)
- Tapered end for a better fit
- Applicator tube provides access to base of cover for easy, mess-free application of gel
- Size: 1.75" – 1" x 48" (4cm-3cm x 122cm)
- Sterile, intended for single-use

Part No.	Kit Content	Sterility	Qty
1-519-2143	Probe cover and 4 tapes	Non-sterile	12
1-519-2140	Probe cover, bite block, sterile gel, and 4 tapes	Sterile	12
1-519-2145	Probe cover, bite block, sterile gel, and 4 tapes	Sterile	12



Protek™ Probe Storage Cover

Product Details

- Protects entire probe including handle after cleaning
- Designated space to apply disinfection label
- Clear polyethylene (not made with natural rubber latex)
- Cover size 4" x 16" (10cmx41cm)
- Non-sterile, intended for single-use

Part No.	Description	Qty
1-519-1160	Dispenser box containing 100 probe covers	4



Protek™ Elastic Band

Product Details

- Use to secure ultrasound probe covers
- Not made with natural rubber latex
- Non-sterile

Part No.	Description	Sterility	Qty
1-519-0126	1.25" x 1/8" Blue Bands	Non-sterile	2150



Protek™ ClingTouch™ Monitor Covers
Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Unrolls for easy application
- Top and bottom adhesive to secure cover
- Clear polyethylene and vinyl (not made with natural rubber latex)
- Sterile, intended for single-use

Part No.	Size	Qty
1-519-0180	17" x 17" (43cm x 43cm)	12
1-519-0181	17" x 12" (43cm x 30cm)	12



Protek™ Ultrasound System Covers
Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedureal cleanup
- Clear polyethylene (not made with natural rubber latex)
- Sterile, intended for single-use

Part No.	Size	Kit Content	Qty
1-519-1008	14.5" x 30" (37cm x 76cm)	CT control rail cover and 6 tapes	50
1-519-1010	6" x 10" (15cm x 25cm)	CT control cradle cover and 3 bands	50
1-519-1011	30" x 30" (76cm x 76cm)	System cover and 2 tapes	24
1-519-1212	12" x 24" (31cm x 61cm)	System cover and 4 tapes	24
1-519-1213	24" x 24" (61cm x 61cm)	System cover and 4 tapes	24
1-519-1624	24" W x 28" H (61cm x 71cm)	System cover w/10" deep pocket and 2 tapes	24
1-519-1628	48" W x 28" H (122cm x 71cm)	System cover w/10" deep pocket and 2 tapes	24
1-519-1217	28" x 17.5" (71cm x 45cm)	Portable system cover w/pocket & adhesive seal	24



Protek™ Tablet Cover
Product Details

- Helps minimize cross-contamination
- Protects equipment
- Self-adhesive strip ensures device is fully enclosed
- Clear polyethylene (not made with natural rubber latex)
- Sterile, intended for single-use

Part No.	Size	Qty
1-519-1028	8" x 10" (20cm x 25cm)	24



Protek™ Smartphone Cover
Product Details

- Helps minimize cross-contamination
- Protects equipment
- Self-adhesive strip ensures device is fully enclosed
- Clear polyethylene (not made with natural rubber latex)
- Sterile, intended for single-use

Part No.	Size	Qty
1-519-1039	4" x 6" (10cm x 15cm)	24



Sanitary MRI Headset Covers

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Elastic opening ensures easy application and secure placement
- Comfortable, soft material (not made with natural rubber latex)
- Non-sterile, intended for single-use

Part No.	Description	Size	Color	Qty
1-719-0601	Sanitary MRI Headset Covers	Large	White	1000



Protek™ PullUp™ Biopsy Instrument Cover

Product Details

- Helps minimize cross-contamination
- Protects equipment
- Reduces post-procedural cleanup
- Patented needle isolator allows needle to pierce through cover while maintaining a sterile barrier
- Removable patented PullUp™ with self-adhesive ensures device is fully enclosed
- Telescopic-fold ensures easy gel and cover application
- Clear polyethylene (not made with natural rubber latex)
- Sterile, intended for single-use

Part No.	Size	Qty
1-519-1003	3.75" x 17" (10cm x 43cm)	24




Protek™ Disinfection Cup

Product Details

- A convenient reservoir for disinfecting ultrasound probes, reusable needle guides, and brackets
- Extended length for maximum probe coverage
- Easy wipe down cleaning
- Slide proof table grips
- Hooks for cord and connection management
- Wall mounting hardware

Part No.	Description	Qty
1-577-6050	Disinfection Cup Kit	1
1-577-6051	Replacement Cups	3



Protek™ Needle Guide Cleaning Brushes

Product Details

- Use with disinfection cup to clean reusable needle guides and brackets
- 7" long (16 Ga)
- Non-sterile

Part No.	Description	Qty/Box
1-235-5105	Cleaning Brushes	5

Notes:

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Notes:

[illegible]

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2023 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1108483 rev 2 20-SEPT-2023 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

USA 616-698-7100
Canada 888-364-7004
www.aspensurgical.com

