


Instructions for Use: REF# HC100 (Adult) and
CHC100 (Child)

1. Remove from packaging and place the positioner the O.R. Table at the cephalad end of the O.R. table.
2. The HC100 and CHC100 is indicated for supine positioning of the head during surgical procedures.
3. Discard after use.



Not made with natural rubber latex



Manufactured For:


PRIME
MEDICAL
1410 E Weisgarber
Knoxville, TN 37909
865-776-8994



QNET BV
Kantstraat 19
NL-5076NP Haaren
The Netherlands