

 symmetry surgical® | **Bovie**®



OBGYN PRACTICE PRODUCTS








THE PROCEDURE SAFETY AND RELIABILITY YOUR PRACTICE DEMANDS

Integrate the Bovie® Specialist | PRO energy-based system with smoke evacuation to perform procedures like LLETZ/LEEP, wart removal and other gynecologic procedures. Treat patients in office and add new revenue streams to your practice.

	 <p>Olsen® Speculum, Graves Vaginal Speculum SERS2 Bovie® Smart Shark Remote Activator Switch SERF Reducer Fitter</p>
<p>Specialist PRO High Energy-based Generator</p> <p>Provide operating room-level technology to your patients in your office with this easy-to-use energy-based unit. Includes:</p> <ul style="list-style-type: none"> • Five modalities with up to 120 watts of power for a variety of tissue effects • Monopolar and bipolar capabilities • Compatibility with most industry-standard accessories • Opportunity for additional revenue streams • Four-year warranty 	<p>Smoke Shark® II Smoke Evacuator</p> <p>Protect yourself, your staff and patients from airborne pathogens and contaminants. Effectively removes hazardous smoke and particles from the operative site, reducing the risk of surgical site infection.</p> <ul style="list-style-type: none"> • Includes 35-hour, long-life filter • Improves visibility of operative site • Immediately evacuates contaminated smoke particles • Operates in conjunction with optional Bovie® Smart Shark Remote Activator Switch for maximum filter life and surgical efficiencies • Compatible with most energy-based generators

		
<p>Disposable LLETZ/LEEP Loop</p> <p>Achieve precise control during biopsies and other excisional procedures</p> <ul style="list-style-type: none"> • 100% safety-tested insulation • Molded, ribbed safety grip • Twenty one sizes available 	<p>Disposable Ball Electrodes</p> <p>Designed to control bleeding, fulgurate and treat diseased tissue safely and efficiently.</p> <ul style="list-style-type: none"> • Includes safety grip insulators • Is available in 5mm and 3mm balls • Is manufactured to the highest standards 	<p>Energy-based Mobile Stand</p> <p>At your fingertips: Keep surgical equipment organized and within reach at all times</p> <ul style="list-style-type: none"> • Mobile stand for A1250S; 18.9" wide x 19.6" deep x 32" high • Locking wheel casters • Adjustable shelving heights

MEDICAL LIGHTS

Bovie® exam lights provide brilliant, color-accurate illumination and clarity. The MI 300 offers two brightness settings and can be easily moved on its mobile base. The MI 750, also available in ceiling-mount and wall-mount options, provides physicians with an economical and energy-efficient LED light for exam and procedure rooms. All Bovie lighting (excluding normal replacement parts) is supported by a three-year limited warranty.



Bovie® MI 150 and MI 300

Bovie® MI 150 and MI 300 exam lights provide brilliant, color accurate illumination and clarity. The MI 300 includes two brightness settings and a mobile base.




Bovie® MI 750

The MI 750 provides physicians with bright, white light for clear and accurate visibility during examinations and procedures. Designed with unmatched lifetime spring-arm technology, 50,000 hour LED life, and 75,000 Lux, this high-quality lighting system easily aligns with any practice budget.

THE NEXT GENERATION OF CLARITY AND COMFORT IN COLOSCOPY EXAMS

The Colpo-Master™ line offers superior, 3D, stereoscopic views with excellent clarity. Designed for maximum efficiency from exam to procedure, this suite of colposcopes features optimized ergonomics and a broad range of magnification.

	
<p>Colpo-Master™ I</p> <p>Easy access to surgical site through swing-arm suspension design.</p>	<p>Colpo-Master™ II</p> <p>Compact center-post mount-style for exams, Colpo-Master™ I and II have 10x more magnification than competitive colposcopes, allowing for enhanced visuals.</p>

Colpo-Master™ I Zoom Suspension Arm Colposcopes & Colpo-Master™ II Center-Post Mount Colposcopes

- 3D stereoscopic viewing
- Zoom magnification range 4-27x with 20x eyepiece
- S-video and digital camera options available for connecting EMR platforms or monitor displays (trinocular head required)
- Dual light-intensity controls (green and white) for increased visibility of tissue variances
- Field-of-view ranges from 58.5mm at 4x to 8.5mm at 27x
- Binocular or camera-ready trinocular heads with 298mm working distance
- Independently adjust focus with dual diopter to eliminate eye strain

Compare What's Best for Your Practice

	LED Illumination	Camera-ready Option	4-27 x Zoom	Swing Arm	7-year Limited Warranty
Colpo-Master™ I	•	•	•	•	•
Colpo-Master™ II	•	•	•		•



SYMMETRY SURGICAL OBGYN PROCEDURAL INSTRUMENTATION

04-3318	
06-2113	
10-1012	

Item Code	Description	Item Code	Description
01-2703	Symmetry® Knife Handle, #3, Standard, 5 in, 127 mm	06-2112	Symmetry® Forceps, Adson Tissue, 1x2 Teeth, 4-3/4 in
01-2804	Symmetry® Knife Handle, #4, Standard, 5-1/4 in, 133 mm	06-2113*	Symmetry® Forceps, Adson Tissue, Serrated, 1x2 Teeth, 4-3/4 in
01-2814	Symmetry® Knife Handle, #4L, Long, 8-3/8 in, 212 mm	06-2512	Symmetry® Forceps, Adson-Brown Tissue, 4-3/4 in (121 mm), 7 x 7 Teeth
02-2011	Symmetry® Scissors, Iris Standard Pattern, Straight, 4-1/2 in	06-2590	Symmetry® Forceps, Callison Tying Tissue, 1.5 mm Teeth, Tying Platform, 1x2 Teeth, 4-7/8 in
02-2111	Symmetry® Scissors, Iris Standard Pattern, Curved, 4-1/2 in	06-8798	Forceps, Heaney Tissue, 2x3 Teeth, 5-1/2 in
02-5002	Symmetry® Scissors, Tenotomy, Straight, Sharp, 4-1/4 in	09-8309	Symmetry® Suction Tube, Frazier, Medium Angle, 9 Fr, 4-1/3 in, 110 mm
02-5003	Symmetry® Scissors, Tenotomy, Curved, Sharp, 4-1/8 in	10-1012*	Symmetry® Needle Holder, Webster, Straight, Blunt, Smooth Jaws, 5-1/4 in, 133 mm
03-0415	Symmetry® Scissors, Operating, Straight, Sharp/Sharp, 6 in	10-1212	Symmetry® Needle Holder, Derf, 5 in, 127 mm
03-0613	Symmetry® Scissors, Operating, Straight, Sharp/Blunt, 5-1/8 in	10-1213	Symmetry® Needle Holder, Halsey, Serrated, 5-3/8 in, 127 mm
03-5017	Symmetry® Scissors, Mayo, Straight, Beveled Blades, 6-3/4 in (173 mm)	10-1813	Symmetry® Needle Holder, Mayo-Hegar, 5-1/4 in, 133 mm
03-5117	Symmetry® Scissors, Mayo, Curved, Beveled Blades, 6-3/4 in	10-1815	Symmetry® Needle Holder, Mayo-Hegar, 6-2/3 in, 160 mm
04-3214	Symmetry® Scissors, Metzenbaum Dissecting, Straight, Standard Pattern, 5-3/4 in	10-3121	Symmetry® Needle Holder, Heaney, 8-1/4 in, 209 mm
04-3314	Symmetry® Scissors, Metzenbaum, Curved, Standard Pattern, 5-3/4 in	14-0609	Symmetry® Scissors, Spencer Stitch, 3-1/2 in
04-3318*	Symmetry® Scissors, Metzenbaum, Curved, Standard Pattern, 7 in (178 mm)	15-3209	Symmetry® Forceps, Hartman Baby Mosquito, Straight, 3-7/8 in
06-0014	Symmetry® Forceps, Dressing, Serrated, 5-1/2 in	15-3309	Symmetry® Forceps, Hartman Baby Mosquito, Curved, 4 in
06-0020	Symmetry® Forceps, Dressing, Serrated, 8 in	15-4412	Symmetry® Forceps, Halsted Mosquito, Straight, 5 in
06-0514	Symmetry® Forceps, Tissue, 1x2 Teeth, 5-3/4 in (140 mm)	15-4512	Symmetry® Forceps, Halsted Mosquito, Curved, 5 in (122 mm)
06-0520	Symmetry® Forceps, Tissue, 1x2 Teeth, 8 in	15-7016	Symmetry® Forceps, Kelly Artery, Straight, 6-1/4 in
06-2012	Symmetry® Forceps, Adson Dressing, Serrated Jaws, 4-3/4 in	15-7114	Symmetry® Forceps, Kelly Artery, Curved, 5-1/2 in (140 mm)

*Image shown above

SYMMETRY SURGICAL OBGYN PROCEDURAL INSTRUMENTATION

52-0610



52-4001




Item Code	Description	Item Code	Description
15-7116	Symmetry® Forceps, Kelly Artery, Curved, 6-1/4 in	32-430	Symmetry® Needle Holder, Crile-Wood, 6 in, 152 mm
15-7316	Symmetry® Forceps, Crile Artery, Curved, 6-1/4 in	32-500	Symmetry® Forceps, Adson Dressing, Tungsten Carbide, 4-3/4 in
16-1016	Symmetry® Forceps, Rochester-Pean Hemostatic, Straight, 6-5/8 in	50-0215	Symmetry® Forceps, Allis Tissue, 4x5 Teeth, 6 in (152 mm)
16-1116	Symmetry® Forceps, Rochester-Pean Hemostatic, Curved, 6-3/5 in	50-0219	Symmetry® Forceps, Allis Tissue, 5x6 Teeth, 7-1/2 in
17-5508	Symmetry® Clamp, Backhaus Towel, 3-1/2 in (89 mm)	50-5540	Symmetry® Retractor, Jackson Vaginal, Small
17-5513	Symmetry® Clamp, Backhaus Towel, 5-1/4 in (133 mm)	50-5541	Symmetry® Retractor, Jackson Vaginal, Medium
19-5024	Symmetry® Forceps, Foerster, Straight, Serrated, 9-3/4 in	50-5542	Symmetry® Retractor, Jackson Vaginal, Large
19-5579	Classic® Clamp, Gemini, 5-1/2 in	52-0610*	Symmetry® Speculum, Graves Vaginal, 1-1/4 in x 4 in Blades, Medium
20-0508	Symmetry® Retractor, Senn, Double-Ended, Sharp, 6-1/4 in	52-2793	Speculum, Auvarad Vaginal, 1-3/4 in x 4 in Blades, 2.5 lbs
20-5001	Symmetry® Hook, Joseph Skin, Single-Prong, 6.0 mm, 6-1/4 in, 158 mm	52-4001*	Symmetry® Sound, Sims Uterine, Graduated In Centimeters, 13-1/4 in, 335 mm
20-5005	Symmetry® Hook, Joseph Skin, Double-Prong, 5 mm Spread, 4.5 mm, 6-1/4 in, 158 mm	52-4426	Symmetry® Forceps, Braun Tenaculum, Straight, 9-1/2 in
22-1418	Symmetry® Retractor, Parker, Large, 7 in, (Set of 2)	54-2555	Symmetry® Scissors, Shortbent Stitch
23-5534	Symmetry® Curette, Sims Uterine, Size 1, Sharp, Malleable Shaft, 11 in	54-4030	Symmetry® Scissors, Operating, Straight, 5-1/2 in
23-7503	Symmetry® Retractor, Deaver, 1 in Wide, Standard Handle, 9-1/4 in	60-1813	Symmetry® Scissors, Stitch, Curved, Delicate, Sharp, Large Ribbon Ring Handles w/Flat Shanks, 4 in
24-2211	Symmetry® Retractor, Weitlaner, Blunt, 2x3 Prong, 4-1/2 in	60-1856	Symmetry® Scissors, Strabismus, Straight, Ribbon Handles, 4-1/4 in
30-5561	Symmetry® Forceps, Bozeman Dressing, Curved, 260 mm	60-1924	Symmetry® Needle Holder, Castroviejo, Straight, w/Lock, 12 mm, 5-1/2 in
30-6556	Symmetry® Forceps, Bishop Harmon Iris, Delicate	60-1925	Symmetry® Needle Holder, Castroviejo, Curved, w/Lock, 12 mm, 5-1/2 in
32-110	Symmetry® Needle Holder, Baumgartner, Tungsten Carbide, 5 in	67-5507	Scissors, Ragnell, Curved, Beveled Blades, Dissecting Tips, 5 in
32-405	Symmetry® Needle Holder, Webster, Tungsten Carbide, Large Ring Handle, 4-3/4 in	86-5013	Symmetry® Scissors, Stevens Tenotomy, Curved, Serrated, Sc, 4 in

*Image shown above



SYMMETRY SURGICAL MINIMALLY INVASIVE SURGERY PROCEDURAL DEVICES



Item Code	Description
LSH-01	Liquid Scope HelpOR™, Liquid Scope WarmOR™, two trocar swabs (5-8mm and 10-12mm), one microfiber cloth
SC-01	Scope CleanOR™, Two trocar swabs (5-8mm and 10-12mm), one microfiber cloth
PS-500-M*	Liquid Scope WarmOR™, Short version (two buttons)
PS-500	Liquid Scope WarmOR™, Long version (three buttons)
PS003	AntifogOR™, 6ml anti-fog solution with radiopaque sponge
SHK-01	Scope HelpOR™, Lap Scope WarmOR™ warming sheath, two trocar swabs (5-8mm and 10-12mm) one microfiber cloth
SW-PA-6V	Lap Scope WarmOR™, Medical grade, reusable 6V AC transformer
RSG-14F	Suture PassOR PRO™, 150mm Suture PassOR, one 10-12mm guide and two trocar swabs (5-8mm and 10-12mm)
RSG-14/18F*	Suture PassOR PRO™, 150mm Suture PassOR, two guides (10-12mm and 15mm) and two trocar swabs (5-8mm and 10-12mm)
RSG-18F-XL	Suture PassOR PRO™, 180mm Suture PassOR, two bariatric (XL) length guides (10-12mm and 15mm) and two trocar swabs (5-8mm and 10-12mm)
PC-120*	Port ClosOR PRO™, 150mm Port ClosOR (Gauge:14GA) and two trocar swabs (5-8mm and 10-12mm)
DUMI-350A*	DUMI ManipulatOR™, 350mm adjustable diagnostic uterine manipulator injector

*Image shown above



3034 Owen Drive | Antioch, TN 37013 | 800.251.3000 (US) | +41 (0) 52 557 59 00
customerservice@symmetrysurgical.com | [symmetrysurgical.com](https://www.symmetrysurgical.com)
©2021 Symmetry Surgical Inc. All rights reserved.
Bovie Brochure OBGYN 06.16.2022 v1