



14410 FluidGard® 160 Anti-Fog Procedure Mask with Extended Shield

Intended Use

Intended to be used in environments where high levels of fluid exposure is a risk.



Complies with ASTM F2100

Fluid Resistance (mmHg):	160
BFE (%) @ 3.0 µm:	≥ 98%
PFE (%) @ 0.1 µm:	≥ 98%
ΔP (mm H ₂ O/cm ²):	< 5
Flame Spread:	Class 1

Key Features

- Complies with ASTM F2100 Level 3 requirements
- Coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Anti-fog foam inhibits fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Made in North America
- Not made with natural rubber latex

FluidGard® 160 Anti-Fog Procedure Mask with Extended Shield

Intended to be used in environments where high levels of fluid exposure is a risk.



Technical Specifications

Part No.	14410
Description	FluidGard® 160 Anti-Fog Procedure Mask with Extended Shield
Color	Blue
Fluid Resistance	160 mm Hg
Filtration - BFE	≥99.8
Filtration - PFE	≥99.0
Delta P	≤3.2
Flammability	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Qty	25/Box, 4 Boxes, 100/Case

Precept® and FluidGard® are a registered trademarks of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and model.

©2021 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1051507 rev 1 18-MAY-2021 ENG – US

Aspen Surgical Products, Inc.
phone +1 616.698.7100 toll-free +1 888.364.7004
aspensurgical.com

