

Liquid Scope HelpOR™

The smart way to achieve and maintain optimal vision




Simple, instant and effective.

The Symmetry Surgical® **Liquid Scope HelpOR™** is a cost-effective and comprehensive scope management system that enhances a surgeon's vision during laparoscopic and robotic surgeries. These scope management systems include a lint-free microfiber cloth to clean the lens, two trocar cleaning swabs to remove debris from the cannula, and the Liquid Scope WarmOR™ in either the two-button size or the more extended, three-button size to address any scope fogging or cleaning issues within seconds.

With one click, the Liquid Scope WarmOR provides heat and reaches peak temperatures of 122°F | 50°C in just five minutes.

Scopes are defogged and cleansed quickly and efficiently without costly capital equipment or electrical connections.



Trocar Cleaning Swabs



Liquid Scope WarmOR™




symmetry surgical® | **Minimally Invasive Surgery**



A SIMPLE AND INSTANT HEAT SOURCE


- Activates quickly and intuitively
- Reaches peak temperature (122°F | 50°C) in 5 minutes
- Retains optimal temperature for up to one hour
- Improves OR safety by eliminating fire, burn, and hot water hazards
- Accommodates all sizes and lengths of endoscopes
- Does not require batteries or electrical connections
- Is an affordable, disposable, and sterile device
- Is latex, DEHP, and alcohol-free.



Fogged view when not using Liquid Scope HelpOR.

CLEAN SCOPES AND TROCARS

- Includes a non-abrasive and soft microfiber cleaning cloth (15x15cm) that is gentle on surfaces without streaking, marring, or scattering
- Includes soft and atraumatic surgical foam tips that remove the blood and other fluids from the trocar and cannula without damaging the trocar valves
- Contains trocar cleaning swabs with radiopaque shafts
- Cleans tight spaces due to rigid yet flexible cleaning swabs



Clear view when not using Liquid Scope HelpOR.

REF / Order #	Description	UOM
LSH-01	Liquid Scope HelpOR™ – Liquid Scope WarmOR™, two trocar swabs (5-8mm and 10-12mm), one microfiber cloth	Box of 10, sterile

Contact a Surgical Specialist to learn more about the Minimally Invasive Surgery portfolio.


symmetry surgical® | Minimally Invasive Surgery
customerservice@symmetrysurgical.com • 800.251.3000 • www.symmetrysurgical.com

08.19.2022 v1