




Scan the code above  
to learn more about  
our Rectal portfolio.

# Rectal Portfolio



symmetry surgical™

# Rectal



Forceps	REC-2
Hooks	REC-2
Ligators	REC-3
Probes	REC-3
Retractors	REC-4
Speculums	REC-7

Questions? Contact us at [customerservice@symmetrysurgical.com](mailto:customerservice@symmetrysurgical.com)  
or learn more at [symmetrysurgical.com](http://symmetrysurgical.com)

 **symmetry surgical™**

# The Surgical Instrument Desk Reference

---

## How to Use Our Catalog

---

The Symmetry Surgical catalog is designed with you in mind, to be your trusted resource for surgical instrumentation. We took great care in creating an easy-to-use, robust catalog and we think you'll be pleased with the design as you become more familiar with our broad portfolio.

## Layout

---

The catalog is organized by surgical specialty *and* product function.

And we know that it is difficult to pick up an instrument and know what specialty it's listed in. To make your job easier, we've sorted the instruments in our portfolio by surgical function, too, including Energy, Laparoscopic, Retractors and Sterilization Solutions.

## U.S. Customer Service

---

Your Symmetry Surgical Customer Service team wants to make ordering surgical instruments the easiest part of your day. We accept orders via phone, fax, email, our website, Global Healthcare Exchange (GHX) and through the mail.

### PHONE

1.800.251.3000

(Available Monday through Friday, 8:30 AM–6:30 PM EST)

### FAX

1.800.342.3272

### EMAIL

[customerservice@symmetrysurgical.com](mailto:customerservice@symmetrysurgical.com)

### ONLINE

[symmetrysurgical.com](http://symmetrysurgical.com)

### MAILING ADDRESS

Symmetry Surgical  
Attn: Customer Service  
3034 Owen Drive  
Antioch, TN 37013

# The Surgical Instrument Desk Reference

## International Customer Service

### PHONE

Symmetry Surgical Switzerland GmbH  
+41 (0) 52 557 59 00  
(Available Monday through Friday, 7:00-16:00 GMT)

Symmetry Surgical, Inc., USA  
1.615.964.5672  
(Available Monday through Friday, 8:30 AM–6:30 PM EST)

### FAX

Symmetry Surgical Switzerland GmbH  
+41 (0) 52 557 59 99

### EMAIL

customerservice@symmetrysurgical.com  
internationalorders@symmetrysurgical.com

### ONLINE

symmetrysurgical.com

### MAILING ADDRESS

Symmetry Surgical Switzerland GmbH  
Bühlstrasse 27  
CH-8200 Schaffhausen  
Switzerland

Symmetry Surgical, Inc.  
3034 Owen Drive  
Antioch, TN 37013

### INTERNATIONAL PAYMENT TERMS/RETURN GOODS POLICY/INSTRUMENT REPAIR & RESTORE PROGRAM/LITERATURE REQUESTS

Please contact your local distributor. Locate a local distributor contact, International Customer Service at +41 (0) 52 557 59 00 or 1.615.964.5672 (US).

Cardiovascular/Thoracic

Energy

Ear, Nose & Throat

General

Laparoscopic

Neurology

Obesity

OB/GYN

Ophthalmology

Orthopedics & Arthroscopy

Pediatrics

Rectal

Retractors

Spine

Urology

Sterilization Solutions

Alphabetical & Numeric Indexes

RECTAL

## Forceps (Biopsy)

### Jackson Rectal Biopsy Forceps

CODE	DESCRIPTION	LENGTH
50-8135	1/8 x 1/4 in (4.0 x 6.0 mm) Bite	14 in 355 mm


### Yeoman Rectal Biopsy Forceps

CODE	DESCRIPTION	LENGTH
50-8093	Square Basket	17 in 430 mm

## Hooks (Crypt)

### Stewart Crypt Hook

CODE	DESCRIPTION	LENGTH
50-9438	Short Style	8 1/2 in 215 mm
50-9439	Long Hook	9 in 230 mm



50-9438

## Ligators (Hemorrhoidal)



### McGivney Hemorrhoidal Ligator Set

CODE	DESCRIPTION	LENGTH
DX-8	Set Includes the Following:	
DX-8-8	Ligator Only	8 in 200 mm
DX-18	Grasping Forceps for DX-18	
DX-2	Ligator Cone (Drum)	

## Probes (Fistula)

### Barr Fistula Probe

CODE	DESCRIPTION	LENGTH
50-9409	Double-Ended	6 in 150 mm



### Buie Fistula Probe

CODE	DESCRIPTION	LENGTH
50-9420	Sterling Silver	8 in 200 mm



### Pratt Fistula Probe

CODE	DESCRIPTION	LENGTH
50-9441	Narrow	8 ¼ in 210 mm
50-9442	Short	8 ¼ in 210 mm
50-9443	Long	8 ¼ in 210 mm



## Probes (Rectal)

### Larry Rectal Probe

CODE	DESCRIPTION	LENGTH
50-9450	Sterling Silver, Malleable	7 in 175 mm




### Pratt Rectal Probe

CODE	DESCRIPTION	LENGTH
50-9410	Sterling Silver, Malleable	8 in 200 mm



## Retractors (Bookwalter®)



### Bookwalter® Rectal Accessory Kit

Kit Includes: 1 Ring and 2 Blades

CODE	DESCRIPTION	DIAMETER	BLADE	DIMENSIONS
50-4546	Kit Includes the Following:			
50-4584	Bookwalter® Round Rings, Small	8 ½ in 216mm		
50-4585	Bookwalter® Rectal Blade, Small		2 7/8 in 73 mm	3 x 7 in 76 x 178 mm
50-4586	Bookwalter® Rectal Blade, Large		3 7/8 in 98 mm	4 x 6 in 101 x 152 mm

Note: Use of this product requires a Bookwalter® Table Post, Bookwalter® Coupler, Bookwalter® Horizontal Bar, and Bookwalter® Ratchet (not included in this accessory set.)

50-4546

## Retractors (Prostatic)



51-9503



51-9504

### Young Prostatic Retractor

CODE	DESCRIPTION	DIMENSIONS		LENGTH	
51-9501	Lateral Notched Blade	$\frac{7}{16} \times 2$ in	11 x 50 mm	8 $\frac{1}{4}$ in	210 mm
51-9503	Notched Blade	1 x 1 in	25 x 25 mm	8 $\frac{1}{4}$ in	210 mm
51-9504	Bifurcated Blade	1 $\frac{3}{4}$ x 1 $\frac{1}{2}$ in	44 x 38 mm	8 $\frac{1}{4}$ in	210 mm



### Barr Anal Retractor

CODE	DESCRIPTION	DIMENSIONS		LENGTH	
50-7325	Self-Retaining with Blades	2 $\frac{3}{4}$ x $\frac{7}{8}$ in	70 x 22 mm	8 in	200 mm



## Retractors (Rectal)



50-9604

### Hill-Ferguson Rectal Retractor


CODE	DIMENSIONS	LENGTH
50-9604	$\frac{7}{8} \times \frac{1}{4}$ in 22 x 6.0 mm	8 in 200 mm
50-9605	$\frac{1}{8} \times 3$ in 3.0 x 76 mm	8½ in 215 mm
50-9606	3 x $1\frac{3}{8}$ in 76 x 35 mm	8½ in 215 mm
50-9607	1½ in 38 mm	8½ in 215 mm



50-9621

### Sawyer Retractor

CODE	DESCRIPTION	BLADE DIMENSIONS	LENGTH
50-9621	Retractor	$1\frac{1}{8} \times 3$ in 28 x 76 mm	11 in 280 mm
50-9622	Retractor	$1\frac{1}{2} \times 3\frac{1}{2}$ in 38 x 89 mm	11 in 280 mm



### Smith-Buie Rectal Retractor

CODE	DESCRIPTION	LENGTH
50-7420	Self-Retaining, with Swivel Blades	4½ in 115 mm

## Speculums (Operating)



### Fansler Operating Speculum

CODE	DESCRIPTION	DIAMETER		LENGTH	
50-9610	Fenestrated Handle, Removable Obturator	1 $\frac{3}{8}$ in	35 mm	8 in	200 mm

## Speculums (Rectal)



### Brinkerhoff Rectal Speculum

CODE	DESCRIPTION	DIMENSIONS		LENGTH	
50-7615	Medium	$\frac{5}{8}$ x 1 in	16 x 25 mm	4 in	100 mm



50-9616

### Chelsea-Eaton Anal Speculum

CODE	DESCRIPTION	DIAMETER	
50-9615	Fenestrated Handle, Removable Obturator	$\frac{7}{8}$ in	22 mm
50-9616	Fenestrated Handle, Removable Obturator	1 $\frac{1}{4}$ in	32 mm

## Speculums (Rectal)



50-9601

### Hinkle-James Rectal Speculum

CODE	DESCRIPTION	DIAMETER	
50-9601	Small	2 $\frac{1}{8}$ in	54 mm
50-9602	Medium	2 $\frac{1}{4}$ in	57 mm



50-7716

### Hirschmann Rectal Speculum

CODE	DESCRIPTION	DIAMETER		LENGTH
50-7716	Small	$\frac{9}{16}$ in	14 mm	2 $\frac{3}{8}$ in 60 mm
50-7719	Medium	$\frac{11}{16}$ in	17 mm	2 $\frac{3}{4}$ in 70 mm



### Pennington Rectal Speculum

CODE	DESCRIPTION	DIAMETER		LENGTH
50-9630	With Obturator, Long Handle, Chrome	1 in	25 mm	5 $\frac{1}{4}$ in 135 mm

## Speculums (Rectal)



Pratt Rectal Speculum

CODE	DESCRIPTION	BLADE (W X D)	LENGTH
50-7320-H	Stainless Steel, with Set Screw	1 x 3 <sup>3</sup> / <sub>8</sub> in 25 x 85 mm	8 <sup>1</sup> / <sub>2</sub> in 215 mm



Sims Rectal Speculum

CODE	DESCRIPTION	BLADE (W X D)	LENGTH
50-7300	With Solid Blade	1/2 x 3 in 13 x 76 mm	6 1/4 in 159 mm
50-7305	With Fenestrated Blade	1/2 x 3 in 13 x 76 mm	6 1/4 in 159 mm