


A close-up photograph of a newborn baby's face, looking slightly to the left. The baby has dark hair and eyes. Two hands are gently holding the baby's head from behind. The baby is wrapped in a white cloth.

Stork® Labor & Delivery  
Product Catalog

**Stork® Labor and Delivery Products**

The birth of a child is one of the riskiest and most costly procedures in a hospital. Sacrificing high quality for cost savings is not an option when liabilities are so high. Our traditional obstetric products have been trusted for well over 100 million births.

All of our brands are produced to high manufacturing standards and exacting quality specifications. Our Labor and Delivery products must pass through our innovative, rigorous laser vision inspection systems to ensure they meet our high quality standards. We have been manufacturing these products and perfecting these processes in the same U.S. facility for decades since their inception. Our brands are built on this experience and our commitment to the highest quality, proven every day, every birth.



**AmniHook® Amniotic Membrane Perforator**

**Product Details**

- 100% laser-vision inspected
- Rounded tip and protected point designed to rupture the amniotic membrane without injury to mother or fetus
- 10" long, with non-glare molded plastic and sides ensure a non-slip grip and precise control
- Disposable, pre-sterilized, individually packaged for convenience

Part No.	Description	Sterility	Qty
9601	Membrane Perforator	Sterile	100 ea/bx

**Double-Grip® Cord Clamp**

**Product Details**

- Safe, disposable, easy-to-apply
- Dual gripping surface prevents slipping during application
- Maintains tight seal and applies constant pressure on any size cord as it shrinks
- Hinge guard confines cord between serrated jaws to help eliminate movement
- No rough edges to snag clothing or get caught on bed linen

Part No.	Description	Sterility	Qty
9423	Cord Clamp	Sterile	125 ea/bx
9411	Cord Clamp	Non-sterile	250 ea/bx



**Cord-Clamp® Clipper**

**Product Details**

- Concealed cutting surface helps eliminate the risk of sharp or pointed objects near baby
- One time single patient use only to prevent risk of cross-contamination
- Designed to work exclusively with our Double-Grip® Umbilical Cord Clamp, our Cord Clamp® Clipper safely cuts through the specially designed clamp hinge for easy removal from baby

Part No.	Description	Sterility	Qty
9441	Single Clipper	Non-sterile	1 ea/bx
9445	Bulk Clipper	Non-sterile	25 ea/bx



**U-Bag® Urine Specimen Collector**

**Product Details**

- Transparent butterfly patch makes application and removal easy and accurate, helping to reduce valuable nursing time
- Compatible and comfortable for male and female anatomy
- Double-chambered construction with no-flow back valve helps keep urine away from skin to prevent irritation and chance of contaminating the urine specimen
- Adhesive patch can be folded in half to send sealed specimen to the lab
- Premie size ideal for NICU use with low birth weight infants
- Designed with an infant's tender skin in mind


**U-Bag® with Porous Cloth Adhesive Single Specimen Style**

Part No.	Description	Sterility	Qty
7540	Premie (10mL)	Sterile	100 ea/bx
7535	Newborn (100 mL)	Sterile	100 ea/bx
7531	Pediatric (200 mL)	Sterile	100 ea/bx

**24-hour Style**

7537	Newborn (100 mL)	Sterile	50 ea/bx
7533	Pediatric (200 mL)	Sterile	50 ea/bx





**Ready-Rolled® FootPrinter**

**Product Details**

- Disposable plastic plate fits easily in the palm of one hand
- 4½" x 2¾" less than ⅛" thick with stackable edges for easy, compact storage
- Non-toxic, black ink is easy to see and wash off of baby
- Great for our Heirloom Quality Birth Certificates, Newborn Identification Forms, or other hospital made mementos
- Designed to guard against over-inking

Part No.	Description	Qty
5402	Plastic FootPrinter	25 ea/pk - 4 pks/bx



**PlastiBell® Circumcision Device**

**Product Details**

- Disposable one piece device plus ligature, clear, plastic construction helps to provide visual access at all times
- Produces a clean line of excision and helps to promote rapid healing
- No post operative care or special dressing required
- Sterile, individually packaged in six sizes and assortment packs
- Eliminates the need for repetitive and costly sterilization required of stainless steel clamps


Part No.	Description	Qty
9210	Hospital Assortment Pack	125 ea/pk
	1.1 cm (25 ea)	
	1.2 cm (35 ea)	
	1.3 cm (45 ea)	
	1.4 cm (10 ea)	
	1.5 cm (5 ea)	
9230	Physicians Office Assortment Pack	25 ea/pk
	1.1 cm (5 ea)	
	1.2 cm (7 ea)	
	1.3 cm (8 ea)	
	1.4 cm (3 ea)	
	1.5 cm (1 ea)	
	1.7 cm (1 ea)	

**Hospital (Bulk) Single Sizes**

Part No.	Description	Qty
9211	1.1 cm	125/ea/bx
9212	1.2 cm	125 ea/bx
9213	1.3 cm	125 ea/bx
9214	1.4 cm	125 ea/bx
9215	1.5 cm	125 ea/bx
9217	1.7 cm	125 ea/bx

**Physician Office Single Sizes**

Part No.	Description	Qty
9231	1.1 cm	25 ea/bx
9232	1.2 cm	25 ea/bx
9233	1.3 cm	25 ea/bx
9234	1.4 cm	25 ea/bx
9235	1.5 cm	25 ea/bx
9237	1.7 cm	25 ea/bx



**Identification Form**

**Product Details**

- Receives and retains every detail of a baby's tiny skin print
- Space for recording pertinent birth information and verifying signatures
- Enough space for two infant footprints and a mother's fingerprint
- Some forms allow for the attachment of infant's ID band
- Printed on high contrast, archival quality paper, perfect for both receiving and preserving infant foot prints.

Part No.	Description	Size	Sterility	Qty
5821	Newborn ID Form	8" x 5"	Non-sterile	50 sheets/pad, 5 pads/box
5823	Newborn ID Form	8" x 11"	Non-sterile	50 sheets/pad, 5 pads/box
5825	Newborn ID Form	8" x 11"	Non-sterile	50 sheets/pad, 5 pads/box
5826	Newborn ID Form	8" x 11"	Non-sterile	50 sheets/pad, 5 pads/box
8403	Patient Name Card	5¾" x ¾"	Non-sterile	208 sheets/pack
8423	Patient Name Card	5¾" x 1"	Non-sterile	144 sheets/pack



Amnihook®, Cord-Clamp®, Double-Grip®, PlastiBell®, Ready-Rolled®, Stork®, and U-Bag® are registered trademarks of Aspen Surgical Products, Inc.

Please note, not all products are available in all regions.

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2022 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1098330 rev 1 19-OCT-2022 ENG – US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

USA 616-698-7100  
Canada 888-364-7004  
[www.aspensurgical.com](http://www.aspensurgical.com)



 aspen surgical®