



15101 Isolation Mask

Intended Use

Intended to be used in environments where exposure to fluids is not a risk.

Key Features

- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning
- Made in North America
- Not made with natural rubber latex

Isolation Mask

Intended to be used in environments where exposure to fluids is not a risk.



Technical Specifications

Part No.	15101
Description	Isolation Mask
Color	Blue
Fluid Resistance	N/A
Filtration - BFE	≥98.7
Filtration - PFE	≥98.9
Delta P	≤2.9
Flammability	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Qty	50/Box, 10 Boxes, 500/Case

Precept® is a registered trademark of Aspen Surgical Products, Inc.

Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2021 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1051511 rev 1 30-MAR-2021 ENG – US

Aspen Surgical Products, Inc.
phone +1 616.698.7100 toll-free +1 888.364.7004
aspensurgical.com

