



High quality realistic ALS training



Available in Light, Medium,
and Dark Skin Tones

SimMan ALS

SimMan ALS provides a mobile, durable solution that will meet the training needs of pre-hospital and in-hospital emergency care providers - from basic assessment to advanced life-support skills. From pre-hospital, on scene assessment and management to definitive care in a hospital, SimMan ALS can fulfill the unique training requirements of emergency healthcare providers.



Multiple Airway Skills/Features

- Controllable open/closed airway; automatically or manually controlled
- Head tilt/Chin lift
- Jaw thrust w/articulated jaw
- Suctioning (Oral & Nasopharyngeal)
- Bag-mask ventilation
- Orotracheal intubation
- Nasotracheal intubation
- Combitube, LMA, I-gel and other airway placement
- Endotracheal tube intubation
- Retrograde intubation
- Fiberoptic intubation
- Transtracheal jet ventilation
- Needle cricothyrotomy
- Surgical cricothyrotomy
- Airway resistance - 3 settings (Open/ Medium/ Closed)
- Right main stem intubation
- Stomach distention
- Airway complications include tongue fallback and tongue edema

Breathing Features

- Simulated spontaneous breathing
- Bilateral and unilateral chest rise and fall
- Normal and abnormal breath sounds
- Lung auscultation sites: posterior and anterior sites
- Oxygen saturation and waveform

Breathing Complications

- Cyanosis
- Pneumothorax w/bilateral needle decompression
- Unilateral & bilateral chest movement
- Chest tube insertion – bilateral

Certifications

- UL, CE, FCC, CSA, HMR

Cardiac Features

- Extensive ECG library
- Heart sounds
- ECG rhythm monitoring (4 wire)
- 12 lead ECG display
- Defibrillation and cardioversion using training pads – ShockLink™ (SimMan ALS w/ ShockLink)
- Defibrillation and cardioversion using live defibrillation pads or paddles (SimMan ALS LiveShock)
- Pacing

Circulation Features

- BP measured manually by auscultation of Korotkoff sounds
- Bilateral carotid, femoral, radial pulses synchronized with ECG
- Pulse strength variable with BP
- Pulse palpation is detected & logged

Vascular Access

- Pre-ported IV access (right arm)
- Intraosseous access (tibia)

Eyes

- Manually set to: Open, closed and partially open
- Set of interchangeable pupils

QCPR

- QCPR feedback and scoring compliant with 2010-2015 Guidelines
- CPR compressions generate palpable pulses, blood pressure wave form, and ECG artifacts
- Realistic compression depth and resistance
- Automatic detection of depth, release and frequency of compressions as well as hand placement

Other Features

- Fully articulated pelvis
- Foley catheterization (without urine output)
- Patient Voice
 - Pre-recorded and custom sounds
- Instructor can simulate patient's voice wirelessly
- Prepared for diagnostic ultrasound training using the Laerdal-SonoSim Ultrasound solution (separate purchase required, for LLEAP operation only).

System Features

- Wireless PC or SimPad PLUS controls simulator remotely
- Control simulations from anywhere on your network
- Manual Mode
 - Precise control "on the fly"
- Automatic Mode
 - Run pre-packaged scenarios
 - Unique, simple controls increase/decrease difficulty & pace
- Integrated video debriefing
- Data logging
- Instructor comments

Patient Monitor (Optional)

- Simulated patient monitor provides concise clinical feedback for physiological parameters
- Simulated parameters for SimPad PLUS operation: ECG, etCO₂, HR, SpO₂, BP RR, and Temperature
- LLEAP operation provides additional simulated parameters including CO₂, CVP, ICP, anesthetic agent, TOF, cardiac output, and more. In addition, the LLEAP monitoring options can display x-ray, custom images, and custom video

Ordering Information

235- 02350 / 02150 / 02250 SimMan ALS (Light / Medium / Dark)

235- 03350 / 03150 / 03250 SimMan ALS LiveShock (Light / Medium / Dark)

Includes: SimMan ALS Manikin, Lithium-ion battery, specially designed clothing, upper teeth (hard), male/female/blank genitalia, BP cuff, set of consumables, and setup guides. In addition, SimMan ALS LiveShock comes with built-in ShockLink unit and special chest skin for defibrillation with real pads or paddles.

Operation and Patient Monitoring Options

- | | |
|-----------|---|
| 204-301xx | SimPad Plus LLEAP |
| 400-10250 | Laptop Instructor PC-Patient Monitor |
| 400-09250 | Tablet Instructor PC-Patient Monitor |
| 400-09550 | Rugged Instructor PC-Patient Monitor |
| 400-29350 | All in One Panel PC Instructor PC-Patient Monitor |

Licenses

- | | |
|-----------|--|
| 204-50150 | LLEAP for SimPad PLUS Software License |
| 400-01050 | LLEAP PC Software License |

Accessories

- | | |
|------------------|--|
| 212-13850/ M / D | Right humeral IO/IV arm
(Light / Medium / Dark) |
| 400-96050 | USB HD Web Cam |
| 212-29650 | Headset and Mic with USB (LLEAP Operation) |
| 200-30950 | SimPad Headset |
| 185-10050 | ShockLink System |
| 390-01350 | Laerdal-SonoSim Ultrasound Solution 2.0 for SimMan (Probe, Software, Tags) |
| 212-07150 | Rugged Carry Case |
| 212-17975 | Soft Carry Case |



Laerdal
helping save lives

For more information, visit www.laerdal.com/SimManALS