

EQUASHIELD®

Female Connector 180 and Spike Adaptor 180



Female Connector 180 and Spike Adaptor 180

Product Overview

EQUASHIELD®'s Spike Adaptor 180 (SA-180) and Female Connector 180 (FC-180) work hand in hand to efficiently prepare and administer multiple chemotherapy drugs to a patient while reducing the number of components needed.

Product Features

/ SA-180 is designed with a pre-bonded EQUASHIELD® Luer Lock Adaptor to a bag spike, allowing a secure and leak-free connection to an IV bag.

/ FC-180 connects to the SA-180 via a pre-bonded Female Connector and contains an IV tubing port, allowing for seamless drug delivery to the patient.

/ Dry disconnections are provided by an innovative elastomeric membrane designed to prevent contamination from drug residue.

Fluid Pathway

SA-180	FC-180
- Polypropylene	- Polypropylene
- Polyisoprene	- Stainless Steel
	- Silicone
	- Polyisoprene

Technical Specification

- **Priming Volume:** Approximately 0.47 ml for SA-180 and 0.09 ml for FC-180
- **Needle Size:** 17G (1.5 mm) for FC-180
- **Needle Length:** 1.2" (31 mm) non-coring pencil point for FC-180
- **Gravity Flow Rate:** A minimum flow rate of 1000 ml of sodium chloride in 10 minutes, according to ISO 8536-4

Online Resources

[Product Page](#)

[How to videos](#)

[EQ Academy](#)

Use With

EQUASHIELD® Female Connector 180
and Spike Adaptor 180 Adaptor Swivel

Ordering Information

Female Connector 180 and Spike Adaptor 180

FC-180

Luer Lock
Y-Site



50



200

SA-180



50



200



[Find a US Sales Representative](#)

[Find a Local Distributor](#)

[Contact EQUASHIELD](#)