

## Hygiena™ Positive Control Kit

The Positive Control Kit validates the efficacy and stability of Hygiena™ ATP test devices. The kit comes with 25 sealed glass vials containing lyophilized adenosine triphosphate (ATP) and sugars to provide predictable results when sampling with test devices.

### Benefits:

- Validates performance of Hygiena test devices
- No lab equipment needed
- No reconstitutions or complicated sample preparation
- Helps to meet GMP and other quality program requirements
- 24-month shelf life at refrigerated temperature (2 – 8 °C)



### Compatible Products

- UltraSnap™ ATP Surface Test
- SuperSnap™ High Sensitivity ATP Surface Test
- SnapShot Universal ATP Surface Test
- SpotCheck Plus Food Safety Residue Test

Catalog No.	Description	Quantity
CK25	ATP Positive Control Kit	1
US2020	UltraSnap™ Surface ATP Test	100
SUS3000	SuperSnap™ High Sensitivity ATP Surface Test	100
See catalog	SnapShot Universal ATP Test	100
GL100	SpotCheck Plus	100



Find support documents, instructional videos, and more at [www.hygiena.com](http://www.hygiena.com)