




AlerTox[®] Sticks

Rapid tests for allergen detection in food, drinks, and surfaces for internal quality control.



Egg



Crustacean



Hazelnut



β-Lactoglobulin



Fish



Almond



Casein



Mustard Seed



Walnut



Total Milk



Soy



Peanut

AlerTox[®] Sticks advantages:

- Quick and Reliable
- Results in 10 minutes
- No Cross-Reactivity
- High Sensitivity

AlerTox[®] Sticks are a practical tool for the control of allergen content in raw materials, final products, and working surfaces.

Industrial application

- Semi-quantification of allergens in food
- Food quality control
- Trace allergen contamination in food, liquid and surfaces
- Implementation and fulfillment of HACCP, ISO 22000, BRC, and IFS compliance



Product Name	Ordering No.	Quantity
AlerTox[®] Sticks Egg	KIT3025	10 tests
AlerTox[®] Sticks β-Lactoglobulin	KIT3018	10 tests
AlerTox[®] Sticks Casein	KIT3021	10 tests
AlerTox[®] Sticks Total Milk	KIT3015	10 tests
AlerTox[®] Sticks Crustacean	KIT3036	10 tests
AlerTox[®] Sticks Fish	KIT3038	10 tests
AlerTox[®] Sticks Mustard Seed	KIT3040	10 tests
AlerTox[®] Sticks Soy	KIT3028	10 tests
AlerTox[®] Sticks Hazelnut	KIT3035	10 tests
AlerTox[®] Sticks Almond	KIT3033	10 tests
AlerTox[®] Sticks Walnut	KIT3092	10 tests
AlerTox[®] Sticks Peanut	KIT3094	10 tests

NOTE: As with any qualitative or semiquantitative test, results should be periodically validated with a quantitative test performed in an outside laboratory.

www.hygiena.com